

# Tooth Brushes



**Prestige Dental Products, Inc**

**[Prestigedentalproducts.com](http://Prestigedentalproducts.com)**

**Tel: 877-772-3888 Fax: 714-632-8821**

# Tooth Brushes



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Minty Finger Tooth Wipes (PacDent)

Textured teeth wipes that can be used anytime, anywhere to quickly clean the plaque off of teeth and freshen breath without the use of water or rinsing. Sugar-free and individually wrapped.



**100 x individually wrapped finger tooth wipes**

54-45491

## Ultra Water Jet Flosser (Waterpik)



- Unique combination of water pressure and pulsation reduces gingivitis, prevents gum disease
- Contemporary compact design
- Includes six tips: two standard, one Pik-Pocket tip, one tongue cleaner, one orthodontic, and one toothbrush

**Each**

46-91281

## Water Pik Flosser (Waterpik)

- Quick, easy, convenient.
- Gentle, yet effective.
- Tip action of 10,000 strokes per minute effectively removes plaque above and below the gum line.
- Single filament easily and comfortably reaches between the teeth.
- Ideal for use with braces, crowns, bridges and other dental work
- Contemporary ergonomic design.
- One-handed flossing.
- Portability on the road.
- Developed by a dentist.
- Flosser is battery operated w/replacement tip cartridge.



**Flosser, Each**

46-91271

**Replacement Tip Cartridge, 30/Pkg**

46-912712

## ClearGrip Adult Toothbrush



- Two-tone soft end rounded bristles with two rows of blue brush use bristles to indicate replacement
- Oval head and narrow neck
- Clear plastic handles with elastomeric grips in assorted colors

**Compact Head, 72/Box**

54-91601

**Full Head, 72/Box**

54-916012

**Full Head w/Angled Handle, 72/Box**

54-9160123

## RotaPower Toothbrush (Gum)



- Compact brush head
- Two circles of bristles that are of different lengths
- Battery operated

### **Benefits:**

- Makes it easier to clean harder-to-reach areas of the mouth
- The different bristles clean both the tooth and along the gum line
- Easy to use and store, requires no charging

**Power Toothbrush, 12/Pkg**

54-11640

**Power Toothbrush, 1/Pkg**

54-11641

## Bright Tooth Brush (TC Dental)



Beautifully designed with an assorted attractive colors. Soft and comfortable bristles provides extra cleaning. Handles are designed with good and firm grip. Each children's toothbrush imprinted with attractive cartoon. Each is individually wrapped in see through plastics.

**For adult, 72 pcs/box**

54-11116

## Child Toothbrush

3-Row, 27 Tufts



**Child Soft, 72/Box**

54-91641

## Interproximal Brushes (Plasdent)



**2 mm-3 mm Tapered Asst Lite Colors 100/Pkg**

54-11114

**2 mm-4 mm Tapered Asst Neon Colors 100/Pkg**

54-11113

**3 mm-7 mm Tapered Asst Neon Colors 100/Pkg**

54-11112



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Disposable Pasted Toothbrush

- 4 colors, popular handle, U S formulated mint paste on soft bristles
- Often called a 'bathroom brush' before seeing dentist/hygienist
- Individually cello wrapped for hygiene & to avoid flat bristles



100/Box

54-91591

### Disposable Toothbrush



- 4 colors, popular handle, U S formulated soft bristles
- Often called a 'bathroom brush' before seeing dentist/hygienist
- Individually cello wrapped for hygiene & to avoid flat bristles

100/Box

54-91581

### Dual-Action Tongue Cleaner (Gum)



Brush and scraper combo for better cleaning and better control of bad breath.

Fresh R Tongue Cleaner, 6/Box

46-11244

### Standard Denture Carebrush (PlasDent)



Firm, Blue Color, 12/box

41-97991

Firm, Magenta Color, 12/box

41-97992

Firm, Assorted Colors, 12/box

41-97993

## Tooth Brushes

### Infant Toothbrush (PlasDent)



Clear Silicone Material Dishwasher Safe, each

54-11111

### Interproximal Brushes (Prestige)

Individually Wrapped, 100/Box 54-91391

### Toothpaste Squeezers (PlasDent)



Asst 6 Bright Colors, 144 pcs/bag 54-11119

### Nylon Clear Toothbrushes



Similar to Oral B

Soft-Youth, 3 Rows 27 Tufts, 72/Box 54-91611

Adult-Soft, 4 Rows 40 Tufts, 72/Box 54-91612

Adult-Soft, 4 Rows 48 Tufts, 72/Box 54-91613

## Prestige Dental Products



Excellent Qualities  
 at Affordable Prices

877-772-3888

www.prestigedentalproducts.com

# Tooth Brushes



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Premium Adult Toothbrushes



Removes plaque from hard-to-reach places

- Color-coded round front tip
- Rubber grips prevent slipping
- Assorted colors
- Compares to Oral B Advantage

**Clear, 72/Box** 54-91621

**Opaque, 72/Box** 54-916212

## Toothbrush Covers (PlasDent)



**Asst Colors 144 pcs/bag** 54-11118

## Smile Tooth Brush (TC Dental)



smile brush

Beautifully designed with an assorted attractive colors. Soft and comfortable bristles provides extra cleaning.

Handles are designed with good and firm grip. Each children's toothbrush imprinted with attractive cartoon.

Each is individually weal wrapped in see through plastics.

**For adult, 72 pcs/box** 54-11117

## Angle Toothbrush (Butler)



Contemporary style, ribbed comfort, longer handle angled at 10 degrees.

**Compact w/34 Tufts, Soft, 12/Pkg** 54-98211

**Full Size w/42 Tufts, Soft, 12/Pkg** 54-98212

## Soft Tooth Brush (TC Dental)



soft brush

Beautifully designed with an assorted attractive colors. Soft and comfortable bristles provides extra cleaning. Handles are designed with good and firm grip. Each children's toothbrush imprinted with attractive cartoon. Each is individually weal wrapped in see through plastics.

**For children, 72 pcs/box** 54-11115

## Supreme Toothbrushes



This massager-style toothbrush offers:

- Tapered head with power tip bristles
- Narrow neck for better visibility and flexibility
- Elastomeric fingers which massage the gums while brushing
- Ribbed tongue cleaner on back of brush head
- Ergonomic handle for better gripping
- Imprinting available

**Assorted Colors, 72/Box** 54-91631

## Orthodontic Toothbrushes



V-Cut, rounded edges, assorted colors, sturdy nylon bristles.

**Soft, Opaque, 6 Dozen** 54-98171





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Toothpaste (Oraline)

Formulas are equivalent to popular branded flavored tooth-pastes.



#### **Fluoride Bubblegum Flavored**

Toothpaste For Kids of All Ages, 0.85oz. 144/Pkg

46-46251

**Fluoride Mint Flavored Toothpaste** For Adults, 0.85oz. 144/Pkg

46-46252

### Premium Denture Carebrush (PlasDent)



Angle Brush &  
 Ergonomic  
 Handle

Soft Muave/White Color, 12/box 41-98001

Soft Tangerine/White Color, 12/box 41-98002

Soft Blue/White Color, 12/box 41-98003

Soft Assorted Colors, 12/box 41-98004

### Aquafresh Tooth Brushes



Infant, Each 54-98192

Kids, X-Soft, Each 54-98193

## Tooth Brushes

### Delicate Post-Surgical Brush (Butler)



Delicate, with ultra gentle 0.004 bristles Ideal for post-surgical cleaning or in the presence of aphthous ulcers, mouth irritations, extractions, implants and grafts.

Brush w /31 Tufts, 12/Pkg

54-98201



# BEST QUALITY DENTAL SUPPLIES



**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

