

Autoclaves & Sterilizers



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Autoclaves & Sterilizers



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Autoclave, Sterilizer Unit (Tuttnauer 2340M)

Tuttnauer 2340M Manual Autoclave, Sterilizer



Manual autoclave is designed to comply any healthcare facility. The bright new easy to use panel was designed with the operator in mind. Our 2340M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certifications

Tuttnauer 2340M 30-11145

Dry Heat Sterilizer (DA)

- Cost efficient
- Sterilize trays with lid on
- Automatic timer
- Even heat, rapid heat-up
- Simple to operate
- Burnout-proof heating element
- Precision thermostat
- No pressure, no moisture
- Chamber sealing door
- Stainless steel, rust resistant



•1-year warranty

- 2-Tray Sterilizer Dimensions, Outside: 15 1/2" W x 10 1/2" H x 9 1/2" D, Inside: 12 1/2" W x 5 1/2" H x 7 1/2" D, Tray: 10" x 6 1/2"
- 3-Tray Sterilizer Dimensions, Outside: 19" W x 14 1/2" H x 11" D, Inside: 15" W x 7 1/2" H x 9" D, Tray: 11 1/4" x 7 1/2"

2 Trays Unit Model 200, 450 Watts 30-92401

3 Trays Unit Model 300, 650 Watts 30-92402



*We Work Hard
So You Don't Have To.*
877-772-3888

Large Capacity Fully Automatic Autoclave/Sterilizer (Tuttnauer 3870EA)



Our largest automatic table top autoclave. Features a closed door active drying system. Perfect for multi-doctor practices, clinics or surgery centers. 15" x 30" chamber accommodates longer instruments or a large number of small instruments. Capacity for 15 full size IMS cassettes. 2 year parts and labor warranty.

Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- Hepa filter 0.2 µm air filter (replaceable)
- Rapid cycles optimize the availability of ready-to-use sterile instruments
- Easy-to-use, single touch control system increases productivity
- 3 Pre-set sterilization programs
- Program and customize cycle parameters to suit your needs and maximize flexibility
- Integrated printer for cycle documentation (optional)
- Interdependent temperature and pressure monitoring
- Cycle information recovery in the case of power failure or cycle interruption
- Insufficient water detector with automatic shut-off

Tuttnauer 3870EA 30-11256

Large Capacity Manual Autoclave/Sterilizer (Tuttnauer 3870M)



Our largest manual table top autoclave. Perfect for multi-doctor practices, clinics or surgery centers. 15" x 30" chamber accommodates longer instruments or a large number of small instruments. Capacity for 15 full size IMS cassettes. 1 year parts and labor warranty.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certifications
- **1 year parts and labor warranty**

Tuttnauer 3870M 30-11257



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Autoclaves & Sterilizers

Rapid COX Dry Heat Sterilizer - 6 minute cycle

- 6 minute cycle.
- Fast turn around
- eliminate corrosion or rusting of instruments
- low maintenance
- No supplies needed

- *one drawer*
- *one removable instrument basket*
- *Basket Removal Tool*
- *Cooling Rack*
- *Instrument Rack*

Features:

- You can adjust the time to: 6 minute cycle time for unwrapped Instrument, or 8 minutes cycle time for hand pieces, or 12 minutes Cycle time for wrapped instrument.
- Contaminated implants or those needing alterations can be re-sterilized Fast.

Unit, 110V

30-92381

Unit, 220V

30-92382

- Avoids added risk and time of making new appointment and re-sedating patient if alteration is needed

- Quick cycle times which reduce needs for increased instrument inventory
- Unit size: 17" x 17" square x 15" high / 432 mm x 432 mm square x 381 mm high

- Instrument basket size: 7.5" x 8.25" x 1.5" high / 229 mm x 216 mm x 38 mm

- Operating Temperature: 375 degrees F (190 degrees C)

ONE YEAR MANUFACTURER WARRANTY



EZ 10 Automatic Autoclave -Sterilizer (Tuttnauer)

Fully Automatic Autoclaves EZ10

With the simplicity of one touch design, all your sterilization and drying needs are fulfilled. The Tuttnauer fully automatic autoclaves fill, sterilize, exhaust and dry at the touch of a button. FORGET IT ONCE YOU SET IT.

Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- High efficiency air pump circulates hot air after the sterilization cycle for fast drying
- Hepa filter 0.2µm air filter (replaceable)
- Pre-programmed cycles for: unwrapped instruments; wrapped/packs; glassware and additional drying cycle
- All program parameters can be customized and stored
- Automatic shut off at the end of both the sterilization and dry cycles
- Audible and visual cycle interruption alert
- Double safety locking device prevents door from opening while chamber is pressurized
- Low water sensor prevents activation of sterilization cycle when there is insufficient water in the chamber or reservoir.
- Constant monitoring of temperature and pressure for sterilization
- Power outage recovery system
- Door design prevents steam from coming in contact with control panel
- Optional internal printer documents: date, time, temperature and pressure
- 510k clearance for dental hand pieces
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- International certifications

•2 year parts and labor warranty



Fully Automatic Kwiklave™ EZ 10K

The Tuttnauer EZ series is also available in a quick cycle model, the Kwiklave (EZ10k). This sterilizer offers all the proven, dependable features of our fully automatic autoclaves, with the added benefit of reducing the total heat up time. This autoclave is the perfect choice for quick instrument turnaround or a high volume office. Speed and reliability are accomplished without sacrificing load size. 2 year parts and labor warranty. Model EZ10k requires a voltage between 220V and 235V for proper operation.

When applicable a buck/boost transformer (0.5KVA) may be required.

Tuttnauer EZ10

30-11205

Tuttnauer EZ10K

30-11206



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Manual Kwiklave Sterilizer (Tuttnauer 2540MK)



The Tuttnauer 2540MK Series is also available in a quick cycle model, the Kwiklave (MK). This sterilizer offers all the proven, dependable features of our manual autoclave, with the

added benefit of reducing the total heat up time. This autoclave is the perfect choice for quick instrument turn-around or a high volume office. Speed and reliability are accomplished without sacrificing load size. 1 year parts and labor warranty. Model 2540MK requires a voltage between 220V and 235V for proper operation.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certification
- **1 year parts and labor warranty**

Tuttnauer 2540MK

30-11274

Tuttnauer Autoclave/Sterilizer (Tuttnauer 2540M)



The 2540M Tuttnauer manual autoclave is designed to compliment any healthcare facility. The bright new easy to use panel was designed with the operator in mind.

Our 2540M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certifications
- **1 year parts and labor warranty**

Tuttnauer 2540M

30-114111

Autoclaves & Sterilizers

Valueklave 1730 Sterilizer (Tuttnauer)



- The Valueklave 1730 Autoclave will satisfy all of your basic sterilization needs at an economical price
- The compact size (**7" chamber**) makes it a perfect primary unit for small offices or a second autoclave for larger practices
- **The fast cycle time** and heated drying delivers a quick turn-around for wrapped or unwrapped instruments
- The simple controls and easy closing door make this unit a breeze to operate
- **High speed** and reliability are built into every Valueklave 1730

Valueklave 1730

30-114141



BEST QUALITY

DENTAL SUPPLIES

Autoclaves & Sterilizers

Statim 5000 G4 Cassette Autoclave (SciCan)



STATIM 5000 G4 Cassette Autoclave of communicating directly with service personnel. The 5000 G4 Cassette Autoclave, through its LCD touch screen control and networking capability, has been designed to provide you with unprecedented ease of use and puts the user in total control.

Features and Benefits:

- Large 3.5" high-resolution touch screen offers a vivid display of messages and current cycle information all with extraordinary clarity
- **Ultra-Fast** 9 minute unwrapped cycle
- 17.5 minute wrapped cycle
- Exceptional instrument turnover
- 15 minutes Dri0Tec Drying Technology
- Reduced Instrument Investment
- Prolongs instrument life

Technical Specifications:

- Unit size: 23.6" x 16.1" x 7.5" / 60 cm x 41 cm x 19 cm
- Cassette Internal dimensions: 15" x 7" x 3" / 38 cm x 18 cm x 7.5 cm
- Weight w/o water: 73lbs/33 Kg
- Power Consumption: 115 V, 50/60 Hz 15A/ 1450 W

ONE YEAR WARRANTY

Unit

30-26011

Prestige Cassette Sterilizer - 2100

It has gained an enviable reputation for reliability, durability and quality. *One fact alone proves that the Prestige Medical 2100 is committed to quality-a lifetime guarantee on the heating element.*

- Simple 1- step operation starts the fully automatic cycle.
- Light sequence indicators show "Power to the unit:" - "Cycle in progress" - "Sterilization cycle complete"
- Electronic detectors turn off power if there is insufficient water charge.

Specifications

- 9 Liter, 2.38 Gal
- Chamber size: 8.26", 210 mm diameter X 9.29", 236 mm height
- Sterilizing Temp: 250F, 121 degrees
- Automatic (electronic) controlled
- Steam Sterilizer Power supply: 110V
- 16 minute cycle



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

- 110V., 1200 Watt
- Vertical Style

Unit

30-26011



2340 M

We Carry Brands Like:

Tuttnauer
Your Sterilization & Infection Control Partners

Buy any EZ Autoclave...



Closed Door System Fully Automatic Autoclaves

- Hepa Filtered Air Drying • Cassette Capacity 3 large + 3 half

EZ Models come with a 2 Year Parts and Labor Warranty

cpac
Equipment Inc.

And More!!!



Call us for Quarterly Promotion

Prestigedentalproducts.com



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

