

Bite Registration



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Bite Registration



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

BiteVance Bite Registration Material (ATW)

BiteVance is a vinyl Polysiloxane Bite Registration Material:

- Optimal hardness, 80+ durometer
- Thixotropic for easy application
- Non-slumping consistency
- Fluffy, mousse-like
- Excellent dimensional stability
- Improves fit & alignment of crowns & bridges
- Used for bite registration, quick die for temporary inlays & onlays, copying transfers, and primary impression with bite tray
- wont compress or flex
- 30 seconds working time, and 50 seconds setting time



Fast Set 2x 50ml Cartridges & 6 mixing tips	16-11129
Fast Set 10 x 50ml Cartridges, No tips	16-12129
Fast Set 60 x 50ml Cart. No tips	16-84651

Colorbite Rock (Zhermack)

Thermochromic VPS for bite registration formulated with pigments that permit the visible control of the various working steps to minimize the margin of error



- **Changes color from purple to pink**
- Thixotropic, does not slump, mousse consistency
- **Rigid set (95 Shore-A)**
- Very precise
- Extra Fast
- **Time in mouth less than 45 seconds**

Purple to Pink To Purple (Purple at room temperature) 12x green mixing tips	16-11132
--	----------

Defend Bite Registration

- optimal hardness, 80+ durometer
- Thixotropic for easy application
- Non-slumping consistency
- Fluffy, mousse-like
- **Excellent** dimensional stability
- Improves fit & alignment of crowns & bridges
- **Used for** bite registration, quick die for temporary inlays & onlays, coping transfers, and primary impression with bite tray
- wont compress or flex



2 x 50ml Cartridges + Tips fast set, work - set time 1 minute	16-11177
Bulk 60 x 50ml Cartridges w/o tips	16-111781
4 x 50ml Cartridges w/o tips	16-111782

Triple Trays (Premier)

Proven Accurate In Over 80 Million Impressions



All Premier Triple Tray designs have very thin, loose webbing which will not distort when the patient closes. Look for Premier on the handle; your assurance of a genuine Triple Tray. The superior performance of Triple Tray, combined with today's very accurate materials, provides you with an ideal impression in an easy technique. Triple Tray has been consistently rated an outstanding product by a leading independent research organization due to "least water absorption and thinnest inner occlusal wafer

Anterior 35/Box	16-11142
Quadrant 40/Box	16-11141
Posterior 48/Box	16-11140
Sideless Posterior 48/Box	16-11143
Assorted Trays 40/Box	16-11144
3/4 Arch Partial 35/Box	16-11145
Ext Quadrant 40/Box	16-11146
Wide Body Posterior 48/Box	16-11147

Bite Registration Accessories (Zhermack)

D2 Dispenser (1:1)	16-11115
Replacement Plunger for Universal Dispenser 50ml	16-11114
D3 Multi-Dispenser (1:1, 2:1, 10:1, 4:1) 1x Applicator / 1x 4:1 plunger/ 1x 1:1/2:1 plunger/ 1x 10:1 plunger/ 1x1: 1/2:1breach/ 1x 4:1/10:1 breach	16-11113
Large mixing tips 48x green mixing tips	16-11112

Bite Registration Trays

The right selection of durable, solid plastic bite trays for different impressions



Bite Trays-Anterior 35 Pcs/Box	16-11120
Bite Trays-Posterior 50 Pcs/Box	16-11119
Bite Trays-Quadrant 35 Pcs/Box	16-11116
Bite Trays-Sideless 50 Pcs/Box	16-11117
Bite Trays-Full Arch 30 Pcs/Box	16-11118



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Blu-Mousse Bite Registration (Parkell)



Superfast 30-Second Set Split Cartridge Contains: 2 cartridges (100ml.) total, mixing tips 37-89901
Split Cartridge Mixing Tips w/Pink Base Pkg. of 30 37-89906

Colorbite D (Zhermack)

Accurate and stable impressions: zero thickness.

• **Changes color from green to yellow**

• Thixotropic, does not slump
• Extremely hard after setting (45 Shore-D)

• Very precise

• **Extra Fast**

• **Time in mouth less than 60 seconds**

Colorbite D, changes color from green to yellow, 2x 50ml cartridges/ 12x green mixing tips



16-11131

Disposable Bite Tray




Anterior Adult 35/Pkg 37-20017
Anterior Jr. 35/Pkg 37-20018
Posterior 50/Pkg 37-20019
Quadrant 35/Pkg 37-20020
Sideless 50/Pkg 37-20021
Assorted 40/Pkg 37-20022

www.PrestigeDentalProducts.com

Bite Registration

Excellent Bite Trays (Plasdent)



Excellent Anterior Trays 40 Pcs/Box 16-11126
Excellent Posterior Trays 50 Pcs/Box 16-11125
Excellent Sideless Posterior Trays 50 Pcs/Box 16-11124
Excellent Quadrant Trays 40 Pcs/Box 16-11123
Excellent Full Arch Trays 40 Pcs/Box 16-11122
Excellent Anterior, Posterior, Sideless Posterior, Quadrant 40pcs/box 10 pcs/each 16-11121

Hygenic Bite Wafers with Foil (Coltene)



U-Shaped Bite Wafers w/ Foil Insert 45-108701
- Light Blue, 24/Pkg

Occlufast CAD (Zhermack)

A-Silicone for bite registration developed to be used with all CAD CAM systems

Advantages:

- Compatible with CAD-CAM systems using optical, laser and contact reading
- Fast: 60 secs permanence in the oral cavity
- Imperceptible consistency for the patient who bites naturally
- After setting occlufast CAD can be finished with a cutter and/or bur
- Product easy to apply in the mouth, both with the bite technique and vestibular side applications

Characteristics:

Thixotropic, does not drip, hard (95 Shore A),Precise, dimensionally stable over time and to disinfection, can be repositioned in the mouth and on the models after setting, light blue color



Call us for Quarterly Promotion

Prestigedentalproducts.com

Bite Registration

•Rigid set (95 Shore-A)

- Very precise
- Extra fast
- Time in mouth less than 60 seconds

Blue 2x 50ml cartridges/12x green mixing tips 16-11133

Occlufast Rock (Zhermack)

A-Silicone for bite registration. Maximum intercuspitation bite registration, centric, protusion or lateral registrations, gnathological evaluations, diagnostic evaluation of cranialcervical-mandibular disorders, orthodontic centric occlusion



Advantages:

- Imperceptible consistency for the patient Easy handling and bubble-free Positioning control
- Accurate and stable impressions: zero thickness
- Easy to trim for excess material removal
- Flow able: it avoids pressure and slipping during positioning
- Accurate: maximum dimensional variation of just -0.05% after 24 hours allows long-term stability
- Non-deformable: it is not affected by any thermal shocks during transport

Characteristics:

Final hardness 95 Shore-A, Fast: 60sec permanence in the oral cavity, Thixotropic purple color

Violet 2x 50ml cartridges/ 12x green mixing tips 16-11134

Opti-Bite Bite Registration (Pac-Dent)



- Fast procedure time
- Generous work time
- Excellent handling and accuracy
- Included ribbon tips ensure easy dispensing

Fast set, 2x 50ml cartridges + 6 Tips 16-98601

Fast set Mint Flavor, 4 50ml cartridges 16-98602



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

Regular set, 2 x 50ml cartridges +6 Tips 16-98603

Super fast, 2 x 50ml cartridges + 6 Tips 16-98604

Teal Mixing Tips, 1:1, 6.5mm, 48/Pkg 16-98605

Regular set, 4 x 50ml Cartridges 16-98606

PolyBite Bite Trays (Dentamerica)

Impression, counter impression, bite registration 3 in one disposable tray



Anterior, 35/Pkg 16-11135


Posterior, 50/Pkg 16-11136

Quadrant, 35/Pkg 16-11137

Full Arch, 30/Pkg 16-11138

Sideless Posterior, 50/Pkg 16-11139

Thermobite Bite Registration Wafers (Dentsply)



- Specially-formulated resin provides highly accurate bite impressions of teeth and horizontal molar area
- Heat according to instructions, bite down, and trim
- Minor perforations allowed, for a correct occlusion with touching teeth surfaces
- Easy to store and transport to labs; resistant to temperature changes during transit

25 Pcs/Pkg 45-14069



*We Work Hard
So You Don't Have
To.*

877-772-3888



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Bite Registration

You'll find Brands like:

Bite Vance ATW



And More!!!



Call us for Quarterly Promotion

Prestigedentalproducts.com



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

