







Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Apex Locators

Endo Combo Package (Parkell)



Includes:

- Formatron Apex Locator

DRAMATRON D


- Digitest II Pulp Vitality Tester


Endo Combo Package - Includes Digitest 29-110301
II Pulp Vitality Tester and Formatron
Apex Locator
Touch-Proof Wire Assembly 29-110302

Apex Locator	
Touch-Proof Wire Assembly	29-110302
Quick Test Probe	29-110303
Replacement Ground Clips	29-110304
Reamer Holder	29-110305
Metal Probes for Pulp Tester	29-110306
Pulp Tester Ground Clip & Lead Wire	29-110307
Precision Probe set for Pulp Testers	29-110308

Foramatron Apex Locator - Parkell

- Precisely pinpoints the apical constriction using high and low-frequency bands
- Ergonomic, Aseptic design
- LED illuminated display for easy viewing
- 4 distinct audio signals reinforce the indicator lights
- Reliable results even in the presence Of all types of fluids.





Formatron D-10 Touch Proof Wire 29-110751 Assembly

Quick Probe 29-110752

I-Root Apex Locator (Meta Biomed)





5th generation Apex Locator by dual frequency type.

Best accuracy in any root canal condition (Dry, Bleeding, Wet, Saline, EDTA, NaOCl or Chlorhexidine etc...) Can measure the apical constriction exactly without mis-measurement. Developed & upgraded based on the technology of e-Magic Finder (EMF-100 Series) Apex Locater.

Measurement by eye and ear at the same time: Available for checking the apex by digital No., graph and beep sound. When the file is closer to the apex, indication sound is intermittent beep sound. When reaches to the apex It provides continuous beep sound and red Apex indication for user's convenient and saving time, it helps to find apex easily and quickly.

Color Display: From graph 2.0 to 1.0, Indicates graph by blue color. From graph 0.9 to 0.5, Indicates graph by green color. From graph 0.4 to -0.5, Indicates graph by red color. From the 0.5mm to -0.5mm. The tooth graphic and Apex indication will be blink with color display.

I-Root A Type Digital Apex Locator,

Rechargeable battery / Without PM 29-11747

Rechargeable Liion Battery, 2200 29-12749 mAh 7.2 Volts, Each

Mark VII Apex Locator (Miltex)

The Mark VII Apex Locator gives dental professionals an exact measurement of the canal length under all conditions, including dry, wet and bleeding. This makes working in the canal easier and more accurate. It precisely measures



the distance from the tip of the endodontic file to the foramen. The Mark VII allows dental professionals to know the depth of the canal by showing when they are 2.0 mm, 1.5 mm, 1.0 mm, 0.5 mm, 0.25 mm away from, and then the apex has been reached.

The Mark VII Apex Locator is easy to use and designed to be placed near the oral cavity, in the operator's field of vision. The apex locator is often clipped to the patient's napkin to ensure that the unit does not accidentally fall. The Mark VII does not require calibration between measuring multiple canals, thereby immediately enabling accurate measurement of another canal.

Features:

- •Can be used in any canal condition
- •Provides Accurate Readings
- •Compact Size, Easy to Use
- •Does not require calibration at the start of a procedure or between canals

Benefits:

•Saves time and canal does not require certain condition

- •Allows dental professional to know where they are in the canal in relation to the apex
- •Can be placed near the oral cavity, in the operator's field of vision
- •Does not require any counter space
- •Reduces procedural and set up time
- •Reduces tray set up time and time it takes to measure multiple canals

Mark VII Unit 29-11982

Root ZX II Apex Locator (J-Morita)



Apex Locator with Low Speed Hand piece

As the world's best selling apex locator, Root ZX II has set the industry standard for accuracy and patient safety.

Invented in 1992, it was the first unit with the ability to work in wet canals - measurement is not affected by the presence or absence of blood, other discharges, or electrolytes. Root ZX® II features patented technology which offers an accuracy rate of 97.5%*. User friendly, its display screen is large and easy to read. The action of the meter in the display corresponds exactly to the tactile sensation of using the file. Other features include: slim, lightweight file holder, no zero-adjustment, automatic calibration, battery power indication, and automatic power off function.

As an optional upgrade, a module enables the unit to function as a low speed handpiece, apex locator, or a combination of both. The module snaps easily onto the back of the unit and offers numerous safety features including Auto Apical Reverse and Auto Apical Slow Down.

APEX LOCATOR FEATURES

- •97.5% accuracy rate*
- •Automatic calibration
- •Ability to work in wet canals
- •Easy to read display screen
- •Upgrade option to Low Speed Handpiece

Long File Holder

A long touch file holder is available to be used with the probe cord. It features a fork-like tip that allows the clinician to get an apex reading by simply touching the file. The file holder is 5 inches long, autoclave safe and ideal for hard-to-reach molar procedures, thus increasing access to the work area.

Indications

Apex location/measurement of root canal length

Canal Measurement Module, Includes:

Probe Cord, File Holders (3), Contrary Electrodes (5), Function Tester, "AA" Alkali Batteries (3)

29-12409

Low Speed Handpiece Module, In-

cludes: Micromotor (TR800), Contra Angle, Handpiece Cord, Handpiece Rest, AR Oil, AC Adapter, Battery (pre-installed) Foot Switch,

29-12408

Plastic Sleeves (20)

Micromotor (TR400)

29-12407

Endodontics

Micromotor (TR800)	29-12405
Contra Head	29-12406
Heandpiece Cord	29-12404
Handpiece Rest	29-12403
Probe Cord (Root ZX II)	29-12402
Probe Cord (Original Root ZX)	29-12401
File Holders (5)	29-12399
Contrary Electrodes (Lip Clips) (5)	29-12400
Long File Holder	29-12398

Root ZX Mini Apex Locator (J-Morita)



The Root ZX Mini offers the same, proven, Root ZX apex location technology which has been independently **evaluated at a 97.5% accuracy rate.** The LCD display screen has maintained excellent readability with a clear, progressive display and high contrast. As with the original model, Root ZX mini offers 3 programmable memory settings. Other features include shock resistance, automatic calibration and power function.

- Compact, conveniently portable design
- Measures 4" x 2 1/3"
- Programmable memory settings
- Easy to read screen with progressive display
- Shock resistance

Mini Unit: Probe Cord, 3 x File Holders, 5 x Contrary Electrodes, 1 x Canal-Length Tester, 3 x DC 1.5 V, AAA size Batteries

RootPro Endo Motor-Apex Locator (Pac-Dent)

CANAL PREPARA-TION & APEX LOCATOR SYTEM

RootProTM Wireless is the most advanced apex locator / endo motor all -in-one system of its kind. With a newly added wireless foot pedal, the RootPro dramatically reduces the clutter and improves your efficiency



The RootProTM endodontic motor with Low-Speed Handpiece is a dual-function device that prepares root canals while displaying accurate measurements of the root canal. The RootProTM offers many automatic safety control options, making it the ideal choice for Dentist or Endodontist. It provides the user with expected predictability in today's dental practice.RootProTM also accurately and precisely locates the position of the file inside a root canal without requiring the user

to set the machine to zero. RootProTM's unique automatic calibration function ensures accuracy and eliminates the effect of changes in temperature, moisture, etc. inside the canal even during treatment. It comes with a wireless foot switch, dramatically reducing clutter and drag during operation, ensuring a user-friendly experience through-

ADVANTAGES

- •Prepares root canals while displaying accurate & precise measurements of the root canal
- •Offers many automatic safety control options: Auto Motor Start/ Stop, Auto Torque Reverse, Auto Apical Reverse
- •Accurately and precisely locates the position of the file inside a root canal
- •No zero-adjustments needed
- •Automatic calibration ensures accuracy & eliminates the effect of changes in temperature and moisture
- •Wireless foot switch for a user-friendly & clutter-free operation

1x Complete RootPro Unit: 1 x Root-

Pro main unit with low-speed handpiece, 1 x probe cord, 3 x file holders, 1 x handpiece rest, 1 x foot switch, 1 x ac adapter, 5 x contrary electrodes, 1 x

canal-length tester, 1 x AR oil 29-25480 1x Low-speed handpiece replacement 29-25481 **Foot Switch** 29-25482 AC Adapter (110V~240V, 50/60Hz) 29-25483 Low-Speed contra-angle handpiece 29-25484 sleeves, 250/box

Endo Handpieces

E&Q Master Obturation Systems (Meta Biomed)



Maximizes convenience and minimizes chair time - Combined filling method of vertical compaction and backfill makes the 3 dimensional cordless obturation system a perfect choice for apical part, accessory and lateral canals.

Full Unit	29-106421
Pen Unit	29-106422
Gun Unit	29-106423
Gutta Percha Bars, 100 Bars/Box	29-106424
Pen Tips, 2 each/pac, (F&FM)	29-106425
Pen Tip, each	29-106426
Gun Needle, 6 each/pack, (23G & 25G)	29-106427
Heat Insulator, 5/Pkg	29-106429
Plunger Head, 2 each/Pkg	29-106430

www.prestigedentalproducts.com



AutoSyringe Irrigation Handpiece



- Go with the flow, controlled flow
- Comfort, Convenience, & Control
- Eliminate hand fatigue associated with the standard syringe irrigation
- Easy to fill, quick

change 20ml canisters

- Safe, controlled flow of irrigating solutions with variable flow rate settings
- Save you MONEY, time, hand fatigue and Waste

Complete Kit: AutoSyringe Handpiece, 3x20ml Canisters, 20 Canister Liners, 3 Autoclavable Sleeves, 100 Plastic Barrier Sleeves & Changing Base

29-108641

Cordless Thermal Condensor Unit (Dentazon)



- lightest cordless warm vertical condensation device
- effectively liquefies and condenses root canal filing materials
- adjustable 3 temperature settings
- color-coded pen tips available in various sizes (XF, F, FM, M, ML)
- excellent apical control over lateral and complicated root canals

1 year warranty

Unit 29-11221

Digitest II Pulp Vitality Tester (Parkell)



Includes: 4 autoclavable probe tips, ground wire, clip, 9-volt battery and durable case.

Metal Probes for pulp test - 2 probes; 1 short 29-110322 & long. They fit the Digitest II & Gentle Pulse

Call us for Quarterly Promotion

Precision Probe Set for Pulp Testers - 1 precision labial and 1 precision lingual probe. They fit the Digitest II and Gentle Pulse

29-110323

Replacement Ground Clips - 2 clips. For use with the Digitest II, Gentle Pulse, and **Foramatron**

29-110324

E&Q Plus Gutta Percha Obturator (Meta Biomed)



Gives the clinician the option to choose from either vertical compaction of warm guttapercha, heat softened gutta-percha injection or

combination of techniques in filling the canals.

- Easy & Quick obturation w/ accurate and predictable filing
- possible best hermetic sealing w/ 3D filing
- convenient and efficient 2-in-1 filing system
- easy to learn and master
- easy touch button allows easy temperature setting

Full Unit

29-100711

Endo Irrigator System (Vista)

- ·Instantaneous.
- fingertip-controlled system •Site-specific pinpoint irrigation is cleaner and more
- efficient
- •Handpiece control allows you to instantly select irrigating solutions with the touch of a button, without leaving the field •Built-in heater allows for the delivery of heated irrigating
- solutions for improved chemical reaction •Improves procedure safety and
- efficiency and significantly decrease procedure time
- •Multiple mounting options
- •OPTIONAL Air Handpiece that includes the Stropko-IrrigatorTM for the delivery of low pressure, moisture free air to your surgical site

Kit Includes:

- •Endo-Irrigator™ II Irrigation & Evacuation System & Power Cord
 - •1 Autoclavable tip manifold
 - •30 assorted Endo-IrrigatorTM irrigating tips
 - •1 Installation Kit
 - •Choice of wall or pole mount

Unit without heat includes: wall or

pole mount, 30 assorted irrigating tips 29-11313

Unit heated 3 bottle includes: wall or

pole mount, 30 assorted irrigating tips 29-11312

1-Button Air/Water Syringe Kit

(Includes 7ft. tubing and metal tip) 29-11314





Prestige Dental Products, Inc. Prestigedentalproducts.com Tel: 877-772-3888

E-Cube Black Endo Rotray System W/Contra Angle (Saeshin)

Compare to "NSK Endo-Mate DT"

Auto reverse feature upon overload sensor detection 9 customizable user programmable preset programs Up to 5 hour batter life back up and electrically powered Constant torque settings for safe control Speed ranges from 125 Rpm - 13,000 Light weight handpiece (3.74 oz) Slim style handpiece (4.25" long) and 3/4" diameter Universal power supply





- OAuto Reverse function
- O Torque control from 0.6 to 5.2 N.cm
- OGear ratio: 16:1 @9 Programmable preset memorizing
- Poot control On/Off switch
- OUp to 5-hour battery life OCompact designed controller including
- various functi O Supply with stable torque
- Operation available with battery or adapter



- (1) AC Adapter / cable (2) Control Unit (2type) (3) On/Off Switch (Option) (4) Handpiece Stand (5) Contra Angle (1EA) (6) Motor (1FA)

Input	100~240V 47/63Hz	Rated input	DC20V 0.5A
	0.9~0.34A	Output	DC4.5V 0.4A
Output	20V 1.42A max	Charging time	5 hours approx.
Dimensions	W 62 X D 118 X H 39mm	Dimensions	W 110 X D 196 X H 139mm
Weight	298g	Weight	582g
Micro Motor	(CUBE EM)	Micro Motor	(CUBE EP)
Rated input	DC4.5V 0.4A	Rated input	DC4.5V 0.4A
Dimensions	¢ 19.8 X L 61.7mm	Dimensions	¢ 20 X L108mm
Weight	116g	Weight	106a
RPM	13,000rpm	RPM	13,000rpm
Contra Angle	Handpiece (ACL(B)-42EM)	Contra Angle	Handpiece (ACL(B)-42EP)
Gear Ratio	16:1	Gear Ratio	16:1
Dimensions	¢ 19.6 X L 95.2mm	Dimensions	¢16.7 X L63mm
Difficitsions		Weight	

E-Cube Black (Auto Reverse) w/Plastic HP & Angle PEDAL ON/OFF	29-11232
E-Cube Blue (Non-Auto Reverse) w/Plastic HP & Angle PEDAL ON/OFF	29-11234
E-Cube Black (Auto Reverse) w/Metal HP & Angle PEDAL ON/OFF	29-11235
E-Cube Black (Auto Reverse) Controller Only	29-11236
E-Cube Blue (Non-Auto Reverse) Controller Only	29-11237
E-Cube Handpiece Type A (On/Off Button) - Plastic	29-11238
E-Cube Handpiece Type B (Non On/Off Button) - Metal	29-11239
E-Cube Charger 110-240V	29-11240
E-Cube On-Off Pedal Only	29-11241
E-Cube Handpiece Stand	29-11242
16:1 Mini Head Contra Angle Plastic	29-11243
16:1 Mini Head Contra Angle Metal	29-11244

Endo A Class Handpiece (Saeyang)

Endo a class



This is a **cordless** endodontic handpiece, which can drive from 120 to 2,000rpm with a torque of 3.0 NCM (16:1). It has 4 program memories and 4 speed ratios (4:1, 10:1, 16:1, and 20:1). Its **VERY** light weight (138g) and portable.

- When the preset torque exceeds the limit; the handpiece stops and immediately runs in reverse automatically
- Its equipped with a quiet motor and variable setting mode
- Contra angle is easy to disconnect and has a **comfortable grip**
- It includes a rechargeable battery with 1 hour and a half operation Made in South Korea

FDA Approved, & 510K approved #K12358

29-101381 **Full Unit**

Call Us Today 833-733-6825



Prestige Dental Products, Inc. Prestigedentalproducts.com Tel: 877-772-3888

Endo Motor S5 Kit (JS-Dental)



Complete Starter Kit Includes: Handpiece, charger, cable, & 2 packages of assorted S5 Rotary Files

29-110941

EndoApex Cordless Obturation System (Dentazon)

Light weight and highly efficient cordless endodontic dual obturation system which offers optimum functionality and reliable apical control along with excellent tactile

- •360° Easy Swivel (no tool required)
- •Complete thermal protection for patient
- •Ultra-fast heating for less chair time
- •Ergonomic and user-friendly operation
- •Maximum portability with stylish design
- •Greatly reduces hand and wrist fatigue
- •Color coded pen tips
- •Delivers reliable and predictable clinical results

1 Year Warranty

2-in-1 Obturation System

29-11315

- 1/4 turn

twisting motion. - Allows a 30degree onequarter turn.

Endodontic C.A. (DA)

Each



- Fully autoclavable 29-92171

29-108651

EndoUltra Cordless Endo Handpiece




- Easy-to Use & cordless handheld device.
- ONLY unit to activate tip frequency to (40,000 Hz) required to create sufficient acoustical streaming and the cavitation necessary to effectively clean, penetrate and remove vapor lock.

Complete Kit Includes: 6 autoclavable NiTi Activator tips which resonate down the entire length of the tip. These unique, multi-use tips will not engage or remove tooth structure. They come in 2 sizes: 15/02 & 25/04 w/

depth markers at 18mm, 19mm, and 20mm. Tip frequency is 40 kHz.

Endomax Cordless Endodontic Handpiece (Dentamerica)



Compare to "NSK Endo-MateTC"

- •Auto-stop and auto-reverse mechanisms for patient safety
- •Heavy-duty high torque motor
- •Reciprocating setting available without use of additional attachment: The file reciprocates in a back and forth motion at 45°.
- •Ergonomic design provides comfort and easy operation
- •Features five torque settings
- •Automatic forward/reverse functions
 - Clockwise rotation
 - •Reverses when torque limit is reached
 - •When the load is removed, rotation continues forward
 - •Auto-Reverse Stop
- •Five rotation speeds available for different applications
- •High capacity lithium-ion battery provides stable and prolonged usage
- •Offers quick and easy coupling
- •Weight: 150g

Endomax Cordless Endodontic

Handpiece Unit 29-11346

Replacement:

Latch Attachment (for Endomax) 29-22346

Replacement:

Battery Pack (for Endomax) 29-22347

Gutta Easy Cordless Thermal Obturation System (Dentazon)

Light weight and highly efficient cordless endodontic obturation system which offers optimum functionality and reliable apical control along with excellent tactile sense.

- •360° Easy Swivel (no tool required)
- •Complete thermal protection for patient
- •Ultra-fast heating for less chair time
- ·Ergonomic and user-friendly operation
- •Maximum portability with stylish design
- •Greatly reduces hand and wrist fatigue
- •Color coded pen tips
- •Delivers reliable and predictable clinical results 1 Year Warranty

Unit 29-11538



Prestige Dental Products, Inc. Prestigedentalproducts.com

Tel: 877-772-3888

Gutta Percha Cutter and Trimmer

In 3 seconds only

This instrument is developed for the easy and safe removal and trimming of excess gutta-percha and thermal dense filling in your endodontic cases. It is composed ofendoremover and trimmer operated by two AAA, NiMH rechargeable batteries, which are included in the package.

ERGÖNOMIC DESIGN

SLIM: To feel comfortable for the treatment in the mouth like dental handpieces

LIGHTWEIGHT: It weighs only 50g together wit 2 cells which leads to easiness in your treatment PRODUCTIVE ENVIRONMENTS

There is no need for alcohol lamp or burner to heat some instrument to remove the excess, also no need for an assistant. **ECONOMICAL:**

If NI-HM batteries are used, you can use this instrument in very powerful and economical ways than disposable alkaline batteries. Normally, it has about 400-500 charging cycles which leads to long lasting usages, also if fully charged, normal usage is between 15-20 times

SAFETY & DURABILITY

This instrument might be semi-permanent only if batteries and tip are refilled, also the patient can not notice the treatment because they can not see compared than the usage of conventional alcohol lamp usage. The patients feels comfortable.

Prior to using this instrument ,please charge and check the batter-

- 1. Push the button only when are ready to operate
- 2. Place the edge of the tip to the canal orifice, and then push the button for 3 or 4 seconds. It produces enough heat to cut the excess within seconds. That means you don't need to hold the button for overlong time.
- 3. After removing the excess ,please pull the instrument carefully not to touch the skin of the patient.
- 4. After it comes cool ,please clean the tip with wet alcohol gauze or sterilize for the next usage.

Cutter & Trimmer Unit: with 4 X

AAA rechargeable Ni-MH batteries, Rapid charger and two cutting and trimmer tips

29-11539

Root ZX Mini Apex Locator (J-Morita)

The Root ZX Mini offers the same, proven, Root ZX apex location technology which has been independently evaluated at a 97.5% accuracy rate. The LCD display screen has maintained excellent readability with a clear, progressive display and high contrast. As with the original model, Root ZX mini offers 3 programmable memory settings. Other features

include shock resistance, automatic calibration and power function.

- Compact, conveniently portable design
- Measures 4" x 2 1/3"
- Programmable memory settings
- Easy to read screen with progressive display
- Shock resistance



Mini Unit: Probe Cord, 3 x File Holders, 5 x Contrary Electrodes, 1 x Canal-Length Tester, 3 x DC 1.5 V, AAA size Batteries

I-Root Apex Locator (Meta Biomed)



5th generation Apex Locator by dual frequency type.

Best accuracy in any root canal condition (Dry, Bleeding, Wet, Saline, EDTA, NaOCl or Chlorhexidine etc...) Can measure the apical constriction exactly without mis-measurement. Developed & upgraded based on the technology of e-Magic Finder (EMF-100 Series) Apex Locater.

Measurement by eye and ear at the same time: Available for checking the apex by digital No., graph and beep sound. When the file is closer to the apex, indication sound is intermittent beep sound. When reaches to the apex It provides continuous beep sound and red Apex indication for user's convenient and saving time, it helps to find apex easily and quickly.

Color Display: From graph 2.0 to 1.0, Indicates graph by blue color. From graph 0.9 to 0.5, Indicates graph by green color. From graph 0.4 to -0.5, Indicates graph by red color. From the 0.5mm to -0.5mm. The tooth graphic and Apex indication will be blink with color display.

I-Root A Type Digital Apex Locator,

Rechargeable battery / Without PM


29-11747

Rechargeable Liion Battery, 2200

29-12749

mAh 7.2 Volts, Each

Integrity Endo Handpiece (Osseo Scientific)





Each 29-92161

Mark VII Apex Locator (Miltex)

The Mark VII Apex Locator gives dental professionals an exact measurement of the canal length under all conditions, including dry, wet and bleeding. This makes working in the canal easier and more accurate. It precisely measures the distance from the tip of the endodontic file to the foramen. The Mark VII



allows dental professionals to know the depth of the canal by showing when they are 2.0 mm, 1.5 mm, 1.0 mm, 0.5 mm, 0.25 mm away from, and then the apex has been reached. The Mark VII Apex Locator is easy to use and designed to be placed near the oral cavity, in



the operator's field of vision. The apex locator is often clipped to the patient's napkin to ensure that the unit does not accidentally fall. The Mark VII does not require calibration between measuring multiple canals, thereby immediately enabling accurate measurement of another canal.

Features:

- •Can be used in any canal condition
- •Provides Accurate Readings
- •Compact Size, Easy to Use
- •Does not require calibration at the start of a procedure or between canals

Benefits:

- •Saves time and canal does not require certain condition
- •Allows dental professional to know where they are in the canal in relation to the apex
- •Can be placed near the oral cavity, in the operator's field of vision
- •Does not require any counter space
- •Reduces procedural and set up time
- •Reduces tray set up time and time it takes to measure multiple canals

Mark VII Unit

29-11982

Root ZX II Apex Locator (J-Morita)



Apex Locator with Low Speed Hand piece

As the world's best selling apex locator, Root ZX II has set the industry standard for accuracy and patient safety.

Invented in 1992, it was the first unit with the ability to work in wet canals - measurement is not affected by the presence or absence of blood, other discharges, or electrolytes. Root ZX® II features patented technology which offers an accuracy rate of 97.5%*. User friendly, its display screen is large and easy to read. The action of the meter in the display corresponds exactly to the tactile sensation of using the file. Other features include: slim, lightweight file holder, no zero-adjustment, automatic calibration, battery power indication, and automatic power off function.

As an optional upgrade, a module enables the unit to function as a low speed handpiece, apex locator, or a combination of both. The module snaps easily onto the back of the unit and offers numerous safety features including Auto Apical Reverse and Auto Apical Slow Down.

APEX LOCATOR FEATURES

- •97.5% accuracy rate*
- Automatic calibration
- •Ability to work in wet canals
- •Easy to read display screen
- •Upgrade option to Low Speed Handpiece

Long File Holder

A long touch file holder is available to be used with the probe cord. It features a fork-like tip that allows the clinician to get an apex reading by simply touching the file. The file holder is 5 inches long, autoclave safe and ideal for hard-to-reach molar procedures, thus increasing access to the work area.

Indications

Apex location/measurement of root canal length

Canal Measurement Module, Includes: Probe Cord, File Holders (3), Contrary Electrodes (5), Function Tester, "AA" Alkali Batteries (3)	29-12409
Low Speed Handpiece Module, Includes: Micromotor (TR800), Contra Angle, Handpiece Cord, Handpiece Rest, AR Oil, AC Adapter, Battery (pre-installed) Foot Switch, Plastic Sleeves (20)	29-12408
Micromotor (TR400)	29-12407
Micromotor (TR800)	29-12405
Contra Head	29-12406
Heandpiece Cord	29-12404
Handpiece Rest	29-12403
Probe Cord (Root ZX II)	29-12402
Probe Cord (Original Root ZX)	29-12401
File Holders (5)	29-12399
Contrary Electrodes (Lip Clips) (5)	29-12400
Long File Holder	29-12398

RootPro CL Endo Cordless Handpiece



It is ergonomically designed with a thin neck and contra angle to provide unobstructed view & superior access in all endo treatments. The simple LCD display and user friendly control streamline clinical work.

- 5 programmable memory
- Auto-reverse function
- Contra angle head can be adjusted 6 positions to suit users needs
- Full charge lasts 80 minutes at rated load
- Auto-off after 10 minutes idle
- -Autoclavable (135°C) contra angle head
- **Speed:** $140 \sim 550$ RPM **Torque:** $0.3 \sim 3.0$ N-cm.

Each 29-108831

RootPro Endo Motor-Apex Locator (Pac-Dent)

CANAL PREPARATION & APEX LOCATOR SYTEM

RootProTM Wireless is the most advanced apex locator / endo motor all-in-one system of its kind. With a newly added wireless foot pedal, the RootPro dramatically reduces the clutter and improves your efficiency



The RootProTM endodontic

motor with Low-Speed Handpiece is a dual-function device that prepares root canals while displaying accurate measurements of the root canal. The RootProTM offers many automatic safety control options, making it the ideal choice for Dentist or Endodontist. It provides the user with expected predictability in today's dental practice.RootProTM also accurately and precisely locates the position of the file inside a root canal without requiring the user



to set the machine to zero. RootProTM's unique automatic calibration function ensures accuracy and eliminates the effect of changes in temperature, moisture, etc. inside the canal even during treatment. It comes with a wireless foot switch, dramatically reducing clutter and drag during operation, ensuring a user-friendly experience throughout.

ADVANTAGES

- •Prepares root canals while displaying accurate & precise measurements of the root canal
- •Offers many automatic safety control options: Auto Motor Start/ Stop, Auto Torque Reverse, Auto Apical Reverse
- •Accurately and precisely locates the position of the file inside a root canal
- •No zero-adjustments needed
- •Automatic calibration ensures accuracy & eliminates the effect of changes in temperature and moisture
- •Wireless foot switch for a user-friendly & clutter-free operation

1x Complete RootPro Unit: 1 x Root-

Pro main unit with low-speed handpiece, 1 x probe cord, 3 x file holders, 1 x handpiece rest, 1 x foot switch, 1 x ac adapter, 5 x contrary electrodes, 1 x canal-length tester, 1 x AR oil

canal-length tester, 1 x AR oil 29-25480 **1x Low-speed handpiece replacement** 29-25481

Foot Switch 29-25482

AC Adapter (110V~240V, 50/60Hz) 29-25483

Low-Speed contra-angle handpiece sleeves, 250/box

Sonic Air MM 1500 Endo System (Medidenta)

All NiTi systems are excellent in creating crown down effects, but they do not irrigate or cleanse lateral canals. These canals are breeding grounds for bacteria, which in efect will increase post-op pain and inevitable retreatment.

29-25484

The Sonic MM1500 will effectively cleanse the lateral canals with its cavitation effect, thus reducing post-op pain and recalls.

Uses:

- Coronal flaring
- Canal shaping
- Solution irrigation
- Cleansing
- Remove an

isthmus

- Cavitate solution

to clean lateral

canals

- Clean buccal and lingual extensions

Requires lubrication after every four autoclaves

- Use only autoclavable oil
- Bur insert / Release: Pull down head
- Connection: 4-hole
- Autoclavable

Accepts: Sonic Files - Rispi Sonic, Shaper Sonic, Trio-Sonic, Retrprep

1 year warranty

MM 1500 Handpiece, each

Therma-Flo Composite Applicator (Vista)

Therma-FloTM Composite Applicator is revolutionizing composite delivery!

Fill bulk restorations in seconds - multiple levels for improved composite handling properties and increased ability to use highly filles composites.



Features:

- •Automated delivery of composite material
- Warms up to four compules for multiple restorations or layering
- •Adjustable heating levels for improved composite handling properties
- •Allows increased ability to use highly filled composites
- •Tip maintains heat level for optimal temperature at point of delivery
- •Cordless, compact, lightweight ergonomic design reduces hand fatigue
- •Compatible with most standard compules

Full Unit: 20 of 2x Step-Down tips, 20 of 4x Step	29-26201
-down tips, Rechargeable handpiece & warming	
base with 3 nose cones	

Step Down tips 2x, 40 Cup/Pkg 29-26202

Step Down tips 4x, 40 Cup/Pkg 29-26203

Therma-Flo Composite Warming Kit (Vista)



For precision placement of highly filled composites...

Therma-FloTM is a warming kit engineered for optimal composite performance.

\Evidence shows that when highly filled composites are heated, they perform with the same clinical ease of flowables!

Vista's new Therma-FloTM Composite Warming Kit offers you an easy way to increase the flowable viscosity of composite material without compromising tensile strength, wear-resistance or aesthetics.

Benefits of HEAT:

- •Increased flowability = Easy extrusion
- •Improved polymerization= Reduced voids
- •Reduced curing time
- •Greater wear resistance

Composite Warming Kit Features:

- •Heats composite up to 155° F
- •Warms multiple compules dome accommodates multiple syringes
- •Compatible with most standard compules and composite syringes •Includes *patent pending* Step Down tips – for precise placement of composite material

110V Kit: (1) Warmer unit (2) Composite syringes (20) 4X Step Down tips (20) 2X Step Down tips

29-26221

220V Kit: (1) Warmer unit (2) Composite syringes (20) 4X Step Down tips (20) 2X Step Down tips

29-26222

Step Down tips 2x, 40 Cup/Pkg 29-26211-1



Step Down tips 4x, 40 cup/Pkg

29-26212-1

Tri Auto ZX Cordless Endodontic Handpiece (J-Morita)

The cordless Tri Auto ZX is the only endodontic handpiece with a built-in apex locator, providing the capability and convenience to electronically monitor the root canal before, during and after instrumentation. With the combined technology and accuracy of the Root ZX apex locator, the Tri Auto ZX can significantly increase accuracy and safety.

The Tri Auto ZX also offers greater control and flexibility with the adjustable torque settings

The choice of automatic or manual mode operations provide versatility.

The 3 automatic functions include:

- Auto Start/Stop
- •Auto Apical Reverse
- Auto

Unit Includes: Handpiece, Contra Head, Charger, Battery, Probe Cord, File Holders (3), Contrary Electrode (5), Function Tester, AR Oil, Practice Block, Instructional Video

Contra Head 29-12482

Handpiece only

 (Battery & Contra Head not included)
 29-12479

 Charger
 29-12481

 Tri Auto ZX Sleeves (500)
 29-12480

Vibring (Dentazon)

•Compatible with 10cc irrigation syringe of most brands

•No special syringe or attachment is required to use

•It is proven to be perfectly safe for every endodontic irrigation procedure and is 100% compatible with every type of irrigation solution



29-12483

Unit 29-89521

Vibringe Sonic Irrigation System (Cavex)

The Vibringe is the first endodontic sonic irrigation system that enables activation of the irrigation solution, acoustic streaming, in the root canal, in only one step.

The activation of the disinfectant enriches and completes your irrigation procedure and improves the success rate of your endodontic treatments, without any hassle!

nd tic

Vibringe® is developed to make the dental practioners life easier. That is why we design easy-to-use products. At Vibringe B.V. we

believe that we should not interrupt into your actual processes and routines too much. Therefore the Vibringe® embraces your actual irrigation procedure, but will improve your success rate by using the patented sonic flow technology which causes acoustic streaming in the root canal.

Superclean: Vibringe® has been clinically proven to realize a more thorough debridement and disruption of the smear layer than conventional hand irrigation. The patented sonic flow technology will realize a full irrigation to the apex and eliminates blockage. The Vibringe® is enabling the solution to reach and disinfect all portions of the canal. One can think of the lateral canals, fins, webs and anastomoses. When using ultrasonic irrigation, the Vibringe® is a perfect complementary device instead of the normal syringe to realize an even higher predictable outcome.

Vibringe® has clinically proven that it removes organic and inorganic debris from hard-to-reach areas *faster* than the conventional syringe.

Vibringe® is shown to be *completely safe* for all endodontic irrigation procedures and is compatible with all irrigation solutions like for example sodium hypochlorite and EDTA.

Vibringe® designed and developed this product in close cooperation with *general practioners and endodontists*. During this extensive research and development stage finally 8 out of 10 *dentists love to work with Vibringe*®. Main reasons for that? It is an easy-to-use device which embraces conventional irrigation techniques and realizes a higher predictable outcome in shorter time!

One Year Warranty

Including Adapter & 15 disposable syringes	29-12534
Disposable 10 ml Luer Lock Sterilized Syringes, 50/Pkg	29-12533
Power Supply	29-12535

Endo Medicaments

Formocresol Bottle


Endodontic medicament for root canal therapy and for pulp devitilization



1oz Bottle 29-16561

Gates Glidden Drills

Gates Glidden Drill (Mani)



Corresion resistance to any kind of sterilization including autoclave, chemical vapor, and super acid water.

Made in Japan

Size . 1,2,3,4,0,0 Assorted 1-0	
Size #1 - 32mm, 6/Box	29-105681
Size #2 - 32mm, 6/Box	29-105682
Size #3 - 32mm, 6/Box	29-105683
Size #4 - 32mm, 6/Box	29-105684
Size #5 - 32mm, 6/Box	29-105685



Size #6 - 32mm, 6/Box	29-105686
Assorted Sizes #1-6 - 32mm, 6/Box	29-105687

Gates Glidden Drills Stainless Steel

•Swiss made latch type instrument for enlarging the coronal portion of the tooth •450-800 rpm

Size 1 ø 0.50, Length Short 28mm	29-11536
Size 2 ø 0.70, Length Short 28mm	29-11524
Size 3 ø 0.90, Length Short 28mm	29-11535
Size 4 ø 1.10, Length Short 28mm	29-11534
Size 5 ø 1.30, Length Short 28mm	29-11532
Size 6 ø 1.50, Length Short 28mm	29-11533
Size 1-6 Asst, Length Short 28mm	29-11531
Size 1 ø 0.50, Length 32mm	29-11530
Size 2 ø 0.70, Length 32mm	29-11529
Size 3 ø 0.90, Length 32mm	29-11528
Size 4 ø 1.10, Length 32mm	29-11527
Size 5 ø 1.30, Length 32mm	29-11526
Size 6 ø 1.50, Length 32mm	29-11525
Size 1-6 Asst, Length 32mm	29-11537

Mity Gates Glidden Nickel-Titanium

- •ISO color coded aluminum handle
- •Instrument used to enlarge the coronal portion of tooth

Size 1, Length 32mm, 6/Pkg	29-12027
Size 2, Length 32mm, 6/Pkg	29-12026
Size 3, Length 32mm, 6/Pkg	29-12025
Size 4, Length 32mm, 6/Pkg	29-12024
Size 5, Length 32mm, 6/Pkg	29-12023
Size 6, Length 32mm, 6/Pkg	29-12022
Size 1-6, Length 32mm, 6/Pkg	29-12021

Best Quality Dental Supplies

Obturation Systems

E&Q Master Obturation Systems (Meta Biomed)



Maximizes convenience and minimizes chair time

- Combined filling method of vertical compaction and backfill makes the 3 dimensional cordless obturation system a perfect choice for apical part, accessory and lateral canals.

Full Unit	29-106421
Pen Unit	29-106422
Gun Unit	29-106423
Gutta Percha Bars, 100 Bars/Box	29-106424
Pen Tips, 2 each/pac, (F&FM)	29-106425
Pen Tip, each	29-106426
Gun Needle, 6 each/pack, (23G & 25G)	29-106427
Heat Insulator, 5/Pkg	29-106429
Plunger Head, 2 each/Pkg	29-106430

Cordless Thermal Condensor Unit (Dentazon)



- lightest cordless warm vertical condensation device
- effectively liquefies and condenses root canal filing materials
- adjustable 3 temperature settings
- color-coded pen tips available in various sizes $(XF, F, FM, \dot{M}, M\dot{L})$
- excellent apical control over lateral and complicated root canals

1 year warranty

Unit 29-11221

EndoApex Cordless Obturation System (Dentazon)

Light weight and highly efficient cordless endodontic dual obturation system which offers optimum functionality and reliable apical control along with excellent tactile sense •360° Easy Swivel (no tool required)

- •Complete thermal protection for patient
- •Ultra-fast heating for less chair time
- •Ergonomic and user-friendly oper-
- •Maximum portability with stylish design
- •Greatly reduces hand and wrist fatigue
- •Color coded pen tips
- •Delivers reliable and predictable clinical results 1 Year Warranty

Call us for Quarterly Promotion

2-in-1 Obturation System 29-11315



Prestige Dental Products, Inc. Prestigedentalproducts.com

Tel: 877-772-3888

Gutta Easy Cordless Thermal Obturation System (Dentazon)

Light weight and highly efficient cordless endodontic obturation system which offers optimum functionality and reliable apical control along with excellent tactile sense.

- •360° Easy Swivel (no tool required)
- •Complete thermal protection for patient
- •Ultra-fast heating for less chair time
- •Ergonomic and user-friendly operation
- •Maximum portability with stylish design
- •Greatly reduces hand and wrist fatigue
- •Color coded pen tips
- •Delivers reliable and predictable clinical results 1 Year Warranty

Unit 29-11538

Therma-Kor Obturators

Made in USA

Therma-Kor is an endodontic obturator comparable to Soft Core

Compare to Thermafil or Densfill Obturators (Dentsply)




#20, 25mm (6/Pk)	39-52412
#25, 25mm (6/Pk)	39-52413
#30, 25mm (6/Pk)	39-52414
#35, 25mm (6/Pk)	39-52415

Hand Files - Stainless Steel

Barbed Broaches, Stainless Steel

ISO color coded handle.



Color Black #6, 0.60, Length 21mm

Asst., Length 21mm

Size XXXX-F, Color Purple #0, 0.20, Length 21mm	29-11167
Size XXX-F, Color White #1, 0.25, Length 21mm	29-11168
Size XX-F, Color Yellow #2, 0.30, Length 21mm	29-11169
Size X-Fine, Color Red #3, 0.35, Length 21mm	29-11170
Size Fine, Color Blue #4, 0.40, Length 21mm	29-11171
Size Medium, Color Green #5, 0.50, Length 21mm	29-11173
Size Coarse	

29-11174

D Finders (Mani) D cr

D-Shaped cross-section has achieved blade past strength, unprecedent-

ed with conventional penetration files.

The D-Shaped cross section allows unprecedented initial access and navigation of calcified and curved canals. Incredibly stiff design & obviates over-cutting

Obviating excessive cutting and intrusion in the calcified root canals or in the curved root canals. Its ISO .02 standard taper allows for a smooth transition of other files to complete instrumentation. D Finders have 1.5 times fracture less resistance than other penetration files


Size 08, 21mm, 6/Pkg	29-109401
Size 08, 25mm, 6/Pkg	29-109402
Size 10, 21mm, 6/Pkg	29-109403
Size 10, 25mm, 6/Pkg	29-109404
Size 12, 21mm, 6/Pkg	29-109405
Size 12, 25mm, 6/Pkg	29-109406
Size 15, 21mm, 6/Pkg	29-109407
Size 15, 25mm, 6/Pkg	29-109408

Endo Micro Pulpa Box (Union Broach)



	Aluminum 30 Holes	
SILVER Color	29-108481	
RED Color	29-108482	,
BLUE Color	29-108483	

Endo Micro Pulpa Round (Miltex)



Finger Pluggers

Size 15, Length 21mm, 6/Pkg	29-11468
Size 20, Length 21mm, 6/Pkg	29-11479
Size 25, Length 21mm, 6/Pkg	29-11480



Size 30, Length 21mm, 6/Pkg	29-11481
Size 35, Length 21mm, 6/Pkg	29-11482
Size 40, Length 21mm, 6/Pkg	29-11483
Size 15-40, Length 21mm, 6/Pkg	29-11485
Size 15, Length 25mm, 6/Pkg	29-11484
Size 20, Length 25mm, 6/Pkg	29-11486
Size 25, Length 25mm, 6/Pkg	29-11478
Size 30, Length 25mm, 6/Pkg	29-11477
Size 35, Length 25mm, 6/Pkg	29-11476
Size 40, Length 25mm, 6/Pkg	29-11467
Size 15-40, Length 25mm, 6/Pkg	29-11469
Size 15, Length 31mm, 6/Pkg	29-11470
Size 20, Length 31mm, 6/Pkg	29-11471
Size 25, Length 31mm, 6/Pkg	29-11472
Size 30, Length 31mm, 6/Pkg	29-11473
Size 35, Length 31mm, 6/Pkg	29-11474
Size 40, Length 31mm, 6/Pkg	29-11475
Size 15-40, Length 31mm, 6/Pkg	29-11487

Finger Spreaders (Mani)

ISO color coded handle.

•Smooth conical shape wit pointed tip. •Designed for compacting of gutta percha and other filling materials	
Size 15, Length 21mm, 6/Pkg	29-11511
Size 20, Length 21mm, 6/Pkg	29-11492
Size 25, Length 21mm, 6/Pkg	29-11491
Size 30, Length 21mm, 6/Pkg	29-11490
Size 35, Length 21mm, 6/Pkg	29-11489
Size 40, Length 21mm, 6/Pkg	29-11488
Size 15-40, Length 21mm, 6/Pkg	29-11513
Size XF, Length 21mm, 6/Pkg	29-11512
Size FF, Length 21mm, 6/Pkg	29-11514
Size MF, Length 21mm, 6/Pkg	29-11515
Size F, Length 21mm, 6/Pkg	29-11516
Size FM, Length 21mm, 6/Pkg	29-11517
Size M, Length 21mm, 6/Pkg	29-11518
Size Asst., Length 21mm, 6/Pkg	29-11519
Size 15, Length 25mm, 6/Pkg	29-11522
Size 20, Length 25mm, 6/Pkg	29-11521

Size 25, Length 25mm, 6/Pkg	29-11493
Size 30, Length 25mm, 6/Pkg	29-11494
Size 35, Length 25mm, 6/Pkg	29-11495
Size 40, Length 25mm, 6/Pkg	29-11510
Size 15-40, Length 25mm, 6/Pkg	29-11509
Size XF, Length 25mm, 6/Pkg	29-11508
Size FF, Length 25mm, 6/Pkg	29-11507
Size MF, Length 25mm, 6/Pkg	29-11506
Size F, Length 25mm, 6/Pkg	29-11505
Size FM, Length 25mm, 6/Pkg	29-11504
Size M, Length 25mm, 6/Pkg	29-11503
Size Asst., Length 25mm, 6/Pkg	29-11502
Size 15, Length 31mm, 6/Pkg	29-11501
Size 20, Length 31mm, 6/Pkg	29-11500
Size 25, Length 31mm, 6/Pkg	29-11499
Size 30, Length 31mm, 6/Pkg	29-11498
Size 35, Length 31mm, 6/Pkg	29-11497
Size 40, Length 31mm, 6/Pkg	29-11496
Size 15-40, Length 31mm, 6/Pkg	29-11520

GPR Gutta Percha Remover (Mani)



- Superior results Fast, easy, & efficient.
- Fast & dynamic removal stress-free and easy to use, w/ torque control and reversal
- No need for solvents
- Max allowable speed 1000+/-500min 0.04 Taper

GPR #70 Size 1S, Green color, 16mm Stainless Steel, 4/Pkg	29-109441
GPR #50 Size 2S, Gold color, 18mm Stainless Steel, 4/Pkg	29-109442
GPR #40 Size 3N, Black color, 21mm Stainless Steel, 4/Pkg	29-109443
GPR #30 Size 4N, Blue color, 21mm Stainless Steel, 4/Pkg	29-109444
GPR Assorted Size #30-#70, 4/Pkg	29-109445

www.PrestigeDentalProducts.com



H-Files Stainless Steel with Endo Sto	n (ATW)
11-Files Staffiess Steel With Endo Sto	D(ATW)
 High Quality Stainless Steel H-Files with Endo Stop FDA and ISO approved Color Coded Handles 	of the state of th
	fi ²
Size 8, Length 21mm, 6/Pkg	29-11636
Size 10, Length 21mm, 6/Pkg	29-11635
Size 15, Length 21mm, 6/Pkg	29-11634
Size 20, Length 21mm, 6/Pkg	29-11633
Size 25, Length 21mm, 6/Pkg	29-11632
Size 30, Length 21mm, 6/Pkg	29-11637
Size 35, Length 21mm, 6/Pkg	29-11638
Size 40, Length 21mm, 6/Pkg	29-11621
Size 15-40, Length 21mm, 6/Pkg	29-11620
Size 45, Length 21mm, 6/Pkg	29-11631
Size 50, Length 21mm, 6/Pkg	29-11630
Size 55, Length 21mm, 6/Pkg	29-11617
Size 60, Length 21mm, 6/Pkg	29-11629
Size 70, Length 21mm, 6/Pkg	29-11639
Size 80, Length 21mm, 6/Pkg	29-11628
Size 45-80, Length 21mm, 6/Pkg	29-11643
Size 8, Length 25mm, 6/Pkg	29-11623
Size 10, Length 25mm, 6/Pkg	29-11622
Size 15, Length 25mm, 6/Pkg	29-11627
Size 20, Length 25mm, 6/Pkg	29-11624
Size 25, Length 25mm, 6/Pkg	29-11625
Size 30, Length 25mm, 6/Pkg	29-11626
Size 35, Length 25mm, 6/Pkg	29-11644
Size 40, Length 25mm, 6/Pkg	29-11645
Size 15-40, Length 25mm, 6/Pkg	29-11660
Size 45, Length 25mm, 6/Pkg	29-11646
Size 50, Length 25mm, 6/Pkg	29-11641
Size 55, Length 25mm, 6/Pkg	29-11640
Size 60, Length 25mm, 6/Pkg	29-11658
Size 70, Length 25mm, 6/Pkg	29-11659
Size 80, Length 25mm, 6/Pkg	29-11647
Size 45-80, Length 25mm, 6/Pkg	29-11648
Size 8, Length 31mm, 6/Pkg	29-11649
Size 10, Length 31mm, 6/Pkg	29-11661
Size 15, Length 31mm, 6/Pkg	29-11662

Size 20, Length 31mm, 6/Pkg	29-11663
Size 25, Length 31mm, 6/Pkg	29-11653
Size 30, Length 31mm, 6/Pkg	29-11652
Size 35, Length 31mm, 6/Pkg	29-11651
Size 40, Length 31mm, 6/Pkg	29-11650
Size 15-40, Length 31mm, 6/Pkg	29-11642
Size 45, Length 31mm, 6/Pkg	29-11664
Size 50, Length 31mm, 6/Pkg	29-11619
Size 55, Length 31mm, 6/Pkg	29-11618
Size 60, Length 31mm, 6/Pkg	29-11655
Size 70, Length 31mm, 6/Pkg	29-11654
Size 80, Length 31mm, 6/Pkg	29-11657
Size 45-80, Length 31mm, 6/Pkg	29-11656

H-Files Stainless Steel (JS Dental)

•ISO color coded plastic handle. •Fast cutting instrument with one continuous flute for fast removal of debris.



Size 08, Length 21mm, 6/Pkg	29-11712
Size 10, Length 21mm, 6/Pkg	29-11699
Size 15, Length 21mm, 6/Pkg	29-11698
Size 20, Length 21mm, 6/Pkg	29-11697
Size 25, Length 21mm, 6/Pkg	29-11696
Size 30, Length 21mm, 6/Pkg	29-11695
Size 35, Length 21mm, 6/Pkg	29-11694
Size 40, Length 21mm, 6/Pkg	29-11693
Size 15-40, Length 21mm, 6/Pkg	29-11710
Size 08, Length 25mm, 6/Pkg	29-11708
Size 10, Length 25mm, 6/Pkg	29-11707
Size 15, Length 25mm, 6/Pkg	29-11706
Size 20, Length 25mm, 6/Pkg	29-11705
Size 25, Length 25mm, 6/Pkg	29-11704
Size 30, Length 25mm, 6/Pkg	29-11703
Size 35, Length 25mm, 6/Pkg	29-11702
Size 40, Length 25mm, 6/Pkg	29-11689
Size 15-40, Length 25mm, 6/Pkg	29-11670





K-File 21mm (Mani)

	K-FILES	9
AAAA(I	SIZE Color	SIZE Color
MANI	※06(000)◎Ⅲ	35 (5)
	08 (00) 💷	40 (6)
MANI, INC. マニー株式会社	10 (0)	45 (7)
	15 (1)	50 (8)
K-FILES (25)	20 (2)	55 (9)
	25 (3)	60 (10)
Length (L): 18,21,25,28,31mm Assorted: 15(1)~40(6) Packed: 6 pcs. to a package 45(7)~80(12) 90(13)~140(18)	30 (4)	70 (11)

Sharp & mirror-like finish edges possible for excellent working characteristics with high quality.

Made in Japan

Size #06 - 21mm , 6/Pkg	29-105691
Size #08 - 21mm , 6/Pkg	29-105692
Size #10 - 21mm , 6/Pkg	29-105693
Size #15 - 21mm , 6/Pkg	29-105694
Size #20 - 21mm , 6/Pkg	29-105695
Size #25 - 21mm , 6/Pkg	29-105696
Size #30 - 21mm , 6/Pkg	29-105697
Size #35 - 21mm, 6/Pkg	29-105698
Size #40 - 21mm, 6/Pkg	29-105699

K-File 25mm (Mani)



Sharp & mirror-like finish edges possible for excellent working characteristics with high quality.

Made in Japan

Size #06 - 25mm , 6/Pkg	29-105701
Size #08 - 25mm , 6/Pkg	29-105702
Size #10 - 25mm , 6/Pkg	29-105703
Size #15 - 25mm , 6/Pkg	29-105704
Size #20 - 25mm , 6/Pkg	29-105705
Size #25 - 25mm , 6/Pkg	29-105706
Size #30 - 25mm , 6/Pkg	29-105707
Size #35 - 25mm , 6/Pkg	29-105708
Size #40 - 25mm , 6/Pkg	29-105709

K-Files Stainless Steel with Endo Stop

- •High quality stainless steel K-Files with Endo Stop
- •FDA and ISO approved
- •Color Coded Handles



Lildodolii	
Size 8, Length 21mm, 6/Pkg	29-11892
Size 10, Length 21mm, 6/Pkg	29-11878
Size 15, Length 21mm, 6/Pkg	29-11879
Size 20, Length 21mm, 6/Pkg	29-11876
Size 25, Length 21mm, 6/Pkg	29-11877
Size 30, Length 21mm, 6/Pkg	29-11875
Size 35, Length 21mm, 6/Pkg	29-11873
Size 40, Length 21mm, 6/Pkg	29-11874
Size 15-40, Length 21mm, 6/Pkg	29-11872
Size 45, Length 21mm, 6/Pkg	29-11870
Size 50, Length 21mm, 6/Pkg	29-11871
Size 55, Length 21mm, 6/Pkg	29-11880
Size 60, Length 21mm, 6/Pkg	29-11881
Size 70, Length 21mm, 6/Pkg	29-11862
Size 80, Length 21mm, 6/Pkg	29-11885
Size 45-80, Length 21mm, 6/Pkg	29-11891
Size 8, Length 25mm, 6/Pkg	29-11890
Size 10, Length 25mm, 6/Pkg	29-11889
Size 15, Length 25mm, 6/Pkg	29-11888
Size 20, Length 25mm, 6/Pkg	29-11887
Size 25, Length 25mm, 6/Pkg	29-11886
Size 30, Length 25mm, 6/Pkg	29-11884
Size 35, Length 25mm, 6/Pkg	29-11883
Size 40, Length 25mm, 6/Pkg	29-11882
Size 15-40, Length 25mm, 6/Pkg	29-11869
Size 45, Length 25mm, 6/Pkg	29-11868
Size 50, Length 25mm, 6/Pkg	29-11855
Size 55, Length 25mm, 6/Pkg	29-11854
Size 60, Length 25mm, 6/Pkg	29-11853
Size 70, Length 25mm, 6/Pkg	29-11852
Size 80, Length 25mm, 6/Pkg	29-11851
Size 45-80, Length 25mm, 6/Pkg	29-11850
Size 8, Length 31mm, 6/Pkg	29-11849
Size 10, Length 31mm, 6/Pkg	29-11848
Size 15, Length 31mm, 6/Pkg	29-11847
Size 20, Length 31mm, 6/Pkg	29-11846
Size 25, Length 31mm, 6/Pkg	29-11856
Size 30, Length 31mm, 6/Pkg	29-11857
Size 35, Length 31mm, 6/Pkg	29-11867
Size 40, Length 31mm, 6/Pkg	29-11866



Size 15-40, Length 31mm, 6/Pkg	29-11865
Size 45, Length 31mm, 6/Pkg	29-11864
Size 50, Length 31mm, 6/Pkg	29-11863
Size 55, Length 31mm, 6/Pkg	29-11861
Size 60, Length 31mm, 6/Pkg	29-11858
Size 70, Length 31mm, 6/Pkg	29-11859
Size 80, Length 31mm, 6/Pkg	29-11860
Size 45-80, Length 31mm, 6/Pkg	29-11845

K-Files Stainless Steel (JS Dental)

•ISO color coded plastic handle •Strong and stable instrument for cleaning and shaping the canal walls





walls	1101	
Size 06, Length 21mm,	6/Pkg	29-11929
Size 08, Length 21mm,	_	29-11928
Size 10, Length 21mm,	6/Pkg	29-11927
Size 15, Length 21mm,	6/Pkg	29-11937
Size 20, Length 21mm,	6/Pkg	29-11926
Size 25, Length 21mm,	6/Pkg	29-11925
Size 30, Length 21mm,	6/Pkg	29-11924
Size 35, Length 21mm,	6/Pkg	29-11923
Size 40, Length 21mm,	6/Pkg	29-11922
Size 45, Length 21mm,	6/Pkg	29-11921
Size 50, Length 21mm,	6/Pkg	29-11920
Size 55, Length 21mm,	6/Pkg	29-11930
Size 60, Length 21mm,	6/Pkg	29-11931
Size 70, Length 21mm,	6/Pkg	29-11932
Size 80, Length 21mm,	6/Pkg	29-11943
Size 15-40, Length 21m	m , 6/Pkg	29-11942
Size 45-80, Length 21m	m , 6/Pkg	29-11941
Size 06, Length 25mm,	6/Pkg	29-11940
Size 08, Length 25mm,	6/Pkg	29-11939
Size 10, Length 25mm,	6/Pkg	29-11938
Size 15, Length 25mm,	6/Pkg	29-11936
Size 20, Length 25mm,	6/Pkg	29-11935
Size 25, Length 25mm,	6/Pkg	29-11934
Size 30, Length 25mm,	6/Pkg	29-11933
Size 35, Length 25mm,	6/Pkg	29-11919
Size 40, Length 25mm,	6/Pkg	29-11918
Size 45, Length 25mm,	6/Pkg	29-11917

Size 50, Length 25mm, 6/Pkg	29-11903
Size 55, Length 25mm, 6/Pkg	29-11902
Size 60, Length 25mm, 6/Pkg	29-11901
Size 70, Length 25mm, 6/Pkg	29-11900
Size 80, Length 25mm, 6/Pkg	29-11899
Size 15-40, Length 25mm, 6/Pkg	29-11897
Size 45-80, Length 25mm, 6/Pkg	29-11898
Size 06, Length 31mm, 6/Pkg	29-11896
Size 08, Length 31mm, 6/Pkg	29-11895
Size 10, Length 31mm, 6/Pkg	29-11894
Size 15, Length 31mm, 6/Pkg	29-11904
Size 20, Length 31mm, 6/Pkg	29-11905
Size 25, Length 31mm, 6/Pkg	29-11906
Size 30, Length 31mm, 6/Pkg	29-11916
Size 35, Length 31mm, 6/Pkg	29-11915
Size 40, Length 31mm, 6/Pkg	29-11914
Size 45, Length 31mm, 6/Pkg	29-11913
Size 50, Length 31mm, 6/Pkg	29-11912
Size 55, Length 31mm, 6/Pkg	29-11911
Size 60, Length 31mm, 6/Pkg	29-11910
Size 70, Length 31mm, 6/Pkg	29-11909
Size 80, Length 31mm, 6/Pkg	29-11908
Size 15-40, Length 31mm, 6/Pkg	29-11907
Size 45-80, Length 31mm, 6/Pkg	29-11893

Mini Clean Gripp (Miltex)



 SILVER Color
 29-108491

 BLUE Color
 29-108492

Mini Endo Organizer Box (Miltex)



SILVER Color, 13 x 5 x 4cm, holds up to 60 Endodontic Finger instruments 29-108501

Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Mity Nickel-Titanium Compactors (JS Dental)

ISO color coded aluminum handle.

Instruments used to compact gutta-percha and other filling materials.

Size 25, Length 21mm, 6/Pkg	29-11992
Size 35, Length 21mm, 6/Pkg	29-11991
Size 45, Length 21mm, 6/Pkg	29-11990
Size 25, Length 25mm, 6/Pkg	29-11989
Size 35, Length 25mm, 6/Pkg	29-11988
Size 45, Length 25mm, 6/Pkg	29-11987

Mity Nickel-Titanium Finger Pluggers (JS Dental)

•ISO color coded plastic handle.

•Smooth conical instrument with flat tip.

•Designed for vertical and lateral compacting of gutta percha and other filling materials

mining materials	
Size 15, Length 21mm, 6/Pkg	29-12141
Size 20, Length 21mm, 6/Pkg	29-12129
Size 25, Length 21mm, 6/Pkg	29-12128
Size 30, Length 21mm, 6/Pkg	29-12127
Size 35, Length 21mm, 6/Pkg	29-12126
Size 40, Length 21mm, 6/Pkg	29-12125
Size 15-40, Length 21mm, 6/Pkg	29-12124
Size 15, Length 25mm, 6/Pkg	29-12123
Size 20, Length 25mm, 6/Pkg	29-12122
Size 25, Length 25mm, 6/Pkg	29-12130
Size 30, Length 25mm, 6/Pkg	29-12131
Size 35, Length 25mm, 6/Pkg	29-12132
Size 40, Length 25mm, 6/Pkg	29-12139
Size 15-40, Length 25mm, 6/Pkg	29-12140

Size 25, Length 31mm, 6/Pkg 29-12136 **Size 30, Length 31mm**, 6/Pkg 29-12135

 Size 30, Length 31mm, 6/Pkg
 29-12135

 Size 35, Length 31mm, 6/Pkg
 29-12134

Size 15, Length 31mm, 6/Pkg

Size 20, Length 31mm, 6/Pkg

Size 40, Length 31mm, 6/Pkg 29-12133 **Size 15-40, Length 31mm**, 6/Pkg 29-12121

Path Finders K-Type-Stainless Steel (JS Dental)



 Size 12, Length 21mm, 6/Pkg
 29-12261

 Size 12, Length 25mm, 6/Pkg
 29-12262

29-12138

Reamers Files-Stainless Steel (JS De	ental)
•ISO color coded plastic	
handle. •Designed for smooth preparation of canal walls	
Size 08, Length 21mm, 6/Pkg	29-12354
Size 10, Length 21mm, 6/Pkg	29-12353
Size 15, Length 21mm, 6/Pkg	29-12337
Size 20, Length 21mm, 6/Pkg	29-12355
Size 25, Length 21mm, 6/Pkg	29-12356
Size 30, Length 21mm, 6/Pkg	29-12357
Size 35, Length 21mm, 6/Pkg	29-12358
Size 40, Length 21mm, 6/Pkg	29-12359
Size 45, Length 21mm, 6/Pkg	29-12360
Size 50, Length 21mm, 6/Pkg	29-12361
Size 55, Length 21mm, 6/Pkg	29-12362
Size 60, Length 21mm, 6/Pkg	29-12352
Size 70, Length 21mm, 6/Pkg	29-12351
Size 80, Length 21mm, 6/Pkg	29-12344
Size 15-40, Length 21mm, 6/Pkg	29-12343
Size 45-80, Length 21mm, 6/Pkg	29-12342
Size 08, Length 25mm, 6/Pkg	29-12341
Size 10, Length 25mm, 6/Pkg	29-12347
Size 15, Length 25mm, 6/Pkg	29-12345
Size 20, Length 25mm, 6/Pkg	29-12346
Size 25, Length 25mm, 6/Pkg	29-12348
Size 30, Length 25mm, 6/Pkg	29-12363
Size 35, Length 25mm, 6/Pkg	29-12350
Size 40, Length 25mm, 6/Pkg	29-12349
Size 45, Length 25mm, 6/Pkg	29-12364
Size 50, Length 25mm, 6/Pkg	29-12340
Size 55, Length 25mm, 6/Pkg	29-12339
Size 60, Length 25mm, 6/Pkg	29-12338
Size 70, Length 25mm, 6/Pkg	29-12336
Size 80, Length 25mm, 6/Pkg	29-12335
Size 15-40, Length 25mm, 6/Pkg	29-12334
Size 45-80, Length 25mm, 6/Pkg	29-12333
Size 08, Length 31mm, 6/Pkg	29-12332
Size 10, Length 31mm, 6/Pkg	29-12331
Size 15, Length 31mm, 6/Pkg	29-12376
Size 20, Length 31mm, 6/Pkg	29-12330



Size 25, Length 31mm, 6/Pkg	29-12375	
Size 30, Length 31mm, 6/Pkg	29-12366	
Size 35, Length 31mm, 6/Pkg	29-12365	
Size 40, Length 31mm, 6/Pkg	29-12367	
Size 45, Length 31mm, 6/Pkg	29-12368	
Size 50, Length 31mm, 6/Pkg	29-12369	
Size 55, Length 31mm, 6/Pkg	29-12371	
Size 60, Length 31mm, 6/Pkg	29-12370	
Size 70, Length 31mm, 6/Pkg	29-12373	
Size 80, Length 31mm, 6/Pkg	29-12372	
Size 15-40, Length 31mm, 6/Pkg	29-12374	
Size 45-80, Length 31mm, 6/Pkg	29-12329	

Reamers PLUS Files-Stainless Steel (JS Dental)

- •ISO color coded plastic handle.
- •Rubber stop and mm markings
 •Designed for smooth preparation of canal walls



Size 08, Length 21mm, 6/Pkg	29-12298
Size 10, Length 21mm, 6/Pkg	29-12311
Size 15, Length 21mm, 6/Pkg	29-12312
Size 20, Length 21mm, 6/Pkg	29-12313
Size 25, Length 21mm, 6/Pkg	29-12314
Size 30, Length 21mm, 6/Pkg	29-12315
Size 35, Length 21mm, 6/Pkg	29-12316
Size 40, Length 21mm, 6/Pkg	29-12317
Size 45, Length 21mm, 6/Pkg	29-12318
Size 50, Length 21mm, 6/Pkg	29-12319
Size 55, Length 21mm, 6/Pkg	29-12320
Size 60, Length 21mm, 6/Pkg	29-12310
Size 70, Length 21mm, 6/Pkg	29-12309
Size 80, Length 21mm, 6/Pkg	29-12299
Size 15-40, Length 21mm, 6/Pkg	29-12300
Size 45-80, Length 21mm, 6/Pkg	29-12301
Size 08, Length 25mm, 6/Pkg	29-12302
Size 10, Length 25mm, 6/Pkg	29-12303
Size 15, Length 25mm, 6/Pkg	29-12304
Size 20, Length 25mm, 6/Pkg	29-12305
Size 25, Length 25mm, 6/Pkg	29-12306
Size 30, Length 25mm, 6/Pkg	29-12307
Size 35, Length 25mm, 6/Pkg	29-12308

l		
	Size 40, Length 25mm, 6/Pkg	29-12321
	Size 45, Length 25mm, 6/Pkg	29-12322
	Size 50, Length 25mm, 6/Pkg	29-12291
	Size 55, Length 25mm, 6/Pkg	29-12290
	Size 60, Length 25mm, 6/Pkg	29-12289
	Size 70, Length 25mm, 6/Pkg	29-12288
	Size 80, Length 25mm, 6/Pkg	29-12287
	Size 15-40, Length 25mm, 6/Pkg	29-12285
	Size 45-80, Length 25mm, 6/Pkg	29-12286
	Size 08, Length 31mm, 6/Pkg	29-12284
	Size 10, Length 31mm, 6/Pkg	29-12283
	Size 15, Length 31mm, 6/Pkg	29-12282
	Size 20, Length 31mm, 6/Pkg	29-12292
	Size 25, Length 31mm, 6/Pkg	29-12293
	Size 30, Length 31mm, 6/Pkg	29-12323
	Size 35, Length 31mm, 6/Pkg	29-12324
	Size 40, Length 31mm, 6/Pkg	29-12325
	Size 45, Length 31mm, 6/Pkg	29-12326
	Size 50, Length 31mm, 6/Pkg	29-12328
	Size 55, Length 31mm, 6/Pkg	29-12327
	Size 60, Length 31mm, 6/Pkg	29-12297
	Size 70, Length 31mm, 6/Pkg	29-12296
	Size 80, Length 31mm, 6/Pkg	29-12295
	Size 15-40, Length 31mm, 6/Pkg	29-12294
	Size 45-80, Length 31mm, 6/Pkg	29-12281

S-Files with mm markings -Stainless Steel (JS Dental)

•ISO color coded plastic handle an d mm marks •Perfect blend of flexibility, durability and cutting ability



Size 08, Length 21mm, 6/Pkg	29-12456
Size 10, Length 21mm, 6/Pkg	29-12443
Size 15, Length 21mm, 6/Pkg	29-12442
Size 20, Length 21mm, 6/Pkg	29-12441
Size 25, Length 21mm, 6/Pkg	29-12440
Size 30, Length 21mm, 6/Pkg	29-12439
Size 35, Length 21mm, 6/Pkg	29-12438
Size 40, Length 21mm, 6/Pkg	29-12437
Size 45, Length 21mm, 6/Pkg	29-12436
Size 50, Length 21mm, 6/Pkg	29-12435



Size 55, Length 21mm, 6/Pkg	29-12444
Size 60, Length 21mm, 6/Pkg	29-12445
Size 70, Length 21mm, 6/Pkg	29-12446
Size 80, Length 21mm, 6/Pkg	29-12455
Size 15-40, Length 21mm, 6/Pkg	29-12453
Size 45-80, Length 21mm, 6/Pkg	29-12454
Size 08, Length 25mm, 6/Pkg	29-12452
Size 10, Length 25mm, 6/Pkg	29-12451
Size 15, Length 25mm, 6/Pkg	29-12450
Size 20, Length 25mm, 6/Pkg	29-12449
Size 25, Length 25mm, 6/Pkg	29-12448
Size 30, Length 25mm, 6/Pkg	29-12447
Size 35, Length 25mm, 6/Pkg	29-12434
Size 40, Length 25mm, 6/Pkg	29-12433
Size 45, Length 25mm, 6/Pkg	29-12432
Size 50, Length 25mm, 6/Pkg	29-12419
Size 55, Length 25mm, 6/Pkg	29-12418
Size 60, Length 25mm, 6/Pkg	29-12417
Size 70, Length 25mm, 6/Pkg	29-12416
Size 80, Length 25mm, 6/Pkg	29-12415
Size 15-40, Length 25mm, 6/Pkg	29-12414
Size 45-80, Length 25mm, 6/Pkg	29-12413
Size 10, Length 28mm, 6/Pkg	29-12412
Size 15, Length 28mm, 6/Pkg	29-12411
Size 20, Length 28mm, 6/Pkg	29-12420
Size 25, Length 28mm, 6/Pkg	29-12421
Size 30, Length 28mm, 6/Pkg	29-12422
Size 35, Length 28mm, 6/Pkg	29-12431
Size 40, Length 28mm, 6/Pkg	29-12430
Size 45, Length 28mm, 6/Pkg	29-12429
Size 50, Length 28mm, 6/Pkg	29-12428
Size 55, Length 28mm, 6/Pkg	29-12427
Size 60, Length 28mm, 6/Pkg	29-12426
Size 70, Length 28mm, 6/Pkg	29-12425
Size 80, Length 28mm, 6/Pkg	29-12424
Size 15-40, Length 28mm, 6/Pkg	29-12410
Size 45-80, Length 28mm, 6/Pkg	29-12423



Endodontics	
S-Finder-Stainless Steel (JS Dental	<u>)</u>
Path finding instrument in S-File configuration	The state of the s
and the same of th	
Size 08, Length 18mm, 6/Pkg	29-12457
Size 08, Length 21mm, 6/Pkg	29-12458
Turbo Files-Stainless Steel (JS Den •ISO color coded plastic handle •Rubber stops and mm markings combination of a K-File and H-File with two •A continuous flutes makes this in- strument effective and stable for fast removal of debris	tal)
Size 08, Length 21mm, 6/Pkg	29-12512
Size 10, Length 21mm, 6/Pkg	29-12520
Size 15, Length 21mm, 6/Pkg	29-12519
Size 20, Length 21mm, 6/Pkg	29-12518

29-12517

29-12516

29-12515

29-12514

29-12513

29-12511

29-12510

29-12521

29-12522

29-12532

29-12531

29-12530

29-12529

29-12528

29-12527

29-12526

29-12525

29-12524

29-12523

29-12509

29-12508

29-12496

Size 25, Length 21mm, 6/Pkg

Size 30, Length 21mm, 6/Pkg

Size 35, Length 21mm, 6/Pkg

Size 40, Length 21mm, 6/Pkg

Size 45, Length 21mm, 6/Pkg

Size 50, Length 21mm, 6/Pkg

Size 55, Length 21mm, 6/Pkg

Size 60, Length 21mm, 6/Pkg

Size 70, Length 21mm, 6/Pkg

Size 80, Length 21mm, 6/Pkg

Size 15-40, Length 21mm, 6/Pkg

Size 45-80, Length 21mm, 6/Pkg

Size 08, Length 25mm, 6/Pkg

Size 10, Length 25mm, 6/Pkg

Size 15, Length 25mm, 6/Pkg

Size 20, Length 25mm, 6/Pkg

Size 25, Length 25mm, 6/Pkg

Size 30, Length 25mm, 6/Pkg

Size 35, Length 25mm, 6/Pkg

Size 40, Length 25mm, 6/Pkg

Size 45, Length 25mm, 6/Pkg

Size 50, Length 25mm, 6/Pkg



Tel: 8//-//2-3888	
Size 55, Length 25mm, 6/Pkg	29-12497
Size 60, Length 25mm, 6/Pkg	29-12498
Size 70, Length 25mm, 6/Pkg	29-12499
Size 80, Length 25mm, 6/Pkg	29-12491
Size 15-40, Length 25mm, 6/Pkg	29-12489
Size 45-80, Length 25mm, 6/Pkg	29-12490
Size 08, Length 31mm, 6/Pkg	29-12488
Size 10, Length 31mm, 6/Pkg	29-12487
Size 15, Length 31mm, 6/Pkg	29-12486
Size 20, Length 31mm, 6/Pkg	29-12495
Size 25, Length 31mm, 6/Pkg	29-12494
Size 30, Length 31mm, 6/Pkg	29-12507
Size 35, Length 31mm, 6/Pkg	29-12506
Size 40, Length 31mm, 6/Pkg	29-12505
Size 45, Length 31mm, 6/Pkg	29-12504
Size 50, Length 31mm, 6/Pkg	29-12503
Size 55, Length 31mm, 6/Pkg	29-12502
Size 60, Length 31mm, 6/Pkg	29-12501
Size 70, Length 31mm, 6/Pkg	29-12500
Size 80, Length 31mm, 6/Pkg	29-12492
Size 15-40, Length 31mm, 6/Pkg	29-12493
Size 45-80, Length 31mm, 6/Pkg	29-12485
Xtreme Turbo Files No Memory (JS	Dental)
•ISO color coded plastic handle	100

Size 15-40, Length 31mm, 6/Pkg	29-12493
Size 45-80, Length 31mm, 6/Pkg	29-12485
Xtreme Turbo Files No Memory	(JS Dental)
•ISO color coded plastic handle •MM markings and rubber stops •A combination instrument for fast removal of debris	
Size 15, Length 21mm, 6/Pkg	29-12591
Size 20, Length 21mm, 6/Pkg	29-12579
Size 25, Length 21mm, 6/Pkg	29-12578
Size 30, Length 21mm, 6/Pkg	29-12577
Size 35, Length 21mm, 6/Pkg	29-12576
Size 40, Length 21mm, 6/Pkg	29-12575
Size 15-40, Length 21mm, 6/Pkg	29-12574
Size 15, Length 25mm, 6/Pkg	29-12573
Size 20, Length 25mm, 6/Pkg	29-12580
Size 25, Length 25mm, 6/Pkg	29-12572
Size 30, Length 25mm, 6/Pkg	29-12581
Size 35, Length 25mm, 6/Pkg	29-12582
Size 40, Length 25mm, 6/Pkg	29-12590

Size 15-40, Length 25mm, 6/Pkg	29-12589
Size 15, Length 31mm, 6/Pkg	29-12587
Size 20, Length 31mm, 6/Pkg	29-12588
Size 25, Length 31mm, 6/Pkg	29-12586
Size 30, Length 31mm, 6/Pkg	29-12585
Size 35, Length 31mm, 6/Pkg	29-12584
Size 40, Length 31mm, 6/Pkg	29-12583
Size 15-40, Length 31mm, 6/Pkg	29-12571

Hand Files-Nickel Titanium

15-40

Eco-NiTi K-File


JS Dental Mtg Inc PO Box 904 Ridgefield, CT 06877 (203) 438-8832

29-11260

29-11271

Eco-NiTi H-File (JS Dental)





•Rubber Stops K-Files and H-•Files of highest

ı		
ı	Size 15, Length 21mm, 6/Pkg	29-11252
ı	Size 20, Length 21mm, 6/Pkg	29-11264
ı	Size 25, Length 21mm, 6/Pkg	29-11265
ı	Size 30, Length 21mm, 6/Pkg	29-11266
ı	Size 35, Length 21mm, 6/Pkg	29-11267
ı	Size 40, Length 21mm, 6/Pkg	29-11268
ı	Size 15-40, Length 21mm, 6/Pkg	29-11269
ı	Size 15, Length 25mm, 6/Pkg	29-11270
ı	Size 20, Length 25mm, 6/Pkg	29-11272
	Size 25, Length 25mm, 6/Pkg	29-11263
ı	Size 30, Length 25mm, 6/Pkg	29-11262
ı	Size 35, Length 25mm, 6/Pkg	29-11261
ı	Size 40, Length 25mm, 6/Pkg	29-11254
ı	Size 15-40, Length 25mm, 6/Pkg	29-11253
ı	Size 15, Length 28mm, 6/Pkg	29-11255
ı	Size 20, Length 28mm, 6/Pkg	29-11256
	Size 25, Length 28mm, 6/Pkg	29-11257
	Size 30, Length 28mm, 6/Pkg	29-11258
	Size 35, Length 28mm, 6/Pkg	29-11259
I		



Eco-NiTi K-File (JS Dental)

•ISO color coded plastic handle.

•Rubber Stop K-Files and H

-Files of highest quality nickel titanium at a very low price



Size 15, Length 21mm, 6/Pkg 29-11273	
Size 20, Length 21mm, 6/Pkg 29-11293	
Size 25, Length 21mm, 6/Pkg 29-11292	
Size 30, Length 21mm, 6/Pkg 29-11291	
Size 35, Length 21mm, 6/Pkg 29-11290	
Size 40, Length 21mm, 6/Pkg 29-11289	
Size 15-40, Length 21mm, 6/Pkg 29-11288	
Size 15, Length 25mm, 6/Pkg 29-11287	
Size 20, Length 25mm, 6/Pkg 29-11286	
Size 25, Length 25mm, 6/Pkg 29-11284	
Size 30, Length 25mm, 6/Pkg 29-11283	
Size 35, Length 25mm, 6/Pkg 29-11282	
Size 40, Length 25mm, 6/Pkg 29-11274	
Size 15-40, Length 25mm, 6/Pkg 29-11275	
Size 15, Length 28mm, 6/Pkg 29-11276	
Size 20, Length 28mm, 6/Pkg 29-11277	
Size 25, Length 28mm, 6/Pkg 29-11278	
Size 30, Length 28mm, 6/Pkg 29-11279	
Size 35, Length 28mm, 6/Pkg 29-11280	
Size 40, Length 28mm, 6/Pkg 29-11281	
Size 15-40, Length 28mm, 6/Pkg 29-11285	

Endo Micro Pulpa Box (Union Broach)



Aluminum 30 Holes

SILVER Color	29-108481
RED Color	29-108482
BLUE Color	29-108483

Size 40, Length 28mm, 6/Pkg

Size 15-40, Length 28mm, 6/Pkg

Endo Micro Pulpa Round (Miltex)

Small aluminum box holds up to 18 endodontic instruments **SILVER** Color 29-108481 **RED** Color 29-108482 **BLUE** Color 29-108483

Flex NiTi K-Files with Endo Stop (ATW)

•FDA and ISO approved •Color Coded Handles •High Quality Nickel

Titanium K-Files with Endo Stop



Size 8, Length 21mm, 6/Pkg	29-11752
Size 10, Length 21mm, 6/Pkg	29-11765
Size 15, Length 21mm, 6/Pkg	29-11766
Size 20, Length 21mm, 6/Pkg	29-11767
Size 25, Length 21mm, 6/Pkg	29-11768
Size 30, Length 21mm, 6/Pkg	29-11769
Size 35, Length 21mm, 6/Pkg	29-11770
Size 40, Length 21mm, 6/Pkg	29-11771
Size 15-40, Length 21mm, 6/Pkg	29-11772
Size 45, Length 21mm, 6/Pkg	29-11773
Size 50, Length 21mm, 6/Pkg	29-11774
Size 55, Length 21mm, 6/Pkg	29-11764
Size 60, Length 21mm, 6/Pkg	29-11763
Size 70, Length 21mm, 6/Pkg	29-11753
Size 80, Length 21mm, 6/Pkg	29-11754
Size 45-80, Length 21mm, 6/Pkg	29-11755
Size 8, Length 25mm, 6/Pkg	29-11756
Size 10, Length 25mm, 6/Pkg	29-11757
Size 15, Length 25mm, 6/Pkg	29-11758
Size 20, Length 25mm, 6/Pkg	29-11761
Size 25, Length 25mm, 6/Pkg	29-11760
Size 30, Length 25mm, 6/Pkg	29-11759
Size 35, Length 25mm, 6/Pkg	29-11762
Size 40, Length 25mm, 6/Pkg	29-11775
Size 15-40, Length 25mm, 6/Pkg	29-11776
Size 45, Length 25mm, 6/Pkg	29-11789
Size 50, Length 25mm, 6/Pkg	29-11790
Size 55, Length 25mm, 6/Pkg	29-11791



C18, Inc.	
Size 60, Length 25mm, 6/Pkg	29-11792
Size 70, Length 25mm, 6/Pkg	29-11793
Size 80, Length 25mm, 6/Pkg	29-11794
Size 45-80, Length 25mm, 6/Pkg	29-11795
Size 8, Length 31mm, 6/Pkg	29-11796
Size 10, Length 31mm, 6/Pkg	29-11751
Size 15, Length 31mm, 6/Pkg	29-11750
Size 20, Length 31mm, 6/Pkg	29-11788
Size 25, Length 31mm, 6/Pkg	29-11787
Size 30, Length 31mm, 6/Pkg	29-11777
Size 35, Length 31mm, 6/Pkg	29-11778
Size 40, Length 31mm, 6/Pkg	29-11780
Size 15-40, Length 31mm, 6/Pkg	29-11779
Size 45, Length 31mm, 6/Pkg	29-11783
Size 50, Length 31mm, 6/Pkg	29-11782
Size 55, Length 31mm, 6/Pkg	29-11781
Size 60, Length 31mm, 6/Pkg	29-11785
Size 70, Length 31mm, 6/Pkg	29-11784
Size 80, Length 31mm, 6/Pkg	29-11749
Size 45-80, Length 31mm, 6/Pkg	29-11786

Hyflex X-File 25mm (Coltene)



- A double-fluted instrument design for enhanced cutting action.
- Made of Nickel titanium for exceptional flexibility

#15-25mm, 6/Pkg 29-109941

Mini Clean Gripp (Miltex)

SILVER Color 29-108491 **BLUE** Color 29-108492





SILVER Color, 13 x 5 x 4cm, holds up to 60 Endodontic Finger instruments

Mity Nickel-Titanium H-Files (JS Dental)

- •ISO color coded plastic handle
- •MM markings and rubber stops.
- •Fast cutting instrument with one continuous flute for fast removal of debris



removal of debris	All Property and the second
Size 10, Length 21mm, 6/Pkg	29-12028
Size 15, Length 21mm, 6/Pkg	29-12040
Size 20, Length 21mm, 6/Pkg	29-12041
Size 25, Length 21mm, 6/Pkg	29-12042
Size 30, Length 21mm, 6/Pkg	29-12043
Size 35, Length 21mm, 6/Pkg	29-12044
Size 40, Length 21mm, 6/Pkg	29-12045
Size 45, Length 21mm, 6/Pkg	29-12046
Size 50, Length 21mm, 6/Pkg	29-12047
Size 55, Length 21mm, 6/Pkg	29-12048
Size 60, Length 21mm, 6/Pkg	29-12039
Size 70, Length 21mm, 6/Pkg	29-12038
Size 80, Length 21mm, 6/Pkg	29-12029
Size 15-40, Length 21mm, 6/P	kg 29-12030
Size 45-80, Length 21mm, 6/Pl	kg 29-12031
Size 10, Length 25mm, 6/Pkg	29-12032
Size 15, Length 25mm, 6/Pkg	29-12033
Size 20, Length 25mm, 6/Pkg	29-12034
Size 25, Length 25mm, 6/Pkg	29-12035
Size 30, Length 25mm, 6/Pkg	29-12036
Size 35, Length 25mm, 6/Pkg	29-12037
Size 40, Length 25mm, 6/Pkg	29-12049
Size 45, Length 25mm, 6/Pkg	29-12050
Size 50, Length 25mm, 6/Pkg	29-12051
Size 60, Length 25mm, 6/Pkg	29-12063
Size 55, Length 25mm, 6/Pkg	29-12064
Size 70, Length 25mm, 6/Pkg	29-12065
Size 80, Length 25mm, 6/Pkg	29-12066
Size 15-40, Length 25mm, 6/P	kg 29-12067
Size 45-80, Length 25mm, 6/P	kg 29-12068
Size 10, Length 31mm, 6/Pkg	29-12069
Size 15, Length 31mm, 6/Pkg	29-12070
Size 20, Length 31mm, 6/Pkg	29-12071
Size 25, Length 31mm, 6/Pkg	29-12062



Size 30, Length 31mm, 6/Pkg	29-12061
Size 35, Length 31mm, 6/Pkg	29-12052
Size 40, Length 31mm, 6/Pkg	29-12053
Size 45, Length 31mm, 6/Pkg	29-12054
Size 50, Length 31mm, 6/Pkg	29-12055
Size 55, Length 31mm, 6/Pkg	29-12056
Size 60, Length 31mm, 6/Pkg	29-12057
Size 70, Length 31mm, 6/Pkg	29-12058
Size 80, Length 31mm, 6/Pkg	29-12059
Size 15-40, Length 31mm, 6/Pkg	29-12060
Size 45-80, Length 31mm, 6/Pkg	29-12072

Mity Nickel-Titanium K-Files (JS Dental)

- •ISO color coded plastic handle.
- •MM markings and rubber stops.
- •Strong and stable instrument for cleaning and shaping the canal walls



Size 08, Length 21mm, 6/Pkg	29-12120
Size 10, Length 21mm, 6/Pkg	29-12107
Size 15, Length 21mm, 6/Pkg	29-12106
Size 20, Length 21mm, 6/Pkg	29-12105
Size 25, Length 21mm, 6/Pkg	29-12104
Size 30, Length 21mm, 6/Pkg	29-12103
Size 35, Length 21mm, 6/Pkg	29-12102
Size 40, Length 21mm, 6/Pkg	29-12101
Size 45, Length 21mm, 6/Pkg	29-12100
Size 50, Length 21mm, 6/Pkg	29-12099
Size 55, Length 21mm, 6/Pkg	29-12098
Size 60, Length 21mm, 6/Pkg	29-12108
Size 70, Length 21mm, 6/Pkg	29-12109
Size 80, Length 21mm, 6/Pkg	29-12119
Size 15-40, Length 21mm, 6/Pkg	29-12118
Size 45-80, Length 21mm, 6/Pkg	29-12117
Size 08, Length 25mm, 6/Pkg	29-12116
Size 10, Length 25mm, 6/Pkg	29-12115
Size 15, Length 25mm, 6/Pkg	29-12114
Size 20, Length 25mm, 6/Pkg	29-12113
Size 25, Length 25mm, 6/Pkg	29-12112
Size 30, Length 25mm, 6/Pkg	29-12111
Size 35, Length 25mm, 6/Pkg	29-12110
Size 40, Length 25mm, 6/Pkg	29-12097

Size 45, Length 25mm, 6/Pkg	29-12096
Size 50, Length 25mm, 6/Pkg	29-12083
Size 55, Length 25mm, 6/Pkg	29-12082
Size 60, Length 25mm, 6/Pkg	29-12081
Size 70, Length 25mm, 6/Pkg	29-12080
Size 80, Length 25mm, 6/Pkg	29-12079
Size 15-40, Length 25mm, 6/Pkg	29-12078
Size 45-80, Length 25mm, 6/Pkg	29-12077
Size 10, Length 31mm, 6/Pkg	29-12076
Size 08, Length 31mm, 6/Pkg	29-12075
Size 15, Length 31mm, 6/Pkg	29-12074
Size 20, Length 31mm, 6/Pkg	29-12084
Size 25, Length 31mm, 6/Pkg	29-12085
Size 30, Length 31mm, 6/Pkg	29-12095
Size 35, Length 31mm, 6/Pkg	29-12094
Size 40, Length 31mm, 6/Pkg	29-12093
Size 45, Length 31mm, 6/Pkg	29-12092
Size 50, Length 31mm, 6/Pkg	29-12091
Size 55, Length 31mm, 6/Pkg	29-12090
Size 60, Length 31mm, 6/Pkg	29-12089
Size 70, Length 31mm, 6/Pkg	29-12088
Size 80, Length 31mm, 6/Pkg	29-12087
Size 15-40, Length 31mm, 6/Pkg	29-12086
Size 45-80, Length 31mm, 6/Pkg	29-12073

Mity Nickel-Titanium Spreaders (JS Dental)

•ISO color coded plastic handle. •Smooth conical instrument with pointed tip.

•Designed for compacting of gutta percha and other filling materials

Size 15, Length 21mm, 6/Pkg	29-12170
Size 20, Length 21mm, 6/Pkg	29-12182
Size 25, Length 21mm, 6/Pkg	29-12183
Size 30, Length 21mm, 6/Pkg	29-12184
Size 35, Length 21mm, 6/Pkg	29-12185
Size 40, Length 21mm, 6/Pkg	29-12186
Size 15-40, Length 21mm, 6/Pkg	29-12187
Size 15, Length 25mm, 6/Pkg	29-12188
Size 20, Length 25mm, 6/Pkg	29-12190
Size 25, Length 25mm, 6/Pkg	29-12181



Size 30, Length 25mm, 6/Pkg	29-12180
Size 35, Length 25mm, 6/Pkg	29-12179
Size 40, Length 25mm, 6/Pkg	29-12171
Size 15-40, Length 25mm, 6/Pkg	29-12172
Size 15, Length 31mm, 6/Pkg	29-12173
Size 20, Length 31mm, 6/Pkg	29-12174
Size 25, Length 31mm, 6/Pkg	29-12175
Size 30, Length 31mm, 6/Pkg	29-12176
Size 35, Length 31mm, 6/Pkg	29-12177
Size 40, Length 31mm, 6/Pkg	29-12178
Size 15-40, Length 31mm, 6/Pkg	29-12189

Mity Nickel-Titanium Turbo Files (JS Dental)

- •ISO color coded plastic handle •MM markings and rubber stops
- •A combination of a K-File and H-File with two continuous flutes makes this instrument effective and



stable for fast removal of debris	
Size 08, Length 21mm, 6/Pkg	29-12205
Size 10, Length 21mm, 6/Pkg	29-12191
Size 15, Length 21mm, 6/Pkg	29-12204
Size 20, Length 21mm, 6/Pkg	29-12206
Size 25, Length 21mm, 6/Pkg	29-12207
Size 30, Length 21mm, 6/Pkg	29-12208
Size 35, Length 21mm, 6/Pkg	29-12209
Size 40, Length 21mm, 6/Pkg	29-12210

Size 30, Length 21 mm, 6/Pkg	29-12208
Size 35, Length 21mm, 6/Pkg	29-12209
Size 40, Length 21mm, 6/Pkg	29-12210
Size 45, Length 21mm, 6/Pkg	29-12211
Size 50, Length 21mm, 6/Pkg	29-12212
Size 55, Length 21mm, 6/Pkg	29-12213
Size 60, Length 21mm, 6/Pkg	29-12203
Size 70, Length 21mm, 6/Pkg	29-12202
Size 80, Length 21mm, 6/Pkg	29-12192
Size 15-40, Length 21mm, 6/Pkg	29-12193
Size 45-80, Length 21mm, 6/Pkg	29-12197
Size 08, Length 25mm, 6/Pkg	29-12196
Size 10, Length 25mm, 6/Pkg	29-12195
Size 15, Length 25mm, 6/Pkg	29-12194
Size 20, Length 25mm, 6/Pkg	29-12198
Size 25, Length 25mm, 6/Pkg	29-12199
Size 30, Length 25mm, 6/Pkg	29-12200

Size 35, Length 25mm, 6/Pkg

Size 40, Length 25mm, 6/Pkg	29-12214
Size 45, Length 25mm, 6/Pkg	29-12215
Size 50, Length 25mm, 6/Pkg	29-12228
Size 55, Length 25mm, 6/Pkg	29-12229
Size 60, Length 25mm, 6/Pkg	29-12230
Size 70, Length 25mm, 6/Pkg	29-12231
Size 80, Length 25mm, 6/Pkg	29-12233
Size 15-40, Length 25mm, 6/Pkg	29-12232
Size 45-80, Length 25mm, 6/Pkg	29-12234
Size 08, Length 31mm, 6/Pkg	29-12235
Size 10, Length 31mm, 6/Pkg	29-12236
Size 15, Length 31mm, 6/Pkg	29-12237
Size 20, Length 31mm, 6/Pkg	29-12227
Size 25, Length 31mm, 6/Pkg	29-12226
Size 30, Length 31mm, 6/Pkg	29-12216
Size 35, Length 31mm, 6/Pkg	29-12217
Size 40, Length 31mm, 6/Pkg	29-12218
Size 45, Length 31mm, 6/Pkg	29-12219
Size 50, Length 31mm, 6/Pkg	29-12220
Size 55, Length 31mm, 6/Pkg	29-12221
Size 60, Length 31mm, 6/Pkg	29-12222
Size 70, Length 31mm, 6/Pkg	29-12223
Size 80, Length 31mm, 6/Pkg	29-12224
Size 15-40, Length 31mm, 6/Pkg	29-12225
Size 45-80, Length 31mm, 6/Pkg	29-12238

Rotary Engine Files

Engine H-Files Stainless Steel (JS Dental)

- •ISO color coded R/A aluminum handle
- •Use in 1/4 hand-piece only
- •Fast cutting instrument with one continuous flute for fast

one continuous flute for fast removal of debris	
Size 10, Length 21mm, 6/Pkg	29-11376
Size 15, Length 21mm, 6/Pkg	29-11360
Size 20, Length 21mm, 6/Pkg	29-11359
Size 25, Length 21mm, 6/Pkg	29-11358
Size 30, Length 21mm, 6/Pkg	29-11357
Size 35, Length 21mm, 6/Pkg	29-11356
Size 40, Length 21mm, 6/Pkg	29-11355
Size 45, Length 21mm, 6/Pkg	29-11354
Size 50, Length 21mm, 6/Pkg	29-11353



CIS, IIIC	
Size 55, Length 21mm, 6/Pkg	29-11352
Size 60, Length 21mm, 6/Pkg	29-11351
Size 70, Length 21mm, 6/Pkg	29-11350
Size 80, Length 21mm, 6/Pkg	29-11349
Size 15-40, Length 21mm, 6/Pkg	29-11348
Size 45-80, Length 21mm, 6/Pkg	29-11361
Size 10, Length 25mm, 6/Pkg	29-11362
Size 15, Length 25mm, 6/Pkg	29-11375
Size 20, Length 25mm, 6/Pkg	29-11374
Size 25, Length 25mm, 6/Pkg	29-11373
Size 30, Length 25mm, 6/Pkg	29-11372
Size 35, Length 25mm, 6/Pkg	29-11371
Size 40, Length 25mm, 6/Pkg	29-11370
Size 45, Length 25mm, 6/Pkg	29-11369
Size 50, Length 25mm, 6/Pkg	29-11368
Size 55, Length 25mm, 6/Pkg	29-11367
Size 60, Length 25mm, 6/Pkg	29-11366
Size 70, Length 25mm, 6/Pkg	29-11365
Size 80, Length 25mm, 6/Pkg	29-11364
Size 15-40, Length 25mm, 6/Pkg	29-11347
Size 45-80, Length 25mm, 6/Pkg	29-11363
E 1 E C. 1 (10 E	

Engine Reamers- Stainless Steel (JS Dental)

- •ISO color coded R/A aluminum handle
- •Use in 1/4 hand-piece
- •Designed for smooth

29-11377
29-11393
29-11394
29-11395
29-11396
29-11397
29-11398
29-11399
29-11406
29-11405
29-11404
29-11403
29-11402

Size 15-40, Length 21mm, 6/Pkg

Size 45-80, Length 21mm, 6/Pkg	29-11392
Size 10, Length 25mm, 6/Pkg	29-11391
Size 15, Length 25mm, 6/Pkg	29-11378
Size 20, Length 25mm, 6/Pkg	29-11379
Size 25, Length 25mm, 6/Pkg	29-11380
Size 30, Length 25mm, 6/Pkg	29-11381
Size 35, Length 25mm, 6/Pkg	29-11382
Size 40, Length 25mm, 6/Pkg	29-11383
Size 45, Length 25mm, 6/Pkg	29-11384
Size 50, Length 25mm, 6/Pkg	29-11385
Size 55, Length 25mm, 6/Pkg	29-11386
Size 60, Length 25mm, 6/Pkg	29-11387
Size 70, Length 25mm, 6/Pkg	29-11388
Size 80, Length 25mm, 6/Pkg	29-11389
Size 15-40, Length 25mm, 6/Pkg	29-11400
Size 45-80, Length 25mm, 6/Pkg	29-11390

Engine	S_Files	Stainless	Steel	CIS D	ental
Lugine	9-1 HC3	Stailless	Sieci	(JO D	entari

•ISO color coded R/A aluminum
handle

- •Rubber stops and mm marks


•Use in 1/4 turn handpiece only •Perfect blend of flexibility, durability and cutting ability	and the state of t
Size 10, Length 21mm, 6/Pkg	29-11407
Size 15, Length 21mm, 6/Pkg	29-11423
Size 20, Length 21mm, 6/Pkg	29-11424
Size 25, Length 21mm, 6/Pkg	29-11425
Size 30, Length 21mm, 6/Pkg	29-11426
Size 35, Length 21mm, 6/Pkg	29-11427
Size 40, Length 21mm, 6/Pkg	29-11428
Size 45, Length 21mm, 6/Pkg	29-11429
Size 50, Length 21mm, 6/Pkg	29-11430
Size 55, Length 21mm, 6/Pkg	29-11431
Size 60, Length 21mm, 6/Pkg	29-11432
Size 70, Length 21mm, 6/Pkg	29-11433
Size 80, Length 21mm, 6/Pkg	29-11434
Size 15-40, Length 21mm, 6/Pkg	29-11435
Size 45-80, Length 21mm, 6/Pkg	29-11422
Size 10, Length 25mm, 6/Pkg	29-11421
Size 15, Length 25mm, 6/Pkg	29-11408
Size 20, Length 25mm, 6/Pkg	29-11409
Size 25, Length 25mm, 6/Pkg	29-11410



Size 30, Length 25mm, 6/Pkg	29-11411
Size 35, Length 25mm, 6/Pkg	29-11412
Size 40, Length 25mm, 6/Pkg	29-11413
Size 45, Length 25mm, 6/Pkg	29-11414
Size 50, Length 25mm, 6/Pkg	29-11415
Size 55, Length 25mm, 6/Pkg	29-11416
Size 60, Length 25mm, 6/Pkg	29-11417
Size 70, Length 25mm, 6/Pkg	29-11418
Size 80, Length 25mm, 6/Pkg	29-11419
Size 15-40, Length 25mm, 6/Pkg	29-11436
Size 45-80, Length 25mm, 6/Pkg	29-11420

Engine T-Files Stainless Steel (JS Dental)

- •ISO color coded R/A aluminum handle
- •Rubber stops and mm markings •A combination of a K-File and
- H-File with two continuous flutes makes this instrument effective and stable for fast removal of debris



debris	
Size 10, Length 21mm, 6/Pkg	29-11437
Size 15, Length 21mm, 6/Pkg	29-11445
Size 20, Length 21mm, 6/Pkg	29-11454
Size 25, Length 21mm, 6/Pkg	29-11455
Size 30, Length 21mm, 6/Pkg	29-11456
Size 35, Length 21mm, 6/Pkg	29-11457
Size 40, Length 21mm, 6/Pkg	29-11458
Size 45, Length 21mm, 6/Pkg	29-11459
Size 50, Length 21mm, 6/Pkg	29-11466
Size 55, Length 21mm, 6/Pkg	29-11465
Size 60, Length 21mm, 6/Pkg	29-11464
Size 70, Length 21mm, 6/Pkg	29-11463
Size 80, Length 21mm, 6/Pkg	29-11462
Size 15-40, Length 21mm, 6/Pkg	29-11461
Size 45-80, Length 21mm, 6/Pkg	29-11446
Size 10, Length 25mm, 6/Pkg	29-11447
Size 15, Length 25mm, 6/Pkg	29-11438
Size 20, Length 25mm, 6/Pkg	29-11439
Size 25, Length 25mm, 6/Pkg	29-11440
Size 30, Length 25mm, 6/Pkg	29-11441
Size 35, Length 25mm, 6/Pkg	29-11442
Size 40, Length 25mm, 6/Pkg	29-11443
Size 45, Length 25mm, 6/Pkg	29-11444

Size 50, Length 25mm, 6/Pkg	29-11453	
Size 55, Length 25mm, 6/Pkg		
Size 60, Length 25mm, 6/Pkg	29-11451	low.
Size 70, Length 25mm, 6/Pkg	29-11450	
Size 80, Length 25mm, 6/Pkg	29-11449	
Size 15-40, Length 25mm, 6/Pkg	29-11460	
Size 45-80, Length 25mm, 6/Pkg	29-11448	
	Size 55, Length 25mm, 6/Pkg Size 60, Length 25mm, 6/Pkg Size 70, Length 25mm, 6/Pkg Size 80, Length 25mm, 6/Pkg Size 15-40, Length 25mm, 6/Pkg	Size 55, Length 25mm, 6/Pkg 29-11452 Size 60, Length 25mm, 6/Pkg 29-11451 Size 70, Length 25mm, 6/Pkg 29-11450 Size 80, Length 25mm, 6/Pkg 29-11449 Size 15-40, Length 25mm, 6/Pkg 29-11460

Hyflex CM NiTi Rotary Files 19mm (Coltene)

CONTROLLED MEMORY EFFECT: HYFLEX CM NiTi

Files have been manufactured utilizing a unique process that controls the material's memory, making the files extremely flexible, with no rebound. This enhancement in NiTi file production technology gives the file the ability to more closely



follow the anatomy of the canal, reducing the risk of the file separation, binding and complications such as ledging, transportation, and perforation.

FATIGUE RESISTANT-LESS BREAKAGE: HYFLEX® CM Niti Files with Controlled memory Effect are up to 385% more

resistant to cyclical fatigue compared to other Niti files, reducing the incidence of file separation.

EASY TO FOLLOW, SAFE FILE SEQUENCE: HYFLEX CM

Niti Files come with a recommended sequence of 3-6 files depending on the size of the root canal, making the procedure straightforward and easy to follow.

08/25, 19mm, 6/Pk

29-11713

Hyflex CM NiTi Rotary Files 21mm (Coltene)

CONTROLLED MEMORY EFFECT:

HYFLEX CM NiTi Files have been manufactured utilizing a unique process that controls the material's memory, making the files extremely flexible, with no re-bound. This enhancement in NiTi file production technology gives the file the ability to more closely follow the anatomy of the canal, reducing the risk of the file separation, binding and complications such as ledging, transportation, and perforation.



FATIGUE RESISTANT-LESS BREAKAGE: HYFLEX® CM

Niti Files with Controlled memory Effect are up to 385% more resistant to cyclical fatigue compared to other Niti files, reducing the incidence of file separation

EASY TO FOLLOW, SAFE FILE SEQUENCE: HYFLEX CM

Niti Files come with a recommended sequence of 3-6 files depending on the size of the root canal, making the procedure straightfor-

06/25, 21mm, 6/Pkg	29-11717
06/30, 21mm, 6/Pkg	29-11716
06/35, 21mm, 6/Pkg	29-11715
06/40, 21mm, 6/Pkg	29-11714



04/15, 21mm, 6/Pkg	29-11727
04/20, 21mm, 6/Pkg	29-11726
04/25, 21mm, 6/Pkg	29-11725
04/30, 21mm, 6/Pkg	29-11724
04/35, 21mm, 6/Pkg	29-11723
04/40, 21mm, 6/Pkg	29-11722
04/45, 21mm, 6/Pkg	29-11721
04/50, 21mm, 6/Pkg	29-11720
04/60, 21mm, 6/Pkg	29-11719
06/20, 21mm, 6/Pkg	29-11718

Hyflex CM NiTi Rotary Files 25mm (Coltene)

300% more resistance to separation HyFlex® CM™ NiTi files with Controlled Memory are up to 300% more resistant to cyclical fatigue compared to other NiTi files* which substantially helps reducing the incidence of file separation. No rebound + Extreme flexibility = Superior Canal Tracking HyFlex® Controlled Memory NiTi files have been manufactured utilizing a unique process



that controls the material's memory, making the files extremely flexible but without the shape memory of other NiTi files. This gives the file the ability to follow the anatomy of the canal very closely, reducing the risk of ledging, transportation or perforation. Regains shape after sterilization = Multi-use The shape and strength of files with straightened spirals can be restored during autoclaving and reused. Files not returning to original shape should be discarded.

04/15 Taper Files, 25mm, 6/Pkg	29-11741
04/20 Taper Files, 25mm, 6/Pkg	29-11740
04/25 Taper Files, 25mm, 6/Pkg	29-11739
04/30 Taper Files, 25mm, 6/Pkg	29-11738
04/35 Taper Files, 25mm, 6/Pkg	29-11737
04/40 Taper Files, 25mm, 6/Pkg	29-11736
04/45 Taper Files, 25mm, 6/Pkg	29-11735
04/50 Taper Files, 25mm, 6/Pkg	29-11734
04/60 Taper Files, 25mm, 6/Pkg	29-11733
06/20 Taper Files, 25mm, 6/Pkg	29-11732
06/25 Taper Files, 25mm, 6/Pkg	29-11731
06/30 Taper Files, 25mm, 6/Pkg	29-11730
06/35 Taper Files, 25mm, 6/Pkg	29-11729
06/40 Taper Files, 25mm, 6/Pkg	29-11728



Lentulo Paste Fillers Stainless Steel (JS Dental)

Latch type instrument designed for delivery of root canal sealers, cements, and medicaments

Size 25, Length 21mm, 4/Pkg	29-11944
Size 30, Length 21mm, 4/Pkg	29-11952
Size 35, Length 21mm, 4/Pkg	29-11951
Size 40, Length 21mm, 4/Pkg	29-11950
Size 25-40, Length 21mm, 4/Pkg	29-11948
Size 25, Length 25mm, 4/Pkg	29-11949
Size 30, Length 25mm, 4/Pkg	29-11947
Size 35, Length 25mm, 4/Pkg	29-11946
Size 40, Length 25mm, 4/Pkg	29-11945
Size 25-40, Length 25mm, 4/Pkg	29-11953

Liberator Nickel Titanium Rotary Files (Miltex)

The liberatorTM rotary nickel-titanium endodontic files are the result of years of research and development as well as feedback from leading clinicians in the field of endodontics. The files incorporate a unique straight flute design and manufacturing process that eliminates the tradi-



tional helical flutes found on virtually all rotary endodontic files. Feel the safety and control - The patented design of the noncutting Roane tip minimizes ledging and transportation while helping to keep liberator files centered in the canal.

Liberator files are designed with straight flutes that will not selfthread, unlike helically fluted files. As proof, a controlled study measure the incidence of self-threading of the Liberator compared to competitive rotary files as shown here.

Feel the superior cutting efficiency - The triangular cross section and lack of radial lands provide sharp cutting edges and reduce torque demand on Liberator files. Because the files operate at high RPM (1,500-200) compared to conventional files (300-500), torque is further reduced according to this basic rule of physics. Efficient debris removal - A vortex created by the high rotational speed effectively removes debris.

Faster, Easy-Use - Liberator files rotating at 1,500-2000RPM cut dentin faster than conventional rotary files at 300RPM. This higher speed also provides excellent tactile feel in canals. Plus, torque control is not required for electric motors. And even air-driven handpieces may be used at 1,500 to 2,000 RPM.

#15 .02 Taper 21mm, 6/Pkg	29-11979
#20 .02 Taper 21mm, 6/Pkg	29-11978
#25 .02 Taper 21mm, 6/Pkg	29-11977
#30 .02 Taper 21mm, 6/Pkg	29-11976
#35 .02 Taper 21mm, 6/Pkg	29-11975



410, 11-	
#40 .02 Taper 21mm, 6/Pkg	29-11974
#15-40 .02 Taper 21mm, 6/Pkg	29-11973
#45 .02 Taper 21mm, 6/Pkg	29-11970
#50 .02 Taper 21mm, 6/Pkg	29-11972
#55 .02 Taper 21mm, 6/Pkg	29-11971
#60 .02 Taper 21mm, 6/Pkg	29-11969
#70 .02 Taper 21mm, 6/Pkg	29-11968
#80 .02 Taper 21mm, 6/Pkg	29-11967
#45-80 .02 Taper 21mm, 6/Pkg	29-11966
#15 .02 Taper 25mm, 6/Pkg	29-11965
#20 .02 Taper 25mm, 6/Pkg	29-11964
#25 .02 Taper 25mm, 6/Pkg	29-11963
#30 .02 Taper 25mm, 6/Pkg	29-11962
#35 .02 Taper 25mm, 6/Pkg	29-11961
#40 .02 Taper 25mm, 6/Pkg	29-11958
#15-40 .02 Taper 25mm, 6/Pkg	29-11959
#45 .02 Taper 25mm, 6/Pkg	29-11960
#50 .02 Taper 25mm, 6/Pkg	29-11957
#55 .02 Taper 25mm, 6/Pkg	29-11956
# 60 .02 Taper 25mm, 6/Pkg	29-11955
#70 .02 Taper 25mm, 6/Pkg	29-11954
#80 .02 Taper 25mm, 6/Pkg	29-11980
#45-80 .02 Taper 25mm, 6/Pkg	29-11981


Mity Nickel Titanium 360 Rotary Files (JS Dental)

- •ISO cc aluminum handle
- •MM markings and stops
- •Full rotation engine driven instrument for preparing a smoother, rounder, more uniform canal

Size 60, Length 21mm, 6/Pkg

Size 70, Length 21mm, 6/Pkg

Size 80, Length 21mm, 6/Pkg



29-12146

29-12145

Size 15-40, Length 21mm, 6/Pkg	29-12143
Size 45-80, Length 21mm, 6/Pkg	29-12155
Size 15, Length 25mm, 6/Pkg	29-12156
Size 20, Length 25mm, 6/Pkg	29-12168
Size 25, Length 25mm, 6/Pkg	29-12167
Size 30, Length 25mm, 6/Pkg	29-12166
Size 35, Length 25mm, 6/Pkg	29-12165
Size 40, Length 25mm, 6/Pkg	29-12164
Size 45, Length 25mm, 6/Pkg	29-12163
Size 50, Length 25mm, 6/Pkg	29-12162
Size 55, Length 25mm, 6/Pkg	29-12161
Size 60, Length 25mm, 6/Pkg	29-12160
Size 70, Length 25mm, 6/Pkg	29-12159
Size 80, Length 25mm, 6/Pkg	29-12158
Size 15-40, Length 25mm, 6/Pkg	29-12157
Size 45-80, Length 25mm, 6/Pkg	29-12142

Mity Nickel Titanium Engine T-Files (JS Dental)

•ISO color coded aluminu	m
handle	

•MM markings and rubber	A Contract of the Contract of	
stops •Engine driven 1/4 turn Mity	A STATE OF THE STA	
Turbo file		
Size 15, Length 21mm, 6/Pkg	29-12020	
Size 20, Length 21mm, 6/Pkg	29-12005	
Size 25, Length 21mm, 6/Pkg	29-12004	
Size 30, Length 21mm, 6/Pkg	29-12003	
Size 35, Length 21mm, 6/Pkg	29-12002	
Size 40, Length 21mm, 6/Pkg	29-12001	
Size 45, Length 21mm, 6/Pkg	29-12000	
Size 50, Length 21mm, 6/Pkg	29-11999	
Size 55, Length 21mm, 6/Pkg	29-11998	
Size 60, Length 21mm, 6/Pkg	29-11997	
Size 70, Length 21mm, 6/Pkg	29-11996	
Size 80, Length 21mm, 6/Pkg	29-11995	
Size 15-40, Length 21mm, 6/1	Pkg 29-11994	
Size 45-80, Length 21mm, 6/1	Pkg 29-12006	
Size 15, Length 25mm, 6/Pkg	29-12007	
Size 20, Length 25mm, 6/Pkg	29-12019	
Size 25, Length 25mm, 6/Pkg	29-12018	
Size 30, Length 25mm, 6/Pkg	29-12017	


29-12016



Size 40, Length 25mm, 6/Pkg	29-12015
Size 45, Length 25mm, 6/Pkg	29-12014
Size 50, Length 25mm, 6/Pkg	29-12013
Size 55, Length 25mm, 6/Pkg	29-12012
Size 60, Length 25mm, 6/Pkg	29-12011
Size 70, Length 25mm, 6/Pkg	29-12010
Size 80, Length 25mm, 6/Pkg	29-12009
Size 15-40, Length 25mm, 6/Pkg	29-12008
Size 45-80, Length 25mm, 6/Pkg	29-11993

Peeso Reamers Stainless Steel Files (JS Dental)

Swiss made latch type instrument for removal of coronal portion of tooth as well as gutta-percha for post placement.



Pro NiTi Rotary Tapered Files - ATW

NiTi Rotary files with good elasticity and appearance




SX , 25mm , 6/Pkg	29-19221
S1, 25mm, 6/Pkg	29-19222
S2, 25mm, 6/Pkg	29-19223
F1, 25mm, 6/Pkg	29-19224
F2, 25mm, 6/Pkg	29-19225

Size 35, Length 25mm, 6/Pkg

F3, 25mm , 6/Pkg	29-19226
SX, 21mm, 6/Pkg	29-19227
S1, 21mm, 6/Pkg	29-19228
S2, 21mm, 6/Pkg	29-19229
F1, 21mm, 6/Pkg	29-19230
F2, 21mm, 6/Pkg	29-19231
F3 , 21mm , 6/Pkg	29-19232

Safet Spring Spiral Fillers Stainless Steel (JS Dental)

Similar to Paste Filler but with a safety spring all the way to the handle for extra flexibility and safer use.



Versa Brush (Vista)

Versa BrushTM is a rotational utility brush which is ideal for a variety of endodontic and restorative procedures



Incredibly versatile, strong and bendable, *Versa Brush*TM *is a*

"must have" for both endodontists and restorative dentists. It is the ideal tool for cleaning endodontic access cavities and removing contaminants left behind, such as sealer components and smear layer.

Used in a latch-grip hand piece at 250rpm or less, the prepared space can be literally scrubbed clean, allowing for more reliable and predictable bonding.

Versa BrushTM is perfect for:

- •Cement removal
- •Post hole cleaning
- •Cleaning coronal access openings
- •Use with EDTA to enhance cleansing and removal of smear layer

40/Pkg 23-26241 **5/Pkg** 23-26242

Xtreme 360° Files No Memory (JS Dental)

•ISO cc aluminum handle

•MM markings and stops

•Rotary engine driven instrument for preparing a smoother, rounder, more uniform canal

Size 15, Length 21mm, 6/Pkg	29-12549
Size 20, Length 21mm, 6/Pkg	29-12537
Size 25, Length 21mm, 6/Pkg	29-12538
Size 30, Length 21mm, 6/Pkg	29-12539



Size 35, Length 21mm, 6/Pkg	29-12540
Size 40, Length 21mm, 6/Pkg	29-12541
Size 15-40, Length 21mm, 6/Pkg	29-12542
Size 15, Length 25mm, 6/Pkg	29-12543
Size 20, Length 25mm, 6/Pkg	29-12544
Size 25, Length 25mm, 6/Pkg	29-12545
Size 30, Length 25mm, 6/Pkg	29-12546
Size 35, Length 25mm, 6/Pkg	29-12548
Size 40, Length 25mm, 6/Pkg	29-12547
Size 15-40, Length 25mm, 6/Pkg	29-12536

Endodontic Stops

Ease Stops (JS Dental)

White dispenser filled with 250 stops of your color choice. Push button and place stop right onto file. Left and right measuring ruler for exact placement of stops.



Color Rust, 250 Stops/Pkg	29-11245
Color White, 250 Stops/Pkg	29-11246
Color Yellow, 250 Stops/Pkg	29-11247
Color Red, 250 Stops/Pkg	29-11248
Color Blue, 250 Stops/Pkg	29-11249
Color Green, 250 Stops/Pkg	29-11250
Color Black, 250 Stops/Pkg	29-11251

Silicone Stops (JS Dental)





Color Rust, Round, 100/Pkg	29-12464
Color Rust, Teardrop, 100/Pkg	29-12465
Color Asst., Round, 100/Pkg	29-12467
Color Asst., Teardrop, 100/Pkg	29-12466

Absorbent Paper Points

Absorbent Paper Points (Meta Biomed)

200 point color coded per Box



Standardized Size 8	29-11112
Standardized Size 10	29-11111

Standardized Size 15	29-11113
Standardized Size 20	29-11114
Standardized Size 25	29-11115
Standardized Size 30	29-11116
Standardized Size 35	29-11117
Standardized Size 40	29-11118
Standardized Size 45	29-11119
Standardized Size 50	29-11120
Standardized Size 55	29-11121
Standardized Size 60	29-11122
Standardized Size 70	29-11123
Standardized Size 80	29-11124
Standardized Size 90	29-11125
Standardized Size 100	29-11126
Standardized Size 110	29-11127
Standardized Size 120	29-11128
Standardized Size 130	29-11129
Standardized Size 140	29-11130
Non-Standardized Size XX-Fine	29-11131
Non-Standardized Size X-Fine	29-11132
Non-Standardized Size Fine	29-11133
Non-Standardized Size Medium	29-11134
Non-Standardized Size Coarse	29-11135
Non-Standardized Size X-Coarse	29-11136

Taper Paper Points - Meta



Special Taper .04 Size 15	29-11137
Special Taper .04 Size 20	29-11138
Special Taper .04 Size 25	29-11139
Special Taper .04 Size 30	29-11140
Special Taper .04 Size 35	29-11141
Special Taper .04 Size 40	29-11142
Special Taper .04 Size 45	29-11143
Special Taper .04 Size 50	29-11144
Special Taper .06 Size 15	29-11149
Special Taper .06 Size 20	29-11150



Special Taper .06 Size 25	29-11151
Special Taper .06 Size 30	29-11152
Special Taper .06 Size 35	29-11153
Special Taper .06 Size 40	29-11154
Special Taper .06 Size 45	29-11155
Special Taper .06 Size 50	29-11156

Gutta Percha Point

GPR Gutta Percha Remover (Mani)



- Superior results Fast, easy, & efficient.
- Fast & dynamic removal stress-free and easy to use, w/ torque control and reversal
- No need for solvents
- Max allowable speed 1000+/-500min 0.04 Taper

Max anowable speed 1000 17 300mm 0.01	ruper
GPR #70 Size 1S, Green color, 16mm Stainless Steel, 4/Pkg	29-109441
GPR #50 Size 2S, Gold color, 18mm Stainless Steel, 4/Pkg	29-109442
GPR #40 Size 3N, Black color, 21mm Stainless Steel, 4/Pkg	29-109443
GPR #30 Size 4N, Blue color, 21mm Stainless Steel, 4/Pkg	29-109444
GPR Assorted Size #30-#70, 4/Pkg	29-109445

Gutta Box (Union Branch)



Color CodedQuantity areLimited

SILVER Color	29-108471
RED Color	29-108472
BLUE Color	29-108473

Gutta Percha Cutter and Trimmer

In 3 seconds only

This instrument is developed for the easy and safe removal and trimming of excess gutta-percha and thermal dense filling in your endodontic cases. It is composed of endo-remover and trimmer operated by two AAA, NiMH rechargeable batteries, which are included in the package.



ERGONOMIC DESIGN

SLIM: To feel comfortable for the treatment in the mouth like dental handpieces

LIGHTWEIGHT: It weighs only 50g together wit 2 cells which leads to easiness in your treatment

PRODUCTIVE ENVIRONMENTS

There is no need for alcohol lamp or burner to heat some instrument to remove the excess, also no need for an assistant.

ECONOMICAL:

If NI-HM batteries are used, you can use this instrument in very powerful and economical ways than disposable alkaline batteries . Normally, it has about 400-500 charging cycles which leads to long lasting usages, also if fully charged, normal usage is between 15-20 times

SAFETY & DURABILITY

This instrument might be semi-permanent only if batteries and tip are refilled, also the patient can not notice the treatment because they can not see compared than the usage of conventional alcohol lamp usage. The patients feels comfortable.

Prior to using this instrument ,please charge and check the batteries.

1. Push the button only when are ready to operate

- 2. Place the edge of the tip to the canal orifice, and then push the button for 3 or 4 seconds. It produces enough heat to cut the excess within seconds. That means you don't need to hold the button for overlong time.
- 3. After removing the excess ,please pull the instrument carefully not to touch the skin of the patient.
- 4. After it comes cool, please clean the tip with wet alcohol gauze or sterilize for the next usage.

Cutter & Trimmer Unit: with 4 X AAA

rechargeable Ni-MH batteries, Rapid

charger and two cutting and trimmer tips 29-11539

Gutta Percha Points (Meta)



Standardized Size 15,	
box of 120 point color coded	29-11542
Standardized Size 20, box of 120 point color coded	29-11543
Standardized Size 25, box of 120 point color coded	29-11544
Standardized Size 30, box of 120 point color coded	29-11545
Standardized Size 35, box of 120 point color coded	29-11546
Standardized Size 40, box of 120 point color coded	29-11547
Standardized Size 45, box of 120 point color coded	29-11548
Standardized Size 50, box of 120 point color coded	29-11549
Standardized Size 55, box of 120 point color coded	29-11550
Standardized Size 60, box of 120 point color coded	29-11551



Standardized Size 70, box of 120 point color coded	29-11552
Standardized Size 80, box of 120 point color coded	29-11553
Non-Standardized Size X-Fine	29-11560
Non-Standardized Size Fine-Fine	29-11561
Non-Standardized Medium-Fine	29-11562
Non-Standardized Fine	29-11563
Non-Standardized Fine-Medium	29-11564
Non-Standardized Medium	29-11565
Non-Standardized Medium-Large	29-11566
Non-Standardized Large	29-11567
Non-Standardized X-Large	29-11568

Protaper Gutta Percha Point (Meta)



- For in sliding Box
- Length market for accurate depth measurement and insertion

Assorted F1, F2, & F3, 60/Box	29-105514
F3 , 60/Box	29-105513
F2 , 60/Box	29-105512
F1 , 60/Box	29-105511

Taper Gutta Percha Point (Meta)



#15 Taper 0.04, 60/Box	29-106391
#20 Taper 0.04, 60/Box	29-106392
#25 Taper 0.04, 60/Box	29-106393
#30 Taper 0.04, 60/Box	29-106394
#35 Taper 0.04, 60/Box	29-106395
#40 Taper 0.04, 60/Box	29-106396
#15-40 Taper 0.04, 60/Box	29-106397
#15 Taper 0.06, 60/Box	29-106398
#20 Taper 0.06, 60/Box	29-106399
#25 Taper 0.06, 60/Box	29-106400
#30 Taper 0.06, 60/Box	29-106401
#35 Taper 0.06, 60/Box	29-106402
#40 Taper 0.06, 60/Box	29-106403


#15-40 Taper 0.06, 60/Box

29-106404

Root Canal Sealers Repairs

Acroseal Root Canal Sealer (Septodont)

This endodontic sealer seals the canal permanently, eliminating the worries about leakage and minimizing the risks of root canal treatment failures





- •Low viscosity allows it to flow into accessory canals for superior fill every time
- •Highly compatible can be used with all gutta-percha obturation techniques, cold or warm
- •Eugenol-free fully compatible with most cements and core materials
- •Automix syringe dispensing systems delivers a homogeneous mix without the mess
- •Radiopaque

Kit: 1 x 8.5g tube of base paste, 1 x

9.5g tube of catalyst paste, 1 crossed-ruled mixing pad

29-94111

Adseal Root Canal Sealer (Meta Biomed)



Resin-based root canal sealer (Compare to Tubli-Seal 10.5q (KRR))

Advantage: Excellent


biocompatibility, Easy to mix

paste-paste, Hermetic sealing ability, Non staining to the teeth, Insoluble in tissue fluids and Good radiopacity

Components: 13.5g dual syringe (Base 9g, Epoxy resin, Calcium phosphate) Catalyst 4.5g (Amines, Bismuth sub carbonate), Spatula,

, Mixing Plate, Direction for use, Dimension : $80\times200\times43$ (mm),







Working Time: 35 Minutes at 23°C (73°F) Setting Time: 45 Minutes at 37°C (99°F)

Film Thickness : 3.3 µm
Biocompatibility

Test M	1ethod	Results
Cytotoxicity (ISO 10993-5)	Agar Diffusion Freshly mixed Set material	mild none
Irritation (ISO 10993-10)	oral	no irritation
Systemic (ISO 10993-11)	acute	safe
Implantation (ISO 10993-6)		acceptable (no edema)

Adseal, 13.5g

29-11161

AH Plus (Dentsply)

Root Canal Sealer Paste-Paste resin system for better handling Very high radiopacity for easy detection



Standard Pkg. of 2 x 4ml

29-11162

AH26 (Dentsply)

Resin based, non-acrylic, eugenol-free sealer Used with gutta percha:

- •Densfil obturators or silver points Radiopaque
- •Working time is 4.5 hours · Set time is 12-48 hours
- Suitable for a cold & warm obturation techniques



Standard Pkg. 8gm powder & 10gm 29-11163 resin

AH 26 Resin (Liquid Only) 10gm tube 29-11164

BioRoot RCS (Septodont)



A bioactive and biocompatible root canal sealer with active biosilicate technology. The mineral based root canal filling material is for permanent canal obturation. BioRoot RCS does not stain teeth and eliminates post-op sensitivity,

Features & Benefits

- Incorporates an easy to use cold

obturation technique, ensuring a long lasting, leak-free seal - High pH (~11) stops bacterial growth & alleviates risk of intracanal reinfection

- Bio compatible & Resin-Free
- Easy to hand-mix & place

29-108091 15gm

Diapex Root Canal Sealer (Diadent)



Pre-mixed calcium hydroxide and iodoform (in a silicone oil base) paste for temporary or permanent root canal treatment. Diapex is highly radiopaque and antibacterial.

2.0g Paste & 20 Tips


29-94131

Tips, 20/Pkg

29-94132

Endomet Root Canal Sealer (Septodont)

- Easy to use
- · Contains thymol iodide, an antiseptic agent
- Radiopaque and non-resorbable
- within the canal
- Stays soft on the mixing pad but rapidly hardens when placed into the tooth canal





catalyst,1 mixing pad

29-94091

Metapex

Metapex (Meta Biomed)

Calcium Hydroxide with **Iodoform** Temporary Root Canal Filling Material with Excellent Radiopacity (Compare with Diapex from Diadent 2.0g paste)

- · Weeping canal
- Formation of a hard tissue barrier(Apexification)
- · Root canal filling for deciduous teeth
- Internal & external root resorption
- Excellent antibacterial effect and radiopacity
- Premixed paste
- Premixed paste
- Disposable tips for prevention of cross-contamination

Metapex, 2.2g 29-11986

Pulp Canal Sealer EWT (Kerr)



- Longer working time 6-8 hours on the pad after mixing. - Ideal for warm gutta percha
- technique.

Q ETA

- Radiopaque - Low solubility in oral fluids

Standard Pkg: 1-10.5gm powder, 1-4ml 29-109931 liquid, 2 dispensing scoop, instruction

www.prestigedentalproducts.com



Root Canal Sealer Silver-Free (Prime Dent)

Excellent sealing capacity



Complete Pkg: 8gm Powder & 10gm resin 29-94101

RSA Roeko Seal Automix (Coltene Whaledent)

A state-of-the-art polydimethylsiloxane based root canal sealer, features the following:

- •Ready for use in a 1:1 automix system
- •Flows into dentine tubules and tiny crevices
- •Sets in both moist and dry conditions
- Insoluable
- •Expands slightly (0.2%)
- •Provides mechanical retention through expansion
- •Permits subsequent revision of the canal
- •Does not contain Eugenol
- •Radiopaque
- Extremely biocompatible
- •Dispenses the sealer homogeneously mixed and free of air

Standard Set: Dual-barrel syringe, 12

flexible mixing tips, 1 cap 29-94121

Refill Mix Tips, 24/Pkg 29-94122

Sealapex Calcium Hydroxide (Kerr)



Calcium Hydroxide Polymeric Root Canal Sealer

12gm Tube Base, 12gm Tube Catalyst

& Mixing Pad

29-16171

Tubli-Seal Endodontic Sealers & Cements (Kerr)

Tubli-Seal is a paste/paste zinc oxide eugenol root canal sealer. Tubli-Seal is light in color, non darkening, and radiopaque.



10gm Tube Base, 1-3.5gm Tube Accelerator & 1 Mixing Pad

29-12484



Call Us Today 833-733-682

Well-Pex Root Canal Filling Material (Vericom)



- Oil-Based Calcium Hydroxide Pase with Iodoform
- Root canal filling material for the infected root canal
- Good radipacity and antibacterial characteristics

- Very stable without any

solidification or separation and highly flow able paste premixed in a


- **EXCELLENT** accessibility to the apex and biocompatibility to the periapical tissue.

2g syringe with disposable tips & **Extension Holder**

29-12986

Chelating Agents

EDTA Solution Aqueous Irrigant (Vista)



•Ideal for root canal preparation & smear layer removal



120cc Refill Bottle 29-87681

3cc Unit Dose.

20 single use patient applications along

with 20 Vista Probe needles, 20/Pkg 29-87682

30ml Bulk Syringe 29-87683

480ml Bottle 29-87684

MD-ChelCream 19% ETDA (Meta Biomed)

Easier Filing with Lubricating

Action (Compare to RC

Prep)

Applications:

- •Root canal preparation with files and reamers
- •Presence of
- calcifications in the canal
- •Narrow canals in posterior teeth
- •Smear layer removal

Advantages

- •Removal of inorganic substances by chelation (EDTA-Ca²)
- •Easier filing with lubricating action
- •Produces the cleanest canal wall after smear layer removal
- •Use MD-ChelCream with discretion to prevent perforation
- •Avoid extrusion of the MD-ChelCream into the apex
- •Keep the syringe capped all the time
- •Store in temperature 10°C~ 25°C (50°F~ 77°F)

2 Sys. 5.4ml 29-11983



Prestige Dental Products, Inc. Prestigedentalproducts.com

Tel: 877-772-3888

MD-Cleanser 17% ETDA (Meta Biomed)

Excellent Smear Layer Removal (Compare to EDTA Solution from PulpDent)

Applications:

- •Presence of calcifications in the canal
- •Narrow canals in posterior teeth
- •Smear layer removal

- Advantages:
 •Removal of inorganic substances by chelation (EDTA-Ca²)
- ·Easier filing with lubricating action
- •Produces the cleanest canal wall after smear laver removal

Components:

•MD-Cleanser: 100mℓ

•1cc syringe

Disposable tips

•Dimension : $50 \times 115 \times 55$ (mm)

•Weight: 138g

100ml w/One Syringe & 3 Tips

29-11984

MD-Clear

Prep-Rite RC (PulpDent)



•17% viscous EDTA gel that is picked up on an endodontic file to facilitate preperation of the root canal

Kit Contains: 4-5g Prep-Rite RC Syringes 29-87691

RC PREP (Premier)





For Chemo-Mechanical **Preparation of Root** Canals

Over 40 years of clinical use and numerous publications demonstrate the clinical efficacy of

RC-Prep for Chemo-Mechanical preparation of Root Canals.

RC-Prep's unique formulation of glycol, urea peroxide and EDTA in a special water soluble base helps remove calcifications and lubricates the canal to permit more efficient instrumentation. Whether you prefer finger or engine driven instruments, RC-Prep allows reamers and files to move easily without binding. RC-Prep is also excellent for use with apex locators as it permits consistently reliable readings.

29-12280 Jar 18gm **Syringe Kit Contains:** 29-70433

5-3cc Syringes & 50 Tips

RC Prep MicroDose (Premier)

Direct Application Into the Canal

RC-PrepTM MicroDose® is convenient, easy to dispense and eliminates the risk of cross-contamination.





Each MicroDose pipette contains the ideal amount of RC-Prep (0.2cc) for more efficient instrumentation during root-canal procedures.

You can apply RC-Prep directly into the pulp chamber or into a dappen dish

After each use, the empty pipette is discarded and the MicroDose syringe can be conveniently autoclaved

Refill (Pkg. of 100) 0.2cc pipettes 29-88632

Well-Prep EDTA Root Canal Cleaner (Vericom)

Compare to: RC Prep (Premier)

EDTA Root Canal Cleaner

Well prep containing ethylenediamine Tetraacetic acid (EDTA) As a chelating agent allows efficient cleaning and easier preparation of the root canal.



Kit: 3 Sys x 6 gm with Tips 29-87701

Irrigating Solutions

CanalClean (Vista)

CanalClearnTM Irrigation Kit

•All your endo solutions in one easy kit

- •Optimizes the time spent on irrigation
- •No mess... Less time
- **5 Easy Steps to Improved Outcomes:**
- 1 Chlor-XTRATM 16ml
- 2 EDTA 6ml
- 3 Chlor-XTRATM 10ml
- 4 Rinse-N-DryTM 6ml
- 5 CHX PlusTM 10ml
- 5 Colored Coded Syringes

Anterior Anatomy, Kits: AA contains: 29-87731 1 (12 cc) Chlor-XTRATM syringe, 1 (3

cc) EDTA syringe, 1 (3 ml) Rinse-N-Dry syringe, 1 (3 ml) CHX-PlusTM syringe, and 4 (30 g) side vent irrigating tips

Complex Anatomy, Kits CA contains: 29-87732

1 (16 cc) Chlor-XTRATM syringe, 1 (6 cc) EDTA syringe, (10 cc) syringe Chlor-Xtra, 1 (6 cc) syringe Rinse and Dry (95% Ethyl Alcohol), 1 (10 cc) syringe CHS-Plus, 5 (30 g) Irrigation needles.

Chlor-Xtra Sodium Hypochlorite (Vista)

The perfect solution for the irrigation, debridement and disinfection of root canals during and after instrumentation

- •4X Faster than ChlorCid®*
- •Dissolve more tissue
- •Improve Disinfection
- •Reduce post-treatment complications

Chlor-XTRATM is a 6% Sodium Hypochlorite formula; enhanced with powerful wetting agents and proprietary





surface modifiers that penetrate into lateral canals and isthmuses by lowering surface tension, The result is a much more powerful oxidizing agent, 2+ times thinner than standard Sodium Hypochlorite and 2.6 times more digestive.

- •2.5 X more wetting ability**
- •2.6 X more digestive**
- •Much more powerful oxidizing agent

 16oz Bottle
 29-87671

 Prefilled Syringe 3cc, 10/Pkg
 29-87612

 Prefilled Syringe 12cc, 10/Pkg
 29-8767123

Sodium Hypochlorite (Vista)

For the debridement and irrigation of root canals during and after instrumentation compare to Ultradent®
ChlorCid®*



3% 16oz (480mL) - 1 Bottle 29-98741 **6% 16oz (480mL) -** 1 Bottle 29-98742

Sodium Hypochlorite 6% (USA)

For the debridement and irrigation of root canals during and after instrumentation. **Gluten Free**



6% Sodium Hypochlorite Solution -

500ml Bottle

6% Enhanced Sodium Hypochlorite Solution - 500ml Bottle

29-109362

29-109361

Irrigating Syringes & Needles

Appli-Vac Bendable Needle Tips (Vista)

1" bendable needle tip Blunt-end, side-port. Universal Luer-style designed. Color -coded for easy ID.



23 Ga, 100/Pkg	29-87751
23 Ga, 20/Pkg	29-87752
27 Ga, 100/Pkg	29-87753
27 Ga, 20/Pkg	29-87754
30 Ga, 100/Pkg	29-87755
30 Ga, 20/Pkg	29-87756

AutoSyringe Irrigation Handpiece



- Go with the flow, controlled flow
- Comfort, Convenience, & Control
- Eliminate hand fatigue associated with the standard syringe irrigation
- Easy to fill, quick

change 20ml canisters

- Safe, controlled flow of irrigating solutions with variable flow rate settings
- Save you MONEY, time, hand fatigue and Waste

Complete Kit: AutoSyringe Handpiece, 3x20ml Canisters, 20 Canister Liners, 3 Autoclavable Sleeves, 100 Plastic Barrier Sleeves & Changing Base

29-108641

Bendable Irrigation Needles (Dentazon)

- •1" multi-purpose bendable irrigation needle tip
- •Excellent flexibility for easy access to any canal
- •Use for root canal irrigation and periodontal pockets
- •Universal Luer-Lock design
- Color coded
- •Non-sterile

 18G: Pink, 1.2 X 25mm, 100/Pkg
 29-11175

 19G: Cream, 1.1 X 25mm, 100/Pkg
 29-11176

 20G: Yellow, 0.9 X 25mm, 100/Pkg
 29-11177

 22G: Black, 0.7 X 25mm, 100/Pkg
 29-11178

 25G: Orange, 0.5 X 25mm, 100/Pkg
 29-11179

 27G: Medium Gray, 0.4 X 25mm, 100/Pkg
 29-11180

Canal Clean Side-Port Irrigation Needles Tips (Dentazon)





Close-end tip and side-port opening

(€₀₄₉₉ ISO9001 ISO13485

Canal Clean irrigation needle can greatly
Increase your patient's safety on root canal and gingival sulci
irrigation treatments. Canal Clean's special design (closed-end tip
& side-port opening) allows you to prevent possible root injury
caused by unexpected forceful injection with conventional
irrigation needle. As a result, secondary inflammation or
flare-up may be also caused by polluted materials. The Canal
Clean is a reliable prevention for such clinical difficulties.

28G: Red, 0.36 X 25mm, Root Canal

Size #40, 25/Pkg 29-11193

30G: Yellow, 0.30 X 25mm, Root Canal

Size #30, 25/Pkg 29-11194



Prestige Dental Products, Inc. Prestigedentalproducts.com Tel: 877-772-3888

CanalClean (Vista)

CanalClearnTM Irrigation Kit

•All your endo solutions in one easy kit

•Optimizes the time spent on irrigation

•No mess... Less time

5 Easy Steps to Improved Outcomes:

- 1 Chlor-XTRATM 16ml
- 2 EDTA 6ml
- 3 Chlor-XTRATM 10ml
- 4 Rinse-N-DryTM 6ml
- 5 CHX PlusTM 10ml
- 5 Colored Coded Syringes

Anterior Anatomy, Kits: AA contains: 29-87731

1 (12 cc) Chlor-XTRATM syringe, 1 (3 cc) EDTA syringe, 1 (3 ml) Rinse-N-Dry syringe, 1 (3 ml) CHX-PlusTM syringe, and 4 (30 g) side vent irrigating tips

Complex Anatomy, Kits CA contains: 29-87732

1 (16 cc) Chlor-XTRATM syringe, 1 (6 cc) EDTA syringe, (10 cc) syringe Chlor-Xtra, 1 (6 cc) syringe Rinse and Dry (95% Ethyl Alcohol), 1 (10 cc) syringe CHS-Plus, 5 (30 g) Irrigation needles.

Endo Irrigation Combo Kit (Pac-Dent)



Monoject style needle tips

 27 Ga, YELLOW, 100/Pkg
 29-87801

 23 Ga, BLUE, 100/Pkg
 29-87802

 30Ga, PURPLE, 100/Pkg
 29-87803

Evacu-N-Irrig Needle Tips (PacDent)



PINK, 100/Pkg	29-87791
GRAY, 100/Pkg	29-87792
YELLOW, 100/Pkg	29-87793
BLACK, 100/Pkg	29-87794
BLUE, 100/Pkg	29-87795


Fill Station (Vista)

The simples, most convenient way to fill your syringes: No mess, No waste, Less time

- •Features syringe-activated valve and low level light indicator
- •Place several units together for convenient filling of multiple solutions
- •Easy to wall mount or use as tabletop modular unit

Available with HEAT for warm delivery of irrigating solutions

Heating hypochlorite solutions greatly increases tissue dissolution; the magnitude of the increase compared with RT [room temperature] varied from 30%



Fill Station TM (no heat)	29-98731
Fill Station TM with heat	29-98732
Empty Bottle	29-98733
Plastic Stopper	29-98734
Replacement Cap	29-98735

Irrigating Syringe W/Luer Lock (TCDental)

Accepts all luer lock based needle. Also applicable for use in surgical site irrigation. Non sterile and disposable. Graduated for 10cc.



Graduated for 10cc, 50/Box

29-87601

Irrigation Needles Tips- Double Side Vented (Wylke)



With this high quality, double side vented, luer lock irrigation needle, you can be sure to get an effective and safe irrigation when performing and doing treatment. Easy and optimized to clean out

the canals with the double side openings, creating a "swirl-effect". Safe with the closed tip, protecting apex.

Grey 27 gauge .040 X 25mm 29-26081 Violet 31 gauge .028 X 25mm 29-26082

Luer Lock Syringes

Disposable, Non-Sterile



12cc, 100/Box 29-98281



6cc , 100/Box	29-98282
3cc , 100/Box	29-98283

Luer-Lock Endo Irrigation Syringes (PacDent)



- •Ideal for surgical site irrigation and compatible with all luer style needles and hubs
- •Disposable, non-sterile.

Caution: Not suitable or intended for use with injection needle tips.

P	
12cc , 100/Box	29-15581
20cc, 50/Box	29-15582
3cc, black, w/easy-twist cap, 100/Box	29-15583
1.2cc, black, 100/Box	29-15584
1.2cc, white, 100/Box	29-15585
3cc, 100/Box	29-15586
6cc, 50/Box	29-15587

Micro Aspirators (Vista)

- •Universal high-speed, low-speed valve fitment available on autoclavable style
- •Provides for fast and efficient micro aspiration
- Virtually eliminates the need for paper points
- •Bendable needle that resists crimping

20g Luer HVE Adap. 20/Pkg 39-24851

Monovac Irrigation Needle Tips

Disposable, Comparable to Monoject® style



29-98341
29-98342
29-98343

www.prestigedentalproducts.com

Monovac Irrigation Syringes with Tips



3cc Syringe, 23 Ga Vented Tip, Monoject Style, Blue, 100Pcs/Box

29-98261

3cc Syringe, 27 Ga Vented Tip,

Blue, 100Pcs/Box 27Ga. Close End, Max-I-ProbeTM Style, 29-98331

PROBEVAC

Monoject Style, Yellow, 100Pcs/Box

29-98262 Yellow, 100Pcs/Box 29-98332

3cc Syringe, 30 Ga Vented Tip,

30Ga. Close End, Max-I-ProbeTM Style, 29-98263

29-98333

Monoject Style, Purple, 100Pcs/Box

Bendable needle 1" monoject style (notched) needle tips, color coded for

Notched Irrigating Needle Tips (PacDent)

Purple, 100Pcs/Box

Probevac Irrigation Pre-Tipped Syringes

23Ga. Close End, Max-I-ProbeTM Style,



Prestige Dental Products, Inc. Prestigedentalproducts.com

Tel: 877-772-3888

Probevac Irrigation Needle Tips

Disposable, Comparable to Max-I-ProbeTM style

3cc Syringes, 23Ga Close End Tip, Max-I-ProbeTM Style, Blue, 100Pcs/Box

29-98251

29-87742

3cc Syringes, 27Ga Close End Tip, Max-I-ProbeTM Style, Yellow, 100Pcs/Box

3cc Syringes, 30Ga Close End Tip, Max-I-

29-98252

29-98253

30 Ga, 100/Pkg

pockets.

Blue 23G, 100Pkg/Tips

Yellow 27G, 100Pkg/Tips

Purple 30G, 100Pkg/Tips

Blue 23 Ga, 100/Pkg

Yellow 27 Ga, 100/Pkg

easy ID

29-87743

with a closed-end

29-26251

29-26252

29-26253

29-87741

ProbeTM Style, Purple, 100Pcs/Box **Prope Irrigation Tips (Vista)**

•1" Bendable needle tip for irrigating subgingival surgical sites and sulcus •Univ Luer-style design. Round, closed end with side-port delivery

Secure Lock Bendable Tips (Vista)



23 Ga, 100/Pkg 29-87771 23 Ga, 20/Pkg 29-87772

27 Ga, 100/Pkg 29-87773

27 Ga, 20/Pkg 29-87774

Bendable

Pre-Tipped Syringes with Probe Irrigating Tips (Vista)

cannula and a side window dispersal point, the OptiProbeTM irrigation probe offers a unique design for unmatched effectiveness, patient comfort, and safety in irrigating root canals and periodontal

OptiProbe Side Open Closed End Needle Tips (PacDent)

Closed end side port needle tip similar to the Max-I-Probe style

3cc Syringes w/23Ga Probe 29-87821

Irrigation Tips, 100/Pkg

3cc Syringes w/27Ga Probe 29-87822

Irrigation Tips, 100/Pkg

BROWN 19 Ga, 80/Pkg

29-87783

YELLOW 20 Ga, 80/Pkg

29-87785

Call us for Quarterly Promotion

•Same design as Appli-VacTM 3/4"

Secure-Lock $^{\text{TM}}$ design to increase safety, ease of use and reduce risk of

cross-contamination.

•Universal luer lock style hub with new

BLACK 22 Ga, 80/Pkg	29-87787
BLUE 25 Ga, 80/Pkg	29-87789

Secure-Lock Irrigating Tips (Vista)

- •Same tip design as Appli-VacTM Irrigating
- •New locking hub design
- •Easier to use, increases safety and reduces cross-contamination
- •For irrigation of canals, pocket and fistulas

BLUE 23 Ga, 80/Pkg	29-87761
YELLOW 27 Ga, 80/Pkg	29-87763
PURPLE 23 Ga, 80/Pkg	29-87765

Syringes with Appli-Vac Irrigation Tips (Vista)

- •Irrigating Needle Tips. 1" bendable needle tip
- •Blunt-end, side-port
- •Universal Luer-style design
- •Color-coded for easy ID
- •Comparable to Monoject® Style




3cc Syringes with 23Ga Appli Vac Irrigation Tips, 100/Pkg	29-87811
3cc Syringes with 27Ga Appli Vac Irrigation Tips, 100/Pkg	29-87812
3cc Syringes with 30Ga Appli Vac Irrigation Tips, 100/Pkg	29-87813
12cc Syringes with 23Ga Appli Vac Irrigation Tips, 100/Pkg	29-87814
12cc Syringes with 27Ga Appli Vac Irrigation Tips, 100/Pkg	29-87815
12cc Syringes with 30Ga Appli Vac Irrigation Tips, 100/Pkg	29-87816

Therma-Flo Step Down Tips (Vista)

Extended Length for precision placement of composite materials.

Step Down tips offer a completely new delivery option for microscope enhanced dentistry and for use with high magnification loupes



compules, our patent pending Step Down tips combine the ultimate in convenience with precise placement of composite materials.

Features:

- •Narrow profile for optimal visualization of work area
- •Extended length for deeper access
- •Narrow profile allows for optimal visualization



Prestige Dental Products, Inc. Prestigedentalproducts.com Tel: 877-772-3888

•Designed for easy extrusion of highly-filled composites

•Available in 2X and 4X reduction for precision placement

•Orifice reduction for greater precision during micro-restorative procedures

Developed for use with the Therma-FloTM Composite Warming Kit, but can be used with any standard compules or compule guns.

Step Down Tips 2x, 40Cup/Pkg 29-26211 Step Down Tips 4x, 40Cup/Pkg 29-26212

Utility Curved Syringes

- Disposable
- •Non-sterile



29-98291 12cc, 50/Box

Utility Curved Syringes (TC Dental)

- •Curved tip provides easy access for maximum canal irrigation
- •Tight fit plunger and see through
- •10 cc and non sterile



29-11748

Irrigating Syringe w/Curve Tips,

10cc 50Pcs/Box

Vpro EndoSafe (Vista)

•Non-clogging VproTM EndoSafeTM is a simple negative pressure, simultaneous irrigation and evacuation device

•When used properly, it

provides safe, easy and cost effective apical irrigation and evacuation.

10/Pkg 29-87721

Calcium Hydroxide Treatment

Calasept Plus Calcium Hydroxide (Wykle)

The most effective calcium hydroxide thanks to its high concentration and optimal calcium release. CALASEPT Plus has a very high content of calci-

um hydroxide (> 41%) giving it a high concentration of hydroxyl ions.

These generate a high pH



value of 12.4 which has a pronounced bactericidal effect. The high concentration of calcium hydroxide in the paste provides for a long lasting release of calcium ions. CALASEPT Plus is a ready-to-use paste in air tight syringes for direct application through the Flexi-needle. CALASEPT Plus is radio-opaque.

Kit: 4 syringes each containing 1,5 ml CA-LASEPT Plus & 20 CALASEPT Plus Flexineedles, luer lock

29-26071

orendo Past

1 syringe containing:

1,5 ml CALASEPT Plus & 5 CALASEPT

Plus Flexi-needles, luer lock

29-2607

Forendo Paste (PulpDent)

Calcium hydroxide paste with iodoform in a silicone oil base

Bactericidal - Bacteriostatic pH > 12 - Radiopaque -Non-Setting

For temporary and intermediate root canal

therapy

Forendo Paste is a strong antibacterial dressing for routine use between office visits to disinfect root canals and for complicated cases when extended calcium hydroxide therapy is indicated.

- Disinfects canals
- •Inhibits bacterial growth
- •Prevents flare-ups
- •Promotes periapical healing
- •Treats infected root canals, abscesses and periapical lesions
- •Treats traumatic injuries, root restoration, fractures and perforations
- •Treats pus, hemorrhage, exudation and weeping canals
- Stimulates apexification and apexogenesis

2.2 gm syringe + 20 applicator tips

29-11523

Irrigation Needles Tips- Double Side Vented (Wylke)



With this high quality, double side vented, luer lock irrigation needle, you can be sure to get an effective and safe irrigation when performing and doing treatment. Easy and optimized to clean out

the canals with the double side openings, creating a "swirl-effect". Safe with the closed tip, protecting apex.

Grey 27 gauge .040 X 25mm 29-26081

Violet 31 gauge .028 X 25mm 29-26082

Metapaste (Meta Biomed)

Calcium Hydroxide with **Barium Sulfate Water**

Soluble Temporary Root Canal Filling Material (Compare to Tempcanal from Pulpdent)

Applications:

- •Weeping canal
- •Formation of a hard tissue barrier(Apexification)
- •Temporary root canal filling
- Internal & external root restoration

Advantages

- •Easy cleaning and removing with good water solubility
- •Excellent antibacterial effect and radiopacity
- •Premixed paste
- •Syringe type for easy delivery of the paste into the canal
- •Disposable tips for prevention of cross-contamination

Metapaste (w/Barium Sulfate & Water

Soluble) 2.2g Syringe with 20 Tips 29-11985



Prestige Dental Products, Inc. Prestigedentalproducts.com Tel: 877-772-3888

Pulpdent Paste (PulpDent)

The Original, Non-Setting, Pre-mixed Calcium Hydroxide Paste

For Root Canal Therapy and Vital Pulp Therapy

Bactericidal - Bacteriostatic - Radiopaque - Non-Setting Internationally acclaimed

as the world's finest treatment for:

- Vital Pulp Therapy
- •Pulpal Healing
- •Dentin Bridge Formation
- •Root Canal Therapy



Pulpdent Paste is the original premixed calcium hydroxide methylcellulose pulpal dressing. It is indicated for direct pulp capping, pulpotomy and pulpal curettage procedures, and it is used routinely as an antimicrobial intracanal dressing between office visits to disinfect root canals and prevent flare-ups.

Pulpdent Paste promotes favorable periapical healing and is the treatment of choice in complicated cases involving infection, root restoration, traumatic injuries and other instances where calcium hydroxide therapy is indicated.

When used for vital pulp therapy, Pulpdent Paste stimulates consistent dentin bridge formation

The new dentin bridge is visible clinically after approximately 2-3

Kit: 3mL syringe + 24 needles (18 gauge

blunt, sterile) 29-12277

3 mL syringe 29-12278

Syringe needles, 24/Pkg 29-12279

PULPDENT

TempCanal

TempCanal (PulpDent)

Temporary Calcium Hydroxide Root Canal Treatment

Paste Bactericidal - Bacteriostatic - Radiopaque

 Non-Setting TempCanal has a pH > 12

and provides the full benefits of calcium hydroxide in root canal therapy.

For routine use as an antimicrobial intracanal dressing between office visits to disinfect the canal and prevent flare-ups.

For complicated cases involving infection, traumatic injury and other instances where calcium hydroxide therapy is indicated.

- **Routine Treatment Between Office Visits:**
- •pH > 12 disinfects canal
- Prevents flare-ups
- •Replaces hazardous formocresol
- •Easily removed with file and irrigation

Treatment of Complicated Cases

- •Promotes favorable periapical healing
- •Treats infection, abscesses and periapical lesions
- •Treats traumatic injuries, avulsed and luxated teeth
- •Treats root resorption, root fractures and perforations
- Stimulates apexification and apexogenesis
- •Treats pus, hemorrhage exudation and weeping canals Preferred Delivery System
- Screw syringe dispensing system
- •Flows through sterile 22- and 25-gauge needles directly into canal
- •No mixing, no condensing
- •Easily removed from canal with file and irrigation
- •Radiopaque



Kit: 3mL syringe + 12 needles 29-12472

Well-Paste Calcium Hydroxide Treatment (Vericom)



- Oil-Based calcium hydroxide paste with lodoform
 Root canal filling material for the infected root canal
- Very stable without any solidification or separation and highly flowable paste premixed in a convenient syringe

- has **EXCELLENT** accessibility to the apex and

biocompatibility to the periapical tissue.

2g Syringe w/ tip & extension holder 29-108101

Endo Boxes & Organizers

Autoclavable Endo Aid (PacDent)

One-finger-held autoclavable organizer gives dental clinicians complete control over endodontic files





•The organizer is made of fully autoclavable material

- •Effectively organizes endodontic files
- •One-finger-held organizer with ruler ensures convenience

Endo aid kit: 1 organizer with integrated ruler, 2 rings to fit all finger sizes and 50 disposable triangular foam refills 26-25591 25 vellow triangular foam refills, 25 teal triangular foam refills 26-25592 50 Triangle foam refills, Gray 26-25593 50 Triangle foam refills, Teal 26-25594 50 Triangle foam refills, Yellow 26-25595 50 Triangle foam, Lavender 26-25596 Endo organizer only 26-25597 50 Rectangular foam refills, Gray 26-25598

Endo Clean Stand & Foam



White Autoclavable 29-11299
Blue Autoclavable 29-11298

Round Endo Foam Blue, Disposable,

48Pcs/Bag 29-11300





Endo Foam Inserts (Plasdent)



White, Disposable, 48Pcs/Bag

29-11309

Endo Micro Pulpa Box (Union Broach)



Aluminum 30 Holes

 SILVER Color
 29-108481

 RED Color
 29-108482

 BLUE Color
 29-108483

Endo Micro Pulpa Round (Miltex)



Small aluminum box holds up to 18 endodontic instruments

 SILVER Color
 29-108481

 RED Color
 29-108482

 BLUE Color
 29-108483

Endo Module (Miltex)



- Small & flat box dedicated to the storage of paper points and Gutta Percha points, but you can also store your endodontic insturments and/or drills.

- Endomodule is made of a box fitted with a 6 compartment tray

 SILVER Color
 27-108511

 RED Color
 27-108512

 BLUE Color
 27-108513

 GREEN Color
 27-108514

Endo Ruler



Stainless Steel, 3" Long 29-107791

Endo Ruler (Pac-Dent)

Autoclavable, stainless steel construction, one-finger-held endo ruler.

- Millimeter scale



1-Endo ruler per box

29-25601

Endo-Cushions hold (IndiGreen)



Disposable Endo-Cushions hold, arrange and clean small instruments during endodontic or implant procedures. Ideal for files, posts, screws, drills and other small items.

Peel-and-stick back secures to patient napkins, tray covers or anywhere handy. Pre-arranging saves set-up time. Dense, grey polyurethane foam.

3"long 50 cushions per bag

29-87921




Endodontic Work Station29-110951File Stand, For Endo Work Station29-110952Foam, For File Stand29-110953Clear File Organizer Set29-110954

Endoring Filecaddy (Jordco)

The Endoring FileCaddy organizer provides a unique bulk storage solution for all endodontic files.

The clinician places multiple files securely into the special foam insert, closes the lid and steam autoclaves the entire FileCaddy organizer—no bag or wrap required. Single handed operation allows easy opening of the protective li



Files are easily identified, selected and transferred to the companion Jordco Endoring II organizer for the current endodontic procedure. The FileCaddy organizer disassembles easily to allow replacement of the foam insert.

Kit: 1 FileCaddy and 12 single use refill
combos29-8987112 FileCaddy single use Blue Foam Refills29-8987212 FileCaddy single use Yellow Foam
Refills29-89873



Endoring II Organizer Gelwells (Jordco)

The Endoring II organizer enables the clinician to directly place,



store, measure and clean endodontic hand and rotary instruments within the operating field

It reduces the risk of crosscontamination and minimizes the passing of sharps between dentists and staff.

Using the Endoring II organizer helps make endodontic procedures safer and faster.

Single Endoring II Premium Kit: 1 En-

doring II assembly (blue), 1 standard metal ruler, 1 e-Foam insert combo pack (48 inserts), 1 pack Gelwells (48 single-use cups) and 1 EndoGel pack (12 1ml MicroDose tubes–USA kits only)

29-89881

E-Foam Blue Pack:

Single package of 48 foam inserts 29-89882

E-Foam Yellow Pack:

Single package of 48 foam inserts 29-89883

Jumbo Steri Box



4 9/16" x 3 3/16" x 2 3/16

 Bur Box, 160 FG Bur
 26-106741

 Endo Box, 144 Endo Instruments
 26-106742

 Bur & Endo Box, 70 Burs, 70 Endo, 4
 26-106743

 Vials
 Vials

Multi Box, Multi-Function 26-106744

Mini Clean Gripp (Miltex)



 SILVER Color
 29-108491

 BLUE Color
 29-108492

Call Us Today 833-733-6825

Mini Endo Organizer Box (Miltex)



SILVER Color, 13 x 5 x 4cm, holds up to 60 Endodontic Finger instruments

29-108501

Endo Burs

Endo Carbide Burs (PacDent)

•151Z can be used for initial endodontic access, sectioning, root seperation, etc. •152EZ, features non -cutting safe-end tip that prevents




penetration of pulp chamber floor, enables easy access to root canals, comparable to Maillefer Endo-Z bur.

151Z 3Pcs/Pkg 19-25042

152EZ 3Pcs/Pkg 19-25043

Endo Diamond Burs (PacDent)

The endo diamond bur's round tip allows access to initial penetration while its



taper shaped shaft flares the pulp chamber. The non-cutting safe tip diamond bur protects chamber floor while removing the pulp chamber roof. The cutting surface is 10mm.

 21mm Total Length, 3Pcs/Pkg
 19-25054


 21mm Total Length, 3pcs/Pkg
 19-25055

 21mm Total Length, 3Pcs/Pkg
 19-25056

 21mm Total Length, 3Pcs/Pkg
 19-25057

Ultrasonic Files

Piezo Endo Chuck tips



E1 120° Angle each 29-11294 (Endo File not included)

E2 90° Angle each 29-11295

(Endo File not included)

Prestige Dental Products, Inc. Prestigedentalproducts.com Tel: 877-772-3888

Ultrasonic NiTi Files

Color Coded used by Ultrasonic scaler units for root canal cleaning




Tip w/excellent tenacity is used for root canal cleaning

File size: 32.5mm±0.5mm Handpiece size: Φ0.68mm~0.71mm K-file taper: 2%

#10, 33mm, 6/Pkg	29-12247
#15, 33mm, 6/Pkg	29-12246
#20, 33mm, 6/Pkg	29-12245
#25, 33mm, 6/Pkg	29-12244
#30, 33mm, 6/Pkg	29-12243
#35, 33mm, 6/Pkg	29-12242
#40, 33mm, 6/Pkg	29-12241

Endodontic Instruments


Broach Removal Pliers (Nordent)



#3 Locking College Plier Locking w/

Grooved Tips (6"/150mm) 40-11293 **Angular** (5 1/2"/ 140mm) 40-11294 **Straight** (5 1/2"/ 140mm) 40-11292

Endo Forceps Stieglitz (Nordent)



45 Degrees 40-11744 **90 Degrees** 40-11743



D/E 5/7, Root Canal Plugger	40-91426
D/E DG16, Explorer	40-91427
D/E Glick 1, Carver	40-91428
S/E D11 T-21, Endo Spreader	40-91429
S/E D11 T-25, Endo Spreader	40-91430


Endodontic Excavator (Nordent)

	A selection of head sizes with Extra-Long	Terminal Shanks
1	#11L	
1	#12L	
	#31L	
	#33L	
#11L Extr	a Long Shank Double End	1

	1
#33L	
	7
#11L Extra Long Shank Double End Spoon DuraLite¬Æ ColorRings,Ñ¢	40-11768
#11L Extra Long Shank Double End Spoon DuraLite® ROUND	40-11754
#11L Extra Long Shank Double End Spoon DuraLite® HEX	40-11753
#11L Extra Long Shank Double End Spoon Standard	40-11751
#12L Extra Long Shank Double End Spoon DuraLite¬Æ ColorRings,Ñ¢	40-11752
#12L Extra Long Shank Double End Spoon DuraLite® ROUND	40-11750
#12L Extra Long Shank Double End Spoon DuraLite® HEX	40-11749
#12L Extra Long Shank Double End Spoon Standard	40-11748
#31L Long Shank Double End Blade DuraLite¬Æ ColorRings,Ñ¢	40-11747
#31L Long Shank Double End Blade DuraLite® HEX	40-11746
#31L Long Shank Double End Blade DuraLite® ROUND	40-11745
#31L Long Shank Double End Blade Standard	40-11755
#31LR Long Shank Double End Spoon DuraLite¬Æ ColorRings,Ñ¢	40-11756
#31LR Long Shank Double End Spoon DuraLite® ROUND	40-11767




1cl. 6//-//2-3888	
#31LR Long Shank Double End Spoon DuraLite® HEX	40-11766
#31LR Long Shank Double End Spoon Standard	40-11765
#32L Long Shank Double End DuraLite¬Æ ColorRings,Ñ¢	40-11763
#32L Long Shank Double End DuraLite® ROUND	40-11762
#32L Long Shank Double End DuraLite® HEX	40-11761
#32L Long Shank Double End Standard	40-11760
#33L Long Shank Double End Spoon DuraLite¬Æ ColorRings,Ñ¢	40-11759
#33L Long Shank Double End Spoon DuraLite® ROUND	40-11758
#33L Long Shank Double End Spoon DuraLite® HEX	40-11757
#33L Long Shank Double End Spoon Standard	40-11764
Explorers (Nordent)	
#DG16	





#5-7 Double End Standard	40-13547
#9-11 Double End DuraLite® Color- Rings™	40-13546
#9-11 Double End DuraLite® ROUND	40-13545
#9-11 Double End Standard	40-13544

Spreader NiTi (Nordent)



4SP Single End Nickle Titanium DuraLite® ColorRings™	40-13543
4SP Single End Nickle Titanium DuraLite® ROUND	40-13530
4SP Single End Nickle Titanium Standard	40-13531
#D11 Single End Nickle Titanium DuraLite® ColorRings TM	40-13532
#D11 Single End Nickle Titanium DuraLite® ROUND	40-13533
#D11 Single End Nickle Titanium Standard	40-13535
#D11T (Thin) Single End Nickle Titanium DuraLite® ColorRings™	40-13534
#D11T (Thin) Single End Nickle Titanium DuraLite® ROUND	40-13536
#D11T (Thin) Single End Nickle Titanium Standard	40-13438
#D11T Single End 25MM Nickle Titanium DuraLite® ColorRings TM	40-13537
#D11T Single End 25MM Nickle Titanium DuraLite® ROUND	40-13539
#D11T Single End 25MM Nickle Titanium Standard	40-13540
#D11T Single End 28MM Nickle Titanium DuraLite® ColorRings™	40-13541
#D11T Single End 28MM Nickle Titanium DuraLite® Standard	40-13529
#D11T Single End 28MM Nickle	40 12542





40-13542

Titanium DuraLite® ROUND

Order Form

Date:												
Ordered By					Delive	er To			Same as Above			
Business Nam	Business Name:					Business Name:						
Address:					Address:							
State/Province	e:				State/Pro	vince:						
Zip/Postal Cod	de:			Zip/Postal			:					
Phone:					Phone:							
Fax:					Fax:							
Contact Name	: :				Contact N	Vame:						
Item Number		D	escr	iption		Qua	antit	y	Unit Price		Amount	
									1			
									 			
	<u> </u>								<u> </u>			
									Sub Total			
Payment				○ Credit (Card				C.O.D. Charge	es		
^ ^ "	~		- \) American Express) MasterCard				Freight Charge	es		
	ne Pay	ment (Bank to Ba	nk)	O ?					Sales Tax			
E-mail					Visa				Total	\perp		
○ Check	ı, D _{ava}	shla Ta		~	_	_						
	A Luj	DIE 10			Number:	-						
					tion Date:	_						
				Carum	older Nam	e:						
Dental Licer	nse #			Year Licensed			7	Exp	piration Date			
DEA#				Year Licensed		,	7	Exp	piration Date			



Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

