

# Handpieces



**Prestige Dental Products, Inc**

**[Prestigedentalproducts.com](http://Prestigedentalproducts.com)**

**Tel: 877-772-3888 Fax: 714-632-8821**

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## High Speed Handpieces - Non Fiber Optic

### Cartridge - Turbine for Handpiece (ND)

#### CARTRIDGE FOR AIRTURBINE HANDPIECE



**CRT-450**  
Cartridge for  
TCP-450M / -450B  
TCP-70QM / -70QB

**CRT-350**  
Cartridge for  
TCP-350M / -350B

**CRT-40P**  
Cartridge for  
TC-40PM / -40PB  
TC-80QM / -80QB

**CRT-35Y**  
Cartridge for  
TC-35YM / -35YB

**CRT-300Z**  
Cartridge for  
TS-300ZT

#### Made in Japan

**Standard Chuck**, for TC-40PM / 40PB  
& TC-80QM / 80QB 35-11176

**Large Size, Standard Chuck**,  
for TC-35YM / 35YB 35-111751

**Push Button Chuck**, for TC-450M /  
450B & TCP-70QM / 70QB 35-111741

**Push Button Chuck**,  
for TC-350M / 350B 35-111731

**Cartridge for TS-300ZT** 35-111732

### Chuck Style Handpiece (Dentex)

Chuck Style Fiber  
Optic SINGLE  
PORT WATER  
SPRAY Hand-  
piece (Dentex),  
Fixed Back End



**Standard Head Non-Optic** 35-111821

**Mini head Non-Optic** 35-111811

### E-Type Straight Handpiece (Nakamura)

#### ES-30A: Straight Handpiece U-nose



**ES-30A** 35-112501

### Evolve High-Speed Handpiece - Neovo

AG Neovo Dental handpieces  
deliver industry-leading cutting efficiency  
and long-lasting precision.



- **Quiet & comfortable** operation, noise is reduced to the minimized level
- ceramic ball bearings, increased durability and reliability are guaranteed by ceramic ball bearings.
- Quattro water spray 4-nozzle water spray ensures efficient cooling and brings comfort to the entire operation
- **Light weighted** design ergonomic yet light weight

body provides ease and comfort throughout each operation  
- **Highly durable** and repeatedly autoclavable at up to 135°C

**Kavo Type Standard**, Head Evolve AG Neovo Highspeed push button swivel type handpiece w/ Quad Water Spray 35-105431

**NSK Type Standard**, Head Evolve AG Neovo Highspeed push button swivel type handpiece w/ Quad Water Spray 35-105432

**Kavo Type Mini**, Head Evolve AG Neovo Highspeed push button swivel type handpiece w/ Quad Water Spray 35-105433

**NSK Mini**, Head Evolve AG Neovo Highspeed push button swivel type handpiece w/ Quad Water Spray 35-105434

**Kavo Type Evolve series swivel AG Neovo Highspeed Quick Connect Swivel** 35-105435

**NSK Type**, Evolve series Swivel AG Neovo Highspeed Quick Connect Swivel 35-105436

### Handpieces with Quick Connect (ND)



#### Made in Japan

**Standard Chuck** 380,000rpm ISO Con 35-82351

**Push-Button (Auto)** Chuck 380,000rpm ISO Con 35-82352

**360° Quick Connect Swivel** for TC-80QM, TCP-70QM 35-82353

### Highspeed Basic Non-Push Button Handpiece (Alfa Air)

**Handpiece Non-Push** 35-117561

**Turbine** 35-117551

**Turbine Deluxe, Made in Japan** 35-11754

### Highspeed Basic Push Button Handpiece (Alfa Air)

**Handpiece Push Button, Made in Japan** 35-117581

**Turbine Deluxe, Made in Japan** 35-117571

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Highspeed Handpiece- Pana-MAX Push Button

**Made in Japan**



Each 35-16531

## Highspeed Handpiece - Push Button (ND)

Push-Button Chuck Type



**TCP-450M (4 holes)**  
**TCP-450B (2 holes)**  
Standard Size Head, 400,000rpm



**TCP-350M (4 holes)**  
**TCP-350B (2 holes)**  
Large Torque Head, 350,000rpm



**Quick-Joint Model**  
**TCP-70QM (4 holes)**  
**TCP-70QB (2 holes)**  
Standard Size Head, 380,000rpm

**Standard Size Head, 400,000rpm,** 35-11764  
4 holes

**Large Torque Head, 350,000rpm,** 35-117651  
4 holes

**Standard Size Head with swivel,** 35-117663  
380,000rpm, 2 holes

## Integrity Leverage Classic Highspeed with Push Button Cap (Osseo Scientific)



**Integrity Leverage Classic Highspeed** 35-117801  
with Push Button Cap

## Integrity Leverage Handpiece with Lever Cap (Osseo Scientific)

Like top brand handpieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best highspeed hand-piece performance and value at a fraction of the cost of the top leading brands.



- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip Grip design

**Handpiece with Lever Cap** 35-117841

## Integrity Push Button Handpiece (Osseo Scientific)

Like top brand handpieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best highspeed handpiece performance and value at a fraction of the cost of the top leading brands.




- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip Grip design


**Highspeed Handpiece Push Button** 35-117891

## Micromotor Handpiece and E-Type Micromotor (Nakamura)

2 LB-25HM



**LB-20EM**  
E-type Micromotor  
Speed Range : 3,000 - 20,000rpm



**LB-25HM**  
Micromotor Handpiece  
Speed Range : 3,000 - 20,000rpm

**LB-30HM** 35-118591

**LB-35EM** 35-118601

**LB-25HM** 35-118611

**LB-20EM** 35-118621

## Micromotor Handpiece Set for Laboratory (Nakamura)

LXC-25FS



**Made in Japan**

**LX-30FS Set Includes:** 35-118641

LX-35V, LB-30HM, AND FC-500

**LX-35FS** 35-118651

**LXC-20FS** 35-118661

**LXC-25FS Set Include:**LXC-12V, 35-118631  
LB-25H, FCC-150 and cradle



# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Portable Micromotor System (Nakamura)

BD-550 (Battery)



BD-20FS



**Set Includes:** BD-550, LB-20EM, 35-119171  
ES-30A, CHR-80A, Cradle

**Set Includes:** BD-550, LB-25HM, 35-119181  
CHR-80A, Cradle

**Battery-Pack Controller** 35-119191

## Prophy Angle U-Type Right Angle HP (Nakamura)

RA-10PS



**RA-10PS** 35-119321

**RA-10PK** 35-119331

**RA-50PS** 35-119341

**RA-50PK** 35-119351

## Push Button Handpiece (Dentex)

Handpiece 4-Hole and  
5-Hole Connector  
Non-Optic (Dentex),  
Push Button, Triple  
Port Water Spray,  
Fixed Back End, Highspeed



**Standard Head Non-Optic** 35-119521

**Mini Head Non-Optic** 35-119531

**Complete Turbine for M3-S** 35-119532

## Push Button Handpiece (Sable)



**ACCESS II:** PB Mini Head, Hand- 35-88661  
piece Only (No Coupler)

**ACCESS II:** PB Mini Head 360/Quick 35-88662  
Connect Coupler Included

**ROTAMAX II:** PB Standard Head, 35-88663  
Handpiece Only (No Coupler)

**ROTAMAX II:** PB Standard Head, 35-88664  
360 Quick Connect Coupler Included

**4 Hole 360 Quick Connect Coupler** 35-88665  
To Fit Access II & Rotamax II

**2 Hole 360 Quick Connect Coupler** 35-88666  
To Fit Access II & Rotamax II

## Push Button Handpiece Fit NSK Swivel Coupler (Dentex)

Push Button Handpiece  
Interchangeable with NSK  
Swivel Coupler (Dentex)  
Swivel Handpiece, Triple  
Port Water Spray



**Standard Head Non-Optic** 35-119581

**Mini Head Non-Optic** 35-119591

## Push Button Handpiece Fit Star Coupler (Dentex)

Push Button Highspeed Handpiece,  
Triple Port Water Spray,  
Swivel Handpiece  
Interchangeable with Star  
Coupler (Dentex)



**Standard Head Non-Optic** 35-119601

**Mini Head Non-Optic** 35-119611

## Push Button Handpiece Interchangeable with Kavo Coupler (Dentex)

Push Button Handpiece  
Interchangeable with  
Kavo Coupler (Dentex),  
Swivel Handpiece



**Standard Head Non-Optic** 35-119541

**Mini Head Non-Optic** 35-119551

## Push Button Handpiece Interchangeable with W&H Coupler (Dentex)

Push Button  
Handpiece  
Interchangeable with  
W&H Coupler  
model RA-24 &  
RA-25 (Dentex)  
Triple Port Water  
Spray, Swivel Handpiece



**Standard Head Non-Optic** 35-119621

**Mini Head Non-Optic** 35-119631



## Handpieces

### Push Button Handpiece M1 (Dentex)

Push Button SINGLE  
PORT WATER SPRAY  
4-Hole and 5-Hole  
Connector Handpiece  
(Dentex), Fixed Back  
End



**Standard Head Non-Optic** 35-119641

**Mini Head Non-Optic** 35-119651

### Push Button Handpiece Midwest Type (Dentex)

Push Button Handpiece  
Interchangeable with  
Midwest/Bien Air  
Coupler (Dentex),  
Swivel Handpiece  
Triple Port Water Spray



**Standard Head Non-Optic** 35-119561

**Mini Head Non-Optic** 35-119571

### Quest Handpieces (Lares)

#### **American Manufactured Highspeeds**


- Quest highsPEEDs feature low-maintenance, long-life, ceramic bearings. The Quest stainless steel exterior and subtle grip design allows for fast, easy cleaning and autoclaving/chemiclaving
- Quest highsPEEDs are priced for value-conscious dentists unwilling to compromise on quality and performance
- The Quest 100 Series features a compact head for easy posterior access
- These handpieces are fast cutting, and have exceptional power and visibility, easy posterior access
- The Quest 300 Series offers a full sized head for powerful torque and smooth, quiet operation
- Quest 300 models feature heavy-duty construction and multi-port water spray



**Push Button Ceramic Bearing: Small** 35-119701  
Head Handpiece, 4 Head, model 110

**Push Button Ceramic Bearing: Full** 35-119691  
Sized Head Handpiece, 4 Head, model 310

### Sheath of E-Type Contra Angle (Nakamura)



**SHS-EC (1:1) Midwest Type** 35-120381

**SHS-ECG (4:1) Midwest Type** 35-120401

**SHS-EGG (16:1) Midwest Type** 35-120391

**SHS-ECN (1:1) NSK or STAR Type** 35-120411

**SHS-EGN (4:1) NSK or STAR Type** 35-120421



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

### SHS-EGGN (16:1) NSK or STAR Type 35-120431

### Quest Push Button for Swivel Handpieces (Lares)



**100 Series**  
Power and Durability in a Compact  
Head

**300 Series**  
Power for Your Toughest Cutting  
Requirements



- Small, compact head for excellent visibility and accessibility
- High speed for fast cutting
- Reduced contra angle for easier posterior access
- Easy, secure push button bur changing
- Chair-side rotor and chuck replacement
- Optional bright twin beam fiber optics
- Optional 360° swivel and quick disconnect
- Autoclavable and chemiclavable
- Low maintenance Ceramic Bearings

- Full-size head for powerful torque
- Heavy duty construction
- Smooth, quiet operation
- Optional 360° swivel and quick disconnect
- Optional bright twin beam fiber optics
- Easy, secure push button bur changing
- Chair-side rotor and chuck replacement
- Autoclavable and chemiclavable
- Low maintenance Ceramic Bearings
- Multiple Port Water Spray

Quest 150 LO PB SW MX	Height 11.9 mm Speed 11320	Diameter 10.2 mm Water Spray Single Port	Weight 50 g Torque 15.3 watts
	450,000 rpm		
Quest 350 LO PB SW MX	Height 14.65 mm Speed 11322	Diameter 13 mm Water Spray Multi-Port	Weight 56 g Torque 20.8 watts
	360,000 rpm		

#### **American Manufactured highsPEEDs**

- Quest highsPEEDs feature low-maintenance, long-life, ceramic bearings. The Quest stainless steel exterior and subtle grip design allows for fast, easy cleaning and autoclaving/chemiclaving
- Quest highsPEEDs are priced for value-conscious dentists unwilling to compromise on quality and performance
- The Quest 100 Series features a compact head for easy posterior access
- These handpieces are fast cutting, and have exceptional power and visibility, easy posterior access
- The Quest 300 Series offers a full sized head for powerful torque and smooth, quiet operation
- Quest 300 models feature heavy-duty construction and multi-port water spray

**130 Quest Non Fiber Optic Small Head HP** 35-119801  
360 Swivel, Push Button, Multiflex Compatible, Ceramic Bearings, w/o Swivel

**330 Quest Non Fiber Optic Standard Head** 35-119791  
HP 360 Swivel, Push Button, Multiflex Compatible, Ceramic Bearings, w/o Swivel

**4 Holes Coupler,** 35-119781  
Multiflex Compatible Coupler

**5 Holes Coupler, Multiflex Coupler** 35-119771



Like us on Facebook

## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Sonic Air MM 1500 Endo System - Medidenta




All NiTi systems are **EXCELLENT** in creating crown down effects, but they do not irrigate or cleanse lateral canals. These canals are breeding grounds for bacteria, which in effect will increase post-op pain and inevitable retreatment.

- Requires lubrication after every four autoclaves
- Use only autoclave oil
- **Accepts:** Sonic Files - Risp Sonic, Shaper Sonic, Trio-Sonic, Retprep

**1 year warranty**

**MM 1500 Handpiece** 329-105421

### Spare Head with Elbow for Star-Titan Type Contra Angle (Nakamura)



MH-30BL Midwest Type	35-120821
MH-30TL Midwest Type	35-120831
MH-20L Midwest Type	35-120841
MH-20ML Midwest Type	35-120851
MH-20FG Midwest Type	35-120861
MH-20AP Midwest Type	35-120771
MH-20PS Midwest Type	35-120761
MH-20PK Midwest Type	35-12067
MH-50PS Midwest Type	35-12068
MH-50PK Midwest Type	35-120691
MH-35BL Midwest Type	35-120701
MHN-30BL NSK or STAR Type	35-120711
MHN-30TL NSK or STAR Type	35-120721
MHN-20L NSK or STAR Type	35-120731
MHN-20ML NSK or STAR Type	35-12074
MHN-20FG NSK or STAR Type	35-12075
MHN-20AP NSK or STAR Type	35-120871
MHN-20PS NSK or STAR Type	35-120811
MHN-20PK NSK or STAR Type	35-120801
MHN-50PS NSK or STAR Type	35-120791

**MHN-50PK NSK or STAR Type** 35-120781

**MHN-35BL NSK or STAR Type** 35-12066

### StarDental 430 SW Handpiece



StarDental 430 Torque has a small head, high torque design that improves overall accessibility and visibility. The high power enables fast and precise removal of tooth surface and amalgam to improve speed and quality of work. The low noise level reduces possibility of auditory damage and can reduce patient anxiety. The coaxial water spray optimizes the direction of the spray and water pressure, cooling the bur and flushing the operatory site

#### **Features & Benefits**

- 100% stainless steel body provides durability
- Electropolished or brushed finish
- Scallop handle eliminates grooves where debris & bacteria are harbored
- Coaxial

**430 SW Handpiece Non-Fiber Optic w/o swivel** 35-105781

**430 SW Lubricated, Non-Fiber Optic w/o swivel** 35-105782

**430 SW LF, Torque w/o swivel** 35-105783

**430 SW Lubricated, Torque w/o swivel** 35-105784

### StarDental 430K Handpiece



Star Dental 430 Torque has a small head, high torque design that improves

overall accessibility and visibility. The high power enables fast and precise removal of tooth surface and amalgam to improve speed and quality of work. The low noise level reduces possibility of auditory damage and can reduce patient anxiety. The coaxial water spray optimizes the direction of the spray and water pressure, cooling the bur and flushing the operatory site

#### **Features & Benefits**

- 100% stainless steel body provides durability
- Electropolished or brushed finish
- Scallop handle eliminates grooves where debris & bacteria are harbored
- Coaxial

**430K, Fixed Back End** 35-110641

**430K Lubricated, Fixed Back End** 35-110642

## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Table-Top Type Controller (Nakamura)



<b>LB-35V</b>	35-121141
<b>LXC-12V</b>	35-121151

### Triple Port Water Spray Handpiece (Dentex)

Chuck Style  
Handpiece 4-Hole and  
5-Hole Connector  
Non-Optic (Dentex),  
Triple Port Water  
Spray, Fixed Back End



<b>Standard Head</b>	35-111831
<b>Mini Head</b>	35-111841
<b>Turbine For M3-S Handpiece</b>	SKU83481

### High Speed Handpieces - Fiber Optic

#### Chuck Style Fiber Optic Handpiece (Dentex)

Chuck Style Fiber  
Optic SINGLE  
PORT WATER  
SPRAY Handpiece  
(Dentex), Fixed  
Back End, 4-Hole  
and 5-Hole Connector



<b>Standard Head Fiber Optic</b>	35-111771
<b>Mini Head Fiber Optic</b>	35-111781

#### Chuck Style Fiber Optic Handpiece with Triple Port Water Spray (Dentex)

Chuck Style  
Handpiece 4-Hole  
and 5-Hole  
Connector Fiber  
Optic (Dentex),  
Triple Port Water  
Spray, Fixed Back End



<b>Standard Head Fiber Optic</b>	35-111801
<b>Mini Head Fiber Optic</b>	35-11179

### Evolve High-Speed Handpiece - Fiber Optic (Neovo)

- Fiber optic  
preparation areas  
are optimally  
illuminated with  
superior  
visibility  
- **Quiet &  
comfortable** operation, noise is reduces to the minimized level  
- Ceramic ball bearings, increased durability and reliability are  
guaranteed by ceramic ball bearings



- Quattro Water Spray 4-Nozzle water spray ensures efficient cooling and brings comfort to the entire operation
- Light-weighted design ergonomic yet light weight body provides ease and comfort throughout each operation
- **High durable and repeatedly autoclavable up to 135°C**

**Standard head,** (work w/ KAVO coupling system) 35-105441

**Standard head,** (work w/ NSK coupling system) 35-105442


**Mini head,** (work w/ KAVO coupling system) 35-105443

**Mini head,** (work w/ NSK coupling system) 35-105444

**Coupling,** (KAVO Compatible) 35-105445

**Coupling,** (NSK Compatible) 35-105446

### Fiber Optic Light Source - Built In



**6 Pin Package Includes:** Control Board, Fiber Optic HP Tubing & Transformer 35-105971

**5 Pin Package Includes:** Control Board, Fiber Optic HP Tubing, Transformer, and Bulb Cartridge 35-105972

**Control Board** (Accept 2 FO Tubing) 35-105973

**Fiber Optic HP Tubing - 6 pin, Grey** 35-105974

**Fiber Optic HP Tubing - 5 pin, Grey** 35-105975

**Fiber Optic Light Bulb Cartridge** 35-105976

**Fiber Optic light source Transformer** 35-105977

### Gemini Fiber Optic Handpiece for Kavo (Johnson-Promident)



*Offering exactly the same  
design and performance of the  
originals for a fraction of the price!*

**Handpiece for KaVo Super Torque** 35-25851  
625/630/640 (large head)

**Handpiece for KaVo Bella Torque** 35-25852  
632/642/643 (small head)



## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Integrity Leverage Fiber Optic Handpiece with Lever Cap (Osseo Scientific)

Like top brand hand pieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best high-speed handpiece performance and value at a fraction of the cost of the top leading brands.

- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip



Each 35-117811

### Integrity Leverage Fiber Optic Handpiece with Push Button Cap (Osseo Scientific)

Like top brand handpieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best high-speed handpiece performance and value at a fraction of the cost of the top leading brands.

- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip Grip design



Each 35-117821

### KAVO Interchangeable Fiber Optic Handpiece (Sable)



**Sable ACCESS S4500KL KAVO Interchangeable: Mini Head (No Coupler)** 35-88651

**Sable ROTAMAX S6500KL KAVO Interchangeable: Torque Head (No Coupler)** 35-88652

**6 Hole KAVO Multiflex LUX Type Coupler** 35-88653

**6 Hole KAVO Multiflex LUX Type Coupler: LED Bulb Long Life** 35-88654

**Sable Light Control Set Complete** 35-88655  
**Includes:** Sable Light Control Board, for 1-2 Handpieces with adjustable time delay & voltage feature Sable Power Transformer 120v, 60Hz, 24VDC 1 Amp 6 Pin, 7 Foot, Straight GRAY Tubing

**Sable Light Control Set Complete (5 Hole)** 35-88656  
**Includes:** Sable Light Control Board, for 1-2 Handpieces with adjustable time delay & voltage feature Sable Power Transformer 120v, 60Hz, 24VDC 1 Amp 5 Pin, 7 Foot, Straight Gray Tubing with Bulb (1000815)

### Push Button Fiber Optic Handpiece (Dentex)

Push Button Fiber Optic SINGLE PORT WATER SPRAY 4-Hole and 5-Hole Connector Handpiece (Dentex) , Fixed Back End



**Standard Head Fiber Optic** 35-119391

**Mini Head Fiber Optic** 35-119381

### Push Button Fiber Optic Handpiece Fit Star Coupler (Dentex)

Push Button Fiber Optic Handpiece Interchangeable with Star Coupler (Dentex), Triple Port Water Spray, Swivel Handpiece



**Standard Head Fiber Optic** 35-119461

**Mini Head Fiber Optic** 35-119471

### Push Button Fiber Optic Handpiece Interchangeable with Kavo Coupler (Dentex)

Push Button Fiber Optic Handpiece Interchangeable with Kavo Coupler (Dentex), Swivel Handpiece



**Standard Head Fiber Optic** 35-119401

**Mini Head Fiber Optic** 35-119411

## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Push Button Fiber Optic Handpiece Interchangeable with Midwest/Bien Air Coupler (Dentex)

Push Button Fiber Optic Handpiece Interchangeable with Midwest/Bien Air Coupler (Dentex), Swivel Handpiece, Triple Port Water Spray



**Standard Head Fiber Optic** 35-119421  
**Mini Head Fiber Optic** 35-119431

### Push Button Fiber Optic Handpiece Interchangeable with NSK Swivel Coupler (Dentex)

Push Button Fiber Optic Handpiece Interchangeable with NSK Swivel Coupler (Dentex), Triple Port Water Spray, Swivel Handpiece



**Standard Head Fiber Optic** 35-119441  
**Mini Head Fiber Optic** 35-119451

### Push Button Fiber Optic Handpiece Interchangeable with W & H Coupler (Dentex)

Push Button Fiber Optic Handpiece Interchangeable with W&H Coupler model RA-24 & RA-25 (Dentex), Triple Port Water Spray, Swivel Handpiece



**Standard Head Fiber Optic** 35-119481  
**Mini Head Fiber Optic** 35-119491

### Push Button Fiber Optic Handpiece with Triple Port Water Spray (Dentex)

Handpiece 4-Hole and 5-Hole Connector, Push Button, Triple Port Water Spray, Fixed Back End FIBER OPTIC



**Standard Head Fiber Optic** 35-119501  
**Mini Head Fiber Optic** 35-119511

### Quest Fiber Optic Push Button for Swivel Handpieces (Lares)



**100 Series**  
Power and Durability in a Compact Head



**300 Series**  
Power for Your Toughest Cutting Requirements



- Small, compact head for excellent visibility and accessibility
- High speed for fast cutting
- Reduced contra angle for easier posterior access
- Easy, secure push button bur changing
- Chair-side rotor and chuck replacement
- Optional bright twin beam fiber optics
- Optional 360° swivel and quick disconnect
- Autoclavable and chemiclavable
- Low maintenance Ceramic Bearings



- Full-size head for powerful torque
- Heavy duty construction
- Smooth, quiet operation
- Optional 360° swivel and quick disconnect
- Optional bright twin beam fiber optics
- Easy, secure push button bur changing
- Chair-side rotor and chuck replacement
- Autoclavable and chemiclavable
- Low maintenance Ceramic Bearings
- Multiple Port Water Spray

Quest 150 LO	Height	Diameter	Weight
PB SW MX	11.9 mm	10.2 mm	50 g
11320	Speed	Water Spray	Torque
	450,000 rpm	Single Port	15.3 watts

Quest 350 LO	Height	Diameter	Weight
PB SW MX	14.65 mm	13 mm	56 g
11322	Speed	Water Spray	Torque
	350,000 rpm	Multi-Port	20.8 watts

**150 Small Head HP 360 Swivel,** 35-119681  
Push Button, Multiflex Compatible, Ceramic Bearings, w/o Swivel

**350 Standard Head HP 360** 35-119661  
**Swivel,** Push Button, Multiflex Compatible, Ceramic Bearings, w/o Swivel

**6 Pin Optic Coupler, Multiflex** 35-119671  
Compatible Coupler w/ Bulb

### Quest Optic Handpieces (Lares)

#### **American-Manufactured Highspeeds**

- Quest high speeds feature low-maintenance, long-life, ceramic bearings
- The Quest stainless steel exterior and subtle grip design allows for fast, easy cleaning and autoclaving/chemiclaving
- Quest high speeds are priced for value-conscious dentists unwilling to compromise on quality and performance
- The Quest 100 Series features a compact head for easy posterior access
- These handpieces are fast cutting, and have exceptional power and visibility, easy posterior access
- The Quest 300 Series offers a full sized head for powerful torque and smooth, quiet operation
- Quest 300 models feature heavy-duty construction and multi-port water spray

**Optic Ceramic Bearing:** 35-119751  
Push Button, Small Head, 5 Holes, model 120

**Optic Ceramic Bearing:** Push 35-119761  
Button, Full Sized Head, 5 Holes, model 320

### TwinPower Turbine High Speed F.O. Handpiece (J-Morita)

J. Morita's TwinPower Turbine now features ceramic ball bearings. Forty percent lighter and 3 times harder than conventional bearings, they offer extended turbine life, reduced operation noise, and less vibration. For maximum torque control, TwinPower's unique double-impeller design provides up to 22 W of power. Other features include mini head option, zero suck back in the air line,



Call us for Quarterly Promotion

## Handpieces

rapid braking within 2 seconds, glass rod optics, and direct connection to popular coupling systems.

The TwinPower operates by air from the drive air nozzles powering the primary impeller. The exhaust air is then directed through fixed fins, which power the secondary impeller. The result is that the handpiece is more powerful and provides constant torque and controlled speed, even under load.

**Fits Morita, Fiber Optic Standard Head,** 35-121421  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits Sirona, Fiber Optic Standard Head,** 35-121431  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits W&H, Fiber Optic Standard Head,** 35-121391  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits Kavo, Fiber Optic Standard Head,** 35-121401  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits NSK, Fiber Optic Standard Head,** 35-121411  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits Morita, Fiber Optic Mini Head,** 35-121381  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits Sirona, Fiber Optic Mini Head,** 35-121371  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits W&H, Fiber Optic Mini Head,** 35-121361  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits Kavo, Fiber Optic Mini Head,** 35-121351  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**Fits NSK, Fiber Optic Mini Head,** 35-121341  
**Includes:** Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool

**TwinPower 4-Hole Coupling CP4-O, w/** 35-121331  
**light (6Pin)**

**TwinPower 4-Hole Coupling CP4-WO,** 35-121321  
**w/ water adjustment and light (6 Pin)**

**TwinPower Turbine 4HEX (Standard** 35-121441  
**Head) capsule cartridge (Ceramic Bearing)**



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**TwinPower 4-Hole Coupling CP5-O, w/** 35-121442  
**light (5 Hole)**

**TwinPower 4-Hole LED Coupling** 35-121443  
**CP4-LD, w/light (6 Pin)**

### Low Speed Handpieces

#### Air Torx E-Type Air Motor (J-Morita)

Air Torx offers efficient, powerful, and constant torque with operational speeds up to 20,000 rpm.

Its innovative fluid dynamics generate about twice as much torque compared to a conventional low



speed motor in the 5,000–10,000 rpm range. Air Torx is comfortable to operate - it is lightweight, perfectly balanced, and delivers smooth power control. Versatile, it can be used for a wide range of tasks such as grinding dentures, adjusting, finishing, and polishing restorations, and for tooth polishing with a prophyl angle. Air Torx is designed for enhanced durability and offers an extended working life. Other features include: forward/reverse drive with continuous speed control, double-lock connection, autoclave safe design, and compatibility with standard ISO-type attachments.


- Powerful torque even at a low speed.
- Long working life
- Perfectly balanced, lightweight design
- Forward/reverse drive with continuous speed control
- Safe and easy double-lock connection
- Autoclavable for complete infection control
- Compatible with all standard ISO-type attachments

**Unit Includes:** Lubricant Spray Nozzle 35-98621

#### Airmotor E-Type (Nakamura)

Max Air Pressure 50-60psi

**Made in Japan**



**Midwest Type, E-Type,** 35-111391  
20,000 rpm 4 Holes

**Midwest Type, E-Type,** 35-111381  
5,000 rpm 4 Holes

**Borden Type, E-Type,** 35-111371  
20,000 rpm 2 Holes

**Borden Type, E-Type,** 35-111361  
5,000 rpm 2 Holes

#### Airmotor Handpiece - ND

- Straight air motor  
low speed dental  
handpiece

- Forward/Reverse

operation with **smooth** start-up

- Accepts all standard contra angle and prophyl angles, including disposables

- 360 degrees swivel nosecone, very strong torque, durable, one





# Handpieces

piece design  
- Autoclavable

**One year Manufacturer Warranty**  
**Made in Japan**

**4 Holes, 20,000 rpm** 35-118561

## Contra Angle Handpiece NRS-E (NSK)

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common handpiece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 degrees C and highly durable.

### **Product Features:**

- 10:1 Reduction Head
- For CA burs(Ø2.35)
- Clean Head System
- Push Button Chuck
- Max. Speed 3,000 min

**NRS-E Contra Angle Handpiece** 35-11186



## Contra Angle Sheath for Midwest Airmotor Type (Nakamura)

### CONTRA ANGLE SHEATH FOR MIDWEST MOTOR



**SHS-MWC Contra Angle Sheath** 35-120651

**EH-50PK Spare Head** 35-120531

**EH-50PS Spare Head** 35-120541

**EH-20L Spare Head** 35-120631

**EH-20ML Spare Head** 35-120621

**EH-20PS Spare Head** 35-120561

**EH-20PK Spare Head** 35-120551

**EH-35BL Spare Head** 35-120571

**EH-30BL Spare Head** 35-120581

**EH-30BLP Spare Head** 35-120591

**EH-30TL Spare Head** 35-120601

**EH-20FG Spare Head** 35-120611

*Visit Our website @*  
[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)



**Prestige Dental Products, Inc.**  
[Prestigedentalproducts.com](http://Prestigedentalproducts.com)  
Tel: 877-772-3888

## Contra Angle with Push Button Head Handpiece FBP-E (NSK)

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common handpiece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 deg. and highly durable.

### **Features:**

- For CA burs(Ø2.35)
- Clean Head System
- Ultra Push Chuck
- Connect directly to all brands of ISO E-type air or electric micromotors 1:1 transmission gear

**1:1 Direct Drive:** For CA burs(Ø2.35) - 35-111881

Clean Head System -Push Button Chuck -  
Ball Bearing -Max. Speed 40,000 min



## E-Type Contra Angle 4:1 (Nakamura)



**EG-30BL Midwest Type** 35-11233

**EG-30TL Midwest Type** 35-112341

**EG-20L Midwest Type** 35-112351

**EG-20ML Midwest Type** 35-112361

**EG-20FG Midwest Type** 35-112371

**EG-20AP Midwest Type** 35-11228

**EF-20APS Midwest Type** 35-11227

**EG-20PS Midwest Type** 35-11218

**EG-20PK Midwest Type** 35-11219

**EG-50PS Midwest Type** 35-11220

**EG-50PK Midwest Type** 35-11221

**EGN-30BL NSK or STAR Type** 35-11223

**EGN-30TL NSK or STAR Type** 35-11222

**EGN-20L NSK or STAR Type** 35-11224

**EGN-20ML NSK or STAR Type** 35-11225

**EGN-20FG NSK or STAR Type** 35-11226


**EGN-20AP NSK or STAR Type** 35-112381

**EGN-20APS NSK or STAR Type** 35-11232

## Handpieces

<b>EGN-20PS NSK or STAR Type</b>	35-11231
<b>EGN-20PK NSK or STAR Type</b>	35-11230
<b>EGN-50PS NSK or STAR Type</b>	35-11229
<b>EGN-50PK NSK or STAR Type</b>	35-11217

### E-Type Contra Angle Handpieces (Nakamura)



**Made in Japan**


<b>EC-30BL Push-Up Latch, Ball-Bearing Head for CA bur, 30,000 rpm</b>	35-11207
<b>EC-30BLP Push Button, Ball-Bearing Head for CA bur, 30,000 rpm</b>	35-11206
<b>EC-20L Standard Latch Head for CA bur, 20,000 rpm</b>	35-11208
<b>ESG-30AR 4:1 Speed Reduction Straight Handpiece</b>	35-112109
<b>ES-30A Straight Handpiece 40,000 rpm</b>	35-11210
<b>ME-80B 2 Holes 5,000 rpm Airmotor</b>	35-11211
<b>ME-80M 4 Holes 5,000 rpm Airmotor</b>	35-112131
<b>ME-20B 2 Holes 20,000 rpm Airmotor</b>	35-11204
<b>EG-20L 4:1 Speed Reduction Standard Latch Head for CA bur</b>	35-11203
<b>EG-20PS 4:1 Speed Reduction Regular Propphy Screw-On Head</b>	35-11193
<b>ME-20M 4 Holes 20,000 rpm Airmotor</b>	35-112151



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### E-Type Durable Handpiece (Dentex)

#### **Durable Stainless Steel**



<b>20,000 RPM Motor</b>	35-112481
<b>1:1 Contra Angle</b>	35-112461
<b>4:1 Reduction Contra Angle</b>	35-112451
<b>Manual Latch Head</b>	35-11242
<b>Push Latch Head</b>	35-11241
<b>Prophy Head Screw Type</b>	35-11240
<b>Quarter Turn Endo Head</b>	35-11239
<b>5,000 RPM Air Motor</b>	35-112471
<b>1:1 Nose Cone</b>	35-112441
<b>4:1 Reduction Nose Cone</b>	35-11243

### E-Type Handpiece Sets (Nakamura)



**Made in Japan**

<b>Set: ES-30A 1:1 Straight Nosecone HP, EC-20L 1:1 Standard Latch Head Contra Angle HP &amp; ME-20M 20,000 rpm Air motor</b>	35-111421
<b>Set: ES-30A 1:1 Straight Nosecone HP, EC-30BL 1:1 Push-Up Latch, Ball-Bearing Head CA HP &amp; ME-20M 20,000 rpm Air motor</b>	35-111411
<b>Set: ES-30A 1:1 Straight Nosecone HP, EC-20L 1:1 Standard Latch Head Contra Angle HP &amp; ME-80M 5000 rpm Air motor</b>	35-111401



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Handpieces

### E-Type Handpiece with Latch Head BB-E (NSK)

E-type External Water Spray – for Varying Treatment Needs

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common handpiece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 deg. and highly durable.

#### **Features:**


- For CA burs(Ø2.35)
- Clean Head System
- Swing latch type
- Connect directly to all brands of ISO E-type air or electric micro motors 1:1 transmission gear
- 40,000 rpm Ball bearing head for smoother and vibration free operation



**NSK BB-E 1:1** 35-112491  
w/Latch Head 40,000 rpm

### E-Type Ultra Light Handpiece (Dentex)

#### **Ultra Light Aluminum**



<b>20,000 RPM Motor</b>	35-112601
<b>1:1 Contra Angle</b>	35-112571
<b>4:1 Reduction Contra Angle</b>	35-112581
<b>Manual Latch Head</b>	35-11255
<b>Push Latch Head</b>	35-11254
<b>Prophy Head Screw Type</b>	35-11253
<b>Quarter Turn Endo Head</b>	35-11251
<b>5,000 RPM Air Motor</b>	35-11256
<b>1:1 Nose Cone</b>	35-11252
<b>4:1 Reduction Nose Cone</b>	35-11259

### Friction Grip Push Button Type Handpiece NFB-E (NSK)

E-type External Water Spray – for Varying Treatment Needs

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common handpiece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 deg. and highly durable.



#### **Features:**

- Standard Head
- Swing Latch type
- For friction grip burs (Ø1.6)
- Max.40,000 min-1
- Connect directly to all brands of ISO E-type air or electric micromotors 1:1 transmission gear

**NFB-E Friction Grip Push Button type** 35-112641  
Handpiece

### Hygiene Prophy Handpiece (ND)



- 360 degree Swivel Hygienist Handpiece
- Compact Size
- Light weight is ergonomically designed to reduce fatigue for long time

use.

- Nearly all the disposable prophy angles can be mounted
- Collet internal mechanical design is wobble free
- 0-5, 000 rpm

**Fully Autoclavable**  
**One year Manufacturer Warranty**  
**Made in Japan**

**Midwest Type 5,000, rpm 4 Holes** 35-117691

### Integrity Low Speed Handpiece Set (Osseo Scientific)

#### **Includes:**

- 20,000 RPM E-Type Motor (NTE-20)
- 1:4 Attachment (NTE-14)
- Standard Latch Head (NTL-01)
- 1:1 Straight Nose Cone (NTS-11)



**Low Speed Handpiece Combination** 35-117881  
**Set includes:** 20,000 RPM E-Type Motor (NTE-20) 1:4 Attachment (NTE-14) Standard Latch Head (NTL-01) 1:1 Straight Nose Cone (NTS-11)

**Low Speed Handpiece Combination** 35-117871  
**Set includes:** 20,000 RPM E-Type Motor (NTE-20) 1:4 Straight Nose Cone (NTS-14) 1:1 Straight Nose Cone (NTS-11)

**Low Speed Handpiece Combination** 35-11786  
**Set includes:** 5,000 RPM E-Type Motor (NTE-5) 1:1 Attachment (NTE-11) Standard Latch Head (NTL-01) 1:1 Straight Nose Cone (NTS-11)

**Low Speed Handpiece Combination** 35-11785  
**Set include:** 20,000 RPM E-Type Motor (NTE-20) 1:4 Attachment (NTE-14) 1:1 Attachment (NTE-11) 2 Standard Latch Heads (2NTL-01)



# Handpieces

## Lynx Low Speed Handpiece (MTI)

### *Air Motor, Low Speed Handpiece System*

- Length: 405"
- Weight: 2.7oz
- Speed: 2,000 to 20,000
- Operating Air Pressure: 35 to 50 psi.
- Air Consumption: 1.4 cfm.
- Connections: 4-, or 2-hole



**2 Hole** 35-118121

**4 Hole** 35-118111

## Micromite Low Speed Handpiece (Micromite-Japan)

Autoclavable, High torque



**5,000RPM 4 Holes** 35-118551

## Micromite Low Speed Handpiece - (ND)

- Straight air motor low speed dental handpiece
- Forward/Reverse operation with **smooth** start-up
- Accepts all standard contra angle and prophy angles, including disposables
- 360 degrees swivel nosecone, very strong torque, durable, one piece design
- Autoclavable



*One year Warranty  
Made in Japan*


**4 holes Midwest Type, 5,000 rpm** 35-105452

## Prophy Angle E-Type 1:1 (Nakamura)

**EC-20PS**  
Prophy Screw-in




**EC-50PS**  
Sealed Prophy Screw-in



**EC-20PK**  
Prophy Snap-on



**EC-50PK**  
Sealed Prophy Snap-on



**Midwest Type, EC-20PS** 35-21546

**Midwest Type, EC-20PK** 35-21547

**Midwest Type, EC-50PS** 35-21548

**Midwest Type, EC-50PK** 35-21549

## Prophy Angle Handpiece Heads (Nakamura)

**AIRMOTOR HANDPIECE FOR HYGIENIST**



**RA-10PS**  
Regular Screw-In  
Prophy Angle

**RA-10PK**  
Regular Snap-On  
Prophy Angle



**RA-50PS**  
Sealed Screw-In  
Prophy Angle

**RA-10PK**  
Sealed Snap-On  
Prophy Angle



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**Regular Screw-In Prophy Angle HP** 35-117692

**Regular Screw-On Prophy Angle HP** 35-117693

**Sealed Screw-In Prophy Angle HP** 35-117694

**Sealed Screw-On Prophy Angle HP** 35-117695

## Push Button Head Contra Angle Handpieces (Nakamura)

30,000rpm, Ball-Bearing  
Push-Button Type Head, For CA  
latch bur



**Made In Japan**

**E-type (ISO Type)** 35-82361  
Connection Contra Angle

**Star Titan-type** Connection Contra Angle 35-82362

**Midwest-type** Connection Contra Angle 35-82363

## Rotary Latch Files Head Handpiece NAC-E316R (NSK)

**E-Type External Water Spray**  
The speed of all Endodontic low speed handpieces is determined by the speed of the micromotor to which they are connected. Micromotors can vary in speed ranges. Micromotor speed should be confirmed before selecting a reduction handpiece.

### **Features:**

- For engine files
- 16:1 Reduction
- 360°rotation
- For Ni-Ti files Standard Head (NAC-Y) Latch Head, (E16R) Contra Angle
- Spray Nozzle (Optional)



**NAC-E16R, Rotary Latch Files Head** 35-120141  
Handpiece

## Slow Speed E Type (Johnson-Promident)



**Classic Slow Speed Set** 35-27241  
(E-Motor + Nose + Contra Angle 30K)

**Classic "E" Type Motor (25,000 RPM)** 35-27242

**Classic "E" Type Motor (5,000 RPM)** 35-27243

## Straight Durable Handpiece (Dentex)

Compact easy-to-use low-speed handpiece for use with doriot style  
contra angles and straight handpiece/  
lab burs. 20,000 rpm

- Straight nose swivels 360°, reducing tubing twist and hand fatigue
- Speed adjusting  
ring controls speed and direction of motor


## Handpieces

- No lubrication required
- Fully autoclavable

**Durable Stainless Steel** 35-120881

### Star Titan Type Low Speed Handpieces (Nakamura)

Superior performance, Very low noise level, High durability, Various system choice, Easy & comfortable use, Autoclavable



**Head Attachment** For Friction Grip bur 20,000rpm 35-92341

**Head Attachment** For CA latch bur 20,000rpm 35-92342

**Head Attachment** For CA latch bur Ball-Bearing 30,000rpm 35-92343

**Head Attachment** For CA latch bur Ball-Bearing Push-Button 30,000rpm 35-92344

**Head Attachment** For Prophy Snap-on 35-92345

**Head Attachment** For Prophy Screw-in 1-64 35-92346

**1:1 Contra Angle Sheath Low Speed Attachment** 35-92347

**4:1 Speed Reduction Torque Multiplier Low Speed Attachment** 35-92348

**Straight Nosecone** Ball-Bearing, U-type Nose 40,000rpm Low Speed Attachment 35-92349

**1:1 Contra Angle Sheath 12-teeth Low Speed Attachment** 35-92350

**20,000rpm Low Speed Air Motor 4-hole Connector** 35-92351

**5,000rpm Low Speed Air Motor 4-hole Connector** 35-92352

### Straight Durable Handpiece M5 (Dentex)

Compact, easy-to-use low-speed handpiece for use with doriot-style contra angles, disposable and metal prophy angles and straight handpiece/lab burs

- Variable speed up to 5,000 rpm

- Straight nose cone swivels 360°, reducing tubing twist and hand fatigue



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888


- Speed-adjusting ring controls speed and direction of motor
- No lubrication required
- Fully autoclavable

**Durable Stainless Steel** 35-120891

### Straight Handpieces E-Type (Nakamura)



ES-30A



ES-35L

**1:1, 40,000 rpm,** 35-120911  
3 Ball Bearings U Type Nose

**1:1, 35,000 rpm,** 35-120901  
2 Ball Bearings For Laboratory

### Straight Handpieces for Midwest Type Airmotor (Nakamura)

**Made in Japan**



MWS-30A  
40,000rpm

**40,000 rpm 3 Ball Bearings** 35-120921

### Straight Nose Cone Handpiece EX5 (NSK)

EX5 E-type External Water Spray – for Varying Treatment Needs

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common handpiece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 deg. and highly durable.

#### **Features:**

- For Prophylaxis Application
- For screw in-cups and brushes
- Connect directly to all brands of ISO E-type air or electric micromotors 4:1 reducer

**NSK EX5, 4:1 reduction Straight Hand-** 35-120931  
**piece**



### Straight Ultra Light Handpiece (Dentex)

Compact easy-to-use low-speed handpiece for use with doriot style contra angles and straight handpiece/lab burs. 20,000 rpm



- Straight nose swivels 360°, reducing tubing twist and hand fatigue
- Speed adjusting ring controls speed and direction of motor
- No lubrication required
- Fully autoclavable

**Ultra Light Aluminum** 35-120951

## Handpieces

### Straight Ultra Light Handpiece M5 (Dentex)

Compact, easy-to-use low-speed handpiece for use with doriot-style contra angles, disposable and metal prophyl angles and straight handpiece/lab burs.



- Variable speed up to 5,000 rpm
- Straight nose cone swivels 360°, reducing tubing twist and hand fatigue
- Speed-adjusting ring controls speed and direction of motor
- No lubrication required
- Fully autoclavable

**Ultra Light Aluminum** 35-120961

### Stripping E-Type Contra Angle HP - ND



- Vertical Reciprocating
- 4:1 SPEED Reduction
- For Eva

Tips  
- Max Speed: 10,000rpm  
- Fit NSK Cable

*Made in Japan*


**Each** 35-105291

### Titan-T Straight Nose Cone Attachment (StarDental)



**Each** 35-108061

### U-Type Contra Angle (Nakamura)



<b>DU-30BL</b>	35-121451
<b>DU-20L</b>	35-121471
<b>DU-20ML</b>	35-121481
<b>DU-20FG</b>	35-121491

Find Us On Social Media



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Surgical Handpieces

### Impact Air 45 Surgical Handpiece (Palisades)



- Designed for use in periodontal and endodontic procedures
- Unique 45° angled head, allowing maximum access and visibility
- available with or without fiberoptics
- Highspeed Air Powered Handpiece
- **Fully Autoclavable/Chemoclavable**

400,000 RPM

**4-Hole Surgical Handpiece** 35-117771

**Fiber Optic 5-Hole, handpiece** 35-117761

**2-Hole, handpiece** 35-11775

**Model 9 Push-Button Conversion** 35-11774

**Kit:** Turbine, head cap, and cap wrench

**Model 9 Push-Button Turbine (Replacement)** 35-11773

**Model 9 Push-Button Chuck** 35-11772

**Model 9 Push-Button Head Cap & Cap Wrench** 35-11771

**IMPACT AIR 45 Autoclavable Lubricant** 35-11770

### Surgical Push Button Handpiece (Pana-MAX)



- **45° Angle**

**Swivel Adapter included** 35-16521

**Standard Connection** 35-16522

### Twin Power Surgical Handpiece (J-Morita)

- **SAFE**, rear-facing exhaust vents

- High torque; up to 20 watts of power

- **QUIET & LONG-LASTING** ceramic ball bearings

**Handpiece 4H-PAR-4HEX-O-45** 35-98821  
(Morita Couplings)

**Handpiece 4H PAR-4HEX-O-KV-45** 35-98822  
(Kavo MULTIflex LUX Coupling)

**Turbine 4HEX-45 Capsule Cartridge** 35-98823  
(Ceramic Bearing)



Call Us Today **833-733-6825**  
Or Visit our website at  
**Prestigedentalproducts.com**



# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Handpiece Controls

### 2 Handpiece Manual Control (Chapman-Huffman)



- **NOTE! Control Requires Filtered/ Regulated Air & Water**
- Manual activation for 2 handpieces
- Adjustable coolant water for one handpiece
- Handpiece tubings (**SPECIFY TYPE BY PART #**)
- Handpiece pressure indicator gauge
- Individual Turbine pressure adjustment
- Instrument holders
- Variable speed disc type foot control
- Asepsis 3-way syringe - pin type
- Control Measures 6" (152mm) wide x 3 1/2" (89mm) high x 1 1/2" (38mm) deep
- Installation Kit

w/ Gray Straight, Handpiece tube	35-111161
w/ Gray Coiled, Handpiece tube	35-111151
w/ Black Straight, Handpiece tube	35-111141
w/ Black Coiled, Handpiece tube	35-111131
w/ Surf Straight, Handpiece tube	35-11112
w/ Surf Coiled, Handpiece tube	35-111171

### 2 Handpiece Manual w/ Coolant Water (Chapman-Huffman)



- **NOTE! Control Requires Filtered/ Regulated Air & Water**
- Manual activation for 2 handpieces
- Adjustable coolant water for both handpieces
- Handpiece tubings (**SPECIFY TYPE BY PART #**)
- Handpiece pressure indicator gauge
- Individual Turbine pressure adjustment
- Instrument holders
- Variable speed disc type foot control
- Asepsis 3-way syringe - pin type
- Control Measures 6" (152mm) wide x 3 1/2" (89mm) high x 1 1/2" (38mm) deep
- Installation Kit

w/ Gray Straight, Handpiece tube	35-111221
w/ Gray Coiled, Handpiece tube	35-111211
w/ Black Straight, Handpiece tube	35-111201
w/ Black Coiled, Handpiece tube	35-111191
w/ Surf Straight, Handpiece tube	35-111181
w/ Surf Coiled, Handpiece tube	35-111231

### 3 Handpiece Automatic Control w/ Floor Valves (Chapman-Huffman)



- **NOTE! Control Requires Filtered/ Regulated Air & Water**
- Automatic activation for 3 handpieces
- Coolant on-off toggle
- Adjustable coolant water for each handpiece
- Asepsis 3-way syringe - pin type
- Individual turbine pressure adjustment

- Master on-off toggle
- Master floor valves less junction box enclosure
- 7' Umbilical w/ air, water, and master switch tubings
- Handpiece holders w/ lock-out lever
- Handpiece pressure indicator gauge
- Variable speed disc type foot control
- Control measures 7 3/4" (197mm) wide x 4" (102mm) High x 2" (51mm) deep
- Handpiece tubings (**SPECIFY TYPE BY PART #**)
- Installation Kit

w/ Floor Valves & Gray Straight, Handpiece tube	35-111281
w/ Floor Valves & Gray Coiled, Handpiece tube	35-11127
w/ Floor Valves & Black Straight, Handpiece tube	35-11126
w/ Floor Valves & Black Coiled, Handpiece tube	35-11125
w/ Floor Valves & Surf Straight, Handpiece tube	35-11124
w/ Floor Valves & Surf Coiled, Handpiece tube	35-111291

### 3 Handpiece Automatic Control without Floor Valves (Chapman-Huffman)



- **NOTE! Control Requires Filtered/ Regulated Air & Water**
- Automatic activation for 3 handpieces
- Coolant on-off toggle
- Adjustable coolant water for each handpiece
- Asepsis 3-way syringe - pin type
- Individual turbine pressure adjustment
- Master on-off toggle
- Master floor valves less junction box enclosure
- 7' Umbilical w/ air, water, and master switch tubings

- Handpiece holders w/ lock-out lever	
- Handpiece pressure indicator gauge	
- Variable speed disc type foot control	
- Control measures 7 3/4" (197mm) wide x 4" (102mm) High x 2" (51mm) deep	
- Handpiece tubings ( <b>SPECIFY TYPE BY PART #</b> )	
- Installation Kit	
w/ Gray Straight, Handpiece tube	35-111341
w/ Gray Coiled, Handpiece tube	35-111331
w/ Black Straight, Handpiece tube	35-111321
w/ Black Coiled, Handpiece tube	35-111311
w/ Surf Straight, Handpiece tube	35-111301
w/ Surf Coiled, Handpiece tube	35-111351




We Work Hard  
So You Don't  
Have To.  
877-772-3888

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Fiber Optic Light Source - Built In



Package includes: H7020+H7045+H7090    Package includes: H7020+H7040+H7070+H7090

**6 Pin Package Includes:** Control Board, Fiber Optic HP Tubing & Transformer    35-105971

**5 Pin Package Includes:** Control Board, Fiber Optic HP Tubing, Transformer, and Bulb Cartridge    35-105972

**Control Board** (Accept 2 FO Tubing)    35-105973


**Fiber Optic HP Tubing - 6 pin, Grey**    35-105974

**Fiber Optic HP Tubing - 5 pin, Grey**    35-105975

**Fiber Optic Light Bulb Cartridge**    35-105976

**Fiber Optic light source Transformer**    35-105977

## Lab Handpiece Control Kit (Chapman-Huffman)



- Handpiece holder and bracket
- Installation instructions

- Control for one dry handpiece
- Air pressure regulator
- Regulator mounting bracket
- Air filter assembly (35 micron)
- Handpiece tubing assembly (4-hole connector w/ 2-hole gray tubing)
- Disc type foot control

**Kit**    35-117981

## Prophy Handpieces

### Hygiene Handpiece (MTI)

- Excellent for hygienic procedures, pins, post excavations and prophys
- This high torque, low speed handpiece is perfectly balanced to reduce hand and arm fatigue
- Fits all straight burs, doriot attachments and disposable prophy heads
- Fixed Forward 5,000 rpm



**Each**    35-89791

## Hygiene Handpiece - Clear Advantage (Johnson-Promident)

### Compare to Midwest RDH

- Same design, performance, and warranty as the original for just half the price
- Lightweight ergonomic design for better positioning and reduced fatigue
- Smooth 360 degree swivel for better comfort and access
- Removable autoclavable color bands for handpiece identification



**Each**    35-117958

## Hygiene Prophy AirLite Handpiece with Friction Grip (MTI)

- Fixed forward 3,000 rpm
- Contoured for comfort
- Super quiet vibration free
- No more torn gloves with our Friction Grip System
- Your choice of attractive colors
- Now with a **2-Year Warranty**



**Blue Color**    35-89801

**Purple Color**    35-89802

**Pink Color**    35-89803

## Hygiene Prophy Handpiece (ND)



- 360 degree Swivel Hygienist Handpiece
- Compact Size
- Light weight is ergonomically designed to reduce fatigue for long time use.

- Nearly all the disposable prophy angles can be mounted
- Collet internal mechanical design is wobble free
- 0-5, 000 rpm

**Fully Autoclavable**  
**One year Manufacturer Warranty**  
**Made in Japan**

**Midwest Type 5,000, rpm 4 Holes**    35-117691

## Hygienist Durable Prophy Handpiece (Dentex)

Durable Stainless Steel M5P Prophy Handpiece (Dentex)

- High Torque Motor
- Designed for use with disposable and metal autoclavable prophy angles
- Quiet, compact and light weight-minimizes hand and wrist fatigue
- Nose Cone swivels 360°
- No Lubrication required
- Fully Autoclavable



**0-5,000 Durable Stainless Steel M5P Hygienist Prophy Handpiece (Dentex)**    35-117671

## Handpieces

### Hygienist Plus Handpiece (Johnson-Promident)



- For disposable prophyl angles **ONLY**
- 5,000 RPM ergonomically designed for

effortless operation

- **ULTRA-Lightweight** for reduced fatigue and easy handling
- 360° swivel for better comfort and access
- Simple push/pull friction grip chuck system

**Hygienist Plus, Each**

35-110981

### iBrite Deplaquer Polisher (Pac-Dent)

A convenient cordless stain remover that cleans better than just brushing and flossing. Its unique fusion of professional technology and ease of use allows for the ultimate removal of tobacco, coffee, tea, wine, and other types of stains even years after they form. The results will be visible immediately after use, and you will be able to smile more confidently. The kit comes with several different types and sizes of cleaning heads to ensure the cleaning cups reach all tooth crevices and gaps, and effectively remove the stains accumulated over time. The included polishing paste contains fluoride for cavity prevention.

#### **Advantages:**

- Cleans better than brushing and flossing alone
- Effective removal of tobacco, coffee, tea, wine, and other types of stains
- Immediate results
- Included paste contains fluoride for cavity prevention

**1x iBrite Deplaquer/Polisher Unit: 1x Handpiece, 1x Remi desensitizing paste pen, 4x AA batteries, 1x Shade guide, 6x Prophyl pastes, 5x Bids with adhesive stickers, 1x Instruction, 6x Prophyl angles (2-Pedo cups, 3-ProAngle Ergo Torque cups, 1-Brush cup)**



### Integrity Hygiene Handpiece (Osseo Scientific)

#### **Lightest Swivel Hygiene Handpiece!**

**Only 2.3oz (60g)**

- Friction grip chucking, just push in or pull out the prophyl head
- Smooth, vibration-free operation up to 4,000 RPM
- Ergonomic design provides perfect balance
- 360° swivel for superb maneuverability
- Compatible with all standard doriot-style disposable metal angles



Brand	Weight (oz)	Length (in)
<b>Prophyhygiene™</b>	<b>2.3</b>	<b>3.9</b>
NSK Hygiene Pro	2.5	4.0
Kavo PROPHYwiz	2.7	4.3
Midwest RDH	3.0	4.4

**Integrity Prophyhygiene**

35-117791



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

### iStar Cordless Prophyl Handpiece (Star Dental)

It is designed to make dental hygiene procedures easier, faster and less stressful. The iStar provides you with the freedom from handpiece hoses and foot pedal. Choose your favorite prophyl angle and go cordless.

- Accepts **ALL** disposable angles
- Cordless
- Disposable Aseptic Sleeves
- Ergonomic, lightweight design produces less stress on hand & Back
- No need to autoclave



**Complete Unit**

35-108891

**Disposable Aseptic Sleeves, 500/Box**

35-108892

### Micromax Cordless Prophylaxis Handpiece (Dentamerica)

- **Cordless** convenience offers complete portability
- **Heavy-duty high torque motor**
- **Forward/reverse** functions
- Reciprocating setting available without use of additional attachment
- **9 rotation speeds** available for different applications
- Ergonomic design provides comfort and easy operation
- **High capacity lithium-ion battery** provides stable and prolonged usage
- Offers quick and easy coupling
- Low battery warning light

#### **Specifications:**

- 100-240V AC 50/60Hz
- Dimensions: L30 x W27 x H197 mm
- Weight: 150g
- Accessories: Charging base, Handpiece, Power Cable



**Unit with Latch Attachment**

35-118541

**Unit with Screw Attachment**

35-118531

**Unit with Snap-On Attachment**

35-118521

**Unit with Straight Attachment**

35-118511

**Latch Attachment**

35-118501

**Screw Attachment**

35-118491

**Snap-On Attachment**

35-118481

**Straight Attachment**

35-118471

### ProMate E Portable Hygiene Handpiece (Pac-Dent)

**A completely portable and rechargeable handpiece.**

This portable hygienist's handpiece is the ultimate in convenience, comfort, and control. No more cords to mess with, lightweight, fully autoclavable, rechargeable, comfortable grip and contour, as well as minimal vibration and chatter. Smoothly control speed with a simple twist to your desired setting anywhere between 2,000 and 20,000 rpm.





# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Advantages:

- Cordless, portable, and rechargeable
- Textured silicon grip reduces vibration and chatter for enhanced comfort & non-slip control
- Fully autoclavable and lightweight
- Consistent, powerful torque delivery
- Variable speed between 2,000 and 20,000 rpm
- Interchangeable autoclavable grips allow color-coding, personalization & easy maintenance
- Works with any standard disposable prophylaxis angle on the market

### 1 X Portable Hygiene Handpiece

with 3 X silicon grips 35-55400

### Dusty Rose silicon grips for hygiene handpiece, 10/pkg

35-55401

### Green silicon grips for portable handpiece, 10/pkg

35-55402

### Blue silicon grips for portable handpiece, 10/pkg

35-55403

### Lavender silicon grips for portable handpiece, 10/pkg

35-55404

## ProMate™ Hygiene Prophylaxis Handpiece (Pac-Dent)

Reduce fatigue and daily micro-traumas on the hand and arm muscles with our new ProMate™ hygiene handpiece. It has a lightweight, ergonomic design that minimizes vibration and chatter as well as a textured silicon grip and contoured shape that provides superior comfort and non-slip control. Maintain the ProMate™ easily by autoclaving the entire handpiece with the grip, or remove and replace the grip for color coding and personalization. Use with your favorite standard disposable prophylaxis angles today and see the difference! The grips are available in blue, green, lavender and dusty rose.



## Advantages:

- **Very economical**, best buy for price & quality
- Soft, non-slip grips help reduce fatigue and daily micro-traumas
- Color-coded grips are ideal for personalization
- **Easy maintenance** with the entire handpiece being autoclavable
- Ergonomic lightweight design ensures superior comfort
- Compared to the best brands on the market
- Works with any standard disposable prophylaxis angle

### 1 X Hygiene Handpiece with 3 X silicon grips

35-55391

### Green silicon grips for hygiene handpiece, 10/pkg

35-55392

### Blue silicon grips for hygiene handpiece, 10/pkg

35-55393

### Lavender silicon grips for hygiene handpiece, 10/pkg

35-55394

### Dusty Rose silicon grips for hygiene handpiece, 10/pkg

35-55395

## Prophy Star 3 Handpieces (Star Dental)



Used for prophylactic procedures. The lube-free, autoclavable handpiece is made of a lightweight anodized aluminum. The longer 360° swivel provides better control of movement. The handpieces operating speeds

are between 0 & 5,000 rpm with an operating pressure between 35 to 40 psi. **Each prophylaxis star includes 4 different grip colors.**

## Features & Benefits

- 4-Hole tubing connection
- Lube Free
- Autoclavable
- Choice of grip colors for color coding

**Complete Set** 35-108901

**Replacement Grips, 4/Pkg** 35-108902

**Real Seal, 4-Line, 3/Pkg** 35-108903

## Electric Motors

### 207S Series (Saeshin)



**207S SET B 106HP 35K Variable** 35-120351

**207S SET B 106HP 45K Variable** 35-120341

**207S SET B 107HP 35K Variable** 35-120361

**207S SET B OPTIMUS HP 35K Variable** 35-120371

**207S SET B OPTIMUS HP 45K Variable** 35-120331

**107S Control box only** 35-120321

### E-File Portable (Saeshin)



**RPM Portable E-File Set in Gray Willow Set** 35-120241

**Ergo Slim Handpiece Only in Gray Willow, Each** 35-120231

**RPM Ergo Slim E-Type Motor in Gray, Each** 35-120251

**RPM On/Off Pedal for Table, Each** 35-120221

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Endo E Class Rotary System (Meta)



**Full Unit, without Handpiece** 29-101391

## Handy I Portable Micromotor (Saeshin)

**Handy I Portable Micromotor Set,** 35-11739  
20K Set

**Handy I Control Box,** 35-11742  
20K Each

**Handy E-Type Handpiece,** 35-11741  
20K RPM Each

**1:1 Straight Nose Cone,** 35-117441  
35K Each

**Handy I Control Box,** 35-11743  
20K with Charger Each

**Handy I Charger for Handy I Con-** 35-117451  
**troller, Each**

**Charger for Handy I Controller,** 35-11740  
230V Each

**Pouch for Handy I Controller Each** 35-11738

## M35 Compact Electric Lab Handpiece System (Buffalo)

**Console:** 5 in. (128 mm)  
**Width x** 6.5 in. ( 166 mm)  
**Depth x** 3.25 in. (82 mm)  
**Height x** 2.9 lb. ( 1.3 kg)  
**Weight**



•High Torque Cutting Power from 0 – 35,000 RPM New high torque handpiece and control system provide excellent cutting power and smooth control in a compact system

•Improved Versatility System is designed for 120V AC operation but can be easily switched to 220V AC via setting on the underside of the Console

•Improved Console and Controls Improved controls provide fully variable speed via console or foot control, simplify operation and protect the controls and handpiece

•Excellent High and Ultra-Low Speed Operation System excels at providing a full range of usable operating speeds from ultra-low to 35,000 RPM

•Heavy Duty Construction Components and system designed to provide reliable and trouble-free operation. **1 year warranty**

### •Sizes:

Handpiece: 6.25 in. (159 mm) Long x 1.1 in. (28.5 mm) Back Diameter .79 in. (20 mm) Front Diameter x 8 oz. (.23 kg) Weight

**120V/220V AC** 35-118131  
(Console+MG Gray HP+FC)

## Magic Wand Lab Handpiece System (Buffalo)

### Features/Benefits:

•Compact electric handpiece system features fully variable cutting speeds from 0 to 35,000 RPM

•Premium black electric handpiece provides smooth, quiet and comfortable operation with a history of long term durability

•Reliable, no-frills, made-in-USA Console

•Fully variable speed conveniently controlled by a heavy-duty, all-metal Variable Speed AC Foot Control

•The convenience and benefits of an electric handpiece without the added extras of higher priced systems


### •Sizes

**Handpiece:** 6 in. (153 mm) Long x 1 1/8 in. (29 mm) Major Diameter, 5/8 in. (16 mm) Nose Diameter x 9.5 oz. (270 gm) Weight

**Console:** 4.25 in. (108 mm) Width x 5.25 in. (133 mm) Depth x 3 1/8 in. (80 mm) Height, 2 lb. 7 oz. (1.1 kg) Weight

**120V AC** (Console + Black HP + FC) 35-118151

**220V AC** (Console + Black HP + FC) 35-118141



## Microlab Brush Motor (Saeshin)

•Feedback control system

•Variable Speed dial system

•Right & Left turning ability

•Hand/Foot speed control system

•Automatic overload protection system

•45,000 rpm high-powered carbon brush motor

•ON/Off and variable speed foot control

•Torque 450 gf.cm

•Input 110 / 230 Volts (50 / 60Hz)

•Output DC 0-32W

•Weight Control box: 3.82 Lb Handpiece: 8.0 oz. (Strong 106)



**MicroLab 350 Control Box,** 35-118411  
Tech2000 HP ON/OFF

**MicroLab 350 Control Box,** 35-118401  
Tech2000 HP Variable

**MicroLab 350 Control Box,** 35-118391  
Tech2000 HP None

**MicroLab 350/450 Control Box Only** 35-118421

**MicroLab 450 Control Box,** 35-118431  
45,000 RPM HP ON/OFF

**MicroLab 450 Control Box,** 35-11844  
45,000 RPM HP Variable

**MicroLab 450 Control Box,** 35-118451  
45,000 RPM HP None

**MicroLab 350/450 CB, E-Type HP, 1:1** 35-118461  
**Straight Nose Cone ON/OFF**

## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**MicroLab 350/450 CB, E-Type HP, 1:1 SNC, 1:1 CA ON/OFF** 35-118381

### Microtorque Motor (Saeshin)



**Microtorque Control Box, Tech2000 Handpiece Variable** 35-118731

**Microtorque Control Box, Slim Nose Handpiece Variable** 35-118741

**Microtorque Control Box, 45,000 RPM Handpiece Variable** 35-118721

**Microtorque Control Box, Tech2000 Handpiece ON/OFF** 35-118711

**Microtorque Control Box, Slim Nose Handpiece ON/OFF** 35-118701

**Microtorque Control Box, 45,000 RPM Handpiece ON/OFF** 35-118691

**35,000 rpm Microtorque Control Box Only (110V)** 35-118681

**35,000 rpm Microtorque Control Box Only (230V)** 35-118761

**45,000 rpm Microtorque Control Box Only (110V)** 35-118751

**45,000 rpm Microtorque Control Box Only (230V)** 35-118671

### Mystisa Pocket Controller Motor (Saeshin)

#### **Includes:**

- MYSTISA POCKET CONTROL BOX
- MYSTISA HANDPIECE
- Length - 5" x Height - 1 1/3" x Width - 2 1/3"
- A/C charger (110v or 230v)
- Leather pouch with belt clip
- Handpiece rubber stand
- Carry Bag
- Under Nail Carbide Bur
- Large Barrel Carbide Bur
- 10 Sanding Bands
- Sanding Band Mandre



**Mystisa Pocket Micromotor Set Complete Set** 35-118851

**Mystisa Handpiece Only Each** 35-118861

**Mystisa Pocket Controller Only , Gray Each** 35-118871

**110V-230V Universal Charger For Mystisa Pocket Cont Each** 35-118881

**Gray Leather Pouch For Mystisa Controller Each** 35-118841

**Mystisa Pocket Controller, Thumb Lab HP** 35-120000

**Mystisa Pocket Controller, Thumb E-Type HP** 35-120002

### Mystisa Pocket Drill (Saeshin)

- Hands Free Portability
- Speed range from 3,000 to 30,000 Rpm
- Runs om AC or battery power
- Runs 5 to 7 hours on a single charge
- Forward/Reverse direction
- 24V/30K Rpm



**Mystisa Table Micromotor Set Complete Set** 35-118891

**Mystisa Handpiece Only Each** 35-118911

**Mystisa Table Control Box Only, Blue Each** 35-118921

**Mystisa Table Control Box Only, Pink Each** 35-118931

**Mystisa Table Control Box Only, White Each** 35-118941

**Charger For Mystisa Table Control Box Each** 35-118901

### Nextgen Brushless Motor (Saeshin)

#### **Nextgen Brushless Motor (Ram)**

- Ergonomic design, Excellent durability, Minimal noise and vibration
- Right & Left turning ability
- Hand/Foot speed control system
- Built-in micro processor
- Self diagnosis, Error display, Work memory storage and handpiece auto recognition function
- Variable speed foot control (Optional)
- Max, 40,000 rpm
- Input 110 / 230 Volts (50 / 60Hz)
- Torque 4 N.cm
- Output 150W
- Weight Control box: 4.2Lb Handpiece: 8.0 oz. (Forte 300II)



**NextGen Micromotor Set with Brushless Handpiece Variable** 35-119001

**NextGen Micromotor Set with E-type Handpiece Variable** 35-118981

**NextGen Control Box Only Each** 35-118971

**NextGen Variable Foot Pedal Only Each** 35-118991



# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**NextGen Brushless Handpiece Only** 35-118961  
Each

**NextGen E-Type Handpiece Only** 35-118951  
Each

## Optimus 35/Optimus 45 HP (Saeshin)



**206 Control Box,** 35-119071  
OPTIMUS 35 HP ON/OFF

**206 Control Box,** 35-119081  
OPTIMUS 35 HP Variable

**206 Control Box,** 35-119091  
OPTIMUS 35 HP None

**206 Control Box,** 35-119101  
OPTIMUS 45 HP ON/OFF

**206 Control Box,** 35-119111  
OPTIMUS 45 HP Variable

**206 Control Box,** 35-119061  
OPTIMUS 45 HP None

## OZ Plus Brushless Motor (Saeshin)

RPM max.

50,000 rpm

Torque 7.8 N.cm

Input 110/ 230

Volts (50/60Hz)

Output 230W

Weight: Control

box: 5.5 lb

Handpiece: 8.9oz

(Forte100)

•Right & Left turning ability

•Hand/Foot speed control system

•Self-diagnosis, Error display, Work memory

•Storage function

•Variable speed foot control (optional)

**OzPlus Set With** 35-119141

Sensor 50K Set B Variable " with Pedal "

**Oz Plus Set With Sensor 50K Set C None "** 35-119131

without Pedal "

**Oz Plus Controller Only (Sensor 50K)** 35-119151

Each

**Oz Plus 50K Handpiece w/Sensor** Each 35-119161

**Variable Pedal Only For Oz Plus** Each 35-119121



## Portable Electric Motor (Saevang)



- Forward/Reverse display **LED**

- RPM Display (2 digit 7 segment)

- **Max Speed:** 30,000 rpm

- Rechargeable battery, when charging green LED blink, when being

Completed green LED on stay  
- LI-ION battery **Included**, 12V, 25W  
- Working **Max 8-4.5 Hours**  
- Over load warning LED  
Adaptor connector included

**Full Unit:** Electric Motor & Straight 35-105111  
Nose HP

**Full Unit:** Electric Motor & Contra 35-105112  
Angle Push Button HP

## Power Brush Type Motor (Saeshin)

- Variable speed dial system
- Right & Left turning ability
- Hand/Foot speed control system
- Automatic overload protection system
- Reset apparatus ability
- ON/Off and variable speed foot control
- Max, 45,000 rpm
- Torque 450 gf.cm
- Input 110 / 230 Volts (50 / 60Hz)
- Output DC 0-32W
- Weight Control box: 4.8Lb Handpiece: 6.7 oz. (strong 110LN)



**Rampower 35 Control Box,** 35-120021  
Tech2000 HP ON/OFF

**Rampower 35 Control Box,** 35-120011  
Tech2000 HP Variable

**Rampower 35 Control Box,** 35-120001  
Tech2000 HP None

**Rampower 35/45 Control Box Only** 35-119991

**Rampower 45 Control Box,** 35-119981  
45,000 RPM HP ON/OFF

**Rampower 45 Control Box,** 35-119971  
45,000 RPM HP Variable

**Ram Power 45 Control Box,** 35-119961  
45,000 RPM HP None

**Ram Power 35 CB with** 35-119951  
Optimus 35K HP ON/OFF

**Ram Power 35 CB with** 35-120061  
Optimus 35K HP Variable

**Ram Power 35 CB with** 35-120031  
Optimus 35K HP None

**Ram Power 45 CB with** 35-120041  
Optimus 45K HP ON/OFF

**Ram Power 45 CB with** 35-120051  
Optimus 45K HP Variable

**Ram Power 45 CB with** 35-119941  
Optimus 45K HP None

## RamPro500N Series (Saeshin)

- Computer controlled microprocessor
- Soft touch keypad control digital speed display
- Hand/foot mode and forward/reverse direction
- 0 - 35,000 or 45,000 RPMs
- Auto cruise control on/off button for continuous power in foot pedal mode



Call us for Quarterly Promotion

Prestigedentalproducts.com

## Handpieces

- Four memory speed reset buttons
- Dual handpiece output jacks
- Safety overload protection
- Forward/Reverse direction
- Hand/Foot selection

**RamPro500N Control Box,**  
Tech2000 HP ON/OFF 35-120081

**RamPro500N Control Box,**  
Tech2000 HP Variable 35-120091

**RamPro500N Control Box,**  
Tech2000 HP None 35-120101

**RamPro500N Control Box Only** 35-120111

**RamPro500N Control Box,**  
45,000 RPM HP ON/OFF 35-120121

**RamPro500N Control Box,**  
45,000 RPM HP Variable 35-120131

**RamPro500N Control Box,**  
45,000 RPM HP None 35-120071

### Rotex 782 Electric Micro Motor (Dentamerica)



**Rotex 782E**

- Ball bearing high torque electric motor
- 1,000 to 35,000 R.P.M. variable speed
- Double insulation heavy duty construction
- Automatic cooling fan for continuous high speed performance
- Forward/Reverse switch Standard package includes Foot Switch, Handpiece, and Handpiece Rest
- Circuit overload protection
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics, and steel
- Unit is available for 100/110V 50/60 Hz, or 220/240V 50 Hz operation

**1 year warranty**

**Unit with Electric Handpiece and**  
Variable Speed Foot Control 35-120211

**Unit with E-Type Straight, Contra An-**  
gle Attachment and Variable Speed  
Foot Control 35-120201

**Handpiece only with cable for Rotex**  
782 35-120191

**Handpiece only with cable for Rotex**  
782E 35-12018

**Straight Nose Cone Attachment, for**  
Rotex 782E 35-12017

**Contra Angle Attachment, For Rotex**  
782E 35-120161

**Replacement Carbon Brushes, 2/Pkg** 35-120151



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Strong 204 Brush Motor (Saeshin)

#### Control Box:

- 35,000 rms
- Torque: 280 gf.cm
- Input 110 - 230 volts
- Compact design



#### Tech2000 Handpiece:

- 35,000 RPM
- Twist type chuck release
- Improved dust seal protection
- Includes 3/32" chuck
- Multi-chuck feature (easy to change chuck)
- Optional chuck sizes available: 1/8", 1/16", or 3MM
- **On/Off Pedal**
- Maintains constant speed

**204 35K Control Box,**  
Tech2000 HP ON/OFF 35-120971

**204 35K Control Box,**  
Tech2000 HP None 35-120981

**204 45K Control Box,**  
45000 RPM HP ON/OFF 35-120991

**204 45K Control Box,**  
45000 RPM HP None 35-121001

### Strong 90 Carbon Brush Motor (Saeshin)

**Basic model, Mini-Compact design,**  
**Less vibration standard carbon**  
**brush motor**

- Variable Speed dial system
- Right & Left turning ability
- ON/Off Foot control
- Max, 35,000 rpm
- Torque 280 gf.cm
- Input 110 / 230 Volts (50 / 60Hz)
- Output DC 0-32W
- Weight: Control box: 2.9 Lb  
Handpiece: 7.1 oz. (Strong 102)



**SS90 Control Box,**  
Tech2000 HP ON/OFF 35-121081

**SS90 Control Box,**  
5000 RPM HP ON/OFF 35-121101

**SS90 Control Box Only** 35-121111

**SS90 On/Off Foot Pedal Only ON/OFF** 35-121091

### Thumb Portable Micromotor (Saeshin)

#### Saeshin Product

Thumb Portable Micromotor  
24V 30,000 Rpm (Ram)

#### Set Consists of:

- Thumb Control Box,
- Charger, Thumb Pouch,
- Either an E-Type Handpiece with a 1:1 Straight Nose Cone or Lab Style Handpiece



- Recharger Battery takes 5 Hours to Fully Charge, **Up to 8 Hours Battery Life**, Variable Speed Control
- Light Weight Handpiece, Chuck Locking, Easy Bit Change, High Torque, Quiet Operations, No Vibration

## Handpieces

<b>Thumb Lab Type Portable Set</b>	35-121201
<b>Thumb E-Type Portable Set With</b> Straight Nose Cone Set	35-121191
<b>Thumb Lab Style Motor Handpiece</b> Only	35-121181
<b>Thumb E-Type Motor Handpiece</b> Only	35-121171
<b>Thumb Control Box Only</b> Each	35-121161
<b>Thumb Charger/Adaptor</b> Each	35-121211
<b>Thumb Pouch</b> Each	35-121231
<b>Thumb Battery</b> Each	35-121221

### V35 Electric Lab Handpiece System 110V (Buffalo)

#### **Features/Benefits:**

- Reliable Gray Electric Handpiece, New design and high-torque motor provide smooth operation and powerful performance with less operator fatigue. Larger bearings, extra-strong chucking mechanism, and proven design promote increased reliability
- Variable Speed Foot Control, Variable speed foot control provides fully variable speed operation from 0 – 35,000 rpm and allows the system to excel at both higher speed and ultra-low speed (0 – 500 rpm) operation
- Dependable V35 Control Console, Four multi-speed options on the control console (14k, 22k, 27k, 35k rpm) simplify console electronics, improve system reliability, and allow ideal speed selection for most cutting needs. Simplified design creates higher power transfer through the console and more torque to the handpiece. Console


**Made in USA!**

<b>120V AC</b> Console + Gray HP + Foot Control	35-12154
<b>220V AC</b> Console + Gray HP + Foot Control	35-12153



### X35 Premium Electric Lab Handpiece System (Buffalo)

- Upgraded Handpiece:** The new black signature handpiece features an improved heavy-duty spindle, upgraded bearings, and more durable square carbon brushes, providing smoother operation and increased reliability
- Improved Console:** The upgraded control console features dual power outlets allowing fast and easy switching between two separate handpieces as desired (one handpiece included with system), saving time-wasting tool changes. Added overload protection helps protect the handpiece from damage during extended heavy grinding
- Superior Control:** The system provides exceptional control at ultra-low Rotational speeds, high torque cutting power, and fully variable speed from 0-35,000 rpm. The X35 comes complete with a fully variable speed foot control for added versatility
- Heavy Duty Construction:** Components and system designed to provide reliable and trouble-free operation.



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**Handpiece:** 6 in. (153 mm) Long x 1 1/8 in. (29 mm) Major Diameter, 5/8 in. (16 mm) Nose Diameter x 9.5 oz. (270 gm) Weight  
**Console:** 7.25 in. (184 mm) Width x 7 in. (177 mm) Depth x 3.5 in. (89 mm) Height x 6.3 lb. (2.9 kg) Weight  
**1 year warranty**

<b>120V AC</b> (Console + Black HP + Foot Control)	35-12156
<b>220V AC</b> (Console + Black HP + Foot Control)	35-12155

### X40 Premium Electric Lab Handpiece System (Buffalo)

#### **New 40,000 RPM High-Torque System Updates Our Classic X35 and Provides Additional Speed and Torque**

- New 40,000 RPM Handpiece 40k blue electric handpiece features an improved heavy-duty spindle, upgraded bearings, and more durable square carbon brushes, providing higher speeds, smoother operation, and increased reliability
- Improved Console The upgraded control console features dual power outlets allowing fast and easy switching between two separate handpieces as desired (one handpiece included with system), saving time-wasting tool changes. Higher voltage output and additional overload control help protect the handpiece from damage during extended heavy grinding
- Superior Control:** The system provides exceptional control at ultra-low rotational speeds, high torque cutting power, and fully variable speed from 0-40,000 rpm. The X40 comes complete with a fully variable speed foot control for added versatility
- Heavy Duty Construction:** Components and system designed to provide reliable and trouble-free operation.

**Handpiece:** 6.25 in. (159 mm) Long x 1 1/8 in. (29 mm) Major Diameter x 5/8 in. (16 mm) Front Diameter x 8.6 oz. (244 gm) Weight  
**Console:** 7.25 in. (184 mm) Width x 7 in. (177 mm) Depth x 3.5 in. (89 mm) Height x 6.3 lb. (2.9 kg) Weight  
**1 year warranty**

<b>120V AC</b> (Console + 40k rpm HP + Foot Control)	35-12158
<b>220V AC</b> (Console + 40k rpm HP + Foot Control)	35-12157
<b>X40 40,000 rpm Blue Electric Handpiece Only- New!</b>	35-63021



### X50 Brushless Electric Lab Handpiece System (Buffalo)

- Ultra-High Torque Cutting Power Brushless high-torque motor and control system provide brute-strength 7.8 N. cm. torque and precise control from 1,000 to 50,000 rpm with no carbon brush maintenance
- New Whisper-Quiet Motor New brushless motor features ultra-quiet and extremely smooth operation. Improved bearings and balance provide highly-





## Handpieces

concentric, vibration-free operation in a sealed, maintenance-free design. Convenient twist lock/unlock chuck provides safe and secure holding of HP-sized instruments

•New Advanced Console Design with Microcomputer Control New intelligent microprocessor control provides fully variable speed via console or foot control. Console provides precise speed control, digital rpm display, and controls for hand/foot and forward/reverse operation.

•Foot Control "Smart Cruise" Feature New smart cruise control circuitry automatically senses when the Foot Control speed is held constant and auto-activates precise rpm control with an audible two beep tone. Feature immediately deactivates with next touch of Foot Control and ready for resetting as desired

•Heavy Duty Construction Components and system designed to provide safe, reliable and trouble-free operation. High impact housing is tough and easy to clean.

•Improved Safety Intelligent safety circuitry is designed to protect the computer controls and the handpiece. The safety-lock on the console speed control dial helps prevent unintentional higher speed operation

**Handpiece:** 6.25 in. (159 mm) Long x 1.125 in. (29 mm) Back Diameter, .75 in. (19 mm) Front Diameter x 8.2 oz. (232 gm) Weight

**Console:** 7 in. (178mm) Width x 9.5 in. (241 mm) Depth x 4.25 in. (108 mm) Height x 7 lb. (3.2 kg) Weight

**1 year warranty**

**120V AC 50,000 RPM Electric Handpiece System** 35-12160

**220V AC 50,000 RPM Electric Handpiece System** 35-12159

### Xi Brush Type Motor (Saeshin)

- Compact / Eegonomic design
- Vertical / Horizontal placement
- Right & Left turning ability
- Hand/Foot speed control system
- Circuit safety system
- Variable Speed dial system (209A)
- ON/Off and variable speed foot control
- Max, 35,000 rpm
- Torque 300 gf.cm
- Input 110 / 230 Volts (50 / 60Hz)
- Output DC 0-32W
- Weight: Control box: 2.6 Lb Handpiece: 8.2 oz. (Strong 110LN)



**XI Control Box, Tech2000 HP ON/OFF** 35-12165

**XI Control Box, Tech2000 HP Variable** 35-12164

**XI Control Box, Tech2000 HP None** 35-12163

**Control Box Only Gray 32V 35/45K** 35-12162

**XI Control Box, 45000 RPM HP ON/OFF** 35-12161

**XI Control Box, 45000 RPM HP Variable** 35-12168

**XI Control Box, 45000 RPM HP None** 35-12167

**XI Control Box, 45000 RPM HP, 2 Collets None** 35-12166



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

## Handpiece Bulbs

### LED Bulbs (Chapman-Huffman)



**LED Bulb** 16-120-00

Fits following Swivel Couplers Chapman-Huffman 20-079-00KAVO MultiflexMK-DENT QC5016, QC5016KKAVO.



**LED Bulb** 16-126-00

Fits KAVO Motors



**LED Bulb** 16-130-00

Fits SIRONA Quick Connectors



**LED Bulb** 16-134-00

Fits W&H Turbines



**LED Bulb** 16-132-00

Fits SIRONA Turbines and SIRONA Motors



**LED Bulb** 16-122-00

Fits BIEN-AIR Motors



**LED Bulb** 16-124-00

Fits BIEN-AIR Quick Connectors



**LED Bulb** 16-128-00

Fits NSK Quick Connector



**LED Bulb** 16-136-00

Fits STAR Quick Connector

**Fits following:** Swivel Couplers 35-118061  
Chapman-Huffman 20-079-00  
KAVO Multiflex MK-DENT  
QC5016, QC5016K KAVO  
cross ref. # 1007.5372

**Fits KAVO Motors** 35-118051

**Fits SIRONA Quick Connectors** 35-118041

**Fits SIRONA Turbines and SIRONA Motors** 35-118031

**Fits W&H Turbines** 35-118021

**Fits BIEN-AIR Motors** 35-118011

**Fits BIEN-AIR Quick Connectors** 35-118001

**Fits NSK Quick Connector** 35-117991

**Fits STAR Quick Connector** 35-118071

## Handpiece Turbines & Couplers - Swivel

### Athena Champion Push Button Turbine (HPP)



**Cap Wrench** 35-112711

**Push Button cap** 35-112721

**Complete Turbine Assembly** 35-112731

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Athena Champion Standard Turbine (HPP)



Cap Wrench	35-11266
Cap	35-112671
Bur Wrench	35-112691
Complete Turbine Assembly	35-112681

## Bien Air Black Pearl Friction Grip Turbine (HPP)



Complete Turbine Assembly	35-112771
---------------------------	-----------

## Bien Air Black Pearl Lab Turbine (HPP)



Bearing Kit	35-112821
Complete Turbine Assembly	35-112831

## Bien Air Black Pearl Push Button Turbine (HPP)



Complete Turbine Assembly	35-113011
---------------------------	-----------

## Bien Air Bora S36-S36L Push Button Turbine (HPP)



Bearing Kit	35-112801
Complete Turbine Assembly	35-112811

## Bien Air Ondine Friction Grip Turbine (HPP)



Bearing Kit	35-112951
Complete Turbine Assembly	35-112971

## Bien Air Ondine Push Button Turbine (HPP)



Bearing Kit	35-112891
Complete Turbine Assembly	35-112881

## Bien Air Prestige Push Button Turbine (HPP)

Bearing Kit	35-97572
Complete Turbine Assembly	35-97571

## Canister Turbine (HPP)



Push Button Turbine	35-11308
Push Button Turbine, Angular H	35-11306
Push Button Turbine, Angular M	35-11307
Push Button Turbine, Phenolic	35-11305
Push Button Turbine, Angular Phenolic H	35-11303
Push Button Turbine, Angular Phenolic M	35-11304
Push Button Turbine, Silver / Maintenance free	35-113101
Push Button Turbine, Maintenance Free	35-113091
Push Button End Cap	35-113111
Cap Wrench	35-11302

## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Canister Turbine Mini (HPP)



<b>Mini Screw Type Turbine H</b>	35-11320
<b>Mini Screw Type Turbine M</b>	35-11319
<b>Mini Screw Type Turbine, Angular H</b>	35-11317
<b>Mini Screw Type Turbine, Angular M</b>	35-11318
<b>Mini Screw Type Turbine, Phenolic H</b>	35-11316
<b>Mini Screw Type Turbine, Phenolic M</b>	35-11314
<b>Mini Screw Type Turbine, Angular Phenolic H</b>	35-11315
<b>Mini Screw Type Turbine, Angular Phenolic M</b>	35-11313
<b>Mini Screw Type Turbine, Silver / Maintenance Free</b>	35-113241
<b>Mini Screw Type Turbine, Maintenance free</b>	35-11321
<b>Mini Bur Wrench</b>	35-113221
<b>Mini Screw Type End Cap</b>	35-113231
<b>Cap Wrench</b>	35-11312

### Canister Turbine Standard (HPP)



<b>Screw Type Turbine H</b>	35-113471
<b>Screw Type Turbine M</b>	35-113451
<b>Screw Type Turbine, Angular H</b>	35-11344
<b>Screw Type Turbine, Angular M</b>	35-11343
<b>Screw Type Turbine, Phenolic H</b>	35-11341
<b>Screw Type Turbine, Phenolic M</b>	35-11342
<b>Screw Type Turbine, Angular Phenolic H</b>	35-11340
<b>Screw Type Turbine, Angular Phenolic M</b>	35-11338
<b>Screw Type Turbine, Silver / Maintenance Free</b>	35-11339

**Screw Type Turbine, Maintenance Free** 35-11336

**Bur Wrench** 35-11337

**Screw Type End Cap** 35-113461

**Cap Wrench** 35-11335

### Cartridge - Turbine for Handpiece (ND)

#### CARTRIDGE FOR AIRTURBINE HANDPIECE



**Made in Japan**

**Standard Chuck, for TC-40PM / 40PB & TC-80QM / 80QB** 35-11176

**Large Size, Standard Chuck, for TC-35YM / 35YB** 35-111751

**Push Button Chuck, for TC-450M / 450B & TCP-70QM / 70QB** 35-111741

**Push Button Chuck, for TC-350M / 350B** 35-111731

**Cartridge for TS-300ZT** 35-111732

### Impact Air Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-97581

**Complete Turbine Assembly with Cap** 35-97582

**Bearing Kit** 35-97583

### Kavo 6000 Series Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-97611

### Kavo 625/630/640 Standard Head Friction Grip Turbine (HPP)





## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**Bearing Kit** 35-114021  
**Complete Turbine Assembly** 35-114031

### Kavo 625/630/640 Standard Head Push Button Turbine (HPP)



**Bearing Kit** 35-97592  
**Complete Turbine Assembly** 35-97591

### Kavo 632/633 Miniature Head Friction Grip Turbine (HPP)

**Bearing Kit** 35-97602  
**Complete Turbine Assembly** 35-97601

### Kavo 632/642/643 Miniature Head Push Button Turbine (HPP)



**Bearing Kit** 35-113651  
**Complete Turbine Assembly** 35-113661

### Kavo 634 Miniature Head Friction Grip Turbine (HPP)



**Bearing Kit** 35-113811  
**Complete Turbine Assembly** 35-113801

### Kavo 634 Miniature Head Push Button Turbine (HPP)



**Bearing Kit** 35-11388  
**Complete Turbine Assembly** 35-113901

### Kavo 635 Mira Lux 3 Push Button Turbine (HPP)



**Push Button Cap** 35-11391  
**Bearing Kit** 35-113941  
**Complete Turbine Assembly** 35-113931  
**Cap Wrench** 35-11392

### Kavo 647/649 Lux 3 Standard Head Push Button Turbine (HPP)



**Push Button Cap** 35-113713  
**Bearing Kit** 35-113711  
**Complete Turbine Assembly** 35-113721  
**Cap Wrench** 35-113712

### KAVO 650 Type Ceramic Turbine (Sable)

**1 year warranty**  
**Made in Canada**



**Each** 35-54251

### Kavo 6500 Series Push Button Turbine (HPP)

**Complete Turbine Assembly** 35-97621

### Kavo 659 Compact Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-114141



[www.PrestigeDentalProducts.com](http://www.PrestigeDentalProducts.com)

Call us for Quarterly Promotion

Prestigedentalproducts.com

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Kavo Back Cap Flat Wrench



Each 35-10018

## Kavo Style Coupler (Dentex)

Non-Optic



Standard Head Non-Optic 35-117941

Mini Head Non-Optic 35-117951

## Kavo Style Fiber Optic Coupler (Dentex)



Standard Head Fiber Optic 35-117961

Mini Head Fiber Optic 35-117971

## Lares 557 Small Head Push Button Turbine (HPP)



Push Button Cap 35-11416

Bearing Kit 35-114171

Complete Turbine Assembly 35-114181

Back Cap Wrench 35-11415

## Lares 757 Large Head Push Button Turbine (HPP)



Complete Turbine Assembly 35-114211

Bearing Kit 35-114201

Back Cap 35-11419

Cap Wrench 35-12419



## Midwest 8000 - 8000i Standard Turbine (HPP)



Complete Turbine Assembly 35-114341

Bearing Kit 35-11432

Bur Wrench 35-114331

Cap Wrench 35-11422

## Midwest Quiet Air Lever & Push Button Turbine (HPP)



Complete Turbine Assembly 35-114451

Complete Turbine Assembly with P.B. Cap 35-11442

Cap Wrench Lever Style 35-11443

Push Button Cap 35-114441

Cap Wrench 35-11436

Bearing Kit 35-11437

## Midwest Quiet Air Standard Turbine (HPP)



Complete Turbine Assembly with Cap 35-114701

Bearing Kit 35-11459

Bur Wrench 35-11468

Cap Wrench 35-11458

## Midwest Stylus Push Button Turbine (HPP)



Complete Turbine Assembly 35-114801

Push Button Cap 35-11477

## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

Cap Wrench	35-11478
Bearing Kit	35-114791

### Midwest Tradition Lever Turbine (HPP)



Complete Turbine Assembly	35-114911
Bearing Kit	35-11482
Cap Wrench Lever Style	35-11481

### Midwest Tradition Standard Turbine (HPP)



Complete Turbine Assembly with Cap	35-115161
Bearing Kit	35-11514
Bur Wrench	35-115151
Cap Wrench	35-11504

### Midwest Tradition XGT Push Button Turbine (HPP)



Complete Turbine Assembly	35-115281
Complete Turbine Assembly with Cap	35-11525
Bearing Kit	35-11526
Push Button Cap	35-115271
Cap Wrench XGT Lever	35-11517

### Ney Hurricane Friction Grip Turbine (HPP)

Complete Turbine Assembly	35-97631
Bearing Kit	35-97632

### NSK 75NL Push Button Small Turbine (Micotech)



Small Size for, NSK, NL75/85 Mini Triple Port 35-109691

### NSK 75NL Push Button Turbine (Micotech)



Complete Turbine Assembly 35-97671

### NSK Mini Phatelus Push Button Turbine (HPP)



Complete Turbine Assembly	35-116011
Bearing Kit	35-116021

### NSK Pana Air Mini Canister (HPP)



Complete Turbine	35-115411
Bearing Kit	35-11537
Bur Wrench	35-11539
Back Cap	35-115401
Back Cap Wrench	35-11529

### NSK Pana Air Mini Head Push Button Canister (HPP)



Complete Turbine	35-113341
Bearing Kit	35-113321




## Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

Back Cap	35-11331
Cap Wrench	35-11325

### NSK Pana Air Standard Head Canister (HPP)



Complete Turbine Assembly	35-11511
Bearing Kit	35-11548
Bur Wrench	35-11549
Back Cap	35-11501
Cap Wrench	35-11543

### NSK Pana Air Standard Head Push Button Canister (HPP)



Complete Turbine Assembly	35-97641
Bearing Kit	35-97642
Cap Wrench	35-97643
Cap	35-97644

### NSK Push and Swivel Push Button Turbine (HPP)



Complete Turbine Assembly	35-11559
Bearing Kit	35-115601

### NSK Push and Swivel Standard Turbine (HPP)

Complete Turbine Assembly	35-97651
Bearing Kit	35-97652

### NSK Style Coupler (Dentex)

Non-Optic



Standard Head Non-Optic	35-119031
Mini Head Non-Optic	35-119021

### NSK Style Fiber Optic Coupler (Dentex)



Standard Head Fiber Optic	35-119041
Mini Head Fiber Optic	35-119051

### NSK Super Grade Push Button Canister (HPP)




Complete Turbine Assembly	35-115711
Bearing Kit	35-11562

### NSK Super Grade Standard Canister (HPP)



Complete Turbine Assembly	35-115811
Bearing Kit	35-115821
Bur Wrench	35-11572

### NSK Torque Phatelus Push Button Turbine (HPP)



Complete Turbine Assembly	35-97661
Bearing Kit	35-97662

### NSK Viper Phatelus Push Button Turbine (HPP)



Complete Turbine Assembly	35-115921
---------------------------	-----------

Call us for Quarterly Promotion

Prestigedentalproducts.com

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**Bearing Kit** 35-11590

## Quest Highspeed Handpiece Couplers (Lares)



**6 Pin Fiber Optic MX Coupler, Multiflex Compatible Coupler (with 0.7 Amp Bulb)** 35-119711

**5 Holes MX Coupler, Multiflex Coupler with Fiber Optic Bundle** 35-119721

**4 Holes Coupler, Multiflex Compatible Coupler** 35-119731

**MX Coupler Fiber Optic (1.5 Amp Bulb) (Lares Power Supply Compatible)** 35-119741

## Star Concentrix Standard Turbine (HPP)



**Complete Turbine Assembly** 35-116211

**Bur Wrench (Hex)** 35-116212

**Cap** 35-116213

## Star Futura I and II Standard Turbine (HPP)



**Complete Turbine Assembly (Square)** 35-116321

**Complete Turbine Assembly (Hex)** 35-116322

**Bur Wrench (Square)** 35-116323

**Bur Wrench (Hex)** 35-116324

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

## Star Ultra Standard Turbine (HPP)



**Complete Turbine Assembly (Square)** 35-116381

**Complete Turbine Assembly (Hex)** 35-116382

**Bur Wrench (Square)** 35-116383

**Bur Wrench (Hex)** 35-116384

## Star Vista 430 Solara Push Button Turbine (HPP)



**Complete Turbine Assembly (Greased Bearing)** 35-116641

**Complete Turbine Assembly (Lube-Free)** 35-116642

**Bearing Kit (Angular Greased Bearings)** 35-116643

**Bearing Kit (Lube-Free)** 35-116644

**Push Button Cap** 35-116645

## Star Vista 430 Standard Turbine (HPP)



**Complete Turbine Assembly (Square)** 35-116551

**Complete Turbine Assembly (Hex)** 35-116552

**Bearing Kit (Angular Greased Bearings)** 35-116553

**Bearing Kit (Lube-Free)** 35-116554

**Bur Wrench (Square)** 35-116555

**Bur Wrench (Hex)** 35-116556

## W&H Adec 890 Series Large Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-117001

## Handpieces

**Complete Turbine Assembly** 35-117002  
(Phenolic Bearing)

**Bearing Kit** 35-117003

### W&H Adec 890 Series Small Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-117071

**Complete Turbine Assembly** 35-117072  
(Phenolic Bearing)

**Bearing Kit** 35-117073

### W&H Adec Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-116721

**Complete Turbine Assembly** 35-116722  
(Phenolic Bearing)

**Bearing Kit** 35-116723

### W&H Adec Synea Large Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-116861

**Complete Turbine Assembly** 35-116862  
(Phenolic Bearing)

**Bearing Kit** 35-116863

### W&H Adec Synea Small Push Button Turbine (HPP)



**Complete Turbine Assembly** 35-116931

**Complete Turbine Assembly** 35-116932  
(Phenolic Bearing)

**Bearing Kit** 35-116933



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

## Handpieces Parts & Accessories

### 2-3 Holes Tubing Adapters (Parts)

**HP Tubing Adaptor, 4-Hole HP To 2/3-Hole Tubing:**



**Adapting:** Adapts

4-Hole HP to

2/3-Hole HP tubing, Provides Drive Air & Coolant Water, Does not provide Coolant Air.

**HP Tubing Adaptor, 4-Hole HP To 3-Hole Tubing:** Adapts

4-Hole HP to 3-Hole HP tubing only, Provides Drive air, water & air coolant

**2/3-Hole HP To 4-Hole Tubing:** Adapts 2 or 3-Hole HP to 4-Hole HP tubing

**HP Tubing Adaptor, 4-Hole HP To 2/3-Hole Tubing** 35-48311

**HP Tubing Adaptor, 4-Hole HP To 3-Hole Tubing** 35-48312

**2/3-Hole HP To 4-Hole Tubing** 35-48313

### 3 Holes Tubing Adapters (Parts)

**3 Hole Borden Adaptor, Plastic:** Lightweight plastic with brass fittings

**3 Hole Borden Adaptor, Metal:** Quality plated brass

**3 Hole Borden Nut & Plastic Connector:** Plated

brass nut and plastic connector, Brass fittings

**3 Hole Borden Nut & Metal Connector:** Plated

brass nut and connector

**3 Hole Borden Nut & Metal Connector:** Plated

brass nut and connector

**3 Hole Borden Nut & Metal Connector:** Plated

brass nut and connector

**3 Hole Borden Adaptor, Plastic** 35-44271

**3 Hole Borden Adaptor, Metal** 35-44272

**3 Hole Borden Nut & Plastic Connector** 35-44273

**3 Hole Borden Nut & Metal Connector** 35-44274

### 4 Hole Tubing Adapters (Parts)



**4 Hole Midwest Adaptor, Plastic** 35-44811

**4 Hole Midwest Adaptor, Metal** 35-44812

**4 Hole Midwest Nut & Plastic Connector** 35-44813

**4 Hole Midwest Nut & Metal Connector** 35-44814



## Handpieces

### 4 Holes Tubing Connectors (Parts)

- Handy for extending HP or Foot Control Tubing
- Plastic body with brass barbs



Each

35-4841

### Bearing Puller



TST710



TST713

Bearing Puller

35-111531

Heavy Duty

35-111521

### Bien Air Cap Wrench (HPP)



Part #: HBA130



Part #: HBA230



Part #: HBP130



Part #: HBP230

Outer Cap Wrench for  
Black Pearl & Bora

35-97781

Inner Cap Wrench for  
Black Pearl & Bora

35-97782

Outer Cap Wrench for Prestige

35-97783

Inner Cap Wrench for Prestige

35-97784

### Bur Changer



Tradition



Pana Air/Naminex/Glow



Quiet Air



Pana Air/Naminex/Glow Mini

Fits: Midwest Tradition, each

35-92231

Fits: Midwest Quiet Air, each

35-92232

Fits: Pana Air/Naminex/Glow, each

35-92233

Fits: Pana Air/Naminex/Glow Mini,  
each

35-92234

### Cap Wrenches



WR100



WR100-1

Universal Back Cap Wrench

35-111481

Back Cap Wrench

35-111471



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

### Fiber Optic Light Source - Built In



H 7070  
Fiber Optic Light Bulb  
Cartridge



H 7040  
Fiber Optic 5-Hole Tubing



H 7020  
Control Board



H 7045  
Fiber Optic 6-Pin Tubing



H 7090  
Fiber Optic Transform



H 7660 - 6 Pin



H 7000 - 5 Hole

Package includes: H7020+H7045+H7090 Package includes: H7020+H7040+H7070+H7090

**6 Pin Package Includes:** Control  
Board, Fiber Optic HP Tubing &  
Transformer

35-105971

**5 Pin Package Includes:** Control  
Board, Fiber Optic HP Tubing,  
Transformer, and Bulb Cartridge

35-105972

**Control Board** (Accept 2 FO Tubing)

35-105973

**Fiber Optic HP Tubing** - 6 pin, Grey

35-105974

**Fiber Optic HP Tubing** - 5 pin, Grey

35-105975

**Fiber Optic Light Bulb Cartridge**

35-105976

**Fiber Optic light source Transformer**

35-105977

### Handpiece Gaskets (HPP)



Part #: HHX4HG  
4 HOLE GASKET



Part #: HHX5HG  
5 HOLE GASKET



Part #: SST004  
WATER LINE GASKET

**4 Hole Gasket, each**

35-117241

**5 Hole Fiber Optic Gasket, each**

35-117231

**2 Hole Gasket, each**

35-11722

**4 Hole Gasket, 5/Pkg**

35-11721

**5 Hole Fiber Optic Gasket, 5/Pkg**

35-11720

**2 Hole Gasket, 5/Pkg**

35-11719

### Handpiece Pull Tester (HPP)

A pull tester is an absolute must when repairing/rebuilding turbines. Our pull testers are simple, but very effective. A precision pin is mounted to a spring loaded shaft that exerts a pull of 8 lbs which more than exceeds the industry standard of 6 lbs.

The standard pull tester is to be used on a turbine that is not in motion.

The dynamic pull tester can be used like the standard, however the test pin is attached that allows you to perform the pull test on a handpiece while it's running.

This is the most accurate way to test the holding power of a spindle.



Part #: THR008



Part #: THR808

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**Dynamic** 35-97682

## Impact Air Cap Wrench (HPP)



**Each** 35-97751

## Kavo Back Cap Flat Wrench



**Each** 35-10018

## Kavo Cap Wrench (HPP)



Part #: HKC130



Part #: HKC230



Part #: HK6230

**Cap Wrench** 35-97721

**Heavy Duty Cap Wrench** 35-97722

**Kavo 635, 647, 649 Cap Wrench** 35-97723

## Kavo Handpiece Spare Parts



TST609



TST657



TST658



TST659

**Kavo Cap Wrench** 35-117931

**Kavo 625/630 Head Expander** 35-117921

**Kavo 634/642 Head Expander** 35-11791

**Kavo 640B Head Expander** 35-11790

## Kavo Slinger and Spray Cover Tools (HPP)



Part #: HK5222



Part #: HK7222



Part #: HKS240

**Kavo 635 Slinger Install Tool** 35-97701

**Kavo 647-649 Slinger Install Tool** 35-97702

**Kavo Dual-Sided Spray Cover Tool** 35-97703

## Lamp Module 6-Pin Style - DCI

- Use with 6-pin current DCI Style power optic tubing
- 3.5VDC, 720 MA
- Fits Current Style DCI
- **RED** Gasket
- Use E-1581 for Kinetics, Proma, Old Style DCI System



**Bulb Cartridge, Power Optic, each** 35-108851

## Lamp Module 5 Hole Style - Kinetics

- Use with 6-pin ISO Style Power-Optic Tubing
- 3.5VDC, 720 MA
- Fits Kinetics, Proma, & previous style DCI
- **BLACK** Gasket
- Use E-1581 (8813) for current DCI systems.



**Bulb Cartridge, for Kinetics St** 35-108861

## Lares Cap Wrench (HPP)



Part #: HL5135



Part #: HL5130

**Heavy Duty Dual Cap Wrench** 35-97741

**Lares 557 Cap Wrench** 35-97742

## Midwest Cap Wrench (HPP)



Part #: HMQ130



Part #: HMQ230



Part #: HMT230



Part #: HMT330



Part #: 753333

**Quiet Air Standard Cap Wrench** 35-97711

**Quiet Air Standard Cap H/D Cap Wrench** 35-97712

**Tradition Standard H/D Cap Wrench** 35-97713

**Multi-Tool Cap Wrench works on:** 35-97714  
XGT, Tradition (L&PB), Quiet Air (L&PB), Stylus PB

**Tradition-L, Quiet Air Lever Cap Wrench (OEM)** 35-97715

## Midwest Quiet Air Bur Changer

- Silver color
- High Quality
- MADE IN THE US



Quiet Air

Call us for Quarterly Promotion


Prestigedentalproducts.com

## Handpieces

Silver Color, each	35-111041
2 Pieces, Silver Color	35-111042
5 Pieces, Silver Color	35-111043

### Midwest Tradition Handpiece Bur Changer

- Silver color
- High Quality
- MADE IN THE US



**Tradition**

Blue Color, each	35-111051
2 Pieces, Blue Color	35-111052
5 Pieces, Blue Color	35-111053

### Midwest Wrenches



**WR-A141**



**WR140**



**WR141**

Midwest Quiet Air Back Cap Wrench	35-118781
Midwest Tradition Bur Wrench	35-118791
Midwest Quiet Air Bur Wrench	35-118771

### NSK Chuck-Type Wrench



Each	35-119011
------	-----------

### Star Cap Wrench (HPP)



Each	35-97761
------	----------

### Turbine Rebuild Kits

Midwest Type Quiet Air, All Models	35-121311
Midwest Type Tradition, All exc. XGT	35-121301
Star 430 Type, All Models	35-121291
Kavo 640/642 Type	35-121281
Kavo 635B Type	35-121271
Kavo 647/649B Type	35-121261
NSK Phatelus/Viper Type	35-121251



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### W&H Adec Cap Wrench (HPP)



Part #: HAS136



Part #: HAS137



Part #: HAS138

Synea TA-96 Cap Wrench	35-97771
Synea TA-97 Cap Wrench	35-97772
Synea TA-98 Cap Wrench	35-97773

### Handpiece Lubricants & Cleaners

### AR Spray Handpiece Lubricant (J Morita)

- Maintenance of Air Turbine
  - Scaler & Motor Handpieces Spray
- Made in Japan**



Spray, 420ml	35-109481
--------------	-----------

### EasyOil Lubrication Adapter (Sable)

The best lubrication adapter. Ever.



2201401  
KaVo  
EasyAdapter



2201402  
W&H  
EasyAdapter



2201403  
Bien Air  
EasyAdapter



2201404  
Sirona  
EasyAdapter



2201405  
NSK  
EasyAdapter



2201406  
Midwest  
EasyAdapter

The universal solution for a trouble-free handpiece care. EasyOil is the universal support that can conveniently gather up to five different adapters to meet all your dental handpiece maintenance needs. An effective handpiece care results in fewer failures and great savings on repair costs. EasyOil allows you to take care of your handpieces in an easy, fast and comfortable way because all you need is already there, readily at hand. With just a few, simple gestures EasyOil will help you make the life of your handpieces longer.

Patented-100% **Made in Italy**

Kavo Easy Adapter	35-98771
W&H Easy Adapter	35-98772
Bien Air Easy Adapter	35-98773
Sirona Easy Adapter	35-98774
NSK Easy Adapter	35-98775
Midwest Easy Adapter	35-98776
Handpiece Maintenance Multi-adapter	35-98777



# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Handpiece Cleaning & Lubrication System (TPC)

### Handpiece Cleaning and Lubrication System

#### Product Features

- Economical, efficient and time-saving
- A fail-safe maintenance solution for your office
- Filtered air used to rotate handpieces for safety
- Care for most handpieces and all handpieces brands
- Disposable mist filtration system
- Solution level indicator
- Detachable door for easy cleaning of unit
- Easy replenishment of solution tank
- Air Pressure gauge
- Oil flow control valve

#### Saves Time

- Simple to Operate
- Optimizes Performance



HG6000  
(3 high speed positions)



HG6005  
• 2 high speed positions  
• 1 low speed E-type motor



HG6025  
• 1 high speed positions  
• 2 low speed E-type motor



**Content:** 3 high speed handpiece positions, Mist filters (pack of 12), Lubrication fluid, Funnel & Air filter 35-12580

**Content:** 2 high speed, 1 low speed E-type Motor positions, Mist filters (pack of 12), Lubrication fluid, Funnel & Air filter 35-12581

**Content:** 1 high speed, 2 low speed E-type Motor positions, Mist filters (pack of 12), Lubrication fluid, Funnel & Air filter 35-12582

**Lubrication Fluid** 35-12583

**Mist Filter 12/Pkg** 35-12584

**Oil absorbent door strips, 20/Pkg** 35-12585

**Kavo Type HP adapter** 35-12586

**NSK Type HP adapter** 35-12587

**Bien Air Type HP adapter** 35-12588

**W&H Type HP adapter** 35-12589

## Handpiece Lubricant (Kavo)



**500ml Bottle** 35-16581

**Nozzle, each** 35-16582

## Handpiece Lubricant Oil (HPP)

General purpose lubricating oil, for use in highspeed and slow speed applications.



Each 35-97731

## Handpiece Lubrication Oil Spray (Yoshida)

Cleaning & Lubrication for push-button ball-bearing latch head Contra Angle  
**Made in Japan**



Each 35-97551

## Handpiece Purge Station with Hepa Filter (Parts)

•Purges excess lubricant and water from handpiece

- Patented Bio-Cleanse HEPA exhaust filter
- Hepa Filter helps protect operator from air-borne pathogens
- Provides for quick and easy handpiece maintenance
- Exhaust reservoir safely collects excess lubricant
- Purges handpiece coolant air and water lines
- Easy view gauge with adjustable air regulator
- Air supply line with 1/4" male quick-disconnect included

- Compact size of a foot control - requires little counter space
- Handpiece adapter required from below

Each 35-117371

**Midwest Adapter** 35-117361

**Star Adapter** 35-11735

**Kavo Adapter** 35-11734

**NSK, Kinetics, Healthco, Viper, Atech Adapter** 35-11733

**Borden Adapter** 35-11732

**Midwest XGT Adapter** 35-11731

**Lares Adapter** 35-11730

**W&H Adapter** 35-11729

**Single Hepa Filter** 35-11728



Like us on Facebook

# Handpieces



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Lubrex Handpieces Cleaning and Lubricating System (Dentamerica)

- Automatic cleaning, purging and lubricating system
- Ergonomic front-loading, easy access design
- Can be used for all high speed, low speed, air scaler, and air motor handpieces
- Up to 4 handpieces can be processed simultaneously
- Quick release air connector
- Automatic fan and filter for vapor removal
- Large reservoir capacity with visible oil level indicator
- Specifications : 100~240V 50/60 Hz AC
- Dimensions : W330 x H430 x D430 mm
- Weight : 13 kg
- Accessories : Main unit, Lubricant, Quick Connector, Hose, Power Cord



**Unit** 35-118081

## Lubricant & Clean Spray

- Effective dental handpieces lubrication
- 2 to 3 seconds twice a day.

**Spray Can 18oz, 550ml** 35-118091

## Lubricant Spray for Kavo (DA)

For Kavo attachments



**12oz Spray Bottle** 35-99031

## Lubrina Automatic Handpiece Lubrication System (J-Morita)

Lubrina cleans and purges handpieces, and a variety of other equipment, via a rapid air blow system. The speed of air flow has been substantially increased resulting in more power. It delivers superior performance in the same amount of time – just 25 seconds per handpiece. Lubrina also offers a long blow out mode, if needed, to help ensure proper maintenance and maximize service life. Rated as highly versatile in 2009 for its setting options, four ports, and compatibility with multiple handpiece and lubricant brands, Morita continues to improve Lubrina's versatility. In addition to 4-hole and 6-pin, it is now compatible with a standard 5-hole connection. Lubrina is also currently the ONLY system on the market that accepts multiple brands of spray lubricant including the use of two different sprays in the same cycle. Powerful, versatile, and easy to operate, Lubrina simplifies the maintenance process and keeps handpieces running smoothly.



## Flexible

Only unit that accepts multiple brands of spray  
Allows two different sprays per cycle  
Accommodates wide variety of attachments & handpiece types  
Short and long mode options

## Saves time & Money

Uses only 1/3 the lubricant (compared to manual lubrication)  
Services up to 4 handpieces in only 100 seconds  
Protects your equipment, helps keep it running smoothly


## Unique

Built-in chuck maintenance system  
Air suction system eliminates mist for clean operation

Lubrina is compatible with the following connections:



6-Pin Connection



4-Hole Connection



5-Hole Connection

**Unit Includes:** AR Spray (420ml), 35-118101  
AR Spray Can Stand, Oil Absorption Pad (5) Oil Absorption Case, Tray, Tray Filter, Air Tube, Power Cord

**Lubrina 4H (Morita Type) Coupling** 35-118102

**Lubrina TR (E-Type Joint) Coupling** 35-118103

**Lubrina Alpha Coupling** 35-118104

**Lubrina Root ZX II Coupling** 35-118105

**Lubrina Tri Auto ZX Coupling** 35-118106

**Lubrina W&H Coupling** 35-118107

**Lubrina Kavo Coupling** 35-118108

**Lubrina Sirona Coupling** 35-118109

## Once-A-Day Handpiece Spray (Palmero)

TopCat Once-A-Day 1 Second Spray, For use in air turbine hand pieces, prophyl angles and contra angles. Can also be used as a lubricant/cleaner for hand pieces. Spray daily. No corrosive (contains hexane). 8.8 oz. aerosol can with extension tube.



**8.8 oz can** 35-99021

## Pana Spray Handpiece Lubricants and Conditioners (NSK)

- Handpiece Lubricants and Conditioners Spray (NSK America)
- PANA Spray Plus is a **specialty formulated general handpiece and air motor lubricant** containing alcohol for better cleaning efficacy.
- Apply after handpiece use to **assure high & low speed handpiece and air motor** performance and durability



**Spray 16oz** 35-94231

## Handpieces

### Spray Lube Surgical Instrument Lubricant (Miltex)

Oil 8oz Anti-Corrosive



8oz oil/bottle

35-99011

### Waterline Treatment

#### BluTab Waterline Maintenance Tablet (Confirm Monitoring)

- Specially designed to be continuously present in your waterlines
- Eliminates the need for weekly shock treatments
- Odorless, tasteless, non-toxic, and safe for equipment
- No need to purge lines at the end of the day



50 tabs/Box for 2 liters of water

35-94152

#### Citrisil Dental Waterline Cleaner (Sterisil)

- Non-toxic, non-corrosive, no taste, no odor
- Just drop tablets into water reservoir and fill with water
- Long lasting
- Eliminates the need to empty bottles and purge and air-dry



**20 Maintenance Tablets & 1 Shock Tablet:** White For 0.7-1 liter bottles Mfg. Part #: C20-W

35-94156

**20 Maintenance Tablets & 1 Shock Tablet:** Blue For 0.7-1 liter bottles Mfg. Part #: C20-B

35-94157

#### DentaPure 365 Day Waterline Treatment Cartridge (Crosstex)



The technology in DentaPure was developed for NASA to ensure water consumed in space was free of bacteria, viruses and any other harmful organisms. It is the

only iodinated water treatment product that is EPA registered as an "Anti-Microbial Waterline Product." DentaPure treats dental unit waterlines from the first day of installation for the life of the cartridge, reducing harmful bacteria exposure to patients and staff. DentaPure is safe for your patients, your staff and the environment. It is the easiest to use and most reliable dental water purifier requiring minimal human interaction.

**"B" Series** - Designed to install into the units independent water supply bottle.



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

DP40B: 60-Day Independent Water Bottle Cartridge  
DP365B: 365-Day Independent Water Bottle Cartridge  
**"M" Series** - Designed to install into the chairs municipal system.  
DP365M: 365-Day Municipal Water Cartridge

**365-Day Independent Cartridge,** 35-108141  
**"B" Series**

**60-Day Independent Bottle Cartridge, "B" Series** 35-108142

**365-Day Municipal Cartridge,** 35-108143  
**"M" Series**

#### Mint-A-Kleen Waterline Disinfectant (Anodia Systems)

- EPA/FDA Registered.
- Fast and ready to use—no mixing or diluting necessary.
- 0.12% Chlorhexidine Gluconate—active ingredient directly contacts waterline walls and bottles for thorough treatment, cleaning, and protection of waterlines and equipment.
- Formulated without harsh chemicals.
- Compatible with all dental units.
- Proven mint-blue, liquid solution—does not corrode, clog, or leave residue in waterlines, bottles, or equipment.



Mint-A-Kleen's weekly treatment prevents bad taste in the water, unpleasant odors, slimy and clogged waterlines, bottles, handpieces and syringes. **Available for sale in ALL 50 States.**

**Case of 10 x 16 oz Bottle**

39-11261

#### Purity Waterline Cleaner (Cory Lab)



The dental waterline cleaner for continuous cleaning of self-contained water systems

**60 Tables/Box**

35-105471

### Variety of Handpieces to choose from:







**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

#### ☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

#### ☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

