

Miscellaneous



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Miscellaneous



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Dapdish Dental Dappen Dishes (PlasDent)



Asst in White, Blue, Yellow and Red 43-11111
200 pcs/box

Asst in White, Blue, Yellow and Red 43-11112
1000 pcs/box

Dappen Dishes (TC Dental)

Dappen Dishes-Glass

Provides an attractive selection of colorful glass dappen dishes. Sturdy and heat resistant. Ideal for mixing acrylics and amalgams, resins sealants and other materials.



Glass, Assorted Colors, 12pcs/box 43-11113

Plastic, Assorted Colors, 200pcs/box 43-11114

Isopropyl Alcohol -70% (DA)



6 x 32 oz. Quart Bottles 43-89051

Nosk Nasal Filter (Meta)

NOSK is a comfortable and lightweight, ultra-thin dual nasal filter that instantly blocks pollen, dust, allergens, and other airborne particles from entering the body.



- Virtually Invisible
- Highly Effective
- Very Convenient

Why use NOSK?

Drug Free: NOSK provides a safe, drug free alternative to allergy or other medicines. No more pills, sprays or side effects, such as drowsiness or dry mouth.

Effective: Reduces the inhalation of pollutants, allergens, dust and other unwanted hazardous substances.

Comfortable: Designed by top Ear, Nose and Throat specialist, NOSK Fits comfortably in your nose, giving you the freedom to talk, and eat while providing unparalleled relief.

Invisible: Virtually invisible to those around you, a thin translucent plastic strip holds NOSK securely in place

10/Pkg 41-98901

Soft Hand Scrub Brush (DA)



- Convenient Grip tight handle makes for easy scrubbing
- The soft, flexible bristles don't irritate skin

Soft, White, 2/Pkg 43-89061

Sticky Applicators (PlasDent)



White, 60 pcs/box 43-11166

Magnifiers & Loupes

LED Headlight System (LumaDent)

Increase efficiency, lighted field wherever you look without shadow. Save your back, the better you see, the faster you can perform

Complete LED Headlight

- Weight: 0.19 ox (5.5 gr). One of the lightest headlamp on the market, it is lighter than a quarter.
- Diameter: 0.64" (16.3 mm), Length: (20mm). Smallest in the market
- Light Output: 40,000 Lux at 14" working distance
- Spot Size: 2.5" at 14", 3.5" at 18"
- Battery Life: 8 hours continuously at maximum intensity, with Low-Battery indicator light
- Battery Charging Time: 3 hours
- Battery size: 3" x 2" x 0.9" (8.8cm x 5.5 cm x 2.3cm)
- Battery Weigh: 5.8 oz (155 gr)

Package includes: LED headlamp, Battery Pack, Charger, Mounts for your loupes/frame, wire accessory)



Comple System: Headlamp, 43-11128
Composite Filter, Mount, Battery Pack, Charger & Carrying Case

Battery Pack: 1 x Litium Ion 43-11129
Battery Pack 1x AC/DC Charger

Composite Resin Orange Filter 43-11130

Battery Pack-Lenser H7, Includes: 43-11131
Battery Pack 500 mA AC/DC
Charger Carrying Case with Swivel
Belt Clip



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Miscellaneous

iZoom Headlight System (DentLight)



Light Intensity	Level 1 = 25,000 lux Level 2 = 44,000 lux Level 3 = 70,000 lux	Level 1 = 16,000 lux Level 2 = 31,000 lux Level 3 = 45,000 lux	Level 1 = 35,000 lux Level 2 = 65,000 lux Level 3 = 110,000 lux
Weight	Lamp = 19 grams Battery = 87 grams	Lamp = 6 grams Battery = 87 grams	Lamp = 24 grams Battery = 87 grams
Size	21 mm (Diameter) x 32mm (Length)	18 mm (Diameter) x 21mm (Length)	28 mm (Diameter)
Beam Size (@12")	1.8 - 3.8 inch	2 inch	1.5 inch
Bulb Life	50,000 hrs	50,000 hrs	50,000 hrs
Replaceable Cord	Yes	No	Yes
Battery Run Time	Level 1 = 10 hrs Level 2 = 5.5 hrs Level 3 = 4 hrs	Level 1 = 12 hrs Level 2 = 7 hrs Level 3 = 5 hrs	Level 1 = 10 hrs Level 2 = 5.5 hrs Level 3 = 4 hrs
Recharge Time	2 hrs	2 hrs	2 hrs
Loupe Adapters	Adapts to all loupes	Adapts to all loupes	Adapts to all loupes
Part Number	8500010	8600110	8700010

Including LED headlamp with power cord, two rechargeable battery packs, universal clip-on adapters (to prescription glasses, safety glass, ant TTL), loupe L-adapter (to TTL), U-adapter (to Flip-up), mounting screws, clip-on battery pouch, universal

43-11117

18 mm (Diameter) x 21mm (Length) / 2 inch

43-11116

28 mm (Diameter) / 1.5 inch

43-11115

Telescopic Loupes (ATW)

Extra wide field of view, large depth of focus. High light efficiency. Advanced mounting system with 3 adjustments for individualizing posture, working angles and achieve perfect eye alignment. Lightweight frames with soft silicone nose pad. Adjustable head strap. Removable side shields. you can put the glass to fit your eye need.



Magnification +250, 17" Working Distance & 4.5" Field of View

43-11167

Magnification +300, 17" Working Distance & 3.5" Field of View

43-11168

Magnification +350, 17" Working Distance & 2.5" Field of View

43-11169

Miscellaneous



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Dental LED & Loupes Light (Miltex)







REF 15051
Flip-Up Dental Composite Filter
(PREVENTS COMPOSITE MATERIALS FROM CURING)

REF 15057
Ion Battery

REF 15056
Replaceable Cable

ASSEMBLY IS EASY FOR MILTEX® LOUPES!






1. Charge battery and place into clip battery pack.

2. Light will mount onto the loupes crossbar; simply push the grooved portion in between the locking knob.

3. Clip battery pack onto your belt and turn on.

4. Position LED Light and its ready for use!

MILTEX® LED LIGHT PART NUMBERS

REF	Description	Mount Type
15050	LED Loupes Light Unit	Miltex® loupes mount
15051	LED Loupes Light Unit	Orascope® loupes mount
15052	LED Loupes Light Unit	Sheervision® loupes mount
15053	LED Loupes Light Unit	Surgitel® loupes mount
15054	LED Loupes Light Unit	Large Universal mount
15055	LED Loupes Light Unit	Heine® loupes mount

NOTE: Not all mounts match all loupes. Contact your Miltex Representative to ensure compatibility.

Highest intensity of any portable dental & surgical Headlight
(80,000 lux @ 14 inches)

- Provides pure white light for true color rendition
- Variable light intensity control allows for professional to adjust light as needed
- Flip up dental composite filter available prevents composite materials from curing
- State-of-the-art lithium batteries (1 included) provides 4.5 Hours of use on a single charge
- Quick charger - 90% charge in 90 minutes reduced charging time
- Portable Headlamp weighs only 0.7 ounces lightweight reduced pressure
- Alternative manufacturer mounts available can fit other loupe brands

Unit: LED Headlight, Battery, Battery Charger, Battery Power Pack (attaches to belt), AC Adapter, Attach mount (Varies per REF) 43-11126

Miltex Loupes mount 43-11125

Orascope Loupes mount 43-11124

Sheervision Loupes mount 43-11123

Surgitel Loupes mount 43-11122

Large Universal mount 43-11121

Heine Loupes mount 43-11120

Flip-up Dental composite Filter 43-11119
(prevents composite materials from curing)

Ion Battery 43-11118

Replaceable Cable 43-11127

Nano Loupe Light (Dentlight)



The most power per pound, and the smallest overall light on the market with the most power at its size of only 6 grams on headlamp and 87 grams on battery pack.

Focused Beam - Unlike other headlights, iZoom has no glares with its unique beam quality giving your patient the most comfort.

Constant Power - With unique heat sink technology and constant current LED driver, Dentlight Nano maintains power throughout the battery run time.

Reliable - With smart battery, cooling design, flexible wire and plug and start power driver, Dentlight Nano works any time you need it.

Smart and Robust Battery - Black anodized aluminum rechargeable battery weighs only 87 grams with seamless backlit membrane buttons that set three power levels and prevents cross-contaminations with easy to disinfect clip-on battery pouch.

Ease of Use and Adaptation - Fit to all major loupes including Orascope, Zeiss, Designs for Vision, Surgitel and iZoom.

Kit Includes: slide-on composite filter, two battery packs, easy-to-clean battery pouch, universal charger, and loupe adapter with articulating arm that fits on all major loupes 43-11157

Headlamp Only 43-11156

Battery Pack 43-11155

Power Adapter 43-11154

Prestige Dental Products



Excellent Qualities
at Affordable Prices

877-772-3888

www.prestigedentalproducts.com



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Miscellaneous

Loupes (Miltex)



KIT CONTENTS

1. MAGNIFICATION LOUPES & FRAME
2. HEAD STRAP
3. FLIP PADDLE
4. LENS COVERS
5. STORAGE CASE
6. CLEANING CLOTH
7. CLEANING SOLUTION
8. SCREWDRIVER

Discipline	Your Height		
	under 5'7"	5'7" to 6'4"	over 6'4"
General Dental	2.5x/3.0x - Regular	2.5x/3.0x - Long	2.5x/3.0x - X Long
Hygiene	2.5x - Regular	2.5x - Long	2.5x - X Long
Endodontics	3.0x/3.5x - Regular	3.0x/3.5x - Long	3.0x/3.5x - X Long
Oral Surgery	3.0x/3.5x - Regular	3.0x/3.5x - Long	3.0x/3.5x - X Long



MILTEX® LOUPE PART NUMBERS

Loupe Description			Miltex® Part Number	Miltex® Part Number
Magnification	Working Distance	Field of View	RED FRAME	BLACK FRAME
2.5x	10"-13" (S)	3.5" (S)	1-5000	1-5011
	12"-15" (R)*	3.8" (R)	1-5001	1-5012
	14"-17" (L)	4.3" (L)	1-5002	1-5013
	16"-20" (XL)	4.6" (XL)	1-5003	1-5014
3.0x	10"-13" (S)	2.8" (S)	1-5004	1-5015
	13"-16" (R)	2.5" (R)	1-5005	1-5016
	16"-21" (L)	3.5" (L)	1-5006	1-5017
	18"-22" (XL)	4.0" (XL)	1-5007	1-5018
3.5x	13"-16" (R)	2.3" (R)	1-5008	1-5019
	16"-20" (L)**	3.0" (L)	1-5009	1-5020
	18"-22" (XL)	3.5" (XL)	1-5010	1-5021

* Most popular loupes for hygiene

** Most popular loupes for endodontic procedures

Based on a report in a 2003 Supplemental Issue of Access on Magnification Use In Dental Hygiene, Magnification loupes extend a dental clinician's career by decreasing the physical strains associated with achieving clear oral cavity visualization while providing superior patient care.

With Set and Forget™ optics, you do not have to be individually measured to get a custom Fit. All Miltex® optics allows you to Set your interpupillary, convergence and declination angle For perfect optical alignment. The optics are quickly and easily locked into your individual setting So your loupes are always perfectly aligned.

Step 1: Choose your magnification: If you are a new loupes user, start with the lowest power you need. Most dental professionals begins using loupes with a magnification power of 2.5x. Remember, as the magnification transmission increases, the field of view, focal range & light transmission decreases.

Step 2: Choose your working distance: All loupes are pre-set to focus at a certain distance and objects will remain in focus within a range of several inches around that pre-set point. Miltex® loupes are available in 4 working distances: Short (10-13"), Regular (12-16"), Long (14-18") and Extra Long (16-20").

The 250 Series (2.5x) is the highest ideal general purpose optic. The 250 series optic will provide a magnified view of approximately 2.5 times the object size. Most first time loupes users choose this magnification.

The 300 Series (3.0x) is an ideal optic for individuals who have worked with a 2.5x magnification power optic and want to step up in power. The 300 series optic provides a view of the entire quadrant of the mouth and has an excellent depth of focus.

The 350 Series (3.5x) is the highest magnification available in a compound lens design. The 350 series optic is an excellent choice for dentists requiring high magnification in lightweight loupes.

2.5x Magnification, Red, 10"-13" Working Distance, 3.5" Field View	43-11152
2.5x Magnification, Black 10"-13" Working Distance, 3.5" Field of View	43-11151
2.5x Magnification, Red 12"-16", Working Distance, 3.8" Field of View	43-11150
2.5x Magnification, Black 12"-16", Working Distance, 3.8" Field of View	43-11149
2.5x Magnification, Red 14"-18", Working Distance, 4.3" Field of View	43-11148
2.5x Magnification, Black 14"-18", Working Distance, 4.3" Field of View	43-11147
2.5x Magnification, Red 16"-20", Working Distance, 4.6" Field of View	43-11146
2.5x Magnification, Black 16"-20", Working Distance, 4.6" Field of View	43-11145
3.0x Magnification, Red 10"-13", Working Distance, 2.8" Field of View	43-11144
3.0x Magnification, Black 10"-13", Working Distance, 2.8" Field of View	43-11143
3.0x Magnification, Red 13"-16", Working Distance, 2.5" Field of View	43-11142



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Miscellaneous

3.0x Magnification, Black 13"-16", Working Distance, 2.5" Field of View	43-11141
3.0x Magnification, Red 16"-21", Working Distance, 3.5" Field of View	43-11140
3.0x Magnification, Black 16"-21", Working Distance, 3.5" Field of View	43-11139
3.0x Magnification, Red 18"-22", Working Distance, 4.0" Field of View	43-11138
3.0x Magnification, Black 18"-22", Working Distance, 4.0" Field of View	43-11137
3.5x Magnification, Red 13"-16", Working Distance, 2.3" Field of View	43-11136
3.5x Magnification, Black 13"-16", Working Distance, 2.3" Field of View	43-11135
3.5x Magnification, Red 16"-20", Working Distance, 3.0" Field of View	43-11134
3.5x Magnification, Black 16"-20", Working Distance, 3.0" Field of View	43-11133
3.5x Magnification, Red 18"-22", Working Distance, 3.5" Field of View	43-11132
3.5x Magnification, Black 18"-22", Working Distance, 3.5" Field of View	43-11153

Storage Boxes

Round Box (PlasDent)



Large Diameter 1 1/2" x Inner Deep 5/8" White 100 pcs/bag with labels	43-11165
Medium Diameter 1 1/2" x Inner Deep 1/2" White 100 pcs/bag with labels	43-11164
Small Diameter 1 1/2" x Inner Deep 5/16" White 100 pcs/bag with labels	43-11163
Medium Diameter 1 1/2" x Inner Deep 1/2" Blue 100 pcs/bag with labels	43-11162
Small Diameter 1 1/2" x Inner Deep 5/16" Red 100 pcs/bag with labels	43-11161
Bulk Large Diameter 1 1/2" x Inner Deep 5/8" White 1000 pcs/case	43-11160
Bulk Medium Diameter 1 1/2" x Inner Deep 1/2" White 1000 pcs/case	43-11159
Bulk Small Diameter 1 1/2" x Inner Deep 5/16" White 1000 pcs/case	43-11158

Dental Supplies



www.prestigedentalproducts.com



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

