

Saliva Ejector System



Prestige Dental Products, Inc

Prestigedentalproducts.com


Tel: 877-772-3888 Fax: 714-632-8821

Saliva Ejector System




Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Flexform Saliva Ejectors (PlasDent)



White with Clear Adapter, 100/Pkg 50-11111

Flexjet Saliva Ejector (PlasDent)



Soft Clear with Blue Tip, 100/Pkg 50-11117

Soft White with Blue Tip, 100/Pkg 50-11116

Soft Clear with White Tip, 100/Pkg 50-11115

Soft White with White Tip, 100/Pkg 50-11114

HVE Adapters (PlasDent)



Saliva Ejector to HVE adapters, 12/Pkg 50-11118

Converts 11mm to 16mm, 12/Pkg 50-11119

HVE Tips Screens (PlasDent)



Clear , 100 pcs/bag 50-11124

HVE Valves (Plasdent)



Lever Beige, each 50-11126

Lever Black, each 50-11127

Rotary Beige, each 50-11128

Rotary Black, each 50-11125

Replacement lever Beige, 5/Pkg 50-11120

Replacement lever Black, 5/Pkg 50-11121

Mirro-Vac Saliva Ejector Mirror (IndiGreen)

Reduce number of hands needed, instruments used and procedure time! Our Mirro-Vac Saliva Ejector Mirror combines evacuator and mirror functions into one efficient instrument. Ideal for sealants, air abrasion, bonding and other dry field procedures. Upper suction inlet relieves tissue grab and ensures anti-fog acrylic mirror stays clear—even under direct exhalation. Size 4 lens with 35° inclination delivers standard mouth mirror feel and function without compromise.




- Easy-lift film protects mirror surface until use.
 - Soft blue color eases patient anxiety.
 - Round, knurled handle improves grip.
- Fully disposable, one-piece construction eliminates need for sterilization. 6" long.

50/pkg

50-87861

Tongue Shield with Saliva Ejector (Dentsply)



Saliva ejector 3/ 44-102891
pkg

Adapter tubes 3/ 44-102892
pkg



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Saliva Ejector - Defend

Disposable saliva ejectors
smooth edges, flexible, comfortable.

Tip designed to ensure highest
suction coupled with maximum
patient comfort.



Made in Italy

Clear with Blue Tip, 100/Pkg	50-16161
White with White Tip, 100/Pkg	50-16162
Clear with Blue Tip, One case 10 x 100/Pkg	50-16163
White with White Tip, One case 10 x 100/Pkg	50-16164

Saliva Ejector Valves (Plasdent)



Beige, each	50-11131
Black, each	50-11132
Lever SE Valve, Beige	50-11135
Lever SE Valve, Black	50-11136
Rotary SE Valve, Beige	50-11137
Rotary SE Valve, Black	50-11138
Snap-On SE tip Beige, 5/Pkg	50-10000
Snap-On SE tip Black, 5/Pkg	50-10001
Clear, 72 pcs/bag One size fits all	50-11133

Saliva Ejector Traps

Disposable traps provides for
replacement metallic traps. No
more cleaning and avoids con-
tamination. Easy to re-
place.



72 pcs/bag	50-11134
-------------------	----------

Saliva Ejector System

Saliva Ejectors

Clear with blue tips, or
white with white tips.
Bends to the right curve for
best comfort to the patient.
Retains the mem-
ory.



White with white tips, 1000/Pkg	50-11140
Clear with white tips, 1000/Pkg	50-11139

Saliva Ejector Levers (Plasdent)



Beige color, 5/Pkg	50-11144
Black color, 5/pkg	50-11146

Saliva Ejectors with Mint Scented (ASA Dental)

Disposable saliva ejectors
smooth edges, flexible,
comfortable.

Tip designed to ensure
highest suction coupled
with maximum patient com-



100/Box	50-11129
1000/Box	50-11130

Snap-On SE Tip (Plasdent)



Black Silicone Autoclavable, each	50-11147
Beige Silicone Autoclavable, each	50-11148
Gray Silicone Autoclavable, each	50-11149

Saliva Ejector System



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Spand-Ezz Saliva Ejector System (Dentsply)

Non-intrusive design works seamlessly with the Spand Ezz cheek retractor.

Creates a dry working field, without being an impediment to the patient or practitioner.



5/Pkg disposable clips/tubes 50-11152

Total-Comfort Aspirators (IndiGreen)



Hard plastic tip aspirators are just that—hard! They pinch lingual surfaces, jab the mandible and aspirate soft tissue. Now there is a truly comfortable solution—so comfortable, it's patented!

Total-Comfort Soft Tip Aspirators feature cushioned, flexible, patient-friendly tips. Tested and proven to be gentle on soft tissue and bone while reducing soft tissue aspiration. Elastic tips are securely bonded for safety. Protected by U.S patent 6,068,477 and other patents pending. 100% non-latex. Made in the USA.

With Blue Tips, 100/Pkg 50-87841

With Clear Tips, 100/Pkg 50-87842

Deluxe Universal Autoclavable OD HVE (Parts)

- Single Piece Contoured Lever/Valve Assembly
- Ultra-smooth Body
- Long Lasting Viton O-rings
- Anodized Aluminum Finish
- Lightweight Swivel QD
- Hangs from Tip in Most Syringe Hangers



With 6-feet of Grey tubing 50-11153

With 6-feet of Sand tubing 50-11154

With 6-feet of Sterling tubing 50-11155

Each 50-11156

90° Autoclavable HVE (Parts)

- Smooth lever action
- Long lasting Viton O-rings
- 90° head for special handling
- Fully autoclavable
- Quick-disconnect



Each 50-11165

Deluxe Autoclavable QD HVE Upper Assembly (Parts)

- Upper assembly only for deluxe QD HVE
- Use while original is being autoclaved



Each 50-11157

Replacement Spool Assembly, deluxe HVE 50-11158

Deluxe Standard Autoclavable OD HVE (Parts)

- Single piece contoured lever/valve assembly
- Hangs from base in most HVE hangers
- Ultra-smooth body
- Long lasting Viton O-Rings
- Anodized aluminum finish
- Lightweight swivel QD



Each 50-11159

With 6-Foot of Grey Tubing 50-11160

With 6-Foot of Sand Tubing 50-11161

With 6-Foot of Sterling Tubing 50-11162

Deluxe Standard Autoclavable QD HVE, Upper Assembly (Parts)

- Upper assembly only for deluxe QD HVE
- Use while original is being autoclaved



Each 50-11163

Replacement spool assembly, deluxe HVE 50-11164

Deluxe Autoclavable QD SE (Parts)

- Single piece contoured lever/valve assembly
- Ultra-smooth body
- Long lasting Viton O-rings
- Anodized aluminum finish
- Lightweight swivel QD
- Fits in most syringe hangers



Each 50-11168

With 6' of tubing 50-11166

Replacement spool 50-11167



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Saliva Ejector System

Contour Vacuum Valve, Slider HVE (Parts)

- Quiet operation
- Smooth external contour without a lever allows one-handed operation
- Quick disconnect and fully autoclavable
- Ready for cleaning and sterilization in seconds



Each 50-11179

Screw-on Tip Saliva Ejector (Parts)

- Secure screw-on style tips
- Choice of Grey or Black tips
- Rubber tip, brass insert
- Autoclavable valve only, remove tip first



SE with Black Tip, each 50-11180

SE with Grey Tip, each 50-11181

SE with Metal QD and Grey Tip, each 50-11182

Standard Rotary HVE (Parts)

- Traditional rotary on/off
- 360° swivel quick-disconnect
- Lightweight comfortable body
- Fits large-bore HVE hangers
- Easy to disassemble for cleaning



Grey HVE, each 50-11183

Black HVE, each 50-11184

Flexjet Multi colored Tips Saliva Ejectors



Clear multi-colored tips, 100/pkg 50-11113

Opaque multi-colored tips, 100/pkg 50-11184

Horizontal Rotary Saliva Ejector (Parts)

- Traditional rotary on/off
- Accepts plastic or metal screen
- Removable push-on tip
- Easy to disassemble for cleaning



Grey SE, each 50-11185

Black SE, each 50-11186

Autoclavable SE No Shut-off (Parts)

- Use for auto-hanger vacuum installations
- Same great features as S-1200, less lever valve



Each 50-11187

SE Adapter for HVE (Parts)

Inserts into HVE to accept SE tips



Each 50-11190

Saliva Ejector Adapter (Parts)

- Allows connection of SE and tubing to an HVE
- Inserts into an HVE tip
- 7/16" OD with a 7/32" barb



Each 50-11191



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

