

Ultrasonic Cleaners



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Ultrasonic Cleaners



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Medflex™ Premium Ultrasonic Cleaners (First Medica)

In order to pursue safe infection control practices and prevent cross contamination the first step is effective cleaning. First Medica's MedFlex™ Premium Ultrasonic Cleaners -provide a crucial step in ridding your instruments of infectious bio-burden. Our high quality stainless steel construction will provide years of reliable service and with the -standard heating element effective cleaning is easy and cost effective. Not only will our Ultrasonic cleaners do an effective job of cleaning, they will look good in your -professional setting!



Stainless steel construction
Acoustic pulse generators produce powerful cavitations for thorough and effective cleaning
Quiet operation
Stainless steel drain for care-free removal of cleaning solution
Scratch and corrosion resistant
Adjustable temperature range for precision cleaning
Adjustable timer
Instrument basket included

3 Liter	55-11122
5 Liter	55-11121
10 Liter	55-11119
Hanging Basket, 3 liter	55-11118
Round Basket	55-11117
Hanging Basket, 5 liter	55-11116
Hanging Basket, 10 liter	55-11113

Ultrasonic Cleaner (ATW)



3L,
with Timer, Basket, Cover w/o Drainage, 120 Watt, Inside Dimension 240×140×90 mm 55-11129

5L,
with Heater, Timer, Basket, Cover & Drainage, 120 Watt, Inside Dimension 240×140×150 mm 55-11130

10L,
with Heater, Timer, Basket, Cover & Drainage, 120 Watt. Inside Dimension 305×245×150 mm 55-11131

Ultrasonic Cleaner (Tuttnauer)

Tuttnauer's Clean & Simple Ultrasonic Cleaning Systems are manufactured with only the highest quality components. These machines provide 60 kHz of quiet cleaning



power which is less aggressive to instruments while still removing all surface contamination. The contemporary design and quiet operation compliment any office environment. Clean & Simple Ultrasonic Systems are available in 2 sizes with a full selection of accessories to meet your specific needs. The 1 gallon system holds a standard selection of instruments or 2 half size cassettes. The 3 gallon holds 3 full size cassettes.

The Tuttnauer CSU1 Table Top Ultrasonic Cleaner

Features:

- High efficiency 60 kHz cavitations for a thorough and quieter cleaning
- Durable vinyl clad steel housing... easy to clean and resists scratching
- Full line of accessories available

•2 1/2, 5 & 10 year Limited Warranty

Unit size: 1 Gallon	55-11132
Unit size: 3 Gallon	55-11134
Unit size: 3 Gallon with Heater	55-11135
Unit size: 1 Gallon with Heater	55-11133
SS Basket for 1 Gallon size	55-11136

Accessory Kit for 1 Gallon size:
includes:positioning cover and 2-600ml beakers with lids and bands 55-11137

SS Basket for 3 Gallons size 55-11138

Accessory Kit for 3 Gallons size:
includes:1/2 size ss basket, positioning cover and 2-600ml beakers with lids a nd bands 55-11139





Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Ultrasonic cleaner Unit (TPC)

- Powerful 40 KHz of cavitations or thorough, quieter cleaning
- Stainless steel cover and tank with welded drain
- Durable vinyl clad steel housing is scratch resistant and easy to clean
- Digital timer for quick and easy process
- Sweep technology ensures exceptional cleaning efficiency



4.0 qt, with basket, Stainless steel lid,
Timer without Drain Kit 55-45681

10.6 qt, with basket, Stainless steel lid,
Timer & Drain Kit 55-45682

Brite Shield (Premier)



Enzymatic action is up to 35% faster than other brands
Chemo-mechanically debrides clogged dental burs, removes alginate and stone from impression trays. Also removes set glass ionomer, carboxylate and zinc phosphate cements. Ideal for removing bio-burden from hand instruments.

A protective shield coats the instrument through a unique passivation process to

restrict rust and corrosion while imparting a lustrous surface to stainless-steel instruments.

800gm/Jar, Container makes up to 34
gallons of finished solution. 55-4921

Efferzyme Ultrasonic Cleaner Tablets (Crosstex)

The first effervescent tablet for all your general purpose instrument and endoscopic cleaning. Fast acting prolease enzymes formulated specifically for cleaning soils from medical instruments.

Just add 2 tablets plus one gallon of water to produce an effective pre-soak or ultrasonic cleaning solution.



52 Tablets/Box 55-24881

Ultrasonic Cleaners

Stain Out (Cory Lab)



Removes tarter and stains from dentures, bridges and all orthodontic appliances.

40 Tabs/ Box 55-105481

Enzymax Powder (Hu-Friedly)

A lighter-weight concentrated powder allows Enzymax to be cost effective while saving valuable storage space. Enzymax powder is specially formulated with a fresh scent.



800gm Makes 55 Gallons 55-24951

EZ-Zyme (Miltex)



- Removes mucous and tissue from hemostats.
- Concentrated, Each 3/4 ounce packet makes 1 gallon of solution.
- Highly effective instrument presoak, ultrasonic and evacuation cleaning solution.

Box of 32 Packets 55-24871

Ultrasonic Cleaner (Yates & Bird)

Cleans any object that can be immersed in liquid. This versatile unit is ideal for both dental office and laboratory use. Compact enough (5" cube overall size) for convenient chair side use, but large enough with its 3 1/2" diameter x 2" deep tank to accommodate full dentures and many other laboratory items. Model 42 generates 42,000 high-intensity cleaning pulses every second that instantly clean any dental diamond, bur, instrument or dental appliance.



Unit Model 42 41-11648

Ultrasonic Cleaners



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Ultrasonic Cleaning Tablets (Defend)



Ultrasonic Enzymatic Tablets

Easy-to-use tablets dissolve in a gallon of water in under a minute to create a powerful, effective cleaning solution. Greatly assists in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents. Contains fast-acting protease enzymes and has a pleasant scent. Also use as an instrument presoak or evacuation system cleaner.

- Safe for amalgam separators
- Biodegradable, neutral pH
- Non-chlorine • Non-corrosive

64 Tablets/Box

55-11142

Citrizyme (Pascal)

Ultra Concentrated Enzyme Cleaner and Instrument Presoak

- Fast-acting dual enzymes remove protein and biofilms quickly
- Leaves no residue
- Premium non-foaming formula
- Powdered concentrate dissolves instantly in warm tap water-no mixing required



900gm, Canister

55-24862

Prozyme Ultrasonic Cleaning Tablets (Cory Lab)



64 Tablets/ Box

55-105641

Clean and Simple Ultra Tablet (Tuttnauer)

Clean & Simple is an ultrasonic / enzymatic solution dispensed in a unique tablet form. Drop 2 tablets into your ultrasonic tank already filled with tap water and the tablets will dissolve in about a minute. This is a non-ionic solution which is non-corrosive, non-toxic, biodegradable and non-caustic.

FEATURES:

- Easy to use tablet form.
- Eliminates handling of heavy 1 gallon containers.
- No bulky storage.
- No pre-mixing required.
- Enzymatic solution effectively breaks down protein.



64 tablets/Box

55-11140

Cetyl-Zyme Pro-Am (Cetylite)



Cetyl-Zyme Pro-Am is a low foaming, bacteriostatic, protease amylase dual enzymatic detergent for cleaning delicate medical and dental instruments/equipment.

It dissolves both protein and starch-containing soils and is useful for cleaning instruments with hard to reach areas such as flexible endoscope channels, box locks

and serrated edges. It may be used as a presoak or in manual, ultrasonic and automated cleaning systems or as an

- One gallon concentrated formula makes 128 gallons of use solution
- Neutral pH Formula
- Non-Corrosive
- Includes no mess easy-to-use 1oz metered pump dispenser
- Pleasant scent
- Use as a presoak or in manual, ultrasonic and automated cleaning systems

**1Gallon Concentrate, 1 oz makes
1gallon of ready to use solution**

55-4931

Enzymax Liquid (Hu-Friedy)

Enzymax is the only specially formulated enzymatic line on the market with dual enzymes and steel guard agents for added corrosion protection & sparkling instruments restoration.



**Box of 40 1/3 oz Packets,
1 Packet makes 1 Gal**

55-24941

**1 Quart Pump Bottle, 1 Packet
makes 1 Gal, 4 Pumps makes 1 Gal**

55-24942

1 Gallon, No Pump

55-24943

Enzymax PAX (Hu-Friedy)

The dissolvable Enzymax PAX are the ultimate in convenience, just drop them in the water. There's no measuring or scooping needed, just use one packet per gallon of warm water.



**Box of 32 powder packets, makes
32 gallons**

55-24961

**Box of 96 powder packets, makes
96 gallons**

55-24962



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

General Purpose Ultrasonic Solutions

Concentrated solution that dilutes 1:10 with water.

Formulated especially to remove debris from a variety of operatory and laboratory items.

1 Gallon	55-89601
4 Gallon	55-89602

Good Vibrations™ (Certol)

Effective Cleaning Action

- Supports ultrasonic cleaning action.
- Removes organic blood debris and inorganic dental materials such as cement.
- Non-Foaming.



Protects Valuable Instruments

- Most effective anti-corrosive agents available.
- Enhances instrument appearance and shine.
- Safe for delicate surgical items and jointed orthodontic pliers.

Economical

- Concentrated Formula.
- 1 ounce makes 1 gallon solution

1 oz. Unit Dose – 24/Box	55-11111
1Gallon with Pump	55-11112

IMS Daily Clean (Hu-Friedy)

Concentrated powder detergent specifically formulated for safe, efficient ultrasonic cleaning and soaking of all dental/surgical instruments and accessories.



1 lb. (Makes 35 gal.)	55-89161
5 lbs. (Makes 160 gal.)	55-89162

Biosonic Enzymatic (Coltene Whaledent)

Enzymatic Ultrasonic Cleaning Concentrate.



8oz. Concentrate	55-24911
-------------------------	----------

Ultrasonic Cleaners

Scooba Ultrasonic Cleaner (Bonart)



By sending ultrasonic pulses through the water for superior cleaning power, the Scooba is ideal for cleaning and removing debris from dental instruments prior to sterilizing.

Include: Main Unit, Cleaning Tray, vented cover, water tube, tube connector

55-24911

ProEZ foam™ (Certol)

Ready-to-Use Foaming Enzymatic Instrument Precleaning Spray

ProEZ foam is an easy to use instrument pre-cleaning spray



Powerful quad enzyme formula quickly dissolves sticky blood soils while unique anti-corrosive system protects and brightens delicate metal surfaces.

ProEZ foam is especially valuable in the dental setting when soiled instruments must be held for extended periods of time. Simply spray cleaning foam on tray of instrument.

24 oz. Spray Bottle	55-11127
1Gallon Refill	55-11128

Biosonic UC-125B Unit (Coltene Whaledent)

This Ultrasonic Cleaner is equipped with programmable LCD display that supports six different languages and features a degas mode, a digital solution timer, easier/faster draining, and customizable operation.



Each	55-105621
-------------	-----------

Ultrasonic Cleaners



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Sonic Cleaner Ultrasonic Cleaner

- Detergent Concentrate.
- Safe, efficient, non-corrosive.
- 5 lbs. makes 160 gallons.



5 Lb. Can

55-89651

Tartar and Stain Remover (DA)



Efficiently removes calculus, tartar, tobacco, and food stains from removable appliances. Use at full strength.

1 Gallon

55-89621

4Gallon

55-89622

Temporary Cement Remover (DA)

Ideal for effectively removing temporary cement from a variety of surfaces. Use at full strength.

1 Gallon

55-89631

4Gallon

55-89632

Biosonic UC-300 Unit (Coltene Whaledent)

The Biosonic UC-300 Ultrasonic Cleaner provides gentle but effective cleaning. The wear-resistant plastic housing which will not stain or dent. The large basket allows you to clean 3 to 4 times as many instruments while the basket dividers facilitate sorting and minimize handling of instruments.



Each

55-109391

Ultrasonic Cleaner - General Purpose (EPR)



Super-concentrate general purpose cleaner designed for most cleaning applications in the dental laboratory.

1 Gallon

55-89641

4Gallon

55-89642

Plaster and Stone Remover (DA)

Removes plaster and stone residue from metal trays and articulators.

1 Gallon

55-89611

4Gallon

55-89612

Dentsonic Ultrasonic Cleaner

Powerful 35KHz of cavitation for thorough, quieter cleaning.
-Stainless steel tank
-Durable plastic housing is scratch resistant and easy to clean



- Digital timer for quick and easy process
- Sweep technology ensures exceptional cleaning efficiency
- Includes basket and lid

0.65 Gallons, heater, digital timer, basket capacitive glass touch control panel.... 55-106571

0.8 Gallons, heater, digital timer, basket durable compact design.... 55-106572

1.6 Gallon, heater, digital timer, basket advanced sweep technology.... 55-106573

2 Gallons, Heater, Digital Timer, Basket capacitive glass touch control panel 55-106574



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Ultrasonic Cleaners

Just to mention a few of our products. To see more visit us online at:
www.prestigedental.com

Best Dental
Supplies!!

Tartar & Stain
Remover

DA



Enzymax
Liquid

Hu-Friedy



Coltene Whaledent

Biosonic Enzymatic



Good Vibrations

Certol



*We Work Hard
So You Don't Have To.*



Call Us at 877-772-3888



www.prestigedental.com

Ultrasonic Cleaner Tuttinauer

Features:

- High efficiency 60 kHz cavitation for a thorough and quieter cleaning
- Durable vinyl clad steel housing... easy to clean and resists scratching
- Full line of accessories available
- 2 1/2, 5 & 10 year Limited Warranty



Follow Us On Social Media!!!

facebook

twitter

LinkedIn

Call us for Quarterly Promotion

Prestigedentalproducts.com



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

