

# Ultrasonic Scalers



**Prestige Dental Products, Inc**

**[Prestigedentalproducts.com](http://Prestigedentalproducts.com)**

**Tel: 877-772-3888 Fax: 714-632-8821**

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Advance Piezo Ultrasonic Scaler Unit 750N (TPC)



### Standard Scaler Features:

- Compact unit housing for easy placement anywhere in the operatory
- Extremely gentle on patient, causes little or no sensitivity
- No heat build up in handpiece tip with extended use
- Autoclavable handpiece
- Extensive power range to allow widest choice of clinical applications
- Fine water spray eliminates excessive overspray onto patient
- Thin sub gingival tips available
- 5 tips included with unit
- 1 year warranty

### Optional Endodontics Unit Features:

- Endodontic file holder available
- Easy and rapid operating for root canal treatment
- Achieve enlargement of canal quickly specifications
- Input rating: 115/230V, 50/60 Hz
- Dimensions: 160mm(W) x 155mm(D) x 70mm(H)
- Weight: 2.1kg
- Frequency: 32 KHz
- 1 year warranty (tips excluded)

<b>Piezo Ultrasonic Scaler &amp; Endodontic Unit</b>	56-17031
<b>Endo files and holder (21mm #25)</b>	56-17032
<b>Piezo scaler tip #1</b>	56-17033
<b>Piezo scaler tip #2</b>	56-17034
<b>Piezo scaler tip #3</b>	56-17035
<b>Piezo scaler tip #4</b>	56-17036
<b>Piezo scaler tip #5</b>	56-17037
<b>Piezo scaler handpiece</b>	56-17038
<b>Tip wrench</b>	56-17039
<b>Piezo scaler O-rings (5 pack)</b>	56-17040
<b>Replacement water line</b>	56-17041
<b>Replacement scaler foot control</b>	56-17042

## ART-MB3 Ultrasonic Scaler Built in type (Bonart)

Power Package includes: PC Board, LCD digital display screen, VR Switch Set, Terminal Block, Magnet Hand Piece w/ Cable Assembly, 2 Inserts. Optional parts includes: Solenoid Valve, Transformer, and Foot Control at additional costs. Available in either 25Khz or 30Khz system.

**Available now in all regions.**

### Warranty Terms:

- **5 year limited** warranty on the internal PC board and electrical components
- **1 year limited warranty** on the hand piece assembly, foot control, power cord, LCD display screen etc.
- 120 days on the expendable part like inserts



<b>Unit, 25 KHZ with 2 Insert</b>	56-11119
<b>Unit, 30 KHZ with 2 Insert</b>	56-11120



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### ART-M1 Ultrasonic Scaler (Bonart)

Powerful easy to use scaling performance that outlasts its competitors. Excellent quality at an affordable price.

Specifications: Operation Frequency: AVAILABLE ONLY IN 25KHZ. Output Power: From 10 to 30 Watts of output power during normal scaling speed. Turbo Function: Increase power during scaling with one touch of a button. C25W to 35W output power. Water Supply: 25-60 PSIG (172-414 Kpa) Optional accessories: 25KHZ inserts ISO Certified; CE Marked; 510K registered. Warranty Terms: 18 month limited warranty on the internal electronic module; 1 year warranty for hand piece, foot control; 120 days warranty for all insert tips.



**Complete Unit with 2 Insert 25KHZ** 56-11112

### Piezo Combo Package (Parkell)

**Combining the power of Piezo technology with an antimicrobial irrigating system**

- TurboPiezo™ Ultrasonic Scaler
- Periosonic® Multi-Fluid irrigator (Model 'W') Three Piezo Tips



**Includes TurboPIEZO scaler unit 110volts w/ 3 tips plus Periosonic Multi-Fluid Irrigator - Model W** 56-110271

### Scaler Sleeves



1 1/2 "Wx7" L stick hold, 500/box 39-105611  
1 5/8 "Wx8" L stick hold, 500/box 39-105612

## Ultrasonic Scalers

### ART-M3 II Ultrasonic Scaler (Bonart)

Has our unique "AGC" (Auto Gain Control) LCD digital display screen that shows scaling mode and power setting. System is available in either 25KHz or 30KHz stand-alone system. Up to a maximum power setting of 50 with auto tune system that locks in the frequency and stabilize power setting. The well-equipped feature called "AGC" (Auto Gain Control) stabilizes the power setting and prevents stalling during scaling operation. System comes with a built in

LCD digital screen that displays the mode and power setting, and allows the user to choose from Scaling, Perio and Endo mode with one turn of the

knob, and an easy to access water control knob at the front control panel. System includes 1 plastic grip internal flow insert and 1 metal grip external flow insert. ISO certified; CE marked; 510K registered. Warranty



**Complete Unit with 2 Insert 30KHZ**  
Inside the box & Another Universal Insert Free 56-11115

**Complete Unit with 2 Insert 25KHZ**  
& Another Universal Insert Free 56-11116

### Biosonic Ultrasonic Scaler System (Coltene Whaledent)

**New Features:**

- Water intensity and power control knobs now rotate in the same direction
- On demand turbo button
- Significantly less expensive than competitors units
- Operates with either 25KHz or 30KHz inserts
- Turbo button for an operator-controlled power boost
- Accepts most competitive magnetostriuctive inserts
- Units small footprint saves space in the operatory
- No bulky external power supply
- 115V Model



**Complete Unit** 56-51247

**One 30k, slim super soft, insert** 56-51248

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## ART-P3 II Piezo Electric Scaler (Bonart)

Up to a maximum power setting of 50 with auto tune system that locks the frequency. The well-equipped feature called "AGC" (Auto Gain Control) stabilizes the power setting and prevents stalling during scaling operation. System comes with a built in LCD digital screen that displays the mode and power setting, allows the user to choose from Scaling, Perio and Endo mode with one turn of the knob, and has an easy to access water control knob at the front control panel. System includes 4 various tips (2 pointed tips & 2 flat tips), 4 individual plastic autoclavable wrenches, 1 aluminum sterilization box. ISO certified; CE marked; 510k registered. Warranty Terms: 5 year limited warranty on the internal PC board and electrical components 2 year limited warranty on detachable piezo hand piece 1 year limited warranty on the hand piece cable, foot control, power cord etc. 120 days on the expendable part like piezo tips.



**Complete Unit with 4 Piezo Tips** 56-11121

## ART-P6 Piezo Electric Scaler (Bonart)

ART-P6 Piezo Electric Compact Sized Scaler

Unit is designed and manufactured to meet the most demanding scaling environment. System is portable and light-weight weighing only 3.52 lbs; it comes with a built in power setting knob (1-20) that allow users to set power output at their desired level of operation. Technical specifications are as follows:



- Power Supply: Available in either 115V or 230V
  - Working Freq: 28KHz to 32KHz
  - Water Supply: 25-60Psi (172-414kPA)
  - Dimension: 7.68"(L)x3.38"(W)x3.27"(H)
  - Weight: 3.52 lbs
  - Handpiece Cable Length: 8.52 feet
  - Footswitch Cable Length: 8.52 feet
  - Operation Environment Temperature: 0-140 degrees F
  - Relative Humidity: 10-90% (non-condensing) system
- includes 4 tips and 1 plastic wrench** in US; 3 tips and 1 plastic wrench in other regions. **5 year power unit warranty**

**Complete Unit with 4 Piezo Tips** 56-11125

## ART-PB3 Piezo Electric Built in Type Scaler (Bonart)



Power Package includes: PC Board, LCD display screen, VR Switch Set, Piezo Titanium Hand Piece w/ Cable Assembly, 4 Tips, 4 wrenches, 1 sterilization box. Optional parts includes: Solenoid Valve, Transformer, Foot Control at additional costs.

### **Warranty Terms:**

**5 year limited warranty** on the internal PC board and electrical components

**2 year limited warranty** on detachable Piezo hand piece


**1 year limited warranty** on the hand piece cable, foot control, power cord etc.

120 days on the expendable part like Piezo tips

**Complete Unit with 4 Piezo Tips** 56-11126

## Ultrasonic NiTi Files

Color Coded used by Ultrasonic scaler units for root canal cleaning



Tip with excellent tenacity is used for root canal cleaning  
File size: 32.5mm±0.5mm  
Handpiece size: Φ0.68mm~0.71mm  
K-file taper: 2%

<b>#10, 33mm, 6/Pkg</b>	29-12247
<b>#15, 33mm, 6/Pkg</b>	29-12246
<b>#20, 33mm, 6/Pkg</b>	29-12245
<b>#25, 33mm, 6/Pkg</b>	29-12244
<b>#30, 33mm, 6/Pkg</b>	29-12243
<b>#35, 33mm, 6/Pkg</b>	29-12242
<b>#40, 33, 6/Pkg</b>	29-12241





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Ultrasonic Scalers

### Manual/Auto - Tune Ultrasonic Scaler (Parkell)



**Includes:** handpiece and autoclavable sheaths

56-110381

### TurboVue Illuminated Magnetostrictive Ultrasonic Scaler (Parkell)



TurboVue Illuminated Magnetostrictive Ultrasonic Scaler is a uniquely engineered device that provides excellent visibility when scaling all areas of the oral cavity. The TurboVue features a light source built into the handpiece allowing a significant amount of light to emit through the 30K, light-transmitting ultrasonic inserts.

#### Features and Benefits:

- Lighted handpiece for better viewing
- External water filter prevents clogs and dripping
- Auto-tuning Technology

This uniquely engineered device provides excellent visibility when scaling all areas of the oral cavity.

TurboVue Scaler: 30K Universal Tip V Insert - 110V scaler w/ attached handpiece, foot controller and water line, In-line water filter w/ spare disk and water connection w/ quick disconnect. Transformer w/ cable	56-110351
30KHz Burnett Power-Tip for TurboVue lighted hand piece	56-110352
30 KHz universal insert for Turbovue lighted hand piece	56-110353
30 KHz straight Perio insert for Turbovue lighted hand piece	56-110354

### Quick Connector for Cavitron Water Line



**Female Male**

**Kit**

Scaler Unit "Male Type"	56-109611
Chair Unit "Female Type"	56-109612
<b>Kit:</b> 1 male & 1 female	56-109613

### Jetpolisher 2000- Air Polisher (Deldent)



Advanced compact, non-clogging airpolishing system for about half the price of standard countertop polisher

**Recommended to clean tooth surfaces prior to bonding procedures, placing sealants, and bleaching.**

- Non-clogging nozzle due to exclusive H.S.T. technology.
- Harmless to titanium implants.
- No installation required, completely portable, and can be moved easily from room to room
- Measure: 3.5" x 3.5" x 7.5"
- Connects directly to existing hand piece tubing (30-35 psi) via 2-3 hole or 4 hole connectors (*kavo, W&H, Sirona connectors*)
- Using existing dental unit's foot pedal
- Lightweight and easy on the wrist
- Fully detachable autoclavable hand piece and water control knob to prevent cross-infection
- Large, easy-view powder reservoir.
  - **FDA cleared for use around implant.**
  - **2 years warranty.**

<b>Midwest Type</b> with 4-Hole connector	56-89691
<b>Kavo Type connector</b>	56-89692
<b>Sirona/Siemens Type connector</b>	56-89693
<b>W&amp;H Type connector</b>	56-89694
<b>JetStream Standard handpiece</b>	56-89695
<b>Booster Handpiece</b> (for heavy stain)	56-89696
<b>Spray Head Tubes</b> (50/Pkg)	56-89697

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## DelSonic 2000 Piezo Scaler (Deldent)



- Sleek, compact, quiet, powerful, and portable.
- Linear tip vibration, for more efficient and more comfortable scaling.
- High frequency 30 kHz design provides greater patient comfort than the competition (less vibration against teeth).
- Cool running, high torque hand piece (minimal tip warmth).
- Fully autoclavable hand piece sleeves, so purchasing additional hand pieces is not necessary.
- Range of tips available for all clinical needs.
- Tips are designed to introduce spray even into deep pockets.
- May be used with any independent fluid supply system dispensing purified water or medicaments.
- 2 year warranty (compare with competitive brands' one year warranty on handpiece).

<b>Full Unit</b>	56-25841
<b>#1 Universal Piezo Tip</b>	56-25842
<b>#2 Round Straight Piezo Tip</b>	56-25843
<b>#2 Round Right Piezo Tip</b>	56-25844
<b>#2 Round Left Piezo Tip</b>	56-25845
<b>#3 Flat Horizontal Piezo Tip</b>	56-25846
<b>#4 Perio Straight Piezo Tip</b>	56-25847
<b>Autoclavable Handpiece Sleeves (set of 3)</b>	56-25848
<b>Tip Changing Tool</b>	56-25849

## JetSonic 2000M Combo - With Cavitron - Type Scaler (Deldent)

Advanced combination scaling and polishing unit with patent H.S.T. air polishing system and advanced scaling design features.

The air polishing portion of the JetSonic 200M has all the features of the JetPolisher 2000 **with the addition of warm water flow.**

The scaler side of this unique unit also offers advanced magneto restrictive technology, including

- Unique sensor technology recognizes and activates the selected hand piece (scaling or polishing)
- Advanced patented digital circuitry **automatically tunes any "stack-type" insert, whether 25K, or 30K.**
- Advanced microchip feedback system maintains optimum scaling power for all settings.
- Power-Surge Mode technology increases power by 15-20% with a "double-click" on the foot switch for removal of tenacious deposits. The pre-selected power setting returns automatically when the foot control is released
- Large, easy-view powder reservoir (no need to stop to add powder).
- Convenient rinse mode available
- Water heating feature
- **2 Years warranty**

**Combo Unit**

56-88681

**Stream Selector 2000 (Fluid supply system):** Unit, 2 Bottles and 2 Quick Disconnect

56-88682

**Bottle for Stream Selector 2000, each**

56-88683





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Piezo Co-Pilot Scaler (Vista)



The piezo **Co-Pilot™** is a beautifully designed and powerful self-contained dual reservoir scaling system. The Co-Pilot™

accepts a wide variety of piezo tips from major manufacturers and can be ordered with either a standard (Acteon\*) or metric (EMS\*) style hand piece.

- Two large, easy-to-fill 350ml reservoir tanks with easy-to-see liquid levels
- Accepts a wide variety of ultrasonic scaling tips for General Scaling, Periodontal, Endodontic and Apical Surgery applications
- Advanced circuitry, touchpad controls, and push-button controlled variable irrigation flow settings
- Compact design for improved productivity and portability (needs NO operatory air/water connections)
- Soft-lit dome and reservoir. Autoclavable handpiece sleeves available in six contemporary colors

#### **System Includes:**

- 1 - Piezo Pilot™ ultrasonic scaler unit,
- 4 - Ultrasonic tips and 4 torque wrenches (Either Hygiene or Endo tip packages)
- 1 - Detachable handpiece (either Universal Acteon\* or Metric EMS\* thread)
- 3 - Color-matched autoclavable handpiece sleeves
- 1 - Quick-disconnect water line and replaceable in-line water filter
- 1 - Omni-directional foot pedal

**Acteon Style Unit, Grey Color** 56-11162

**Acteon Style Unit, Purple Color** 56-11163

**Acteon Style Unit, Blue Color** 56-11164

**Acteon Style Unit, Yellow Color** 56-11165

**Acteon Style Unit, Blue/Green Color** 56-11166

**Acteon Style Unit, Pink Color** 56-11167

**EMS Style Unit, Grey Color** 56-11168

**EMS Style Unit, Purple Color** 56-11169

**EMS Style Unit, Blue Color** 56-11170

**EMS Style Unit, Yellow Color** 56-11171

**EMS Style Unit, Blue/Green Color** 56-11172

**EMS Style Unit, Pink Color** 56-11173

## Ultrasonic Scalers

### Piezo Handpiece (Bonart)



Piezo handpiece truly defines the innovative modern piezo technology developed by Bonart. The outer casing of the hand piece is made of a specialized material that can withstand steam sterilization of 135 degrees Celcius. Its newly improved transducer design offers a more powerful scaling, endo and perio vibrating frequency with less pain and great ease. It comes with a protective cap to be installed at the bottom portion of the hand piece during autoclave session to prolong the life of the entire hand piece. ISO certified; CE Pending. 510k registered!. 2 year limited warranty.

**Each**

**56-11184**

### ProSonic Air Scaler (Pac-Dent)



The ProSonic™ air scaler is a pneumatically driven device that generates an elliptical oscillation of the working tips in the sound range (< 6.5 kHz) with longitudinal and transverse movements. It is as powerful as it gets without being an Ultra Sonic, and has been clinically proven to be efficient and gentle with a wide range of uses, including: Hygiene, Periodontics, Minimally invasive caries therapy, Cavity shaping and finishing, Endodontics and Cementation. It features a controlled water spray that prevents the "drowning effect" often experienced by patients during deep scaling. The scaler is entirely autoclavable / chemically, and the tips can also be dry heat sterilized. To further safeguard the operator and patient from cross contamination, it also comes with replaceable tactile silicone grips which also relieves finger fatigue. The ProSonic™ has a slim and lightweight design that permits complete scaling even in tight interproximal areas. Packaged with 5 scaling tips - Universal, Sickle, Perio, Perio Left and Perio Right

- Sonic technology effectively removes calculus deposits and dental stains.
- Connects to existing air driven tubings and foot control
- Lightweight - helps to relieve wrist and finger fatigue
- Can be safely used on patients with pacemakers
- Fully autoclavable
- Ergonomic & color coded silicone grips
- Selection of tips for different scaling needs

#### **ProSonic™ Air Scaler Unit: 1x**

ProSonic air scaler, 1 x universal scaling tip, 1 x sickle scaling tip, 1 x perio scaling tip

**56-55401**

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Piezo Pilot Scaler Unit (Vista)



Introducing the Piezo Pilot™ a New Generation of ultrasonic scalers. Featuring state-of-the-art circuitry and a compact modern design, the Pilot is a lightweight powerhouse. The soft-lit dome and removable / autoclavable hand piece sleeves are available in six contemporary colors to adapt to any operatory decor. The unit also features a magnetic handpiece holder that can be positioned on either side of the Pilot or onto any magnetic surface.

A complement to any office decor, the Piezo Pilot™ is a powerful, lightweight, and beautifully designed piezo scaler that allows the hygienist to work quickly and effortlessly. The Pilot's state-of-the-art circuitry offers a wide range of amplitude settings for regular hygiene procedures and all ultrasonic endodontic procedures, including certain bone surgery applications. It comes with a choice of either four ultrasonic hygiene or four ultrasonic endodontic tips, your choice of either a Universal (P5) or Metric (EMS) handpiece, and is compatible with tips from a variety of major tip manufacturers, including Vista's full line of piezo tips.

**Acteon Style Unit**, Grey Color with Hygiene 8 Tips Package SU1, SU6, HC, PTU & 8 Torque Wrenches + Interchangeable colored Domes & Autoclavable hand piece sleeves

56-11186

**EMS Style Unit**, Grey Color with Hygiene 8 Tips Package SU1, SU6, HC, PTU & 8 Torque Wrenches + Interchangeable colored Domes & Autoclavable hand piece sleeves

56-11198

## Autoscaler (Southeast Inst)

### TRUSTED QUALITY & PERFORMANCE

Special Periodontal & Routine Scaling

Detachable handpiece

Assembly

Precision water control

**Made in USA**

**Best in Class warranty**

- 6 Years for circuit

boarsa

- 2 Years for all other components

The Autoscaler is the premier in Ultrasonic

Accepts the corresponding frequency: **DENTSPLY® CAVI-**

**TRON®** Brand ( TFI, P, SLI and FSI inserts

The Autoscaler allows for feather touch removal of stain and calculus



Autoscaler 25 KHz Model 56-110551

Autoscaler 30 KHz Model 56-110552

Handpiece Assembly for Autoscaler 25 or 30 K 56-110553

Insert 25K DF10 56-110554

Insert 30K DF10 56-110555

## Integra Self - Contained Ultrasonic Scaler System (Parkell)



The Integra features an expanded low-power range for patient comfort during subgingival scaling.

Designed for today's infection-control concerns, it features an easy-to-clean waterproof membrane control panel, autoclavable handpiece sheaths and a built-in safety tum-off feature.

Complete Unit

56-110581





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Scalex 800 Ultrasonic Scaler (Dentamerica)



The **scalex 800 Ultrasonic Scaler** is **25,000** scaling strokes per second, solid state power unit.

- Streamlined compact size and slim light weight handpiece
- Convenient foot control

- No power loss at maximum pressure
- Unit is available in 100/110V 50-/60 Hz, or 220/240V 50 Hz operation
- Dimensions 22L x 19W x 7.5H cm, weight 3.6kg

**Unit 25K, Tips not Included** 56-11267


**F4 Scaler Universal, External Water Tip 25K each** 56-11268

**F1 Scaler Universal, Internal Water Tip 25K each** 56-11269

### Scalex 880 Piezo Scaler (Dentamerica)

The Scalex 880 Piezo Scaler has State of the art piezo-electric technology is self-tuning and load adjusting for maximum smooth operation

- Streamlined compact size and slim lightweight handpiece.
- Innovative water-cooled circuitry designed for greater stability during operation
- Extremely gentle on the patient
- The Scalex 880 Piezo Scaler has Available a wide variety of high quality tips to correspond with specific clinical applications
- Three universal tips included. **1 year warranty**



**Unit with 3 Universal A Tips** 56-11276

**Piezo Handpiece only** 56-11277

**Piezo Scaler Tip (Universal A) each** 56-11278

**Piezo Scaler Tip (Universal B) each** 56-11279

**Piezo Scaler Tip (Chisel) each** 56-11280

**Piezo Scaler Tip (Perio) each** 56-11281

**Endo Chuck each** 56-11282

## Ultrasonic Scalers

### Scalex 850 Self Contained Ultrasonic Scaler (Dentamerica)

Scalex™ 850 offers self-contained solution tank for total portability, adding therapeutic medicaments



The **Scalex 850 Self-Contained Ultrasonic Scaler** has all the cleaning power

you need with performance guaranteed for one full year. Self-contained solution tank allows for convenient portability and the addition of medicated solutions to scaler therapy. Automatic tuning activates inserts in a magnetostrictive cleaning pattern for feather-touch removal of calculus and stains under all loading conditions—no power loss under pressure. Each system includes base unit with graduated 1000ml self-contained solution tank, handpiece with extra-long 7 ft. cord and one 30kHz insert. 110V. **1 year manufacturer's warranty.**

**Unit, 30K with 1 Tip** 56-11273


**F6 Scaler Universal External Water Tip 30K each** 56-11274

**F3 Scaler Universal Internal Water Tip 30K each** 56-11275

**Water Tank** 56-12275

### Scalex 830 Ultrasonic Scaler (Dentamerica)

- 30,000 scaling strokes per second, solid state power unit.
- Streamlined, compact size and slim lightweight hand piece
- Innovative water-cooled circuitry designed for greater stability during operation
- No power loss at maximum pressure
- Available selection of a wide variety of high quality inserts t correspond to specific clinical applications
- Unit for 25K scaling strokes is available
- Specifications : 110/220V 50/60Hz AC
- Dimensions : L140 x W140 x H50 mm
- Weight : 1 kg
- Accessories : Control Box, Hand piece, Foot Control, Water Connector, Tube, Power Cable



**Unit 30K with 1 Tip** 56-11272

**F6 Scaler Universal External Water Tip 30K each** 56-11271

**F3 Scaler Universal Internal Water Tip 30K each** 56-11270

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Scalex 890 Self-Contained Piezo Scaler (Dentamerica)

Scalex 890 offers self-contained solution tank for total portability, adding therapeutic medicaments



### The **Scalex 890 Self-Contained Piezo Scaler**

has all the cleaning power you need with performance guaranteed for one full year. Self-contained solution tank allows for convenient portability and the addition of medicated solutions to scaler therapy. State of the art Piezo Electric Self-tuning And load adjusting technology for smooth and effortless operation.

Innovative water cooled circuitry designed for greater stability during operation.

<b>Unit with 3 Universal A Tips</b>	56-11289
<b>Piezo Handpiece only</b>	56-11288
<b>Piezo Scaler Tip (Universal A) each</b>	56-11287
<b>Piezo Scaler Tip (Universal B) each</b>	56-11286
<b>Piezo Scaler Tip (Chisel) each</b>	56-11285
<b>Piezo Scaler Tip (Perio) each</b>	56-11284
<b>Endo Chuck each</b>	56-11283

## Sonic Air Scaler (TPC)

### Features:

- Includes 3 tips AS1, AS2, AS3
- Tip wrench included
- Light weight
- Portable
- Saves time
- Autoclaveable
- Available in 2 or 4 hole

### Package includes the following:

- AS5000 Air sonic scaler handpiece
- ASI Universal Tip
- AS2 Sickle Tip
- AS3 Perio Tip
- A776 Tip Wrench
- **6Month warranty**



<b>4 Hole Kit</b>	56-17051
-------------------	----------

<b>2 Hole Kit</b>	56-17052
<b>ASI Universal Tip</b>	56-17053
<b>AS2 Sickle Tip</b>	56-17054
<b>AS3 Perio Tip</b>	56-17055

## Ultrasonic Inserts

### After Five Plus Ultrasonic Inserts (Hu-Friedy)

#### Sub-gingival

- External fixed water delivery
- Thin, probe-sized design
- Provides smooth, deep access into periodontal pockets

- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use

AFTER FIVE PLUS - For Subgingival



<b>Left, 25 kHz, Teal Handles, each</b>	56-90975
<b>Straight, 25 kHz, Blue Handle, each</b>	56-90976
<b>Right, 25 kHz, Red Handle, each</b>	56-90977
<b>Set, 25 kHz, each</b>	56-90978
<b>Left, 30 kHz, Teal Handles, each</b>	56-90979
<b>Straight, 30 kHz, Dark Blue Handle, each</b>	56-90980
<b>Right, 30 kHz, Red Handle, each</b>	56-90981
<b>Set, 30 kHz, each</b>	56-90982

### After Five Ultrasonic Inserts (Hu-Friedy)

#### Sub-gingival

- External water delivery
- Provides smooth, deep access into periodontal pockets
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use
- Thinner Tip, Slender Taper with Metal Handles

HU-FRIEDY

AFTER FIVE - For Subgingival



<b>Straight, 25 kHz, each</b>	56-90983
<b>Set of 3:</b>	
Left, Straight, Right, 25 kHz, each	56-90984
<b>Straight, 30 kHz, each</b>	56-90985
<b>Set of 3:</b>	
Left, Straight, Right, 30 kHz, each	56-90986
<b>Left, 3. kHz, each</b>	56-90987



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Biosonic Ultrasonic Inserts (Whaledent)

BioSonic Inserts - variety of most popular tip designs in both 25 kHz and 30 kHz

Better performance reliability and longevity through patented manufacturing process  
Ergonomically designed grip for operator's comfort and reduced hand stress  
Maximum comfort in popular tip designs



<b>Biosonic Ultrasonic Scaler Tip Slim #10</b> (25kHz)	56-21110
<b>Biosonic Ultrasonic Scaler Tip Slim #10</b> (30kHz)	56-21111
<b>Biosonic Ultrasonic Scaler Tri-Bend #1000</b> (25kHz)	56-21112
<b>Biosonic Ultrasonic Scaler Tri-Bend #1000</b> (30kHz)	56-21113
<b>Biosonic Ultrasonic Scaler Tip Opti-Mist</b> (25kHz)	56-21114
<b>Biosonic Ultrasonic Scaler Tip Opti-Mist</b> (30kHz)	56-21115
<b>Biosonic Ultrasonic Scaler Tip Slim OptiMist #10</b> (25kHz)	56-21116
<b>Biosonic Ultrasonic Scaler Tip Slim OptiMist #10</b> (30kHz)	56-21117
<b>Biosonic Scaler O-Rings</b> Pkg of 12	56-21118

### Satin Swivell Ultrasonic Inserts (Hu-Friedy)

- Exclusive 360° rotation
- Color-coded silicone grips
- Large handle diameter reduces finger pinching
- Internal water delivery
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use



<b>10-25K, Satin Swivel Insert,</b> each, UI25SS10	56-90950
<b>100-25K, Satin Swivel Insert,</b> each, UI25SS100	56-90951
<b>1000-25K, Satin Swivel Insert,</b> each UI25SS1000	56-90952
<b>10-30K, Satin Swivel Insert,</b> each UI30SS10	56-90947
<b>100-30K, Satin Swivel Insert,</b> each UI30SS100	56-90948
<b>1000-30K, Satin Swivel Insert,</b> each, UI30SS1000	56-90949

## Ultrasonic Scalers

### ART-3 Piezo - Implant motor combo (Bonart)

Piezo/Implant motor combo. The Most Trusted, Dependable and Exceptional Technology



PiezoART units have become the most wanted systems in the world.

6 tips included

**Complete Unit:** L.E.D. hand piece 36-110201 & low speed optical hand piece included

### Streamline Ultrasonic Inserts (Hu-Friedy)



- Internal water delivery
- The #3 Beavertail design is excellent for removal of gross calculus and orthodontic cement
- Use #10 Universal for the removal of all types of debris
- Use #100 for removing light deposits
- Use #1000 for removing tenacious stains and calculus
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use.

<b>Lavender, 25 kHz, #10 Universal</b>	56-90953
<b>Black, 25 kHz, #100 Thin</b>	56-90954
<b>Orange, 25 kHz, #1000 Triple Bend</b>	56-90955
<b>Yellow, 25 kHz, #3 Beavertail</b>	56-90956
<b>Lavender, 30 kHz, #10 Universal</b>	56-90957
<b>Black, 30 kHz, #100 Thin</b>	56-90958
<b>Orange, 30 kHz, #1000 Triple Bend</b>	56-90959
<b>Yellow, 30 kHz, #3 Beavertail</b>	56-90960

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Swivel Direct Flow Ultrasonic Inserts (Hu-Friedy)

- Patented 360° rotation for improved access and adaptability
- Exceptional water delivery focuses the water flow directly to the tip for improved visibility and increased patient comfort
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use



Univ #10, Swivel, 25K, each	56-90931
Univ #10, Swivel, 30K, each	56-90932
Thin #100, Swivel, 25K, each	56-90933
Thin #100, Swivel, 30K, each	56-90934
After Five, Left, 25K, each	56-90935
After Five, Straight, 25K, each	56-90936
After Five Set, 25K, 3/Pkg	56-90937
After Five, Left, 30K, each	56-90938
After Five, Right, 30K, each	56-90939
After Five, Straight, 30K, each	56-90940
After Five 30K, 3/Pkg	56-90941
#1000 Triple Bend, 25K, each	56-90942
#1000 Triple Bend, 30K, each	56-90943

## Endo Chuck Peizo Tips (ATW)




E1 120 Degree Angle, each (Endo File Not Included)	29-11294
E2 90 Degree Angle, each (Endo File Not Included)	29-11295

## Ultrasonic Scaler ART Inserts (Bonart)




SPA Series Insert  
Internal water flow



P-10 Series Insert  
Universal external water flow



P-100L Insert  
Left angle slim external water flow



P-100R Insert  
Right angle slim external water flow



P-100 Series Insert  
Straight angle slim external water flow



IF-50 Streamline Insert  
Universal internal water flow



IF-5 Beavertail Insert  
Beavertail internal water flow



IF-100 Series Insert  
Slim internal water flow

### SPA Series– Internal water flow insert:

Bonart Medical Technology, Inc. holds the patent to the SPA series of internal water flow inserts. The SPA inserts spray water close to the tip utilizing a silicone sleeve that pushes the water from the back of the tip to the front. This feature allows for a more focused water spray. **Autoclavable, Made in USA**

### P-10 Series– Universal external water flow insert

The Bonart universal external water flow insert series (P-10) is a straight angled universal insert with an external water tubing that sits above the tip. The P-10 series is designed to be compatible with most major brand magnet scalers. **Autoclavable, Made in USA**

### P-100L - Left angle slim external water flow insert

The Bonart P-100L insert is a left-angled insert with a slim tip. Water flows from an external water tube that sits above the tip. The P-100L insert is designed to be compatible with most major brand magnet scalers. **Autoclavable, Made in USA**





Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### P-100R– Right angle slim external water flow inserts

The Bonart P-100R insert is a right-angled insert with a slim tip. Water flows from an external water tube that sits above the tip. The P-100L insert is designed to be compatible with most major brand magnet scalers. **Autoclavable, Made in USA**

**P-100 Series– Straight angle slim external water flow insert**  
The P-100 series is a straight-angled insert with a slim tip. The slim tip is best used for perio scaling of calculus on sensitive areas of the mouth. Water flows from an external water tube that sits above the tip. The P-100 insert is designed to be compatible with most major brand magnet scalers. **Autoclavable, Made in USA**

### IF-50 Streamline Insert

The IF-50 Streamline #10 ultrasonic insert has a universal sized tip. Water flow comes from the back of the tip. This tip is designed to be compatible with most major brand magnet scalers. **Autoclavable, Made in USA**

### IF-5 Beavertail Insert

The IF-5 Beavertail insert has a spatula-style tip for removing hidden calculus more easily. Water flows from the back of the tip. This tip is designed to be compatible with most major brand magnet scalers. **Autoclavable, Made in USA**

### IF-100 Slim Series Insert

The FS slim series insert has a slim tip for better perio scaling of calculus on sensitive areas of the mouth. Water flow comes from the back of the tip. This tip is designed to be compatible with most major brand magnet scalers. **Autoclavable, Made in USA**

Internal water flow insert 30KH, each	56-11306
Internal water flow insert 25KH, each	56-11305
Universal External water flow insert, 30KH each	56-11304
Universal External water flow insert, 25KH each	56-11303
Left angle slim external water flow insert 30KH, each	56-11302
Left angle slim external water flow insert 25KH, each	56-11301
Right angle slim external water flow insert 30KH, each	56-11300
Right angle slim external water flow insert 25KH, each	56-11299
Straight angle slim external water flow insert 30KH each	56-11298
Straight angle slim external water flow insert 25KH each	56-11295
Universal internal water flow insert 30KH, each	56-11297
Universal internal water flow insert 25KH, each	56-11294

## Ultrasonic Scalers

<b>Slim internal water flow insert</b> 30KH, each	56-11293
<b>Slim internal water flow insert</b> 25KH, each	56-11292
<b>Beavertail internal water flow insert</b> 30KH, each	56-11296
<b>Beavertail internal water flow insert</b> 25KH, each	56-11307

### Ultrasonic Scaler Inserts (Dentamerica)

#### ※ULTRASONIC INSERTS

Reorder No.80005  
Universal Insert 25K  
(External Water)



Reorder No.83005  
Universal Insert 30K  
(External Water)



Reorder No.80006  
Universal Insert 25K  
(Internal Water)



Reorder No.83006  
Universal Insert 30K  
(Internal Water)



Reorder No.80007  
Chisel Type Insert 25K  
(External Water)



Reorder No.83007  
Chisel Type Insert 30K  
(External Water)



Reorder No.80008  
Chisel Type Insert 25K  
(Internal Water)



Reorder No.83008  
Chisel Type Insert 30K  
(Internal Water)



Reorder No.80011  
Endodontic Insert 25K  
(Canalsonic Tip)



#### SPECIFICATIONS :

- \* INPUT : 115V 60Hz or 230V 50Hz 3A
- \* OUTPUT : AC 24V 2.4A
- \* FREQUENCY : 25,000 or 30,000 Hz/sec.
- \* DIMENSION : W140xD140xH45mm
- \* WEIGHT : 1000g

F4 External Water, 25K Universal Insert	56-11313
F1 Internal Water, 25K Universal Insert	56-11312
External Water, 25K Chisel Type Insert	56-11311
Internal Water, 25K Chisel Type Insert	56-11310
F6 External Water, 30K Universal Insert	56-11309
F3 Internal Water, 30K Universal Insert	56-11308

### Ultrasonic Scaler - Turbo Sensor (Parkell)



Complete Unit	56-108261
---------------	-----------



Find us on:  
**facebook®**

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Piezo Periodontic Tips



P1

Each

56-11133

## O-Rings For ART Ultrasonic Insert




5/Pkg

56-106431

## Piezo Bone Surgery Tips (Vista)

- Acteon\* Style
- Mectron\* Style



BT1 (Acteon Style)	56-11160
BT3 (Acteon Style)	56-11161
ST1 (Acteon Style)	56-11159
ST2 (Acteon Style)	56-11158
ST3 (Acteon Style)	56-11157
ST4 (Acteon Style)	56-11156
LT1 (Acteon Style)	56-11155
BT1 (Mectron Style)	56-11154
BT3 (Mectron Style)	56-11153
ST1 (Mectron Style)	56-11152
ST2 (Mectron Style)	56-11151
ST3 (Mectron Style)	56-11150
ST4 (Mectron Style)	56-11149
LT1 (Mectron Style)	56-11148

## Piezo Endodontic Tips (Vista)

**GE20 20mm Tip from Bend Angle** Same geometry as the GE20D without the diamond coating. This tip reacts less aggressively and is not recommended for cements.

**GE20D Diamond Coated 20mm Tip from Bend Angle** Used in a variety of endodontic procedures, including the

removal of dental pulp, re-treatment, retrieval of separated instruments, identifying canal openings and the removal of cements. Vista Dental recommends usage of irrigation as a cooling agent when removing crowns and posts to avoid injury.



### **4 Acteon Style Tips**

(GE20, GE20D, GE2D-2 & GE2D-2) 56-11174

**Piezo (Acteon, Vista P5, Tulsa, Spartan, NSK, Hu - Friedy) Style Tips (GE20) each**

56-11175

**Piezo (Acteon, Vista P5, Tulsa, Spartan, NSK, Hu - Friedy) Style Tips (GE20D) each**

56-11176

**Piezo Acteon Style Tips (GE2D-2) each**

56-11177

**Piezo Acteon Style Tips (GE2D-2) each**

56-11178

### **4 EMS Style Tips**

(GE20, GE20D, GE2D-2 & GE2D-2) 56-11179

**Piezo EMS Style Tips (GE20) each**

56-11180

**Piezo EMS Style Tips (GE20D) each**

56-11181

**Piezo EMS Style Tips (GE2D-2) each**

56-11182

**Piezo EMS Style Tips (GE2D-2) each**

56-11183

**Piezo (Amdent\* (US30), Piezo Advantage) Style Tips (GE20) each**

56-12183

**Piezo (Amdent (US30), Piezo Advantage) Style Tips (GE20D) each**

56-12184

## Piezo Hygiene Tips (Vista)

**SU1 Light to Medium - Heavy Cases of Supra-Gingival Calculus,**

Utilizing the side radius of the initial 2-3mm of the top of the tip in a sweeping motion on the tooth

enamel provides a smooth scaled surface with minimal need for secondary hand instrumentation.

**SU6 Designed for Subgingival Scaling** Long and thin design allows for a less evasive entry into subgingival pockets resulting in minimal patient sensitivity and reduced bleeding





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

**HC Universal Curette Design for Sub and Supragingival Scaling** Designed for deep pocket access and powerful enough for medium to heavy calculus removal. Use the edge of the tip to guide the motion for smooth end results.

**PTU Subgingival Pocket Scaling** With a design to simulate the natural curvature of a tooth, this universal tip provides for easy access into pockets with minimal patient discomfort. Extremely thin design helps address patient sensitivity and allows tip to create a smooth scaled surface.

**4 Tips (SU1, SU6, HC & PTU) 56-11185**


### Piezo Scaler Tips (Dentamerica)

#### ※PIEZO SCALER TIPS



**Endo Chuck, each 56-11214**  
**Piezo Scaler Tip (Universal A) each 56-11213**  
**Piezo Scaler Tip (Universal B) each 56-11212**  
**Piezo Scaler Tip (Chisel) each 56-11211**  
**Piezo Scaler Tip (Perio) each 56-11210**

### Piezo Ultrasonic Tip Kits (Vista)



Tip kits offer great value and provide the clinician with the tips most likely used for a specific procedure in one complete package in both Acteon\* Style and EMS\* Style. Vista Dental's piezo tips kits include:

## Ultrasonic Scalers

- General Scaling (SU1, SU2, SU3, SU4, SU5, SU6)
- Apical Surgery (S70D, S90D, SPL, SPR, SA, SAD)
- Gracey Mini Curette (HC-1, HC-2, HC-3, HC-4, HC-7, HC-8)
- Periodontal (DST, HC, PTU, PTL, PTR, SU6)
- Retreatment (GE20, GE20D, GE40D, PRT, CT4, VPro™)
- Restorative (GEBD-1, GEBD-2, GEFD-1, GEFD-2)
- Bone Surgery (BT1, BT3, ST1, ST2, ST3, ST4, LT1)
- Hygienist Kit: (DST, HC, PTU, PTL, PTR, SU6)

**EMS\* Style:** Mini Piezon\*, Scalex 880\*, Bonart\*, Super-Son EIE\*, Master 400\*, Parkell (Piezo)\*, Dent America\*

**Amdent\* Style:** Piezo +

**Satelec\* Style:** Vista P5\*, NSK\*, ENAC\*, Sybron Endo\*, Obtura Spartan\*, ADEC\*, J.Morita\*, Dentsply/Tulsa\*, Amdent\*

**Cavity Prep Kit (Acteon Style) 56-11264**

**General Scaling Tip Kit:** SU1, SU2, SU3, SU4, SU5, SU6 one each Plus 6 Torque Wrenches (Acteon Style) 56-11263

**Perio Scaling Kit (Acteon Style) 56-11262**

**Hygienist Tip Kits:** DST, HC, PTU, PTL, PTR, SU6 one each Plus 6 Torque Wrenches (Acteon Style) 56-11261

**Retreatment Tip Kit:** GE20, GE20D, GE40D, PRT, CT4, VPro™ one each Plus 2 Wrenches (Acteon Style) 56-11260

**Restorative Tip Kit:** GEBD-1, GEBD-2, GEFD-1, GEFD-2 one each Plus 2 Wrenches (Acteon Style) 56-11259

**Apical Surgery Tip Kit:** S70D, S90D, SPL, SPR, SA, SAD one each (Acteon Style) 56-11258

**Bone Surgery Tip Kit:** BT1, BT3, ST1, ST2, ST3, ST4, LT1 one each (Acteon Style) 56-11257

**Gracey Mini Curette Tip Kit:** HC-1, HC-2, HC-3, HC-4, HC-7, HC-8 one each (Acteon Style) 56-11256

**General Scaling Tip Kit:** SU1, SU2, SU3, SU4, SU5, SU6 one each (Amdent Style) 56-11255

**Perio Scaling Kit (Amdent Style) 56-11254**

**Cavity Prep Kit (EMS Style) 56-11253**

**Bone Surgery Tip Kit:** BT1, BT3, ST1, ST2, ST3, ST4, LT1 one each (EMS Style) 56-11252

**Gracey Mini Curette Tip Kit:** HC-1, HC-2, HC-3, HC-4, HC-7, HC-8 one each (EMS Style) 56-11251

*[www.prestigedental.com](http://www.prestigedental.com)*



# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Bone Surgery Tip Kit:

BT1, BT3, ST1, ST2, ST3, ST4, LT1  
one each (Mectron Style) 56-11250

**General Scaling Tip Kit:** SU1, SU2,  
SU3, SU4, SU5, SU6 one each Plus 6  
Torque Wrenches (EMS Style) 56-11249

**Perio Scaling Kit (EMS Style)** 56-11248

**Hygienist Tip Kit:** DST, HC, PTU,  
PTL, PTR, SU6 one each Plus 6  
Torque Wrenches (EMS Style) 56-11247

**Retreatment Tip Kit:** GE20, GE20D,  
GE40D, PRT, CT4, VPro™ one each  
Plus 2 Wrenches (EMS Style) 56-11246

**Restorative Tip Kit:** GEBD-1, GEBD-  
2, GEFD-1, GEFD-2 one each Plus 2  
Wrenches (EMS Style) 56-11245


**Apical Surgery Tip Kit:** S70D, S90D,  
SPL, SPR, SA, SAD one each (EMS  
Style) 56-11244

**Bone Surgery Tip Kit:** BT1, BT3,  
ST1, ST2, ST3, ST4, LT1 one each  
(Piezotome Style) 56-11265

## Ultrasonic Piezo Tips (Vista)

Hygiene/Perio  
**SU1 Light to  
Medium -  
Heavy Cases of  
Supra-Gingival  
Calculus**

Utilizing the side  
radius of the ini-  
tial 2-3mm of the  
top of the tip in a  
sweeping motion  
on the tooth  
enamel provides a  
smooth scaled  
surface with mini-  
mal need for sec-  
ondary hand in-  
strumentation.



**SU2 Scaling of High Density Supra-Gingival Calculus Deposits**

Utilizing the flat edge of the top of the tip and the side of the tip perpendicular to the tooth surface, use an up-and-down motion on the tooth surface to shear off calculus.

**SU3 Removal of Stains (Tobacco, Coffee & Other Similar Stains)**

Utilizing both the round edge of the top of the tip and the side of the tip on the tooth surface, use an up and down motion on the tooth surface to shear off stains.

## SU4 General Scaling of Sensitive Areas

Less abrasive than the SU1 - SU3 tips, this tip is ideal for use in the periodontal pockets and with other sensitive locations.

## SU5 Designed for Easy Interproximal Scaling

The thin design of the SU5 tip allows for entry into most interproximal spaces.

## SU6 Designed for Subgingival Scaling

Long and thin design allows for a less evasive entry into subgingival pockets resulting in minimal patient sensitivity and reduced bleeding.

**EMS\* Style:** Mini Piezon\*, Scalex 880\*, Bonart\*, Super-Son EIE\*, Master 400\*, Parkell (Piezo)\*, Dent America\*

**Amdent\* Style:** Piezo +

**Satelec\* Style:** Vista P5\*, NSK\*, ENAC\*, Sybron Endo\*, Obtura Spartan\*, ADEC\*, J.Morita\*, Dentsply/Tulsa\*, Amdent\*

## SU1 - Advantage / Amdent™

Thread Style, each 56-11328

SU1 - P5 Thread Style, each 56-11327

SU1 - EMS™ Thread Style, each 56-11325

SU2 - P5 Thread Style, each 56-11323

## SU2 - Advantage / Amdent™

Thread Style, each 56-11324

SU2 - EMS™ Thread Style, each 56-11322

SU3 - P5 Thread Style, each 56-11321

## SU3 - Advantage / Amdent™

Thread Style, each 56-11320

SU3 - EMS™ Thread Style, each 56-11319

SU4 - P5 Thread Style, each 56-11318

## SU4 - Advantage / Amdent™

Thread Style, each 56-11317

SU4 - EMS™ Thread Style, each 56-11316

SU5 - P5 Thread Style, each 56-11315

## SU5 - Advantage / Amdent™

Thread Style, each 56-11314

SU5 - EMS™ Thread Style, each 56-11329

SU6 - P5 Thread Style, each 56-11330

## SU6 - Advantage / Amdent™

Thread Style, each 56-11331

SU6 - EMS™ Thread Style, each 56-11332






Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Ultrasonic Scalers

### Piezo Scaler Tips (Bonart)

**ART Piezo Tip: BE-1 / BS-1:** is a pointed universal tip best used for general deposit removal in all quadrants.

**ART Piezo Tip: BE-2 / BS-2:** is a flat crescent tip best used for removing resin-free ortho cements. This tip may also be used for supragingival calculus removal.



ART Piezo Tip: BE-1 / BS-1  
Universal pointed tip



ART Piezo Tip: BE-2 / BS-2  
Flat chisel tip



ART Piezo Tip: BE-3 / BS-3  
Rounded flat tip

**ART Piezo Tip: BE-3 / BS-3:** is a rounded flat tip. The rounded edge is designed for gentle treatment.

**ART Piezo Tip: BS-PFR:** is a right-angled tip designed for perio procedures and reaching difficult-to-reach spots near the posterior teeth.



ART Piezo Tip: BS-PFR  
Right-angled perio tip



ART Piezo Tip: BS-P  
Slim perio tip



ART Piezo Tip: BS-PFL  
Left-angled perio tip

**ART Piezo Tip: BS-P:** is a slim perio tip designed for supra and subgingival scaling.



ART Piezo Tip: BS-PFU  
Straight-angled perio tip



ART Piezo Tip: BS-10P  
45° angled perio tip



ART Piezo Tip: BS-H2L  
Left-angled curved perio tip

**ART Piezo Tip: BS-PFL:** is a left-angled tip designed for perio procedures and reaching difficult-to-reach spots near the posterior teeth.

**ART Piezo Tip: BS-PFU:** is a straight-angled universal tip specially designed for perio procedures and reaching difficult-to-reach spots near the posterior teeth.



ART Piezo Tip: BS-H2R  
Right-angled curved perio tip



ART Piezo Tip: BS-10Z  
Flat headed perio tip

**ART Piezo Tip: BS-10P:** is a 45° angled tip recommended for procedures on shallow pockets.

**ART Piezo Tip: BS-H2L:** is a left-angled curved tip specially designed for perio functions and reaching difficult-to-reach spots near the posterior teeth. This tip is recommended for treatment on premolars, molars, furcation, and narrow spaces.

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**ART Piezo Tip: BS-H2R:** is a right-angled curved tip specially designed for perio functions and reaching difficult-to-reach spots near the posterior teeth. This tip is recommended for treatment on premolars, molars, furcation, and narrow spaces.

**ART Piezo Tip: BS-10Z:** is a wide straight-angled tip specially designed for perio functions and reaching difficult-to-reach spots near the posterior teeth. This tip is recommended for work on medium pockets less than 4mm.

**BE tips** are designed to be **EMS** compatible like **ART-P3, ART-SP2**

**BS tips** are designed to be compatible with **Satelec Acteon** like **ART-P3II, ART-P4, ART-P6**

<b>BE-1 Universal Pointed Piezo Tip for EMS type, each</b>	56-11230
<b>BS-1 Universal Pointed Piezo Tip for Satelec Acteon type, each</b>	56-11231
<b>BE-2 Flat chisel Piezo Tip for EMS type, each</b>	56-11223
<b>BS-2 Flat chisel Piezo Tip for Satelec Acteon type, each</b>	56-11229
<b>BE-3 Rounded Flat Piezo Tip for EMS type, each</b>	56-11228
<b>BS-3 Rounded Flat Piezo Tip for Satelec Acteon type, each</b>	56-11227
<b>BS-PFR Right-Angled perio Tip for Satelec Acteon type, each</b>	56-11226
<b>BS-P Slim Perio Tip for Satelec Acteon type, each</b>	56-11225
<b>BS-PFL Left-Angled Perio Tip for Satelec Acteon type, each</b>	56-11224
<b>BS-PFU Straight-Angled Perio Tip for Satelec Acteon type, each</b>	56-11222
<b>BS-10P 45° Angled Perio Tip for Satelec Acteon type, each</b>	56-11221
<b>BS-H2L Left-Angled Curved Perio Tip for Satelec Acteon type, each</b>	56-11220
<b>BS-H2R Right-Angled Curved Perio Tip for Satelec Acteon type, each</b>	56-11219
<b>BS-10Z Flat Headed Perio Tip for Satelec Acteon type, each</b>	56-11218

## Piezo Tip Wrenches (Vista)

**Torque - EMS\* Style**



**Universal - Acteon\* EMS\* - Amdent\* Styles**



**Torque - Acteon\* Style**

### Tip Wrenches - 3 Styles

- Torque - Acteon\* Style (Frosted Beige)
- Torque - EMS\* Style (Blue)
- Universal - Acteon\* EMS\*, Amdent\* Style

**Universal - Acteon, EMS, Amdent** 56-11216

**Torque - Acteon (Frosted Beige)** 56-11215

**Torque - EMS (Blue)** 56-11217

## Sonic Air Scalers (MTI)



### ScaleMaster

- The perfect choice for air-driven scalers.
- Gentle, yet powerful to remove stains and calculus.
- Attach to any 2-or 4-hole connections.

### Black color With:

4 Holes & Three differently shaped scaling tips and wrench 56-89811

### Blue color With:

4 Holes & Three differently shaped scaling tips and wrench 56-89812

### Black color With:

2 Holes & Three differently shaped scaling tips and wrench 56-89815

### Blue color With:

2 Holes & Three differently shaped scaling tips and wrench 56-89816




Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Sonic Scalers Tips (MTI)



<b>C-Perio Tip</b>	56-89821
<b>B-Sickle Tip</b>	56-89822
<b>A-Universal Tip</b>	56-89823
<b>Sonic Scaler Tip Wrench</b>	56-89824

### Star Titan Scalers Type Tips (MTI)



<b>Perio Scaling Tips</b>	56-89831
<b>Sickle Scaling Tips</b>	56-89832
<b>Universal Scaling Tips</b>	56-89833
<b>Star Type Tip Wrench</b>	56-89834

### Foot Control Assembly Cavitron Select SPS BOOST Type (Chapman-Huffman)



- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot
- Same light foot pressure to operate as original equipment switch
- Round pin male plug as illustrated with 8 foot cable
- Two position switch
- Fits DENTSPLY CAVITRON Select SPS BOOST
- DENTSPLY cross ref. # 565162003

<b>Fits Dentsply Cavitron Select SPS Boost</b>	56-88091
--	----------

## Ultrasonic Scalers

### Perio Flow Ultrasonic Adapter (Vista)



- Transforms Piezo ultrasonic handpiece into a versatile subgingival pocket irrigator!
- Attaches any style of luer-lock irrigation tip to your piezo ultrasonic scaler
- Autoclavable
- Available in three styles to fit all Vista Piezo Ultrasonic units, Satelec and Tulsa P5 units, EMS and Amdent™ Piezo ultrasonic systems
- Great for subgingival irrigation after scaling

<b>Perio-Flow™ - P5 / Satelec Style Adapter</b>	56-11138
---	----------

<b>Perio-Flow™ - Piezo Advantage™ / Amdent Adapter</b>	56-11137
--	----------

<b>Perio-Flow™ - EMS Style Adapter</b>	56-11139
--	----------

### O-Ring for Cavitron Insert Replacement (Plasdent)



<b>Black For Cavitron Insert, 12/Pkg</b>	56-11132
--	----------

<b>Green For Cavitron Insert, 12/Pkg</b>	56-11131
--	----------

### O-Ring Ultrasonic Insert (Hu-Friedy)

The O-Ring acts as a water seal when the insert is placed into the open end of the handpiece.



<b>30K, 12/Pkg</b>	56-90961
--------------------	----------

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Advance Piezo Ultrasonic Scaler Unit 850LED (TPC)

Compact unit housing for carry placement anywhere in the operator

No heat build up in hand piece tip

Autoclavable  
hand piece

Bottle water  
or external  
water source

Easy and  
rapid operat-  
ing for root  
canal treat-  
ment



Achieve enhargement of canal quickly

15 Watts, 24 - 28 KHz and water pressure 30-145 PSI

EMS Piezo Tips Type

1 year manufacturer warranty

Piezo sacler, hand piece, (G1, G2, G4, P4, P5 Piezo tips), foot control, tip wrench, water bottle and endo file holder + files

**Endo files and holder** (21mm #25) 56-109802

**G1** Piezo scaler tip 56-109803

**G2** Piezo scaler tip 56-109804

**G4** Piezo scaler tip 56-109805

**P3** Piezo scaler tip 56-109806

**P4** Piezo scaler tip 56-109807

**P5** Slim Piezo scaler tip 56-109808

Piezo scaler handpiece w/ LED 56-109809

**Tip wrench** 56-109810

Replacement o-rings **5/Pkg** 56-109811

Replacement water line 56-109812

Replacement scaler foot control 56-109813

## MicroProphy II (Danville)

It removes superficial stains rapidly, even in areas a prophy cup won't reach. It cleans interproximally and around orthodontic fixtures and removes stains in darkened grooves to diagnose lesions and is ideal for pit & fissure preparation. Nozzle rotates 360° and is autoclavable.



**Micro Prophy II** Air polisher unit 56-109761

**Nozzle** .023 56-109762

**Nozzle** .026 56-109763

**Nozzle** .023 stainless steel 56-109764

**Nozzle** .026 Carbide 56-109765

**Rubber Cups** 25/Pk 56-109766

## Oral Surgery Combo Package (Parkell)

**Essential oral surgery and prophylaxis devices for your operator.**

- Sensimatic™  
700SE Electrosurge  
-TurboSENSOR™  
Ultrasonic Scaler  
- Three-Insert Perio  
Pak- your choice of  
25K or 30K inserts.



**w/ 25K inserts**-includes sensimatic 700SE electrosurge, turboSENSOR scaler, 25K, 3 insert Perio Pak 27-110261

**w/ 30K inserts....** 27-110262





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Ultrasonic Scalers

### Perio Combo Package (Parkell)

An indispensable perio-dontal kit.

- TurboSENSOR™  
Magnetostrictive  
Scaler
- Periosonic® Multi-  
Fluid Irrigator  
Model "W"

Three-Insert Perio Pak  
 (your choice of 25K or  
 30K inserts)



w/ **25K inserts**-includes 56-110281  
 turboSENSOR, Periosonic Multi-  
 Fluid Irrigator model "W", and 3  
 insert Perio Pak (25K)

w/ **30K inserts**.... 56-110282

### PowerMAX Ultrasonic scaler (TPC)



#S1025  
25K Universal Scaler Insert  
(metal handle)

#S1030  
30K Universal Scaler Insert  
(metal handle)



#S4025  
25K Universal Slim Scaler Insert  
(plastic handle)

#S4030  
30K Universal Slim Scaler Insert  
(plastic handle)



#S2025  
25K Universal Slim Scaler Insert  
(metal handle)

#S2030  
30K Universal Slim Scaler Insert  
(metal handle)



#S3025  
25K Universal Scaler Insert  
(plastic handle)

#S3030  
30K Universal Scaler Insert  
(plastic handle)

#### **Features:**

- Available in 25Khz
- Precision automatic tuning circuitry
- Linear power control for exact power setting  
 Accepts Cavitron®, P, SLI, TFI, and FSI scaler inserts

- Compact and lightweight
- Auto-ranging 110-240 VAC input
- Includes one ultrasonic scaler insert and foot control
- 1 year warranty

#### **Specifications:**

- Weight: 6 lbs
- Dimensions: 8"(W) x 8"(D) x 3"(H)
- Output power to tip: Operatory adjust approximately 10-40 watts
- Handpiece cable: 7.5'
- Power supply to unit: 110-240 VAC Input, 50/60 Hz
- Water flow and supply pressure: 18 cc to 30 cc per minute flow from insert tip
- Supply pressure to scaler is 25-60 PSIG city water or bottle system

**With 25K universal scaler insert & Foot control** 56-45872

**25K universal insert**—external water flow 56-45873

**30K universal insert**- external water flow 56-45874

**25K universal slim insert**- external water flow 56-45875

**30K universal slim insert**- external water flow 56-45876

**25K universal insert**- plastic handle- internal water flow 56-45877

**30K universal insert**- plastic handle- internal water flow 56-45878

**25K universal slim insert**- plastic handle- internal water flow 56-45879

**30K universal slim insert**- plastic handle- internal water flow 56-45880

**25K Replacement hand piece w/ Tubing** 56-45881

**Power Max replacement foot control** 56-45882

**Replacement water line 6' long with male quick disconnect (with shut off)** 56-45883



*We Work Hard  
 So You Don't Have  
 To.*

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## TurboVue Combo Package (Parkell)

**Increase visibility. Reduce strain.**

- TurboVue® Illuminated Ultrasonic Scaler
- Burnett Power-Tip V light-transmitting insert
- Periosonic® Multi-Fluid Irrigator (Model W)



<b>TurboVue Combo Package - TurboVue Illuminated Scaler, Burnett Power-Tip V 30K light-transmitting insert and Periosonic® Multi-Fluid Irrigator (Model W)</b>	56-110291
<b>TurboVue Burnett Power-Tip V (30KHz)</b>	56-110292
<b>TurboVue Straight Perio Tip V (30KHz)</b>	56-110293
<b>TurboVue Universal Tip V (30KHz)</b>	56-110294

## Ultrasonic Scaler Tips Holder



<b>No 1 each w/o scaler tips</b>	56-106271
<b>No 2 each w/o scaler tips</b>	56-106272

## Piezo Cavity Preparation Scaler Tips - SATELEC, DTE & NSK (Yamm)




<b>SB1-SAT Cavity Preparation, each</b>	56-109071
<b>SB2-SAT Cavity Preparation, each</b>	56-109072
<b>SB3-SAT Cavity Preparation, each</b>	56-109073
<b>SBL-SAT Cavity Preparation, each</b>	56-109074

**SBR-SAT Cavity Preparation, each** 56-109075

## Piezo Crown Removal Scaler Tips - SATELEC, DTE & NSK (Yamm)

Compatible with

- VARIOS Piezo scaler
- HENRY SCHEIN Piezo scaler ACCLEAN PRO
- GNATUS Piezo scaler
- HU-FRIEDY Piezo scaler
- DENTAMERICA Piezo scaler
- BONART Piezo scaler
- ROLLENCE Piezo scaler



<b>C1 Crown Removal Tip, each</b>	56-109081
<b>C2 Crown Removal Tip, each</b>	56-109082

## Piezo Endo Scaler Tips - SATELEC, DTE & NSK Type (Yamm)



Compatible with

- VARIOS Piezo scaler
- HENRY SCHEIN Piezo scaler ACCLEAN PRO
- GNATUS Piezo scaler
- HU-FRIEDY Piezo scaler
- DENTAMERICA Piezo scaler
- BONART Piezo scaler
- ROLLENCE Piezo scaler

<b>E3 Endo Tip, each</b>	56-109061
<b>E4 Endo Tip, each</b>	56-109062
<b>E5 Endo Tip, each</b>	56-109063
<b>E4D Endo Tip, each</b>	56-109064
<b>E5D Endo Tip, each</b>	56-109065

*[www.prestigedental.com](http://www.prestigedental.com)*

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Piezo Endodontic Tips - Satrlec Type



**Fit to:** DTE "WOODPECKER" / BONART ART-P3II, ART-P4, ART-P6 / GNATUS / **SATELEC** / ROLLENCE-S / HUFRIEDY / NSK VARIOS / DENTAMERICA / HENRY SCHEIN Acclean Pro

**ED3DT with diamond & with Titanium Plated** to remove the calcification and bad filling material of pulp cavity, overhang pulp chamber.

**ED4DT with diamond & with Titanium Plated** to enlargement of root canal wall. (Tip end with diamond of 23mm length)

**ED5DT with diamond & with Titanium Plated** to enlargement of root canal wall. (Tip end with diamond of 13mm length)

**ED3DT with Diamond & Titanium Plated, each** 56-106311

**ED4DT w/ Diamond end for 23mm length & Titanium Plated, each** 56-106312

**ED5DT with Diamond end 13mm length & Titanium Plated, each** 56-106313

## Piezo Handpiece (Bonart)



Piezo handpiece w/o tip

56-11184


## Piezo Perio Scaler Tips - SATELEC, DTE & NSK Type (Yamm)



**P4 Perio Tip, each** 56-109051

**P4D Perio Tip, each** 56-109052

## Piezo Periodontic Tips - Satelec type



**Fit to:** DTE "WOODPECKER" / BONART ART-P3II, ART-P4, ART-P6 / GNATUS / **SATELEC** / ROLLENCE-S / HUFRIEDY / NSK VARIOS / DENTAMERICA / HENRY SCHEIN Acclean Pro

**PD1** Pointed to remove subgingival calculus

**PD3** to remove subgingival and deep perio calculus.

**PD3D with diamond** for perio treatment to Level off the surface of endo during the periodontal flap surgery.

**PD3DT with diamond & with Titanium Plated** for perio treatment to Level off the surface of endo during the periodontal flap surgery

**PD1, each** 56-106301

**PD3, each** 56-106302

**PD3D with Diamond, each** 56-106303

**PD3DT with Diamond & Titanium Plated, each** 56-106304


**Torque Wrench** 56-106305



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Ultrasonic Scalers

### Piezo Scaler Tips



**Fit to:** DTE "WOODPECKER" / BONART ART-P3II, ART-P4, ART-P6 / GNATUS / **SATELEC** / ROLLENCE-S / HU-FRIEDY / NSK VARIOS / DENTAMERICA / HENRY SCHEIN Acclean Pro

**GD1** to remove calculus and bacterial plasue on supragingival, neck of teeth and hard calculus on the teeth point

**GD3** to remove the calculus and bacterial plasue on subgingival interdental calculus.

**GD4** to remove all supragingival calculus and bacterial plasue.

**GD1**, each 56-106291

**GD3**, each 56-106292

**GD4**, each 56-106293

### Piezo Scaling Scaler Tips - SATELEC, DTE & NSK Type (Yamm)



**G2 Scaling Tip**, each 56-109041

**G6 Scaling Tip**, each 56-109042

### Piezo Tip Wrench



**Metal wrench**, each 56-106171

### Piezo Tips Wrench - Satelec Type



**Torque wrench**, each 56-110591

### Piezo Ultrasonic Gracey Mini Curette Tips (Vista)



The long shank and slender design allow these tips to be used in places that are out of reach with regular sized tips. For instance, when furcations are present or when teeth overlap (crowding) because the bite was never corrected with braces.

- Longer shank than regular curettes for

better access into deeper pockets.



# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

- Shorter blades (on average 1.2 mm) permit better scaling in NARROW pockets, reduce tissue damage and offer great adaptation to the root surface.
- Slender blades aid in tactile sensitivity and enhance deep pocket access and adaptation to the root surfaces.
- Tips will fit: EMS\* STYLE: Mini Piezon\*, Scalex 880\*, Bonart\*, SuperSon EIE\*, Master 400\*, Parkel (Piezo)\*, Dent America\* - Acteon\* STYLE: Vista P5\*, NSK\*, ENAC\*, Sybron Endo\*, Obtura Spartan\*, ADEC\*, J.Morita\*, Dentsply/Tulsa\*, Amdent\*

<b>HC-Mini-1</b> (Acteon Style)	56-11243
<b>HC-Mini-2</b> (Acteon Style)	56-11242
<b>HC-Mini-3</b> (Acteon Style)	56-11241
<b>HC-Mini-4</b> (Acteon Style)	56-11240
<b>HC-Mini-7</b> (Acteon Style)	56-11239
<b>HC-Mini-8</b> (Acteon Style)	56-11238
<b>HC-Mini-1</b> (EMS Style)	56-11237
<b>HC-Mini-2</b> (EMS Style)	56-11236
<b>HC-Mini-3</b> (EMS Style)	56-11232
<b>HC-Mini-4</b> (EMS Style)	56-11233
<b>HC-Mini-7</b> (EMS Style)	56-11234
<b>HC-Mini-8</b> (EMS Style)	56-11235



Order Online at  
PrestigeDental.com





**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

