

Vacuum Forming



Prestige Dental Products, Inc


Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Vacuum Forming

Vacuum Forming Machines

Air Vac XQ Vacuum Former (Keystone)



The Keystone Air Vac XQ Dental Vacuum Forming Machine is a work horse Dental vacuum forming machine. It's price, features and reliability make it one of the best dental vacuum formers on the market. The Air Vac XQ has greater accuracy and definition when processing bleaching trays, sports mouth guards, implant stints, temporary splints, custom impression trays, denture bases, copings, temporary crowns and bridges and orthodontic appliances.

- The first and only patented 'Xtra Quiet', American made vacuum forming machine in the industry
- Utilizes compressed air to create all types of dental appliances in less time, with precise adaptation and superior accuracy
- **5 year warranty**
- UL, CUL, and CE rated

Air Vac XQ Vacuum Former 110V 57-11111

Air Vac XQ Vacuum Former 220V 57-11112

Sta-Vac Premium Electric Vacuum Forming Machine (Buffalo)

Features/Benefits:

Versatile New ConvertTable Combination Vacuum Table and Bead Basket unique new table design can be used in quick and easy flat table mode or flipped over and used in the versatile bead basket mode for infinite model positioning and improved vacuum flow

High Performance Electric Vacuum Forming System Dual-Stage electric vacuum motor provides powerful suction and highly accurate adaptations— no compressed air connections required.

Accurate Frame and Alignment System spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table

Heavy Duty Construction Solid Construction (7"W x 16" D x 11" H, 32 lbs) for maximum durability and long life

Convenient and Safe Design Unique housing design prevents motor contamination, shields heater, and promotes ergonomic operation



Sta-Vac Vacuum Former, 120V
AC, Low Profile, High Powered 57-11469

Sta-Vac Vacuum Former, 220V
AC, Low Profile, High Powered (CE) 57-11470



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Sta-Vac II Premium Electric Vacuum Forming Machine (Buffalo)

Features/Benefits:

Versatile New Convert-Table Combination Vacuum Table and Bead Basket unique new table design can be used in quick and easy flat table mode or flipped over and used in the versatile bead basket mode for infinite model positioning and improved vacuum flow

High Performance Electric Vacuum Forming System Dual-Stage electric vacuum motor provides

powerful suction and highly accurate adaptations— no compressed air connections required.

Accurate Frame and Alignment System spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table

Heavy Duty Construction Solid Construction (7"W x 16" D x 11" H, 32 lbs) for maximum durability and long life

Convenient and Safe Design Unique housing design prevents motor contamination, shields heater, and promotes ergonomic operation

Sta-Vac II Vacuum Former, 120V
AC, w/2 ConvertTable Tables 57-11467

Sta-Vac II Vacuum Former, 220V
AC, w/2 ConvertTable Tables 57-11466

Replacement ConvertTable Table for Sta-Vac II Vacuum Former 57-11465

Replacement Steel Beads for Sta-Vac II (and Accu-Vac), 3 lb. Pkg. 57-11468



Econo-Vac Vacuum Forming Machine Heavy Duty (Buffalo)

Dual Stage vacuum motor
Square center post
Self-leveling frame
Heavy-duty construction
Frame (7 1/2" wide x 10 1/2" deep x 15" high, 32 lbs)



Econo-Vac Vacuum Former,
120V AC, High Powered 57-11332

Econo-Vac Vacuum Former,
220V AC, High Powered (CE) 57-11333



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

The Machine III Vacuum Former (Keystone)

This is the most popular Dental Vacuum Forming machine in the dental industry. Its price, features and reliability make it the first choice in dental vacuum formers. More power means greater accuracy and definition when processing bleaching trays, sports mouth guards, implant stints, temporary splints, custom impression trays, denture bases, copings, temporary crowns and bridges and orthodontic appliances.

- Uses 5" x 5" resin sheets
- 22 lbs



The Machine III 110V 57-11471

The Machine III
 220V CE Euro Cord 57-11472

Tray-Vac Vacuum Forming Machine (Buffalo)



Excellent Performance Heavy duty electric vacuum motor provides powerful and immediate suction, creating accurate adaptations, no compressed air connections required.

Gasket-Free Easy Frame System spring loaded, self-leveling, gasket-free frame with clever quick-release "fingernail-saving" latch securely handles a variety of 5" x 5" sheet materials and maintains excellent alignment with the flat vacuum table

Solid Construction Solid construction (7 1/2" W x 9 1/2" D x 15" H, 17lbs) for improved durability, long life and virtually no maintenance in a lighter-weight design

Convenient and Safe Design heater hood rotates behind the center post for safety and convenience

Wide Range of Available Accessories Please see the complete line of Buffalo Dental vacuum forming equipment and supplies including a wide range of Vacuum Forming Sheets, Block-out Compound, Silicone Spray, Tray Handles, Knives and Burs, the ThermoKnife Thermal Cutting Knife, and PolyShoK Sports Mouth guard Materials.

Tray-Vac Vacuum Former, 120V
 AC, Light Weight, Small, Powerful 57-11479

Tray-Vac Vacuum Former, 220V
 AC, Light Weight, Small Powerful (CE) 57-11480

Vacuum Forming

Vacuum Forming Materials

Base Plate Material (Keystone)

- Three thicknesses
 (.060, .080 & .100)
- Pure, high impact styrene 5" x 5" sheets
 - Low melt temperature used for denture base or wax bite plates
 - Packages of 100 or 625 sheets



Also available in round

.060 Base Plate 5" x 5" 100/Pkg	57-11117
.080 Base Plate 5" x 5" 100/Pkg	57-11116
.100 Base Plate 5" x 5" 100/Pkg	57-11115
.060 Base Plate 5" x 5" 625/Pkg	57-11114
.080 Base Plate 5" x 5" 625/Pkg	57-11113
.100 Base Plate 5" x 5" 625/Pkg	57-11118

Coping Material (Keystone)

- Excellent for temporary bridges
- Get good Definition
 - 5" x 5" sheets
 - Three thicknesses
 (.020, .030 & .040)



Also available in round

.020 Coping 5" x 5" 50/Pkg	57-11123
.030 Coping 5" x 5" 50/Pkg	57-11122
.040 Coping 5" x 5" 50/Pkg	57-11121
.020 Coping 5" x 5" 625/Pkg	57-11120
.030 Coping 5" x 5" 625/Pkg	57-11119
.040 Coping 5" x 5" 625/Pkg	57-11127

Crown & Bridge Material (Keystone)

- Thin, clear rigid material
- 5" x 5" sheets
- Acrylic does not bond to material
- For temporary bridges
- Packages of 50 or 625 sheets



Also available in round

.020 Crown & Bridge 5" x 5" 50/Pkg	57-11128
.020 Crown & Bridge 5" x 5" 625/Pkg	57-11129

Vacuum Forming



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

Duraclear Sheets Round (Dentsply)

Best material for combination of strength, durability, and clarity. Unsurpassed clarity when formed. Can be cleaned with almost any method without staining or wear. Extremely durable, will not crack, bend out of shape, or break. Will bond to acrylic, making it easy to modify for other appliance uses. Specific uses include invisible retainers, minor tooth movement, temporary bridges, aligners, and bite planes.



.030" 125mm 20 sheets/box	57-11139
.030" 120mm 20 sheets/box	57-11138
.040" 125mm 20 sheets/box	57-11137
.040" 120mm 20 sheets/box	57-11136
.060" 125mm 10 sheets/box	57-11135
.060" 120mm 10 sheets/box	57-11134
.080" 125mm 10 sheets/box	57-11133
.080" 120mm 10 sheets/box	57-11132

Duraclear Sheets Square (Dentsply)

Best material for combination of strength, durability, and clarity. Unsurpassed clarity when formed. Can be cleaned with almost any method without staining or wear. Extremely durable, will not crack, bend out of shape, or break. Will bond to acrylic, making it easy to modify for other appliance uses. Specific uses include invisible retainers, minor tooth movement, temporary bridges, aligners, and bite planes.



.030" 5" 20 sheets/box	57-11147
.030" 125mm 20 sheets/box	57-11146
.040" 5" 20 sheets/box	57-11145
.040" 125mm 20 sheets/box	57-11144
.060" 5" 10 sheets/box	57-11143
.060" 125mm 10 sheets/box	57-11142
.080" 5" 10 sheets/box	57-11141
.080" 125mm 10 sheets/box	57-11140

Duraforce Sheets (Dentsply)

Superior conformity to the model when heated, little or no shrink. Cooling spray is not required.



Great long term wear. Strong and durable, will not crack or break. Strong enough to create divots and initiate tooth movement if needed. Can be easily cleaned without heavily staining the material. Good clarity. Not bondable to acrylic. Specific uses include final retainers, bruxism splints, temporary bridges, and minor tooth movement.

Square .030 Box of 100 sheets	57-11155
Square .040 Box of 100 sheets	57-11153
Square .060 Box of 50 sheets	57-11151
Square .030 Bag of 10 sheets	57-11149
Square .040 Bag of 10 sheets	57-11156
Square .060 Bag of 5 sheets	57-11158

Duraguard Sheets Round (Dentsply)



Specially formulated EVA heats and forms easily, with great conformity to the model. Specific uses include customer mouth guards, soft splints, and whitening trays. Material forms soft and flexible, with good durability. 25 sheets/box

Aqua .040" 5"	57-11206
Black .040" 5"	57-11223
Blue .040" 5"	57-11224
Blue Light .040" 5"	57-11177
Clear .040" 5"	57-11176
Green .040" 5"	57-11178
Orange .040" 5"	57-11179
Purple .040" 5"	57-11180
Red .040" 5"	57-11182
Silver .040" 5"	57-11181
White .040" 5"	57-11183
Yellow .040" 5"	57-11184
Multi Color .060" 5"	57-11165
Aqua .080" 5"	57-11167
Black .080" 5"	57-11168
Blue .080" 5"	57-11169
Blue Light .080" 5"	57-11170
Clear .080" 5"	57-11172
Green .080" 5"	57-11171
Orange .080" 5"	57-11173
Purple .080" 5"	57-11190
Red .080" 5"	57-11191
Silver .080" 5"	57-11192
White .080" 5"	57-11209



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Vacuum Forming

Yellow .080" 5"	57-11211
Multi Color .080" 5"	57-11210
Aqua .120" 5"	57-11212
Black .120" 5"	57-11213
Blue .120" 5"	57-11214
Blue Light .120" 5"	57-11215
Clear .120" 5"	57-11216
Green .120" 5"	57-11217
Orange .120" 5"	57-11218
Purple .120" 5"	57-11219
Red .120" 5"	57-11220
Silver .120" 5"	57-11221
White .120" 5"	57-11222
Yellow .120" 5"	57-11208
Multi Color .120" 5"	57-11207
Aqua .160" 5"	57-11193
Black .160" 5"	57-11194
Blue .160" 5"	57-11195
Blue Light .160" 5"	57-11196
Clear .160" 5"	57-11197
Green .160" 5"	57-11199
Orange .160" 5"	57-11198
Purple .160" 5"	57-11200
Red .160" 5"	57-11201
Silver .160" 5"	57-11202
White .160" 5"	57-11204
Yellow .160" 5"	57-11203
Multi Color .160" 5"	57-11205

Youth Mouthguard (Keystone)



Also available in round

- 5" x 5", square sheets, .120 thick
- Soft EVA, for energy absorption
- Great for little kids in sports

.120 Soft EVA 5" x 5", 25/Pk	57-11494
.120 Red 5" x 5", 25/Pk	57-11493
.120 Blue 5" x 5", 25/Pk	57-11492
.120 Green 5" x 5", 25/Pk	57-11491
.120 Yellow 5" x 5", 25/Pk	57-11490
.120 Black 5" x 5", 25/Pk	57-11489
.120 White 5" x 5", 25/Pk	57-11488
.120 Soft EVA 5" x 5", 300/Pk	57-11487
.120 Soft EVA 5" x 5", 300/Pk	57-11487
.120 Red 5" x 5", 300/Pk	57-11486
.120 Blue 5" x 5", 300/Pk	57-11485
.120 Green 5" x 5", 300/Pk	57-11484
.120 Yellow 5" x 5", 300/Pk	57-11483
.120 Black 5" x 5", 300/Pk	57-11495

The Machine III Vacuum Former (Keystone)



Pro-Form Dual Laminates-E Gasket (Keystone)



To See More Products Visit
www.prestigedental.com

Vacuum Forming

Prestige Dental Products



**Excellent Qualities
at Affordable Prices**

877-772-3888

www.prestigedentalproducts.com

Pro-Form Mouthguard Resin Sheets (Keystone)

Soft, for energy
absorption

- 5" x 5" sheets, .160 thick
- Mouthguards & Niteguards
- Packages of 25 or 300 sheets



Also available in round

.160 Clear 5" x 5" 25/Pkg	57-11351
.160 Red 5" x 5" 25/Pkg	57-11349
.160 Blue 5" x 5" 25/Pkg	57-11350
.160 Green 5" x 5" 25/Pkg	57-11348
.160 Yellow 5" x 5" 25/Pkg	57-11347
.160 Black 5" x 5" 25/Pkg	57-11346
.160 Orange 5" x 5" 25/Pkg	57-11345
.160 Eva White 5" x 5" 25/Pkg	57-11344
.160 Purple 5" x 5" 25/Pkg	57-11343
.160 Clear 5" x 5" 300/Pkg	57-11342
.160 Red 5" x 5" 300/Pkg	57-11341
.160 Blue 5" x 5" 300/Pkg	57-11340
.160 Green 5" x 5" 300/Pkg	57-11339
.160 Yellow 5" x 5" 300/Pkg	57-11338
.160 Black 5" x 5" 300/Pkg	57-11337
.160 Orange 5" x 5" 300/Pkg	57-11336
.160 White 5" x 5" 300/Pkg	57-11335
.160 Purple 5" x 5" 300/Pkg	57-11334

Mouthguard Straps (Keystone)

Keystone has developed an easy-to-attach helmet strap used in the production of custom-made mouthguards. These straps have a safety release feature that allows the mouthguard to



Prestige Dental Products, Inc.

Prestigedentalproducts.com


Tel: 877-772-3888

break away under 12 pounds of pressure. Available in black and clear. Sold in packages of 12 and 25.

Black 12/Pkg	57-11354
Black 25/Pkg	57-11353
Clear 12/Pkg	57-11352
Clear 25/Pkg	57-11355

Duraguard Sheets Square (Dentsply)

Specially formulated EVA heats and forms easily, with great conformity to the model. Specific uses include customer mouthguards, soft splints, and whitening trays. Material forms soft and flexible, with good durability.



Aqua .040" 5", 25 sheets/box	57-11329
Black .040" 5", 25 sheets/box	57-11277
Blue .040" 5", 25 sheets/box	57-11279
Blue Light .040" 5", 25 sheets/box	57-11280
Clear .040" 5", 25 sheets/box	57-11282
Green .040" 5", 25 sheets/box	57-11281
Orange .040" 5", 25 sheets/box	57-11283
Purple .040" 5", 25 sheets/box	57-11284
Red .040" 5", 25 sheets/box	57-11285
Silver .040" 5", 25 sheets/box	57-11286
White .040" 5", 25 sheets/box	57-11287
Yellow .040" 5", 25 sheets/box	57-11288
Multi Color .040" 5", 25 sheets/box	57-11289
Aqua .060" 5", 25 sheets/box	57-11291
Black .060" 5", 25 sheets/box	57-11290
Blue .060" 5", 25 sheets/box	57-11278
Blue Light .060" 5", 25 sheets/box	57-11276
Clear .060" 5", 25 sheets/box	57-11328
Green .060" 5", 25 sheets/box	57-11327
Orange .060" 5", 25 sheets/box	57-11318
Purple .060" 5", 25 sheets/box	57-11266
Red .060" 5", 25 sheets/box	57-11265
Silver .060" 5", 25 sheets/box	57-11268
White .060" 5", 25 sheets/box	57-11267
Yellow .060" 5", 25 sheets/box	57-11269
Multi Color .060" 5", 25 sheets/box	57-11270
Aqua .080" 5", 25 sheets/box	57-11271
Black .080" 5", 25 sheets/box	57-11272
Blue .080" 5", 25 sheets/box	57-11273
Blue Light .080" 5", 25 sheets/box	57-11274
Clear .080" 5", 25 sheets/box	57-11275




Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Green .080" 5" , 25 sheets/box	57-11292
Orange .080" 5" , 25 sheets/box	57-11293
Purple .080" 5" , 25 sheets/box	57-11394
Red .080" 5" , 25 sheets/box	57-11311
Silver .080" 5" , 25 sheets/box	57-11312
White .080" 5" , 25 sheets/box	57-11313
Yellow .080" 5" , 25 sheets/box	57-11314
Multi Color .080" 5" , 25 sheets/box	57-11315
Aqua .120" 5" , 25 sheets/box	57-11316
Black .120" 5" , 25 sheets/box	57-11317
Blue .120" 5" , 25 sheets/box	57-11319
Blue Light .120" 5" , 25 sheets/box	57-11320
Clear .120" 5" , 25 sheets/box	57-11321
Green .120" 5" , 25 sheets/box	57-11322
Orange .120" 5" , 25 sheets/box	57-11323
Purple .120" 5" , 25 sheets/box	57-11324
Red .120" 5" , 25 sheets/box	57-11325
Silver .120" 5" , 25 sheets/box	57-11309
White .120" 5" , 25 sheets/box	57-11310
Yellow .120" 5" , 25 sheets/box	57-11295
Multi Color .120" 5" , 25 sheets/box	57-11296
Aqua .160" 5" , 25 sheets/box	57-11297
Black .160" 5" , 25 sheets/box	57-11298
Blue .160" 5" , 25 sheets/box	57-11299
Blue Light .160" 5" , 25 sheets/box	57-11300
Clear .160" 5" , 25 sheets/box	57-11301
Green .160" 5" , 25 sheets/box	57-11302
Orange .160" 5" , 25 sheets/box	57-11303
Purple .160" 5" , 25 sheets/box	57-11304
Red .160" 5" , 25 sheets/box	57-11305
Silver .160" 5" , 25 sheets/box	57-11306
White .160" 5" , 25 sheets/box	57-11307
Yellow .160" 5" , 25 sheets/box	57-11308
Multi Color .160" 5" , 25 sheets/box	57-11326

Pro-Flex & Dual Combo (Keystone)

An excellent way to make niteguards. First, vacuum form a dual laminate onto the patient's cast. This gives a hard bite plate with a built-in soft liner. Then, modify the occlusal surface by adding the semi-flexible Pro-Flex NFC acrylic resin. The final result is a hard bite plate modified to the patient's needs that is slightly flexible with a soft inner liner for patient comfort. Two great products combined into one terrific final appliance.



Dual Laminate Kit	57-11364
Clear NFC 160gm Powder w 4oz Liquid Kit	57-11365

Vacuum Forming

Pro-Form Fun Mouthguards (Keystone)

Pro-form has added some new designs to its extensive line of laminated mouthguards. They are all made using the same laminating procedure to insure the same Pro-form quality.




Also available in round

Blue Polkadot 5" x 5" , 12/Pk	57-11391
White Polkadot 5" x 5" , 12/Pk	57-11390
Red Polkadot 5" x 5" , 12/Pk	57-11389
Fang 5" x 5" , 12/Pk	57-11388
Blood 5" x 5" , 12/Pk	57-11392

Pro-Form Dual Laminates - E- Gasket (Keystone)

The Pro-form Dual Laminate is a unique combination of two materials formed into a single laminated sheet. It was made especially for use in bruxing cases. The soft liner assures patient comfort and the hard surface assures long-lasting wear. Also, acrylic will bond to the hard surface to allow the dental technician to modify the occlusal surface. Available in 3 different thicknesses (2, 3, & 5 mm). 12 laminated sheets per package.



Also available in round

2mm E-Gasket 5" x 5" , 12/Pk	57-11130
3mm Dual Laminates 5" x 5" , 12/Pk	57-11131
5mm Anti Snore Device 5" x 5" , 12/Pk	57-12132

Pro-Form Niteguards (Keystone)

This laminate or layered soft/soft material was made specifically for nighttime bruxers. The layering gives the thermoplastic a more abrasion-resistant characteristic that will withstand grinding better than traditional soft sheet materials. Available in packages of 12 laminates.



Also available in round

Niteguards 5" x 5" , 12/Pk	57-11363
-----------------------------------	----------

Vacuum Forming



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Pro-Form Tie-Dye Mouthguards (Keystone)



Also available in round

Pro-form has added some new Tie Dye designs to its extensive line of laminated mouthguards. They are all made using the same laminating procedure to insure the same Pro-form quality.

Wilderness 5" x 5", 6/Pkg	57-11360
Midnight 5" x 5", 6/Pkg	57-11359
Shamrock 5" x 5", 6/Pkg	57-11358
Flamingo 5" x 5", 6/Pkg	57-11357
Rainbow 5" x 5", 6/Pkg	57-11356
All In 5" x 5", 6/Pkg	57-11362
Patriot 5" x 5", 6/Pkg	57-11361

Pro-Form Mouthguard Laminates (Keystone)

Pro-form Mouth-guard Laminate is a strong two layer material made of high-energy absorbing EVA vinyl's. 0.150" thick in dual-color sheets. Dual-color Mouth guard Laminate in assorted two-color combos. 5"x 5", 12/Pk.



Also available in round

Blue/Green	57-11386
Black/White	57-11385
Blue/Red	57-11384
Black/Yellow	57-11382
Black/Red	57-11383
Purple/Yellow	57-11381
Maroon/Yellow	57-11380
Blue/Yellow	57-11379
Red/White	57-11378
Orange/White	57-11377

Blue/White	57-11376
Orange/Blue	57-11375
Green/Yellow	57-11374
Black/Orange	57-11367
Purple/White	57-11368
Green/White	57-11369
Red/Yellow	57-11370
Maroon/White	57-11371
Red/Clear	57-11372
Black/Blue	57-11373
Orange/Clear	57-11366
Assorted 2 Color	57-11387
Clear 5" x 5", 12/Pk	57-11403
Orange 5" x 5", 12/Pk	57-11402
Yellow 5" x 5", 12/Pk	57-11401
Maroon 5" x 5", 12/Pk	57-11400
White 5" x 5", 12/Pk	57-11399
Green 5" x 5", 12/Pk	57-11398
Black 5" x 5", 12/Pk	57-11397
Red 5" x 5", 12/Pk	57-11396
Blue 5" x 5", 12/Pk	57-11395
Purple 5" x 5", 12/Pk	57-11394
Pink 5" x 5", 12/Pk	57-11394
Assorted Solid 5" x 5", 12/Pk	57-11404
Red/White/Blue	57-12388
Canada Red/White/Red	57-12389

Pro-Form Mouthguard Tri-color Laminates (Keystone)

Pro-form Mouth-guard Laminate is a strong two layer material made of high-energy absorbing EVA vinyl's. 0.150" thick in tri-color sheets.



Also available in round

Red/White/Black 5" x 5", 12/Pk	57-11432
Red/White/Blue 5" x 5", 12/Pk	57-11431
White/Green/Black 5" x 5", 12/Pk	57-11430
Red/White/Red Canada 5" x 5", 12/Pk	57-11429
Blue/Yellow/Blue 5" x 5", 12/Pk	57-11428
Red/Black/Red 5" x 5", 12/Pk	57-11427



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Black/White/Green 5" x 5", 12/Pk	57-11426
Green/Clear/Yellow 5" x 5", 12/Pk	57-11425
Orange/White/Orange 5" x 5", 12/Pk	57-11421
Orange/Blue/Orange 5" x 5", 12/Pk	57-11422
Purple/Blue/Orange 5" x 5", 12/Pk	57-11423
Red/White/Red 5" x 5", 12/Pk	57-11424
Red/Yellow/Blue 5" x 5", 12/Pk	57-11420
Blue/White/Yellow 5" x 5", 12/Pk	57-11419
Black/Blue/Red 5" x 5", 12/Pk	57-11418
Yellow/Blue/Yellow 5" x 5", 12/Pk	57-11417
Orange/Black/Orange 5" x 5", 12/Pk	57-11416
Purple/White/Purple 5" x 5", 12/Pk	57-11415
Yellow/Red/Yellow 5" x 5", 12/Pk	57-11414
Green/Black/Green 5" x 5", 12/Pk	57-11413
Blue/White/Blue 5" x 5", 12/Pk	57-11412
Black/White/Yellow 5" x 5", 12/Pk	57-11411
Purple/Black/Purple 5" x 5", 12/Pk	57-11410
Black/Yellow/Green 5" x 5", 12/Pk	57-11409
Blue/Black/Yellow 5" x 5", 12/Pk	57-11408
Clear/Green/Clear 5" x 5", 12/Pk	57-11407
Clear/Red/Clear 5" x 5", 12/Pk	57-11406
Black/Green/Yellow 5" x 5", 12/Pk	57-11405
Assorted 3 Color 5" x 5", 12/Pk	57-11433

Retainer Material (Keystone)

- Stiff, clear material
- 5" x 5", sheets
- Acrylic sticks to material
- Invisible orthodontic retainers
- Packages of 50 or 625 sheets



Also available in round

.030 Super Clear Retainer 5'x 5", 25/Pk	57-11435
030 Super Clear Retainer 5'x 5", 50/Pk	57-11434
030 Super Clear Retainer 5'x 5", 625/Pk	57-11436

Vacuum Forming

Soft EVA Sheets (Keystone)

- Soft, clear material
- 5" x 5" sheets
- Easily formed and trimmed
- Bleaching trays, soft brushing
- Custom fluoride trays, model duplication
- Packages of 25 or 300 sheets



Also available in round

.040 Soft EVA 5" x 5", 25/Pk	57-11443
.060 Soft EVA 5" x 5", 25/Pk	57-11442
.080 Soft EVA 5" x 5", 25/Pk	57-11441
.120 Soft EVA 5" x 5", 25/Pk	57-11440
.040 Soft EVA 5" x 5", 300/Pk	57-11439
.060 Soft EVA 5" x 5", 300/Pk	57-11438
.080 Soft EVA 5" x 5", 300/Pk	57-11437
.120 Soft EVA 5" x 5", 300/Pk	57-11444

Splint Material (Keystone)

- Stiff, clear material
- 5" x 5", sheets
- Acrylic sticks to material
- Packages of 25, 50 or 625 sheets
- For hard splints, thin splits, surgical trays, orthodontic retainers and bleaching trays



Also available in round


1 .020 Splint 5" x 5", 50/Pk	57-11464
3 .040 Splint 5" x 5", 50/Pk	57-11463
4 .060 Splint 5" x 5", 50/Pk	57-11461
5 .080 Splint 5" x 5", 50/Pk	57-11460
6 .100 Splint 5" x 5", 50/Pk	57-11459

Vacuum Forming



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Tuff Form Sheets (Dentsply)



Patented arch form indentation. Provides additional material surface area, which eliminates thinning along the incisal edge, and provides a retainer with uniform thickness throughout the appliance. Will bond to acrylic, making it useful for creating modified appliances. Superior conformity to the model when heated. Material has little or no shrink. Material forms to the contours of the model, and cooling spray is not required. A very secure fit to the model is assured. Great clarity when formed, creating a true invisible retainer for patients. Will not get brittle and crack. Material forms thin, making it comfortable for the patient, and increasing patient compliance. Specific uses include invisible retainers, minor tooth movement, temporary bridges, aligners, and bite planes.

Clear Square .040", Box of 50 sheets 57-11481

Clear Round .040", Box of 50 sheets 57-11482

Tray Material (Keystone)



Also available in round

- Two thickness
- Pure, high impact styrene
- 5" x 5", sheets
- Custom impression trays
- Packages of 100 or 625 sheets
- Spacer foam included

.125 Tray Material 5" x 5", 100/Pk 57-11477

.150 Tray Material 5" x 5", 100/Pk 57-11476

Specialty Breast Cancer Pink Mouthguard (Keystone)



Specialty breast cancer pink mouthguard

Breast Cancer Pink Mouthguard

5" x 5", 12/Pk

57-11446

Spacer Foam (Keystone)



- 5" x 5" square sheets
- Use as a spacer material
- Makes cast easier to remove
- Packages of 25
- One size works for all materials
- Reusable

Foam Spacer 5" x 5", 25/Pk

57-11445

Thermoform Sheet Intro Kit (Dentsply-Glenroe)



Special introductory kit featuring 10 sheets each of Duraclear, Duraforce, Tuff Form, and Duralight. 40 sheets/box

Square

57-11474

Round

57-11473



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

